

Cleanup of Hazardous Waste Site Samples by Solid Phase Extraction

ENVI-Florisil solid phase extraction tubes meet the requirements for pesticide recovery from water, sediment, and soil samples taken from hazardous waste sites, as described in US EPA Contract Laboratory Program methodology. Mean recovery of nine pesticides was 82% or better, well within the specifications of the CLP method. A procedure for recovering the pesticides is described.

Key Words:

- pesticides • solid phase extraction • CLP methods

The US Environmental Protection Agency (EPA) Contract Laboratory Program (CLP) method for monitoring pesticides in water, sediment, and soil samples from hazardous waste sites specifies use of Florisil® solid phase extraction cartridges with stainless steel or Teflon® frits for sample cleanup (1). This mandatory cleanup procedure significantly reduces matrix interferences caused by polar compounds.

ENVI™-Florisil solid phase extraction tubes fulfill the pesticide recovery requirements of the EPA CLP pesticide method. These 6mL polypropylene tubes with stainless steel frits contain 0.5g or 1.0g of Florisil packing. The bed size used depends on the final extract volume required for injection. Manual injection requires a 1.0mL final extract volume and a 0.5g tube. If your autosampler operates reliably with 1.0mL of final extract, use a 0.5g tube, otherwise, prepare 2.0mL of sample extract using a 1g tube.

Table 1. Pesticide Recovery After Extraction on ENVI-Florisil SPE Tubes

Compound	Concentration (ng/mL)	Recovery (%)	
		1.0g Bed	0.5g Bed
α-BHC	25	81.8	90.0
γ-BHC	25	87.3	92.4
Heptachlor	25	90.6	89.2
Endosulfan I	25	91.8	101.6
Dieldrin	50	94.9	92.4
Endrin	50	108.2	93.7
4,4'-DDD	50	104.9	96.1
4,4'-DDT	50	95.6	102.0
Methoxychlor	250	110.0	91.7

To demonstrate the suitability of ENVI-Florisil tubes for the CLP pesticide method, we measured recoveries of the nine pesticides specified in the CLP method, for use as a cartridge performance check (Table 1). Samples of a standard containing the nine

pesticides were prepared and concentrated as specified in the CLP method. The ENVI-Florisil tubes were conditioned with 5mL of 90:10 (v/v) hexane:acetone, then 0.5mL of the pesticide mix was applied. The sample was eluted with 9.0mL of the hexane:acetone solution at a slow, dropwise rate, followed by 2mL of hexane. The eluant was concentrated to 1.0mL. GC analysis was performed on an SPB™-608 capillary column (30m x 0.53mm ID 1.0µm film).

Table 1 summarizes recovery data for tubes containing 0.5g and 1.0g beds of a typical lot of Florisil packing at mid-point analyte concentration levels. Mean recoveries were also determined for five lots of 0.5g tubes using the same procedure. Both studies confirmed that ENVI-Florisil tubes satisfy the CLP pesticide method requirement that target compounds be recovered at 80 to 120%.

ENVI-Florisil tubes are specifically designed and tested to ensure reduced matrix interferences and high analyte recovery for the EPA CLP pesticide method. We also offer a complete selection of products — including solid phase extraction tubes, columns, and standards — for other EPA methods. See our latest catalog for more information.

Ordering Information:

Description	Cat. No.
ENVI-Florisil tubes, 6mL, stainless steel frits	
1.0g pk. of 30	57053
0.5g pk. of 30	57046

Pesticide Individual Standard Mix A

5µg/mL each: α-BHC, γ-BHC, heptachlor, endosulfan I;
10µg/mL each: dieldrin, endrin, 4,4'-DDD, 4,4'-DDT;
50µg/mL methoxychlor. In 1mL hexane.

48796

SPB-608 Fused Silica Capillary Column

30m x 0.53mm ID, 1.0µm film

made-to-order

Reference

1. US EPA Contract Laboratory Program Statement of Work for Organic Analysis, Multimedia Multiconcentration, Doc. No. OLM01.0.

Trademarks

ChromFax, ENVI, SPB — Sigma-Aldrich Co.

Florisil — U.S. Silica Co.

Teflon — E. I. du Pont de Nemours & Co., Inc.

Fused silica columns manufactured under HP US Pat. No. 4,293,415.

Free Technical Literature From Supelco

Supelco produces a wide variety of technical literature, all of which is free. You may request this literature several different ways:

- Use the business reply card located at the back of our catalog.
- Call our Ordering and Customer Service department (800-247-6628 or 814-359-3441).
- Use our ChromFax™ service.

How to Use ChromFax

ChromFax is Supelco's automated digital data system, which allows you to instantly obtain technical literature any time of the day or night. Supelco technical bulletins, application notes, lab hints, and product specification sheets can be sent to your Fax machine.

1. To use ChromFax you must dial from a touch-tone phone. If you do not have a touch-tone phone, call our Ordering and Customer Service department and we will fax the publication to you.
2. Have the phone number of your Fax machine ready, and dial **800-652-8069 or 814-359-5748**. Remember you must have the *complete* Fax number (011 + country/city code + Fax number) when you are calling from outside the US or Canada, otherwise the document cannot be transmitted.
3. When ChromFax answers, follow the simple recorded instructions. Select the **ChromFax number** of the document you want. Be certain you order documents by their ChromFax numbers and NOT their literature numbers.
4. Enter your Fax number and hang up. The requested information will be sent to you within minutes. To report a problem with the ChromFax system contact our Technical Service department.

Technical Literature Available on ChromFax

The listing of documents available on ChromFax is updated every month. We encourage you to request a list of available ChromFax documents each time you use the ChromFax system.

Contact our Technical Service Department
(phone 800-359-3041 or 814-359-3041,
FAX 800-359-3044 or 814-359-5468)
for expert answers to your questions.

Note 28

For more information, or current prices, contact your nearest Supelco subsidiary listed below. To obtain further contact information, visit our website (www.sigma-aldrich.com), see the Supelco catalog, or contact Supelco, Bellefonte, PA 16823-0048 USA.

ARGENTINA · Sigma-Aldrich de Argentina, S.A. · Buenos Aires 1119 **AUSTRALIA** · Sigma-Aldrich Pty. Ltd. · Castle Hill NSW 2154 **AUSTRIA** · Sigma-Aldrich Handels GmbH · A-1110 Wien
BELGIUM · Sigma-Aldrich N.V./S.A. · B-2880 Bornem **BRAZIL** · Sigma-Aldrich Quimica Brasil Ltda. · 01239-010 São Paulo, SP **CANADA** · Sigma-Aldrich Canada, Ltd. · 2149 Winston Park Dr., Oakville, ON L6H 6J8
CZECH REPUBLIC · Sigma-Aldrich s.r.o. · 186 00 Praha 8 **DENMARK** · Sigma-Aldrich Denmark A/S · DK-2665 Vallensbaek Strand **FINLAND** · Sigma-Aldrich Finland/YA-Kemia Oy · FIN-00700 Helsinki
FRANCE · Sigma-Aldrich Chimie · 38297 Saint-Quentin-Fallavier Cedex **GERMANY** · Sigma-Aldrich Chemie GmbH · D-82041 Deisenhofen **GREECE** · Sigma-Aldrich (o.m.) Ltd. · Ilioupoli 16346, Athens
HUNGARY · Sigma-Aldrich Kft. · H-1067 Budapest **INDIA** · Sigma-Aldrich Co. · Bangalore 560 048 **IRELAND** · Sigma-Aldrich Ireland Ltd. · Dublin 24 **ISRAEL** · Sigma Israel Chemicals Ltd. · Rehovot 76100
ITALY · Sigma-Aldrich s.r.l. · 20151 Milano **JAPAN** · Sigma-Aldrich Japan K.K. · Chuo-ku, Tokyo 103 **KOREA** · Sigma-Aldrich Korea · Seoul **MALAYSIA** · Sigma-Aldrich (M) Sdn. Bhd. · Selangor
MEXICO · Sigma-Aldrich Quimica S.A. de C.V. · 50200 Toluca **NETHERLANDS** · Sigma-Aldrich Chemie BV · 3330 AA Zwijndrecht **NORWAY** · Sigma-Aldrich Norway · Torshov · N-0401 Oslo
POLAND · Sigma-Aldrich Sp. z o.o. · 61-663 Poznań **PORTUGAL** · Sigma-Aldrich Quimica, S.A. · Sintra 2710 **RUSSIA** · Sigma-Aldrich Russia · Moscow 103062 **SINGAPORE** · Sigma-Aldrich Pte. Ltd.
SOUTH AFRICA · Sigma-Aldrich (pty) Ltd. · Jet Park 1459 **SPAIN** · Sigma-Aldrich Quimica, S.A. · 28100 Alcobendas, Madrid **SWEDEN** · Sigma-Aldrich Sweden AB · 135 70 Stockholm
SWITZERLAND · Supelco · CH-9471 Buchs **UNITED KINGDOM** · Sigma-Aldrich Company Ltd. · Poole, Dorset BH12 4QH
UNITED STATES · Supelco · Supelco Park · Bellefonte, PA 16823-0048 · Phone 800-247-6628 or 814-359-3441 · Fax 800-447-3044 or 814-359-3044 · email: supelco@sial.com

H

Supelco is a member of the Sigma-Aldrich family. Supelco products are sold through Sigma-Aldrich, Inc. Sigma-Aldrich warrants that its products conform to the information contained in this and other Sigma-Aldrich publications. Purchaser must determine the suitability of the product for a particular use. Additional terms and conditions may apply. Please see the reverse side of the invoice or packing slip.