

Determination of Bis(2-chloroisopropyl) ether Using Chiral Capillary GC

Commercially available bis(2-chloroisopropyl) ether is only about 50% pure. The existing commercial sources of this material contain three optical isomers and an achiral contaminant, 1,1'-oxybis[3-chloropropane]. Under US EPA method conditions, the contaminant and the isomers elute closely and are typically not separated. Detailed chromatography, using an α -DEX 120 chiral capillary GC column, reveals the contamination now known to be present in commercially available bis(2-chloroisopropyl) ether.

Key Words:

- bis(2-chloroisopropyl) ether
- chiral capillary GC
- optical isomers

Supelco regularly receives comments about chemical standards, ranging from suggestions for new formulations to questions about the purity of material used in manufacturing. Over an extended period, several questions about the purity of the bis(2-chloroisopropyl) ether (BCIPE) used in our manufacturing process have been received. BCIPE, a byproduct in the manufacture of bis(chloromethyl) ether, is a toxic compound that is included as a constituent in a number of US EPA methods. Analysis is typically by capillary GC using an SPB™-5 (5% diphenyl/95% dimethylpolysiloxane) column.

Concerns about the purity of this raw material were not resolved through routine techniques, such as GC/MS or GC/IR. Available NMR data on this material indicate the probable existence of contaminant, but the contaminant is not verifiable by comparative technique. An additional technique, differential scanning calorimetry, failed to indicate the presence of contaminants in BCIPE.

Structurally, BCIPE contains two stereocenters and is optically active (Figure A). In molecules containing more than one stereocenter, the total number of stereoisomers will not exceed 2^n , where n is the number of stereocenters. Thus, four enantiomeric forms of BCIPE are theoretically expected, but close examination reveals the existence of a plane of symmetry within the molecule. This symmetry means that one optical isomer pair exists in the meso form and exhibits achiral behavior. Although separation of optical isomers is not expected under normal GC conditions, an α -DEX™ 120 chiral capillary GC column will resolve the BCIPE optical isomers (Figure B).

The α -DEX 120 capillary column contains 20% (w/w) permethylated α -cyclodextrin in a polysiloxane co-phase — offering a unique selectivity capable of resolving non-meso

Figure A. Bis(2-chloroisopropyl) ether (BCIPE) Isomer Structures






Figure B. Total Ion Chromatogram of BCIPE Isomers on an α -DEX 120 Capillary Column

Column: α -DEX 120, 30m x 0.25mm ID, 0.25 μ m film
 Cat. No.: 24310
 Oven: 70°C to 200°C at 2°C/min
 Carrier: helium, 20cm/sec
 Det.: Mass Spec (Scan 18-500 amu)
 Inj.: Direct, 1 μ L



stereoisomeric forms. When a BCIPE sample was analyzed on this column, three peaks were observed instead of the anticipated two (Figure B). A mass spectrum of the peaks confirmed two BCIPE isomers, and a single achiral contaminant. The contaminating compound, 1,1'-oxybis[3-chloropropane] is very similar in fragmentation pattern to BCIPE — so similar that, unless the analytical search is specific, it is unlikely that this contaminant will be discerned (Figure C).

Figure C. Fragmentation Patterns of BCIPE and 1,1'-Oxybis[3-chloropropane]



One generally accepted route for synthesizing commercial BCIPE involves the condensation reaction of 1-chloro-2-propanol. Supelco learned that this starting material is only commercially available in a technical grade containing about 30% of an alcohol impurity. It is this impurity which forms the contaminant, 1,1'-oxybis[3-chloropropane]. This contaminant was observed in all commercially available BCIPE samples reviewed. The BCIPE studied was found to be only about 50% pure on the basis of comparative peak areas. (This value is estimated because an authentic reference BCIPE of 98% or better purity is currently unavailable.)

Supelco is presently exploring synthetic routes to prepare BCIPE free of 1,1'-oxybis[3-chloropropane]. Until this work is completed, we will determine the purity of every lot of BCIPE used in our manufacturing, and report this information on the Certificate of Analysis accompanying each product. While Supelco typically corrects for raw material purity during the gravimetric preparation of standards, it was decided not to do this for BCIPE. Because commercially available BCIPE is approximately 50% pure, and the impurities are not easily resolved, the apparent amount of BCIPE in a purity corrected solution would be double that normally expected in a prepared reference solution.

If BCIPE is a constituent in your analytical work, the results you are obtaining are probably high by about 50%. Extra consideration should thus be given the data you report about this compound. If it is critical that you obtain an accurate measurement of the BCIPE in your sample, and you are not using a standard obtained from Supelco, quantitation of your material using the Supelco α -DEX 120 capillary column is highly recommended.

Ordering Information:

α-DEX 120 Fused Silica Capillary GC Column	
30m x 0.25mm ID, 0.25 μ m film	24310
Bis(2-chloroisopropyl)ether, neat	
100mg	48498
Bis(2-chloroisopropyl)ether	
5000 μ g/mL in methanol	40040-U

DEX and SPB are trademarks of Supelco, Inc.

Contact our Technical Service Department (phone 800-359-3041 or 814-359-3041, FAX 814-359-5468) for expert answers to your questions.

For more information, or current prices, contact your nearest Supelco subsidiary listed below. To obtain further contact information, visit our website (www.sigma-aldrich.com), see the Supelco catalog, or contact Supelco, Bellefonte, PA 16823-0048 USA.

ARGENTINA · Sigma-Aldrich de Argentina, S.A. · Buenos Aires 1119 AUSTRALIA · Sigma-Aldrich Pty. Ltd. · Castle Hill NSW 2154 AUSTRIA · Sigma-Aldrich Handels GmbH · A-1110 Wien
 BELGIUM · Sigma-Aldrich N.V./S.A. · B-2880 Bornem BRAZIL · Sigma-Aldrich Quimica Brasil Ltda. · 01239-010 São Paulo, SP CANADA · Sigma-Aldrich Canada, Ltd. · 2149 Winston Park Dr., Oakville, ON L6H 6J8
 CZECH REPUBLIC · Sigma-Aldrich s.r.o. · 186 00 Praha 8 DENMARK · Sigma-Aldrich Denmark A/S · DK-2665 Vallensbaek Strand FINLAND · Sigma-Aldrich Finland/YA-Kemia Oy · FIN-00700 Helsinki
 FRANCE · Sigma-Aldrich Chimie · 38297 Saint-Quentin-Fallavier Cedex GERMANY · Sigma-Aldrich Chemie GmbH · D-82041 Deisenhofen GREECE · Sigma-Aldrich (o.m.) Ltd. · Ilioupoli 16346, Athens
 HUNGARY · Sigma-Aldrich Kft. · H-1067 Budapest INDIA · Sigma-Aldrich Co. · Bangalore 560 048 IRELAND · Sigma-Aldrich Ireland Ltd. · Dublin 24 ISRAEL · Sigma Israel Chemicals Ltd. · Rehovot 76100
 ITALY · Sigma-Aldrich s.r.l. · 20151 Milano JAPAN · Sigma-Aldrich Japan K.K. · Chuo-ku, Tokyo 103 KOREA · Sigma-Aldrich Korea · Seoul MALAYSIA · Sigma-Aldrich (M) Sdn. Bhd. · Selangor
 MEXICO · Sigma-Aldrich Quimica S.A. de C.V. · 50200 Toluca NETHERLANDS · Sigma-Aldrich Chemie BV · 3330 AA Zwijndrecht NORWAY · Sigma-Aldrich Norway · Torshov · N-0401 Oslo
 POLAND · Sigma-Aldrich Sp. z o.o. · 61-663 Poznań PORTUGAL · Sigma-Aldrich Quimica, S.A. · Sintra 2710 RUSSIA · Sigma-Aldrich Russia · Moscow 103062 SINGAPORE · Sigma-Aldrich Pte. Ltd.
 SOUTH AFRICA · Sigma-Aldrich (pty) Ltd. · Jet Park 1459 SPAIN · Sigma-Aldrich Quimica, S.A. · 28100 Alcobendas, Madrid SWEDEN · Sigma-Aldrich Sweden AB · 135 70 Stockholm
 SWITZERLAND · Supelco · CH-9471 Buchs UNITED KINGDOM · Sigma-Aldrich Company Ltd. · Poole, Dorset BH12 4QH
 UNITED STATES · Supelco · Supelco Park · Bellefonte, PA 16823-0048 · Phone 800-247-6628 or 814-359-3441 · Fax 800-447-3044 or 814-359-3044 · email: supelco@sial.com

Supelco is a member of the Sigma-Aldrich family. Supelco products are sold through Sigma-Aldrich, Inc. Sigma-Aldrich warrants that its products conform to the information contained in this and other Sigma-Aldrich publications. Purchaser must determine the suitability of the product for a particular use. Additional terms and conditions may apply. Please see the reverse side of the invoice or packing slip.