

CARRYOVER MITIGATION USING NEEDLE WASH SOLVENT CHEMISTRY AND AUTOSAMPLER FEATURES OF A UPLC-MS SYSTEM

Waters

THE SCIENCE OF WHAT'S POSSIBLE.™

Chris DesJardins, Zhimin Li, Patricia McConville
Waters Corporation, Milford, Massachusetts

INTRODUCTION

The greater sensitivity provided by mass detection is driving a need for improved carryover mitigation in LC applications. Carryover has many sources, but the most common source of carryover is related to the autosampler and occurs when sample components adhere to or absorb onto the outside of the needle. A key attribute of modern autosamplers is the ability to effectively remove sample from all active surfaces of the autosampler. The design of the sample flow path in the autosampler, features controlling needle washing, and the chemistry of the wash solvent all play a role in eliminating detectable carryover.

In this study, an ACQUITY UPLC H-Class PLUS system configured with an ACQUITY QDa Mass Detector is utilized to examine the impact of wash solvent chemistry and autosampler features for controlling carryover. Granisetron HCl, an antiemetic and anti-nausea agent commonly used in cancer therapy¹, was selected as the model compound for this study. The study is designed to quantify carryover to 0.0002% or 2pg on column, to this end the ACQUITY Diverter Valve, installed on the ACQUITY QDa, was programmed to divert the flow to waste when the highly concentrated granisetron challenge sample was injected. This precaution was taken to avoid saturation of the QDa mass detector which would interfere with accurate quantitation.

METHODS

UPLC System: ACQUITY UPLC H-Class PLUS System (176810005), ACQUITY QDa detector and a single column heater with the Active Pre-heating (CH-A).

Column: UPLC HSS T3, 1.8 µL 3.0 x 5.0 mm (186004679)

Column temp.: 35 °C

Sample temp.: 15 °C

Injection volume: 1 µL

Flow rate: 0.9 mL/min

Needle wash solvent:
Water:acetonitrile, 90:10 and 50:50
Water:methanol, 90:10 and 50:50
Acetonitrile, 100%
Methanol, 100%

Needle wash mode :
Default (6 seconds post-injection)
6 seconds pre & post-injection
12 seconds pre- and post- injection

Mobile Phase A: 0.1% Formic acid in water

Mobile Phase B: 0.1% Formic acid in acetonitrile

Gradient: Isocratic
80:20 (mobile phase A : mobile phase B)

Run time: 3 minutes,

SIR: 313.1 m/z

Data Management:
Chromatography Data Software
Empower 3 FR3 Hotfix 1

SAMPLES

Granisetron HCl (Catalog#: PHR1616) was purchased from MilliporSigma.

All solutions were prepared in 85:15, water:acetonitrile (diluent). The granisetron challenge solution was prepared at 1 mg/mL. This solution was then serially diluted to 0.0002, 0.0005, and 0.0010% of the challenge solution (2, 5 and 10 ng/mL respectively) which were used to generate the calibration curve for quantitation. These calibrators represent 2, 5 and 10 pg on column. The concentration of the calibrators was selected based on the mass detection linear dynamic range. Sample diluent was used for pre-standard and post-challenge blank sample injections.

RESULTS & DISCUSSION

STUDY DESIGN AND QUANTIFICATION OF CARRYOVER

Granisetron HCl was analyzed on an ACQUITY UPLC H-Class PLUS System configured with an ACQUITY QDa w/ ACQUITY Diverter Valve². Three replicate injections were performed for each wash solvent and needle wash mode tested.

The method selected to evaluate carryover employed a highly concentrated challenge solution. The ACQUITY Diverter Valve, installed on the ACQUITY QDa, was programmed to divert the flow to waste when the highly concentrated granisetron challenge sample was injected. This precaution was taken to avoid saturation of the QDa mass detector which would interfere with accurate quantitation.

A three-point standard calibration curve was generated to quantify carryover of the highly concentrated challenge sample. The calibration curve represents 0.0002% to 0.001% of the challenge sample. The concentration of the calibrators was selected based on the mass detection linear dynamic range. The sequence for this study was as follows: pre-standard blank, standard solutions (levels 1, 2 & 3), challenge solution, followed by 5 post-challenge blanks (Figure 1).

% Carryover was calculated using the following equation:

$$\% \text{ Carryover} = \left(\frac{\text{calculated concentration in post-challenge blank}}{\text{concentration of the challenge sample}} \right) * 100$$

The function of needle wash solvent is to solubilize the compound(s) of interest that remain on the sample needle after injection. Therefore, the composition of an optimal wash solvent (Figures 2 & 3) and needle wash mode (Figure 4) is application and compound specific.^{3,4} Important considerations when formulating wash solvent include the choice and ratio of organic to aqueous solvent to the wash solvent.^{5,6} Granisetron HCl is readily soluble in water⁷ therefore varying concentrations of (aq) and (org) in the needle wash solvent were investigated to demonstrate the impact of wash solvent composition.




Figure 1: Mass Detected Calibrators. The mass detected chromatographic results for the pre-standard blank and three calibrators. A flat baseline was established before the calibrators and challenge solution. The s/n for the calibrators was ≥ 9.




Figure 2: Carryover, Acetonitrile-Based Solvents. Carryover results of triplicate injections for granisetron using acetonitrile-based solvents as the needle rinse solvent. Needle wash program: Default, 6 sec. post-injection.




Figure 3: Carryover, Methanol-Based Solvents. Carryover results of triplicate injections for granisetron using methanol-based solvents as the needle rinse solvent. Needle wash program: Default, 6 sec. post-injection.




Figure 4: Needle wash Mode Impact on Carryover. Chromatographic results of observed carryover with varying needle wash modes using 100% ACN as the needle wash solvent. The extra needle rinse time improved carryover 3-fold.

CONCLUSION

Optimization of the needle wash solvent and wash settings are important aspects of LC method development. The needle wash solvent should be able to solubilize the components of the sample. A common starting place for many reversed phase applications is a mixture of water and the strong solvent used in the gradient. Investigation of the wash settings available on the autosampler should be included when optimizing needle wash during LC method development.

The efficient design of the ACQUITY UPLC H-Class PLUS FTN easily managed carryover from granisetron HCl providing greater flexibility to manage carryover and increased reliability of quantitation for LC applications with minimal impact to cycle time.

References

1. National Center for Biotechnology Information. PubChem Compound Database; CID=5284566, <https://pubchem.ncbi.nlm.nih.gov/compound/5284566> (accessed Nov. 29, 2018).
2. ACQUITY Diverter Valve Configuration [715005336 Rev. A](https://www.waters.com/715005336_Rev_A). January 2017.
3. Dolan, J. Autosampler carryover. LCGC Europe. Volume 19, Issue 10, pg 522-529. <http://www.chromatographyonline.com/autosampler-carryover-3?id=&sk=&date=&pageID=4>.
4. Dolan, J., Cooley, L. Reproducibility and Carryover-A Case Study. LCGC Volume 19, Number 3 pg 290-296. March 2001. <http://www.chromatographyonline.com/reproducibility-and-carryover-case-study-0>.
5. ACQUITY UPLC Sample Manager-Flow Through Needle PLUS Series Overview and Maintenance Guide [715005708 Rev. A](https://www.waters.com/715005708_Rev_A). March, 2018.
6. ACQUITY UPLC H-Class, H-Class Bio, and I-Class Series Solvent Considerations [715005742 Rev. B](https://www.waters.com/715005742_Rev_B). June 2018.
7. O'Neil, M.J., (Ed.) (2006). The Merck Index: An Encyclopedia of Chemicals, Drugs, and Biologicals (14th ed.). NJ: Merck. pp. 782-783.