

LIPIDOMIC PROFILING OF BLADDER CANCER PATIENTS USING A COMPACT LC-*oa*-TOF

Lisa Reid¹, Emmanuelle Claude¹, Adam King¹, Robert Plumb², Giorgis Isaac² and Lauren Mullin²
1: Waters Corporation, Wilmslow, UK; 2: Waters Corporation, Milford, MA, US

INTRODUCTION

- HRMS provides critical information for various Omics studies, affording identification of unique and reliable biomarkers of disease.
- Here, we explore the deployment of an easy-to-use, smart, compact benchtop LC-*oa*-TOF MS coupled with advanced statistical processing for a typical lipidomics workflow.
- This feasibility study was performed on plasma samples collected from bladder cancer patients, alongside healthy controls.



Progenesis[®] QI



Sample Preparation
100µL plasma from healthy control and bladder cancer were protein precipitated with 400µL IPA.

UPLC
System: Acquity UPLC I-Class
Column: UPLC CSH C₁₈ (100 x 2.1 mm, 1.7µm)
MPA: 60:40 ACN:water, 0.1% formic acid and 10mM ammonium formate
MPB: 90:10 IPA:ACN, 0.1% formic acid

MS
MS System: RDa
Mass Range: 50-2,000 m/z
Ionization Polarity: +/-
Informatics
Instrument Control: UNIFI 1.9.4
Processing: Progenesis QI with EZinfo with HMDB Searching

CONCLUSIONS

- LC-*oa*-TOF MS provides suitable dynamic range for labelled PC, LPC and SM lipids.
- Markers of significant association with bladder cancer patient samples were tentatively identified as a variety of lipid classes.
- This work has demonstrated a facile workflow for lipidomic profiling utilising a compact LC-*oa*-TOF MS.

Linearity of MS response for lipid classes afforded by compact LC-*oa*-TOF MS.




Figure 1: Using the Avanti SPLASH LipidoMIX, MS response is plotted from 40-20,000 ng/mL for 15:0-18:1(d7) PC; 4-40,000 ng/mL for d18:1-18:1(d9) SM; and 0.625-6250 ng/mL for 18:1(d7) Lyso PC. Mass errors were largely < +/- 5ppm for all injections.

Using advanced statistical processing, tentative identifications were proposed based on accurate mass as possible lipids of significance between bladder cancer and control groups.




Figure 2: Principal component analysis (PCA) from control and bladder cancer patient samples, and S-plot from which markers of interest are selected and then searched against HMDB tentative identification.