

# Quantification of Oligonucleotide Therapeutics in Plasma Using a Generic Kit-Based Approach

Margot Lee and Mary Trudeau

Waters Corporation, Milford, MA USA

CONTACT INFORMATION: [Margot\\_Lee@Waters.com](mailto:Margot_Lee@Waters.com) [Mary\\_Trudeau@Waters.com](mailto:Mary_Trudeau@Waters.com)

## PURPOSE

Developing LC-MS oligonucleotide bioanalytical sample preparation workflows are complex and laborious with no universal extraction method, often requiring skilled scientists and long method development times to achieve high extraction efficiency for highly sensitive, accurate and robust assays. This work leverages a generic kit-based approach with automated sample preparation and extraction to quantify multiple oligonucleotide therapeutics from plasma biomatrices, with minimal to no sample extraction method development while achieving high oligonucleotide recoveries (>75%), with excellent batch, day and user reproducibility.

## OBJECTIVE(S)

- Highlight a standardized, streamlined approach to oligonucleotide bioanalytical sample preparation using a prototype sample preparation and extraction kit.
- Demonstrate accurate and reproducible bioanalytical quantification of therapeutic oligonucleotides using the Prototype AX-SPE Microplate Kit automated on the Andrew+™ Pipetting Robot.

## METHOD(S)

### Oligonucleotides

- Oligonucleotide deoxythymidine nucleotide (dT) mix, containing a mix of 15-35 nucleotides, Waters Milford Mass (P/N 186004135).
- Gene expression modulators (GEM91 and GEM132), a 25-mer phosphorothioated antisense oligonucleotide and 20-mer phosphorothioated antisense oligonucleotide with 2' methoxy caps (Nitto Denko Avecia, Milford, MA)
- N-Acetylgalactosamine (GalNAc) conjugated -siRNA oligonucleotide (Gifted by Alnylam Pharmaceutical)
- Single-stranded DNA (ssDNA), 20-mer oligonucleotide, Waters Milford Mass (P/N 186009451)
- A lipid-conjugated 16-mer antisense oligonucleotide with 5' palmitate modification, a phosphorothioated backbone and terminal methoxy ethyl modifications (BioSearch Technologies, Lystrup, Denmark)

### LC-MS Analysis

LC-MS/MS analysis was performed using a Waters Xevo® TQ-XS tandem quadrupole MS (ESI-) and chromatographic separation using an ACQUITY I-Class PLUS UPLC® system and ACQUITY PREMIER Oligonucleotide BEH C<sub>18</sub>, 1.7 μm, 2.1 x 50 mm column (Waters, P/N 186009452). A flow rate of 0.6 mL/min and shallow gradient was employed using 1% HFIP (Hexafluoro-2-propanol) 0.1% DIPEA (N,N-Diisopropylethylamine) in H<sub>2</sub>O and acetonitrile. Total analysis time was 5 minutes. Injection volumes of the extracted sample were 1-30 μL.

### Sample Preparation and Extraction

Oligonucleotide concentrations in plasma were between 0.025-1, μg/mL and/or 0.01-1 pmol/μL. Starting sample volumes used for Prototype AX-SPE Microplate Kit development and evaluation were between 12.5-300 μL, with 100 μL being the optimal starting volume.




Figure 1. Graphical representation of the oligonucleotide sample preparation protocol using the oligonucleotide Prototype AX-SPE Microplate Kit.




Figure 2. Andrew+ Pipetting Robot, its OneLab deck Layout, dominos and pipettes used for plasma sample preparation extraction using the Prototype AX-SPE Microplate kit.

## RESULT(S)

Oligonucleotide	OST 15-ODT	OST 20-ODT	OST 25-ODT	OST 30-ODT	OST 35-ODT	GEM 91	GEM 132	GalNAc	Lipid Conj
Plasma Recovery	83.3	89.0	93.6	88.7	87.7	86.1	93.0	95.6	74.6
RSD	6.2	1.9	4.5	8.9	6.7	5.4	4.1	4.9	4.8
% Matrix Effects	-0.7	-14.5	-4	-12.6	1.4	-0.4	-10.3	-6.2	-35.2

Table 1. Prototype AX-SPE Microplate Kit extraction performance (no internal standard correction) demonstrating high plasma recovery (1 hour digestion @55°), low matrix effects, and intra-assay RSDs ≤15 % for a diversity of oligonucleotides.




Figure 3. Prototype AX-SPE Microplate Kit extraction performance (no internal standard correction) demonstrating high plasma recovery (1 hour digestion @55°), and inter-batch and intra-assay RSDs ≤15 % for a diversity of oligonucleotides.

Oligonucleotide Prototype Sample Prep Kit Reproducibility									
User-to-User % Plasma Recovery Variation									
	OST 15-ODT	OST 20-ODT	OST 25-ODT	OST 30-ODT	OST 35-ODT	GEM 91	GEM 132	GalNAc	Lipid Conj
User-User Day 1	1.7	0.1	6.5	-0.6	7.7	5.3	4.0	-0.2	12.6
User-User Day 2	3.0	7.4	13.6	-5.5	1.8	1.8	-0.2	-1.9	4.0
Day-to-Day % Plasma Recovery Variation									
Day-to-Day User 1	1.6	-3.8	-0.3	-8.4	8.3	0.4	-5.2	-3.6	-15.3
Day-to-Day User 2	2.9	3.4	6.8	-13.3	2.3	-3.0	-9.4	-5.2	-23.9

Table 2. Inter-day/inter-user Prototype AX-SPE Microplate Kit extraction performance (no internal standard correction) demonstrating ≤15 % difference in oligonucleotide plasma recoveries (N=4) from 100 μL starting sample volume across 2 users and 2 days using the standard starting protocol (1 hour digestion @55°).




Figure 4. Prototype AX-SPE Microplate Kit flexibility demonstrating linear LC-MS response for the GEM91 with extracted plasma volumes 12.5-300 μL (A) and for LC injection volumes from 1-30 μLs (B).

to download a copy of this poster, Visit [www.waters.com/posters](http://www.waters.com/posters)




Figure 5. Prototype AX-SPE Microplate Kit Andrew+ Pipetting Robot automated extraction performance (no internal standard correction) demonstrating high oligonucleotide plasma recovery, reproducibility, and comparability to manual preparation (1 hour digestion @55°).

Sensitive, Linear, Accurate & Precise					
Calibration Curve Statistics					
Analyte	Range	Weighting	Linear Regression	% Accuracy Range	% CV Range
GEM91	0.250-1000 ng/mL	1/x	>0.99	85.4-114.7	2.01-11.4
GalNAc				85.2-114.4	2.01-13.44
GEM132				85.9-119.2*	1.97-9.67
ss DNA (20-mer)	0.5-1000 ng/mL			85.7-112.4	0.99-13.87

\*% Accuracy of 119.2 for LOQ - Acceptable per Bioanalytical method validation guidelines

QC Statistics					
Analyte	QC Level	Expected concentration (ng/mL)	Mean observed concentration (ng/mL) (N=3)	Mean % Accuracy (N=3)	Mean % CV (N=3)
GEM91	LQC	0.75	0.74	98.17	6.42
GalNAc			0.69	92.30	2.90
GEM132			0.78	104.07	5.13
ss DNA (20-mer)	MQC	50	0.69	92.63	8.69
GEM91			52.97	105.95	2.82
GalNAc			49.79	99.56	6.42
GEM132	HQC	750	51.72	103.43	7.41
ss DNA (20-mer)			55.15	110.29	0.99
GEM91			756.55	100.87	13.61
GalNAc			733.97	97.87	13.44
GEM132			748.28	99.8	6.77
ss DNA (20-mer)			763.63	101.8	4.66

Figure 6. Prototype AX-SPE Microplate Kit Andrew+ Pipetting Robot automated extraction performance (no internal standard correction) demonstrating linear, accurate and precise quantitation.

## CONCLUSION(S)

This application successfully highlights a standardized, streamlined kit-based approach to oligonucleotide bioanalytical sample preparation. Using the Prototype AX-SPE Microplate Kit with optimized protocol, and pre-measured, QC-verified, lot-traceable detergent-free reagents, high analytical performance was achieved. Percent plasma recoveries of > 75%, with excellent inter and intra-kit, day-to-day and user-to-user reproducibility (RSDs ≤ 15%) for a diverse range of plasma volumes and oligonucleotide therapeutics and including unmodified, highly modified and both GalNAc and Lipid conjugated modalities was achieved. Automated sample preparation with Andrew+ Pipetting Robot, ensured linear, accurate, robust and reproducible quantitation performance with % accuracies within 15% of expected values and RSDs ≤15 %.

Waters Technology Library Literature Codes  
720008066EN  
720008067EN  
720008068EN

Waters  
THE SCIENCE OF WHAT'S POSSIBLE.™