



Agilent Prep LC Columns for Small Molecules and Biomolecules

**MAINTAIN RAPID, RELIABLE  
SEPARATIONS AS YOU SCALE-UP**



**Agilent Technologies**

## FLEXIBLE, COST-EFFECTIVE OPTIONS FOR SCALING AND PREPARATIVE SEPARATIONS

Scaling up from analytical HPLC to preparative HPLC can be costly and time-consuming. To avoid wasting supplies—and valuable samples—you must reliably predict the separation you will achieve when you scale to manufacturing levels.

Agilent can help you rise to the challenge—whether you are scaling a routine analytical method, or maintaining precise separations throughout every phase of production. Our wide range of preparative particle size columns and phases are designed for high loadability—and employ similar column packing materials as their analytical counterparts. Bulk material is also available in select materials, and can be ordered through our Custom Ordering Process at [www.agilent.com/chem/customlc](http://www.agilent.com/chem/customlc)

### Inside: our latest column technologies for increasing product yield and system throughput

**Agilent Prep LC columns for small molecules** ..... Page 3

A cost-effective first choice for purifying mg to g quantities. Available in both C18 and bare silica

**Additional Agilent Prep columns for small molecules**..... Page 4

- ZORBAX PrepHT columns: up to 7  $\mu\text{m}$  preparative scale
- Pursuit/Pursuit XRs and Polaris columns: scale current methods to 10  $\mu\text{m}$

**Agilent Prep LC columns for biomolecules**..... Page 5

- ZORBAX PrepHT columns: high purity, recovery, and throughput
- PLRP-S columns: thermally and chemically stable
- PL-SAX and PL-SCX columns: stability over a range of linear velocities

**Ordering information** ..... Page 8



Agilent Prep columns for HPLC let you purify microgram to gram quantities of sample with consistent results—making method transfer simple and predictable.

## PURIFY LARGE AMOUNTS OF TARGET COMPOUNDS

### Agilent Prep LC columns: A reliable first choice for scale-up

These highly loadable columns let you purify mg to g quantities of product, and are available in a variety of formats to accommodate most preparative samples. Other advantages include:

- Rapid method development, plus easy scalability from 4.6 to 50 mm
- Fast purification: 21.2 mm id cartridges increase throughput
- Combined high surface area and bonding density maximize loading
- Exceptional column stability and loadability up to pH 10
- Columns, guards, and bulk media are available in both C18 and bare silica



Agilent Prep LC column

#### Specifications


Bonded Phase	Pore Size	Surface Area	Temp. Limits	pH Range	Endcapped	Carbon Load
C18	100Å	400 m <sup>2</sup> /g	60 °C*	2.0-10.0	Yes	24%
Silica	100Å	400 m <sup>2</sup> /g	**	1.0-8.0	Yes	N/A

Specifications represent typical values only.

\*Temperature limits are 60 °C (up to pH 8) and 40 °C (from pH 8-10).

\*\*Temperature limits for bare silica are determined by the pH of the mobile phase.

#### Agilent Prep-C18 column has higher sample loading at pH 2 and 10



##### Conditions pH 2

Column: 4.6 x 150 mm, 5 µm  
 Mobile phase: 0.1% TFA in Water/ACN (ratio was adjusted at 41%-45%B to have same k)  
 Flow: 1.0 mL/min  
 Injection: 20 µL  
 Sample: Oxybutynin in DMSO

##### Conditions pH 10

Column: 4.6 x 150 mm, 5 µm  
 Mobile phase: 10 mM ammonia in Water/ACN (ratio was adjusted at 80%-85%B to have same k)  
 Flow: 1.0 mL/min  
 Injection: 20 µL  
 Sample: Oxybutynin in DMSO

Sample loading comparison of oxybutynin (a basic compound) at pH 2 and 10 in DMSO, using an Agilent Prep-C18 column and two non-Agilent columns.

Get fast, reproducible results with exceptional scalability.  
 Learn more at [www.agilent.com/chem/prep](http://www.agilent.com/chem/prep)

## Perform high-resolution scale-up with Agilent ZORBAX PrepHT columns

ZORBAX PrepHT columns are designed for rapid analytical to preparative scale-up *without* loss of resolution. For complex separations on larger columns (21.2 mm id; 150 mm length and longer), the 7 µm particle size balances high efficiency and high loadability. You also get the benefits of:

- Easy scale-up with ZORBAX phases
- Fast preparative separations—up to 2000 mg
- High efficiency and yield with 5 µm to 7 µm particles
- Easy installation: finger-tight connections effective up to 5000 psi/350 bar
- Analytical phase selectivity for your prep separations

Bulk material is available for ZORBAX 7 µm material, and can be ordered through our Custom Ordering Process at [www.agilent.com/chem/customlc](http://www.agilent.com/chem/customlc)

### ZORBAX PrepHT Column Specifications

Bonded Phase	Pore Size	Surface Area	pH Range	Endcapped	Carbon Load
Eclipse XDB-C18	80Å	180 m <sup>2</sup> /g	2.0-9.0 *	Double	10%
Eclipse XDB-C8	80Å	180 m <sup>2</sup> /g	2.0-9.0 *	Double	7.6%
SB-C18	80Å	180 m <sup>2</sup> /g	0.8-8.0 *	No	10%
SB-C8	80Å	180 m <sup>2</sup> /g	1.0-8.0 *	No	5.5%
SB-Phenyl	80Å	180 m <sup>2</sup> /g	1.0-8.0 *	No	5.5%
SB-CN	80Å	180 m <sup>2</sup> /g	1.0-8.0 *	No	0.04%
SB-Aq	80Å	180 m <sup>2</sup> /g	1.0-8.0 *	No	proprietary
Bonus-RP	80Å	180 m <sup>2</sup> /g	2.0-9.0	Triple	9.5%
Extend-C18	80Å	180 m <sup>2</sup> /g	2.0-11.5 *	Double	12.5%
Rx-SIL	80Å	180 m <sup>2</sup> /g	0-8.0	No	
Rx-C18	80Å	180 m <sup>2</sup> /g	2.0-8.0 *	No	12%
CN	70Å	300 m <sup>2</sup> /g	2.0-7.0	Yes	7%
NH2	70Å	300 m <sup>2</sup> /g	2.0-7.0	Yes	4%
SIL	70Å	300 m <sup>2</sup> /g	0-8.0	No	

\* At pH 6-9, highest column stability for all silica-based columns is obtained by operating at temperatures <40 °C and using lower buffer concentrations in the range of 0.01-0.02 M.



## Agilent Pursuit, Pursuit XRs, and Polaris columns help you maintain linearity

Pursuit and Polaris prep columns are available in the same phases as our Pursuit and Polaris analytical columns, with diameters up to 50 mm. Choose from chemistries of 3 µm or 5 µm in particle sizes, offering scalability up to 10 µm particle size.

- **Pursuit and Pursuit XRs prep columns** are prep-scalable for Pursuit and Pursuit XRs phases, and feature high surface-area silica. Choose from particle sizes to 10 µm, and column diameters up to 50 mm.
- **Polaris columns** are prep-scalable for Polaris phases, and are available in diameters of 10.0 mm and 21.2 mm, with particles up to 10 µm.

### Pursuit Column Specifications

Bonded Phase	Pore Size	Surface Area	pH Range	Endcapped	Carbon Load
Pursuit C18	200Å	200 m <sup>2</sup> /g	2 - 9	Yes	12.9%
Pursuit C8	200Å	200 m <sup>2</sup> /g	2 - 9	Yes	7.4%
Pursuit Diphenyl	200Å	200 m <sup>2</sup> /g	2 - 8	Yes	7.3%
Pursuit PFP	200Å	200 m <sup>2</sup> /g	2 - 9	Yes	6.3%

### Pursuit XRs Column Specifications

Bonded Phase	Pore Size	Surface Area	pH Range	Endcapped	Carbon Load
Pursuit XRs C18	100Å	440 m <sup>2</sup> /g	2 - 9	Yes	22%
Pursuit XRs C8	100Å	440 m <sup>2</sup> /g	2 - 9	Yes	15%
Pursuit XRs Diphenyl	100Å	440 m <sup>2</sup> /g	2 - 8	Yes	14.6%
Pursuit XRs Si	100Å	440 m <sup>2</sup> /g	2 - 8	Yes	N/A

### Polaris Column Specifications

Bonded Phase	Pore Size	Surface Area	Pore Volume	Endcapped	Carbon Load
Polaris C18-A	180Å	200 m <sup>2</sup> /g	1.1 cm <sup>3</sup> /g	Yes	13.8%
Polaris C18-Ether	180Å	200 m <sup>2</sup> /g	1.1 cm <sup>3</sup> /g	Yes	12.1%
Polaris Amide C18	180Å	200 m <sup>2</sup> /g	1.1 cm <sup>3</sup> /g	Yes	15%
Polaris Si-A	180Å	200 m <sup>2</sup> /g	1.1 cm <sup>3</sup> /g	N/A	N/A
Polaris C8-A	180Å	200 m <sup>2</sup> /g	1.1 cm <sup>3</sup> /g	Yes	7.4%
Polaris C8-Ether	180Å	200 m <sup>2</sup> /g	1.1 cm <sup>3</sup> /g	Yes	12.1%
Polaris NH2	180Å	200 m <sup>2</sup> /g	1.1 cm <sup>3</sup> /g	Amide	5.5%

# PRECISELY SCALE SEPARATIONS OF PROTEINS, PEPTIDES, AND OTHER LARGE MOLECULES


## Achieve high purity and recovery for biomolecules with Agilent ZORBAX PrepHT columns

ZORBAX PrepHT columns are designed for rapid analytical to preparative scale-up without loss of resolution. Benefits include:

- Higher throughput: Large column diameters and mechanically stronger ZORBAX particles allow for flow rates up to 100 mL/min
- Packed with 5  $\mu\text{m}$  to 7  $\mu\text{m}$  particles for high efficiency and yield
- Available in a variety of bonded phases—StableBond 300 $\text{\AA}$ , C18, C8, C3, and CN—for optimized resolution and loadability under any conditions

For complex separations on larger columns (21.2 mm id; 150 mm length and longer), we recommend the 7  $\mu\text{m}$  particle size to balance high efficiency and high loadability.

### SEC fraction of proteins with reversed-phase confirmation of purity



**Peak-based fractionation.** Fractions were re-analyzed using a reversed-phase C8 column. This confirmed the exact fractionation procedure using peak-based fraction trigger mode.

#### Conditions

Columns: Agilent Bio SEC-3, 300 $\text{\AA}$ , 7.8 x 300 mm, 3  $\mu\text{m}$   
 Agilent ZORBAX 300SB-C8, 4.6 x 50 mm, 8  $\mu\text{m}$   
 Detector: Agilent 1260 Infinity Bio-inert Quaternary LC System

#### SEC and fraction collection

Flow Rate: 1 mL/min  
 Gradient: Isocratic  
 Injection Volume: 30  $\mu\text{L}$   
 Thermostat autosampler and FC: 8  $^{\circ}\text{C}$   
 Temperature TCC: RT  
 DAD: 280 nm/4 nm  
 Ref.: OFF  
 Peak Width: >0.05 minutes  
 (1.0 second response time)(5 Hz)

#### Re-analysis – reversed-phase C8

Flow Rate: 1 mL/min  
 Gradient: 0 min 5% B, 95% C, 10 min 95% B, 5% C, Run time: 10 min, Stop time: 10 min  
 Injection Volume: 100, 50, and 10  $\mu\text{L}$   
 Thermostat autosampler and FC: 8  $^{\circ}\text{C}$   
 Temperature TCC: 70  $^{\circ}\text{C}$   
 DAD: 280 nm/4 nm  
 Ref.: OFF  
 Peak Width: >0.05 minute  
 (1.0 second response time) (5 Hz)

Get fast, reproducible results with exceptional scalability.  
 Learn more at [www.agilent.com/chem/prep](http://www.agilent.com/chem/prep)

## Agilent PLRP-S columns: High stability under demanding conditions

PLRP-S columns give you additional options for optimizing capacity and resolution – two key parameters for maximizing purification throughput. The thermal and chemical stability of PLRP-S makes it ideal for purifications in which sample preparation, compound elution, and column regeneration are performed under extreme conditions. Other advantages include:


- Less method development time: Particle sizes range from 3  $\mu\text{m}$  to 50  $\mu\text{m}$  for scale-up from the  $\mu\text{g}/\text{mg}$  discovery stage to multi-kg cGMP production
- The 600 L batch size for the production of PLRP-S media enables multiple preparative columns to be packed from the same batch of media, eliminating the need to validate multiple batches which saves time and increases productivity
- Excellent chemical stability, up to 1 M NaOH, permits sanitation and regeneration for longer column lifetime

The PLRP-S media, consisting of rigid poly(styrene/divinylbenzene) particles, are available in a range of pore sizes for small molecule, synthetic biomolecule, and macromolecule purification.

### Specifications

pH Range	Buffer Content	Organic Modifier	Temperature Limits	Maximum Pressure
1-14	Unlimited	1-100%	200 °C	5-8 $\mu\text{m}$ : 3000 psi (210 bar) 3 $\mu\text{m}$ : 4000 psi (300 bar)

### Purification of a 25-mer trityl-off oligonucleotide and analytical quantitation of the fraction using Agilent PLRP-S 100Å, 4.6 x 50 mm



#### Conditions

Column:	Agilent PLRP-S 100Å PL1512-1300 4.6 x 50 mm, 3 $\mu\text{m}$
Mobile Phase:	A: 100 mM Triethylammonium acetate (TEAA) B: 100 mM TEAA in 25:75 Acetonitrile: water
Flow Rate:	1 mL/min
Gradient:	25% B 0 min, 35% B 2 min, 45% B 22.5 min, 45% B 23 min, 25% B 23.05 min, 25% B 26 min
Temperature:	80 °C

The full length oligo(n) can be purified from the failure sequences, including n-1.

## Choose Agilent PL-SAX and PL-SCX columns for prep-to-process applications

These robust ion-exchange columns can be used over a wide range of linear velocities, with fast loading of dilute solutions and wash cycles. Both the polymeric PL-SAX and PL-SCX materials are chemically and thermally stable over a range of HPLC conditions. Thermal stability lets you use denaturing conditions and stabilizing/solubilizing agents to purify target compounds, such as those encountered when purifying synthetic oligonucleotides with self-complementary sequences.


You can also count on:

- Excellent reproducibility and long column life: Our dynamic axial compression (DAC) column hardware packing is straightforward and highly efficient
- Strong, hydrophilic ion-exchange functionalities are covalently linked to a chemically stable polymeric particle facilitating purifications over a wider pH range
- Rapid HPLC flow rates and equilibration for reduced purification cycle times
- Faster sanitation and clean-up through increased column stability

The 1000Å pore size is ideal for high-capacity purifications, while the 4000Å gigaporous particles with improved mass transfer are best for large biomolecules and high-speed, high-resolution purifications.

### Specifications

	PL-SAX	PL-SCX
<b>Matrix</b>	Fully polymeric	Fully polymeric
<b>Pore Sizes</b>	1000Å, 4000Å	1000Å, 4000Å
<b>Particle Sizes</b>	10 µm, 30 µm	10 µm, 30 µm
<b>Bead Form</b>	Rigid spherical	Rigid spherical
<b>Functionality</b>	Quaternary amine	Sulfonic acid
<b>Pressure Stability</b>	3000 psi	3000 psi
<b>Temperature Stability</b>	80 °C	80 °C
<b>pH Range</b>	1 - 14	1 - 14
<b>Eluent Compatibility</b>	All anion-exchange buffers	All cation-exchange buffers
<b>Packed Bed Density</b>	0.39 g/mL	0.39 g/mL



Agilent offers instrumentation, columns, software, and services to help you meet all of your purification challenges.

To learn more about performing fast, reproducible biomolecule purification, visit [www.agilent.com/chem/prep](http://www.agilent.com/chem/prep)

## Ordering Information

# COLUMNS FOR SMALL MOLECULES

### Agilent Prep LC Columns

Size (mm)	Particle Size (µm)	C18	Silica	Size (mm)	Particle Size (µm)	C18	Silica
4.6 x 250	10	440910-902	440910-901	21.2 x 50	10	446910-702	446910-701
4.6 x 150	10	443910-902	443910-901	21.2 x 150	5	443905-702	443905-701
4.6 x 100	10	449910-902		21.2 x 100	5	449905-702	449905-701
4.6 x 250	5	440905-902	440905-901	21.2 x 50	5	446905-702	446905-701
4.6 x 150	5	443905-902	443905-901	Guard Cartridge Hardware		420420-901	420420-901
4.6 x 100	5	449905-902	449905-901	Guard Cartridges 2/pk		420212-902	420212-901
4.6 x 50	5	446905-902	446905-901	30.0 x 250	10	410910-302	410910-301
10.0 x 250	10	440910-802	440910-801	30.0 x 150	10	413910-302	413910-301
10.0 x 150	10	440905-802	440905-801	30.0 x 100	10	419910-302	419910-301
10.0 x 100	10	443905-802	443905-801	30.0 x 100	5	449905-302	449905-301
10.0 x 50	5	446905-802	446905-801	30.0 x 50	5	446905-302	446905-301
Guard Cartridge Hardware		420420-901	420420-901	50.0 x 250	10	410910-502	410910-501
Guard Cartridges 2/pk		420212-902	420212-901	50.0 x 150	10	413910-502	413910-501
21.2 x 250	10	410910-702	410910-701	50.0 x 100	10	419910-502	419910-501
21.2 x 150	10	413910-702	413910-701	50.0 x 100	5	449905-502	449905-501

### Agilent ZORBAX PrepHT Cartridge Columns

Size (mm)	Particle Size (µm)	Eclipse XDB-C18	Eclipse XDB-C8	SB-C18	SB-C8
21.2 x 250	7	977250-102	977250-106	877250-102	877250-106
21.2 x 150	7	977150-102	977150-106	877150-102	877150-106
21.2 x 150	5	970150-902	970150-906	870150-902	870150-906
21.2 x 100	5	970100-902	970100-906	870100-902	870100-906
21.2 x 50	5	970050-902	970050-906	870050-902	870050-906
PrepHT Guard Cartridge, 17 x 7.5	5	820212-925	820212-926	820212-920	820212-915
Guard Cartridge Holder		820444-901	820444-901	820444-901	820444-901
PrepHT Endfittings, 2/pk		820400-901	820400-901	820400-901	820400-901

Size (mm)	Particle Size (µm)	SB-CN	SB-Phenyl	SB-Aq	Bonus-RP
21.2 x 250	7	877250-105	877250-112	877250-114	878250-101
21.2 x 150	7			877150-114	878150-101
21.2 x 150	5			870150-914	868150-901
21.2 x 100	5			870100-914	868100-901
21.2 x 50	5			870050-914	868050-901
PrepHT Guard Cartridge, 17 x 7.5	5	820212-915	820212-915	820212-933	820212-928
Guard Cartridge Holder		820444-901	820444-901	820444-901	820444-901
PrepHT Endfittings, 2/pk		820400-901	820400-901	820400-901	820400-901

Size (mm)	Particle Size (µm)	Extend-C18	Rx-SIL	Rx-C18
21.2 x 250	7		877250-101	877967-102
21.2 x 150	7			
21.2 x 150	5	770150-902		
21.2 x 100	5	770100-902		
21.2 x 50	5	770050-902		
PrepHT Guard Cartridge, 17 x 7.5	5	820212-930	820212-919	820212-914
Guard Cartridge Holder		820444-901	820444-901	820444-901
PrepHT Endfittings, 2/pk		820400-901	820400-901	820400-901

Size (mm)	Particle Size (µm)	ODS (C18)	C8	CN	NH2	SIL
21.2 x 250	7	877952-102	877952-106	877952-105	877952-108	877952-101
PrepHT Endfittings, 2/pk		820400-901	820400-901	820400-901	820400-901	820400-901



### Agilent Pursuit Prep Columns

Size (mm)	Particle Size (µm)	C18	C8	Diphenyl	PFP
10.0 x 250	5	A3000250X100	A3030250X100	A3040250X100	A3050250X100
10.0 x 250	10	A6002250X100	A3032250X100		
10.0 mm Guard		A3000050G100	A3032050G100	A3000050G100	
21.2 x 250	5	A3000250X212			
21.2 x 250	10	A6002250X212	A3032250X212		
21.2 mm Guard		A3000030G212	A3000030G212	A3000030G212	

### Agilent Pursuit XRs Prep Columns

Size (mm)	Particle Size (µm)	XRs C18	XRs C8	XRs Diphenyl	XRs Diol	XRs Si
21.2 x 250	10	A6002250X212				A6004250X212
21.2 x 250	5	A6000250X212		A6020250X212	A6100250X212	
21.2 x 150	5	A6000150X212	A6010150X212			
21.2 x 100	5	A6000100X212	A6010100X212	A6020100X212		
21.2 x 50	5	A6000050X212				
30.0 x 250	10	A6002250X300				A6004250X300
30.0 x 150	10	A6002150X300				
30.0 x 250	5	A6000250X300	A6010250X300			
30.0 x 150	5	A6000150X300				
30.0 x 100	5	A6000100X300		A6020100X300		
50.0 x 250	10	A6002250X500		A6022250X500		A6004250X500

### Agilent Polaris Prep Columns

Size (mm)	Particle Size (µm)	C18-A	C18-Ether	Amide C18	Si-A
10.0 x 250	5	A2000250X100	A2020250X100	A2006250X100	
10.0 x 250	10			A2008250X100	
21.2 x 250	5	A2000250X212	A2030250X212		A2003250X212
21.2 x 250	10	A2002250X212			A2004250X212

Size (mm)	Particle Size (µm)	C8-A	C8-Ether	NH2
10.0 x 250	5	A2010250X100	A2030250X100	A2013250X100
21.2 x 250	5	A2010250X212		A2013250X212

To place your order now, visit [www.agilent.com/chem/prep](http://www.agilent.com/chem/prep)

## Ordering Information

# COLUMNS FOR BIOMOLECULES

### Agilent ZORBAX PrepHT 300 StableBond Columns

Size (mm)	Particle Size (µm)	300SB-C18	300SB-C8-8	300SB-CN	300SB-C3
<b>PrepHT Cartridge Columns (require endfittings kit 820400-901)</b>					
21.2 x 250	7	897250-102	897250-106	897250-105	897250-109
21.2 x 150	7	897150-102	897150-106		897150-109
21.2 x 150	5	895150-902	895150-906		895150-909
21.2 x 100	5	895100-902	895100-906		895100-909
21.2 x 50	5	895050-902	895050-906		895050-909
PrepHT Endfittings, 2/pk		820400-901	820400-901	820400-901	820400-901
PrepHT Guard Cartridge, 17 x 7.5	5	820212-921	820212-918	820212-924	820212-924
Guard Cartridge Hardware		820444-901	820444-901	820444-901	820444-901

### Agilent PLRP-S Columns

Size (mm)	Particle Size (µm)	PLRP-S 100Å	PLRP-S 300Å	PLRP-S 1000Å	PLRP-S 4000Å
100 x 300	30			PL1812-3102	PL1812-3103
100 x 300	15-20	PL1812-6200	PL1812-6201	880995-902	880995-906
100 x 300	10-15	PL1812-6400	PL1812-6401	883995-902	883995-906
100 x 300	10	PL1812-6100	PL1812-6101	860950-902	860950-906
100 x 300	8	PL1812-6800	PL1812-6801	863973-902	863973-906
50 x 300	8	PL1712-6800	PL1712-6801	861973-902	861973-906
50 x 150	30			PL1712-3702	PL1712-3703
50 x 150	15-20	PL1712-3200	PL1712-3201	863974-302	863974-306
50 x 150	10-15	PL1712-3400	PL1712-3401		861973-306
50 x 150	10	PL1712-3100	PL1712-3101	PL1712-3102	PL1712-3103
50 x 150	8	PL1712-3800	PL1712-3801	883750-902	883750-906
25 x 300	15-20	PL1212-6200	PL1212-6201		863750-906
25 x 300	10-15	PL1212-6400	PL1212-3401	861772-902	861775-906
25 x 300	10	PL1212-6100	PL1212-6101	865750-902	865750-906
25 x 300	8	PL1212-6800	PL1212-6801	861630-902	
25 x 150	30			PL1212-3702	PL1212-3703
25 x 150	10	PL1212-3100	PL1212-3101	PL1712-3102	PL1712-3103
25 x 150	8	PL1212-3800	PL1212-3801	5185-5920	5185-5920
25 x 50	10			PL1212-1102	PL1212-1103

### Agilent PLRP-S Bulk Media


Particle Size (µm)	Unit	PLRP-S 100Å	PLRP-S 300Å	PLRP-S 1000Å	PLRP-S 4000Å
50	1 kg	PL1412-6K00	PL1412-6K01	PL1412-6K02	
	100 g	PL1412-4K00	PL1412-4K01	PL1412-4K02	
30	1 kg			PL1412-6702	PL1412-6703
	100 g			PL1412-4702	PL1412-4703
15-20	1 kg	PL1412-6200	PL1412-6201	861973-906	
	100 g	PL1412-4200	PL1412-4201		
10-15	1 kg	PL1412-6400	PL1412-6401		
	100 g	PL1412-4400	PL1412-4401		
10	1 kg	PL1412-6100	PL1412-6101	PL1412-6102	PL1412-6103
	100 g	PL1412-4100	PL1412-4101	PL1412-4102	PL1412-4103
8	1 kg	PL1412-6800	PL1412-6801		

### Agilent PL-SAX and PL-SCX Columns

Size (mm)	Particle Size (µm)	PL-SAX 1000Å	PL-SAX 4000Å	PL-SCX 1000Å	PL-SCX 4000Å
100 x 300	30	PL1851-3102	PL1851-3103	PL1845-3102	PL1845-3103
100 x 300	10	PL1851-2102	PL1851-2103	PL1845-2102	PL1845-2103
50 x 150	30	PL1751-3702	PL1751-3703	PL1745-3702	PL1745-3703
50 x 150	10	PL1751-3102	PL1751-3103	PL1745-3102	PL1745-3103
25 x 150	30	PL1251-3702	PL1251-3703	PL1245-3702	PL1245-3703
25 x 150	10	PL1251-3102	PL1251-3103	PL1245-3102	PL1245-3103
25 x 50	10	PL1251-1102	PL1251-1103	PL1245-1102	PL1245-1103
7.5 x 150	8	PL1151-3802	PL1151-3803		
7.5 x 50	8	PL1151-1802	PL1151-1803	PL1145-1802	PL1145-1803

### Agilent PL-SAX and PL-SCX Bulk Media

Particle Size (µm)	Unit	PL-SAX 1000Å	PL-SAX 4000Å	PL-SCX 1000Å	PL-SCX 4000Å
30	1 kg	PL1451-6702	PL1451-6703	PL1445-6702	PL1445-6703
	100 g	PL1451-4702	PL1451-4703	PL1445-4702	PL1445-4703
10	1 kg	PL1451-6102	PL1451-6103	PL1445-6102	PL1445-6103
	100 g	PL1451-4102	PL1451-4103	PL1445-4102	PL1445-4103



#### Quickly find the best columns for your application

Pressed for time? Agilent's FREE LC Columns and Sample Prep Navigator lets you:

- Explore columns based on method parameters, compound, and method
- Find sample prep guidance and method optimization tips
- Get selection assistance from Agilent chromatography experts
- Obtain technical support contact information

Get started now at [www.agilent.com/chem/navigator](http://www.agilent.com/chem/navigator)

To place your order now, visit [www.agilent.com/chem/prep](http://www.agilent.com/chem/prep)

## For more information

To learn more about Agilent prep LC columns for small molecules and biomolecules, visit

[www.agilent.com/chem/prep](http://www.agilent.com/chem/prep)

Find an Agilent customer center in your country:

[www.agilent.com/chem/contactus](http://www.agilent.com/chem/contactus)

U.S. and Canada:

**1-800-227-9770**

[agilent\\_inquiries@agilent.com](mailto:agilent_inquiries@agilent.com)

Europe:

[info\\_agilent@agilent.com](mailto:info_agilent@agilent.com)

Asia Pacific:

[inquiry\\_Isca@agilent.com](mailto:inquiry_Isca@agilent.com)



## Broad chemistries put you in control of your analyses

Agilent manufactures columns and media that suit nearly every technique for small molecules and biomolecules – allowing you to scale-up from prep to manufacturing levels. What's more, you can be confident that Agilent's meticulous end-to-end production oversight ensures the highest column consistency and performance.

## Purify your samples with maximum efficiency

For isolation and purification of your samples with highest purity and recovery, Agilent InfinityLab LC Purification Solutions offer high-performance instrumentation, columns, software, and services for analytical to preparative-scale workflows. A comprehensive and scalable portfolio based on a single platform gives you the choice to tailor a system to meet your laboratory's current and future needs.

Learn more: [www.agilent.com/chem/liveprepc](http://www.agilent.com/chem/liveprepc)



This information is subject to change without notice.

© Agilent Technologies, Inc. 2017  
Printed in the USA June 9, 2017  
5991-3155EN



**Agilent Technologies**