

## CHEMICAL & ENERGY ANALYSIS

# SPECIATION OF CHLORINATED HYDROCARBONS IN REFORMATE USING THE AGILENT 7200 GC/Q-TOF

Solutions for Your Analytical Business  
Markets and Applications Programs

### Authors

*Frank David, Bart Tienpont  
and Pat Sandra*

Research Institute for Chromatography,  
Kortrijk, Belgium

*Fabien Chainet, Michel Thomas  
and Eric Lemaire*

IFP Energies Nouvelles, Solaize, France

*Sofia Nieto*

Agilent Technologies, Santa Clara, USA



### ABSTRACT

The presence of chlorine in petroleum products can be harmful to equipment and refining processes. In addition to total chlorine determination, hyphenated techniques are needed to identify and quantitate the different chlorinated species. Speciation of chlorinated hydrocarbons in reformat was performed using the Agilent 7200 GC/Q-TOF high resolution mass spectrometer. Using extracted ion chromatograms at accurate mass, excellent selectivity and sensitivity were achieved, enabling the detection and quantitation of chlorinated hydrocarbons in the very complex reformat matrix. These compounds were determined at concentrations below 2 pg on-column, corresponding to 0.5 ppm (mg/kg) in the reformat matrix.

### INTRODUCTION

The presence of chlorine in the petroleum industry is of considerable importance because of the potential formation of hydrogen chloride in liquid or gas streams [1]. This chloride-containing gas can deactivate downstream catalysts and cause undesired reactions. Chlorine is a poison recovered on Pd/Al<sub>2</sub>O<sub>3</sub> catalysts during selective hydrogenation [2]. The presence of HCl in a hydrocarbon matrix leads to the formation of chlorinated hydrocarbons and can also promote the polymerization of olefins to produce “green” oils, mainly containing C6-C18 hydrocarbons [3]. Even if the concentration of HCl is relatively low, it can still interfere with operating processes that use hydrogen and can also cause corrosion problems in equipment such as pipes, valves and compressors.

In addition, HCl is considered a hazardous material and must be eliminated before release to the environment [1]. Recently, a review on total chlorine determination was published by Doyle et al. [4]. Chlorine can be measured using X-ray fluorescence (XRF) with a quantitation limit (LOQ) of 4 mg/kg in petroleum products and derivatives [5], by ICP-OES [4] and also by ICP/MS [6]. Micro-coulometric methods are currently the most widely used techniques for the determination of total chlorine, using standard test methods ASTM D4929 [7] or NF EN 14077, with an LOQ of around 1 mg/kg (ppm).



However, these methods determine the total chlorine concentration rather than identifying specific organochlorine compounds. In the petroleum industry, speciation of chlorine is necessary as these compounds need to be identified and quantitated, in order to develop systems for trapping these molecules. To our knowledge, there are no speciation methods for chlorine in petroleum products available in the literature.

Typically, chlorinated organic compounds are analyzed by gas chromatography (GC) and liquid chromatography (LC) coupled to mass spectrometry or specific detectors such as electron capture detector (ECD), atomic emission detector (AED) and inductively coupled plasma mass spectrometry (ICP/MS). GC-ECD lacks specificity and the ECD detector response is strongly dependent on the organic molecule [8, 9]. Chlorobenzene, for instance, gives a very low response in ECD. GC-AED can also be used, but chlorine detection is not very sensitive and the co-elution of hydrocarbons at very high levels (3-4 orders of magnitude higher than the chlorine compounds) exceeds the selectivity of the detector [10]. Despite spectral interferences at  $m/z$  35 and low ionization potential for chlorine, GC-ICP/MS can be performed, as shown by Peters et al. [11], for the detection of three compounds (dichloromethane, trichloroethane and trichloroethylene) in pentane with LODs around 2 mg/L. However, GC-ICP/MS configurations are not widely available.

Due to the complexity of reformat samples, containing about 250 hydrocarbon species, and the chlorine concentration at low (or sub-) ppm level, selective and sensitive techniques must be used. Compound specific analysis can be performed using two-dimensional GC (heart-cutting) [12] or GC-MS/MS in multi-reaction monitoring (MRM) mode [13]. However, these techniques can only be applied to selected target compounds that are available as standards and for which retention time and mass spectra (MRM transitions) can be determined. In organochlorine speciation, the availability of a screening (non-target) method that allows for the detection and quantitation of chlorinated hydrocarbons (alkanes and aromatics) in gasoline, naphtha or reformat samples at low ppm levels would offer interesting possibilities for the petrochemical industry.

In this application note, the use of the Agilent 7200 GC/Q-TOF high resolution mass spectrometer for the analysis of chlorinated hydrocarbons in reformat is described. Ion extraction at accurate masses is used for detection and quantitation of the target compounds.

## **EXPERIMENTAL**

**Samples:** A reformat sample from IFP Energies Nouvelles (Solaize, France) was used as the test matrix. Basic properties of this reformat are listed in Table 1. A stock solution of 10 chlorinated hydrocarbons was prepared in heptane at a concentration of 100 mg/kg (ppm). The reference compounds are listed in Table 2. A diluted stock solution was also prepared in heptane at 10 ppm. Spiked matrix solutions at concentrations between 0.5 and 20 ppm (per chlorinated compound) were prepared by spiking the reformat test matrix with the stock solution.

### **GC/Q-TOF instrument:**

Analyses were performed on an Agilent 7890B GC equipped with a split/splitless inlet. Injection was performed using an Agilent 7693 ALS. The GC was coupled to an Agilent 7200 Q-TOF mass spectrometer.

### **Analytical conditions:**

0.2  $\mu$ L of sample was injected in split mode (1/50) at 250°C. Separation was performed on a 100 m x 0.25 mm ID x 0.5  $\mu$ m HP-PONA column (19091Z-530). The carrier gas was helium at 1.5 mL/min constant flow. The oven was programmed from 35°C (5 min) at 4°C/min to 250°C (10 min). The transfer line was set at 280°C.

The 7200 Q-TOF was operated in MS mode using electron impact (EI) ionization. The source temperature was 230°C. The scan rate was 5 Hz in HR (high resolution) mode and the mass range was 50 to 300 amu. No internal mass referencing was used, but the instrument was mass calibrated before each run using the keyword command (MassCal) in the sequence.

## RESULTS AND DISCUSSION

The total ion chromatogram obtained for the analysis of the diluted calibration stock (10 ppm) is shown in Figure 1. Taking into account the injection volume (0.2  $\mu\text{L}$ ) and the split ratio (1/50), the injected amount corresponds to 40 pg per compound on-column. Seven compounds were easily detected. However, compounds 2, 3 and 4 were masked by the heptane solvent peak. For the detected compounds, accurate mass spectra were obtained and searched against the standard NIST library. The chloroalkanes from chloropentane to chlorononane, all have the ion at  $m/z$  91.031 in common.




Fig 1: Base Peak Chromatogram of the analysis of a test solution containing chlorinated hydrocarbons in heptane (10 ng/ $\mu\text{L}$  concentration, 40 pg on-column).

Using the formula generator option, this fragment ion was correctly identified as  $[\text{C}_4\text{H}_8\text{Cl}]^+$ , as illustrated by the high resolution mass spectrum of chloroheptane (compound 7) in Figure 2. At  $m/z$  93.028 the chlorine isotope peak is detected. For these ions, the mass error was found to be lower than 4 ppm.




Fig. 2: High resolution accurate mass spectrum of peak 7 (1-chloroheptane)

Using extracted ion chromatograms at accurate mass, very high selectivity can be obtained. This is illustrated in Figure 3. In Figure 3A, the total ion chromatograms are compared for a reformate sample (top) and the same reformate sample spiked at 1 ppm level with the chloroalkanes (bottom). These profiles show the high complexity of the reformate sample, mainly containing aromatic hydrocarbons (toluene, xylenes, C3-, C4-, C5- benzenes) and paraffins. If the ion chromatograms are extracted using a window of  $\pm 0.5$  amu, as is done on low resolution instruments such as a single quadrupole or a triple quadrupole instruments, the chloroalkane compounds cannot be selectively detected in the spiked sample.




Fig. 3A: Comparison of TIC obtained for reformat sample (top) and reformat sample spiked at 1 ppm level with chlorinated hydrocarbons (bottom).




Fig. 3B: Comparison of EIC obtained for reformat sample (top) and reformat sample spiked at 1 ppm level with chlorinated hydrocarbons (bottom). EIC performed at unit mass resolution (91.0309 +/- 0.5 amu)




Fig. 3C: Comparison of EIC obtained for reformat sample (top) and reformat sample spiked at 1 ppm level with chlorinated hydrocarbons (bottom). EIC performed at high mass resolution (91.0309 +/- 20 ppm)

This is illustrated in Figure 3B, showing the EIC at  $m/z$  91.0309  $\pm$  0.5 amu. No difference can be observed between the sample (top) and the spiked sample (bottom). The main peaks detected in these EICs are the aromatic hydrocarbons that give the well known tropylium ion at  $m/z$  91 in their mass spectra. This fragment corresponds to  $[C_7H_7]^+$  and the exact mass is 91.0542. The first 4 peaks correspond to toluene (19.5 min), ethylbenzene, m- and p- xylene and o-xylene (24 - 26 min). These compounds all have the tropylium ion as the most abundant ion in their mass spectra and low resolution MS is unable to differentiate the tropylium ion from the  $[C_4H_8Cl]^+$  ion. If ions are extracted at a mass accuracy of  $\pm$  20 ppm, the target compounds can be selectively detected, as illustrated in Figure 3C. In the unspiked sample, no chloroalkanes were detected (top chromatogram), while compounds 6, 7, 9 and 10 were detected in the spiked sample.

The EIC signals at accurate mass for a standard (10 ppm of compounds in heptane, corresponding to 40 pg on-column) and for a spiked matrix sample (10 ppm compounds in matrix) are compared in Figure 4. The peak areas are quite similar, indicating that no quenching or ion suppression occurs in the matrix. Only slight differences are noted in retention time and peak width, caused by matrix effects on the chromatography (solvent effects).

Using this method, all chlorinated compounds were detected at low ppm or even sub-ppm levels in the matrix (low ppb levels on column). To take maximum advantage of the selectivity, the ions chosen for EIC extraction all contained chlorine (see Table 2). Even though in some cases chlorine-containing ions were of relatively low abundance in a compound spectrum, it allowed efficient detection of the compounds in the very complex matrix. This is illustrated in Figure 5 for chlorobenzene using the ion at 112.007 and in Figure 6 for benzylchloride using the ion at  $m/z$  126.023. Ion chromatograms were extracted at  $\pm$  20 ppm  $m/z$  window. The high resolution mass spectra are shown in these figures and the mass accuracy was calculated using the formula generation option. Chlorobenzene showed a mass error of -3.03 ppm and benzylchloride a mass error of -3.17 ppm. It should be emphasized that these mass errors were obtained in the complex reformat matrix.




Fig. 4: Comparison of EIC obtained for a 10 ppm standard solution (top) and reformat sample spiked at 10 ppm level with chlorinated hydrocarbons (bottom). EIC performed at high mass resolution (91.0309  $\pm$  20 ppm)




Fig. 5: Detection of chlorobenzene in reformat (spiked at 10 ppm level) - EIC @112.0074 +/- 20 ppm  
Mass spectrum is given in insert




Fig. 6: Detection of benzylchloride in reformat (spiked at 10 ppm level) - EIC @126.0231 +/- 20 ppm  
Mass spectrum is given in insert

PROPERTIES	METHODS	VALUES
Density (g/cm <sup>3</sup> )	NF EN ISO 12185	0.8316
Boiling range (°C)	ASTM D2887	82-197
RON	IFPEN 9302	101
Paraffins (%wt)		5.0
Isoparaffins (%wt)		14.7
Naphthenes (%wt)		1.1
Aromatics (%wt)		78.4
Olefins (%wt)		0.8
Chlorine (mg/kg)	NF EN 14077	< 2

Table 1

PEAK	TR (MIN)	SOLUTE	FORMULA	MW	EXACT MASS	FRAGMENT
1	9.949	2-chloro-2-methylpropane	C <sub>4</sub> H <sub>9</sub> Cl	92	77.0158	C <sub>3</sub> H <sub>6</sub> Cl
2	16.463	1-chloro-2-methylbutane	C <sub>5</sub> H <sub>11</sub> Cl	106	62.9996	C <sub>2</sub> H <sub>4</sub> Cl
3	16.747	2-chloro-pentane	C <sub>5</sub> H <sub>11</sub> Cl	106	62.9996	C <sub>2</sub> H <sub>4</sub> Cl
4	18.784	1-chloropentane	C <sub>5</sub> H <sub>11</sub> Cl	106	62.9996	C <sub>2</sub> H <sub>4</sub> Cl
5	23.390	chlorobenzene	C <sub>6</sub> H <sub>5</sub> Cl	112	112.0074	C <sub>6</sub> H <sub>5</sub> Cl
6	23.955	1-chlorohexane	C <sub>6</sub> H <sub>13</sub> Cl	120	91.0309	C <sub>4</sub> H <sub>8</sub> Cl
7	28.834	1-chloroheptane	C <sub>7</sub> H <sub>15</sub> Cl	134	91.0309	C <sub>4</sub> H <sub>8</sub> Cl
8	30.921	benzylchloride	C <sub>7</sub> H <sub>7</sub> Cl	126	126.0231	C <sub>7</sub> H <sub>7</sub> Cl
9	33.334	1-chlorooctane	C <sub>8</sub> H <sub>17</sub> Cl	148	91.0309	C <sub>4</sub> H <sub>8</sub> Cl
10	37.484	1-chlorononane	C <sub>9</sub> H <sub>19</sub> Cl	162	91.0309	C <sub>4</sub> H <sub>8</sub> Cl

Table 2

## CONCLUSION

Sub-ppm chlorinated hydrocarbons can be analyzed in reformat samples using GC/Q-TOF high resolution MS analysis. The high selectivity offered by accurate mass ion extraction allows the target compounds to be detected and quantified in this extremely complex matrix. Sensitivity is better than 2 pg on column, corresponding to 0.5 ppm (mg/kg) in the reformat matrix.

## REFERENCES

- [1] O.K. Karan, M.A. Ay, K. Karhaman, A. Selmen, "Combating green oil formation in a CCR reformer", PTQ 2013, Q3, 29.
- [2] B. Didillon, J. Cosyns, C. Cameron, D. Uzio, P. Sarrazin, J.P. Boitiaux, "Industrial evaluation of selective hydrogenation catalyst poisoning", Stud. Surf. Sci. Catal., 1997, 111, 447.
- [3] J. Cosyns, O. Ducreux, Q. Debuisschert, and F. Lepeltier, Removal of chlorinated compounds in hydrocarbon cuts, (6-6-2012)
- [4] A. Doyle, A. Saavedra, M.L.B. Tristão, L.A.N. Mendes, R.Q. Aucélio, "Spectrometric methods for the determination of chlorine in crude oil and petroleum derivatives — A review", Spectrochim. Acta B, 2013, 86, 102.
- [5] ASTM D6443, Standard Test Method for Determination of Calcium, Chlorine, Copper, Magnesium, Phosphorus, Sulfur, and Zinc in Unused Lubricating Oils and Additives by Wavelength Dispersive X-ray Fluorescence Spectrometry (Mathematical Correction Procedure) (2004).
- [6] F.G. Antes, M.F.P. dos Santos, R.C.L. Guimarães, J.N.G. Paniz, E.M.M. Flores, V.L. Dressler, "Heavy crude oil sample preparation by pyrohydrolysis for further chlorine determination", Anal. Methods, 2011, 3, 288.
- [7] ASTM D4929, Standard Test Methods for Determination of Organic Chloride Content in Crude Oil (2007).
- [8] A. Schwarz, K.G. Heumann, "Two-dimensional on-line detection of brominated and iodinated volatile organic compounds by ECD and ICP-MS after GC separation", Anal. Bioanal. Chem., 2002, 374, 212.
- [9] S. Pedersen-Bjergaard, S.I. Semb, J. Vedde, E.M. Brevik, T. Greibrokk, "Comparison of GC-ECD, GC-MS and GC-AED for the determination of polychlorinated biphenyls in highly contaminated marine sediments", Chromatographia, 1996, 43, 44.
- [10] C. Brede, S. Pedersen-Bjergaard, "State-of-the art of selective detection and identification of I-, Br-, Chlorine-, and F-containing compounds in gas chromatography and liquid chromatography", J. Chromatogr. A, 2004, 1050, 45.
- [11] G.R. Peters, D. Beauchemin, "Versatile interface for gas chromatographic detection or solution nebulization analysis by inductively coupled plasma mass spectrometry: preliminary results", J. Anal. At. Spectrom., 1992, 7, 965.
- [12] F. Chainet, C.P. Lienemann, J. Ponthus, M. Courtiade, O.F.X. Donard, "Development of heart-cutting multidimensional gas chromatography coupled to time of flight mass spectrometry for silicon speciation at trace levels in gasoline samples", J. Chromatogr. A, 2012, 1264, 80.
- [13] N.S. Chary, S. Herrera, M.J. Gómez, A.R. Fernández-Alba, "Parts per trillion level determination of endocrine-disrupting chlorinated compounds in river water and wastewater effluent by stir-bar-sorptive extraction followed by gas chromatography-triple quadrupole mass spectrometry", Anal. Bioanal. Chem., 2012, 404, 1993.



**Solutions for your analytical business**

Markets & Applications Programs

[www.solutions-to-win.com](http://www.solutions-to-win.com)

The Measure of Confidence



**Agilent Technologies**

Agilent Products are for Research Use Only.  
Not for use in diagnostic procedures.  
Information, descriptions and specifications in this publication are subject to change without notice.

© Agilent Technologies, Inc. 2015  
Published in USA, March 26, 2015  
5991-5684EN