



The Best of Both Worlds: LC/Q-TOF as a Method to Detect a Targeted List of 35 Drugs and Metabolites in Urine with Retrospective Data Mining Capabilities

Application Note

Clinical Research

Authors

Natalie N. Rasmussen¹, Julie Cichelli³,
Carrie Adler³, Adam Barker^{1,2}, and
Frederick G. Strathmann^{1,2}

¹ ARUP Institute for Clinical and
Experimental Pathology
Salt Lake City, UT

² Department of Pathology
University of Utah,
Salt Lake City, UT

³ Agilent Technologies, Inc.
Santa Clara, CA USA

Abstract

Compliance testing continues to be an essential component of the drug testing market. As an alternative to conventional immunoassay screens with reflex to confirmation, we developed a dual identification strategy for drug screening and confirmation with one injection for clinical research purposes. Making use of the quadrupole-time-of-flight mass spectrometer's (Q-TOF) quadrupole, we combined fragmentation pattern with accurate mass, retention time, and isotope fidelity to identify analytes of interest with increased accuracy.



Agilent Technologies

Introduction

High-resolution accurate mass liquid chromatography mass spectrometry is ideally suited for rapid analysis of common drugs of abuse and other illicit drugs. A highly sensitive and specific method has been developed for the qualitative determination of 35 drugs and metabolites in urine. The target compounds were detected and confirmed at their respective cutoffs using an Agilent LC/Q-TOF system. An Agilent 6550 LC/Q-TOF coupled to an Agilent 1290 Infinity UHPLC stack provided analyte analysis in a fast 3.3 minute method. This method was developed for common drugs of abuse in urine to screen and confirm target compounds in a single run, with the ability to retrospectively review scan data for possible unknowns at a later date.

Compound identification is based on accurate mass, retention times, isotopic spacing and abundance, MS/MS spectral library matching, and overall score. MassHunter Qualitative Analysis uses peak picking algorithms to generate an overall score based on how the observed identification criteria correspond with theoretical and empirical values. Semiquantitation of compounds is determined with a one-point calibration forced through the origin based on the respective drug cutoffs designated by a previously published method¹.

Internal standards

- 6am-d6
- α OH Alprazolam-d5
- Diazepam-d5
- Meperidine-d4
- Methamphetamine-d5
- Morphine-d3
- Norfentanyl-d5

Experimental

LC Configuration and parameters

Configuration		
Agilent 1290 Infinity II high speed pump (G7120A)		
Agilent 1290 Infinity II multisampler (G7167B)		
Agilent 1290 Infinity II multicolumn thermostat (G7116B)		
Analytical column	Agilent InfinityLab Poroshell 120 EC-C18, 2.1 × 50 mm, 2.7 μ m LC column (699775-902)	
Column temperature	60 °C	
Injection volume	2 μ L	
Mobile phase A	5 mM ammonium formate in water, pH 3.5	
Mobile phase B	Optima grade methanol	
Flow rate	1 mL/min	
Gradient	Time (min)	%B
	0.0	10
	0.5	10
	3.0	95
	3.3	95

Q-TOF Mass spectrometer configuration and parameters

Configuration	
Agilent 6550 Quadrupole Time of Flight Mass Spectrometer	
Electrospray, Positive ionization	
AutoMSMS mode, using a preferred list	
Drying gas temperature	250 °C
Drying gas flow	15 L/min
Nebulizer pressure	30 psi
Sheath gas temperature	400 °C
Sheath gas flow	12 L/min
Nozzle voltage	0 V
Capillary voltage	3,500 V
Fragmentor voltage	125 V
Skimmer voltage	65 V
Octopole RF	750 V
Mass range	MS 100–1,000 m/z
	MS/MS 50–500 m/z
Acquisition rate	MS 8 spectra/s; MS/MS 4 spectra/s
Collision energies	10, 20, and 40 V
Detector rate	Extended dynamic range, high sensitivity slicer position

Chemicals and reagents

All of the drug and metabolite compounds, as well as the deuterated internal standards, were purchased from Cerilliant (Round Rock, Texas). Methanol was purchased from ThermoFisher Scientific (Waltham, MA), and β -glucuronidase was purchased from IMCS (Columbia, SC).

Sample preparation

Samples, single point calibrator, and QC materials in matrix were prepared using the following procedure. First, 200 μ L of sample was added to a collection plate. Next, 180 μ L of sodium bicarbonate buffer containing seven deuterated markers, 34 μ L resorufin glucuronide solution, and 16 μ L of β -glucuronidase were added and the sample was gently mixed. Following a 15-minute room temperature incubation, 300 μ L of the sample mixture was transferred to a Biotage Isolute SLE+ 96-well plate (Charlotte, NC, USA) for further purification. Positive pressure was applied and the sample was allowed to absorb for 5 minutes. Samples were eluted with 0.8 mL ethyl acetate, dried under nitrogen, then resuspended using a 90:10 blend of mobile phase A and B.

Data analysis

System control and data acquisition were performed by Agilent MassHunter Quadrupole Time of Flight Acquisition Software (B.06) with Qualitative Analysis (B.07 SP2) Software. Data were analyzed using Find by Formula and library search software functions. Acceptance criteria included a match to retention time (within 0.015 minutes), isotope spacing and abundance, accurate mass (within 5 ppm), MS/MS spectral library matching, and overall score (>85).

Drug and metabolite cutoffs (ng/mL)

Compound	Cutoff limit
6-Monoacetyl-morphine	10
7-Aminoclonazepam	40
<i>alpha</i> -Hydroxyalprazolam	20
Alprazolam	40
Amphetamine	100
Buprenorphine	5
Clonazepam	20
Codeine	40
Diazepam	50
Fentanyl	2
Hydrocodone	40
Hydromorphone	40
Lorazepam	60
MDA (3,4-methylenedioxyamphetamine)	200
MDEA (3,4-methylenedioxyethamphetamine)	200
MDMA (methylenedioxymethamphetamine)	200
Methamphetamine	400
Midazolam	20
Morphine	20
Norbuprenorphine	20
Nordiazepam	50
Norfentanyl	2
Norhydrocodone	100
Normeperidine	50
Noroxycodone	100
Noroxymorphone	100
Oxazepam	50
Oxycodone	40
Oxymorphone	40
Phentermine	100
Ritalinic acid	100
Tapentadol	100
Tapentadol-O-sulfate	200
Temazepam	50
Zolpidem	20

Results and Discussion

De-identified samples (n = 160), previously analyzed with an Agilent 6230 Time of Flight mass spectrometer (TOF) were re-analyzed using the LC/Q-TOF system. Samples were prepared by adding sodium bicarbonate buffer with internal standard mix, resorufin glucuronide, and β -glucuronidase. Resorufin glucuronide acted as a hydrolysis control. After a short incubation period, the sample mixture was transferred to a supported liquid extraction plate for further purification. Samples were eluted with ethyl acetate, nitrogen dried, and resuspended in mobile phase. Samples were then analyzed using the LC/Q-TOF system and a reversed-phase C18 column. A single-point calibrator at the compound cutoffs was analyzed along with quality control samples.

The Agilent 6550 Q-TOF and 6230 TOF results showed excellent correlation within the 160-sample data set. Most discrepancies, such as finding additional metabolites, were due to the higher analytical sensitivity achieved with the 6550 Q-TOF instrument. In four cases, zolpidem and 7-aminoclonazepam were detected by the 6550 QTOF system that were not originally reported. These samples had concentrations below the cutoff defined with the 6230 TOF assay.

Five times within the data set, the 6550 Q-TOF system identified highly abundant drugs based on known parent drug metabolites. Alternatively, these samples could have been diluted for confirmatory identification of the parent drug.

Within the original 160-sample data set, nine samples required additional confirmation using a liquid chromatography triple quadrupole mass spectrometer (LC/MS/MS) system. These nine samples either had an unusual drug metabolite pattern or suboptimal quality metrics required for acceptance with the 6230 TOF assay. In contrast, the 6550 Q-TOF was able to identify and obtain matching fragmentation spectra for all of the analytes present in those nine samples without any additional testing.

Conclusions

The presented LC/Q-TOF method is highly specific and analytically sensitive, providing a fast and accurate way of identifying analytes. It eliminates the need for LC/MS/MS confirmation of suboptimal results because of the ability to fragment and match to a library. The method uses the ion transmission abilities of a TOF with the enhanced specificity of high-mass-accuracy detection of fragments, providing better specificity in comparison to conventional LC/TOF and LC/MS/MS methods. This method can be applied to a targeted screen, while the simultaneous collection of TOF data provides a mechanism for retrospective data analysis.

Reference

1. McMillin, G. A.; *et al.* A hybrid approach to urine drug testing using high-resolution mass spectrometry and select immunoassays. *Am. J. Clin. Pathol.* **2015 Feb**; *143*(2), 234-40.

www.agilent.com/chem

For Research Use Only. Not for use in diagnostic procedures.

This information is subject to change without notice.

© Agilent Technologies, Inc., 2017
Published in the USA, March 28, 2017
5991-7955EN



Agilent Technologies