

Rapid and Sensitive Analysis of Taurine in Human Serum and Beagle Plasma using Liquid Chromatography-Tandem Mass Spectrometry

Priyanka Chitranshi, Evelyn H. Wang, and Christopher T. Gilles, Shimadzu Scientific Instruments, Columbia, MD, USA 21046

Objective

This work focuses on the development and evaluation of an LC-MS/MS method for the identification of taurine in biological samples such as human serum and beagle plasma.

Introduction

Taurine is a non-essential aminosulfonic acid, often administered externally as a supplement in commercial products such as infant formulas, energy drinks, and contact lens solutions. Typically, chromatographic techniques such as (LC-UV, GC and ion chromatography) and capillary electrophoresis are used for the analysis of taurine. However, the above methods may require extensive sample preparation, and/or lack in sensitivity and specificity required for biological fluids. We report here a sensitive and specific method for the direct determination of taurine in human serum and beagle plasma using liquid chromatography-tandem mass spectrometry (LC-MS/MS). The method was developed and optimized on a Shimadzu LCMS-8050 triple quadrupole mass spectrometer coupled to a Nexera UHPLC liquid chromatograph.

Experimental




Table 1. Instrument parameters used for the analysis of taurine in serum and beagle plasma.

	Nexera UHPLC	LCMS-8050
Column:	PFPP column (150 × 2.1 mm; 5 μm)	Nebulizing Gas: 2 L/min
Mobile phase:	A: 0.03 % Formic acid in water; B: 0.03 % Formic acid in methanol	Heating Gas: 10 L/min
Flow Rate:	0.3 mL/min	Interface Temp.: 300 °C
Oven Temp.:	25°C	DL Temp.: 250 °C
Injection volume	Method 1: 0.5 μL Method 2: 1.0 μL	Heat Block Temp.: 400 °C
Internal Standard	Method 1: Sulfanilic acid Method 2: Aspartic acid	Drying Gas: 10 L/min

Compounds	Mode	Precursor ion (m/z)	Product ion (m/z)	CE (V)
Taurine	ESI -	124.9	79.9	22
Sulfanilic acid (IS)	ESI -	172.0	79.95	27
Aspartic acid (IS)	ESI -	132.1	88	13

Results and Discussion




Figure 1. Representative MRM chromatogram for taurine (0.05 μg/mL in serum) using sulfanilic acid as internal standard

- Calibration curve was generated for taurine in serum and beagle plasma.
- Calibrant range in human serum was 0.025-50 μg/mL and 0.5-60 μg/mL in beagle plasma.
- The curves demonstrated excellent linearity ($r \geq 0.999$).
- No carry over was observed.
- % accuracy of QC injected at low (0.3 μg/mL), medium (1 μg/mL), and high (25 μg/mL) concentration was calculated in serum between 95.2-98.8%.

Table 2. Detection of taurine in unknown beagle plasma samples.

Unknown #	Average (μg/mL)	SD (±)
1	7.4	0.2
2	13.4	0.3
3	15.7	0.8
4	21.1	1.2
5	17.7	0.5
6	36.5	1.2
7	10.2	0.8
8	9.7	0.9
9	12.5	0.1
10	14.8	0.9
11	7.5	0.2
12	16.4	0.0

- The LC-MS/MS method was used for the analysis of taurine in unknown beagle plasma samples.
- The internal standard was modified from sulfanilic acid to aspartic acid to aid counteract matrix effects observed with beagle plasma not previously observed with human serum.
- Percent accuracy of low, mid, and high level QC samples was between 97.6-103.4% with aspartic acid as internal standard.
- Data for 12 unknown beagle plasma samples is shown in table 2. Excellent precision was observed with standard deviation in 0.0-1.2 μg/mL range in beagle plasma.

Conclusions

- A simple, rapid, and sensitive quantitative MRM method was developed for the analysis of taurine in human serum and beagle plasma.
- The method was linear in 0.025-50 μg/mL range with an r value >0.999 . The method demonstrated excellent accuracy (95.2-98.8%) with matrix matched QC samples.
- The method was successfully applied for the quantification of taurine in 12 unknown beagle plasma samples with minor modifications.

References

1. Yu et. al. J. Sep. Sci. 2016, 36, 3837-3844.
2. Chitranshi, P. and Gilles, C.T. Analysis of taurine in beagle plasma using LCMS-8050. SSI-LCMS-105 Application News.

Authors would like to thank Ms. Kristin Campbell Wisconsin Veterinary Diagnostic Laboratory, Madison WI for generously providing beagle plasma samples.
For Research Use Only. Not for use in diagnostic procedures.