

Multi-stage Mass Spectrometry up to MS⁴ on a QToF System

Valerie Koppen,¹ Russell Mortishire-Smith,² Filip Cuyckens¹

¹Janssen R&D, Beerse, Belgium

²Waters Corporation, Wilmslow, United Kingdom

APPLICATION BENEFITS

This application note illustrates how multi-stage mass spectra up to MS⁴ can be generated on a SYNAPT® G2 HDMS® System in a single experiment. This experiment avoids the need to trigger on, isolate, and fragment specific ions, as is done in directed MSⁿ methods. Moreover, the duty cycle remains constant, differing from the performance of ion-trap instruments, where the duty cycle for each transition decreases as the number of monitored ions/transitions increases.

WATERS SOLUTION

SYNAPT G2 HDMS System

KEY WORDS

T-Wave™ ion mobility, IMS, MS³, MS⁴
structure elucidation

INTRODUCTION

The integration of ion mobility separations (IMS) with high-resolution tandem mass spectrometry offers multiple and distinct capabilities and advantages. First, as an orthogonal dimension of separation, IMS can enhance selectivity in MS and MS^E data sets.¹⁻³ Second, IMS can serve as a powerful tool for structure elucidation, as indicated by its ability to characterize positional isomers of metabolites from generated collision cross-section measurements.⁴⁻⁶ Finally, when combined with tandem-MS, IMS enables multi-stage mass analysis.

Here we describe generating multi-stage mass analyses up to MS⁴ using a Waters® SYNAPT G2 High Definition® MS (HDMS) System. Multi-stage mass analysis is defined in this application note as follows: MS² = 1st generation product ion spectra, MS³ = 2nd generation product ion spectra, and MS⁴ = 3rd generation product ion spectra. Midazolam N-glucuronide was chosen as the test compound and the resultant spectra compared with those obtained on an LTQ-Orbitrap mass spectrometer (Thermo Fisher Scientific).

EXPERIMENTAL

Sample description

Midazolam N-glucuronide, purchased from Toronto Research Chemicals, was dissolved in methanol, to a final concentration of 20 ng/ μ L. The solution was infused directly into the mass spectrometer at 10 μ L/min via the instrument's internal syringe pump.

MS conditions

MS system:	SYNAPT G2 HDMS
Ionization mode:	ESI+
Capillary voltage:	3.0 kV
Cone voltage:	70 V
Extraction cone:	3 V
Trap collision energy (CE):	35 V
Transfer collision energy (CE):	20 V or 40 V
IM gas:	N ₂
Wave velocity:	1100 m/s
Wave height:	40 V
IM gas pressure:	3.11 mbar




Figure 1. Midazolam N-glucuronide and tentative structure of the main product ions.




Figure 2. SYNAPT G2 HDMS and generation of MS³ and MS⁴ pathways.

RESULTS AND DISCUSSION

The geometry of the SYNAPT G2 System appears in Figure 2 together with the mechanisms by which multistage MS spectra can be obtained. Briefly, MS³-like spectra can be generated in two ways. They can be generated by in-source fragmentation, followed by selection of the first-generation product ion in the quadrupole and fragmentation in the trap or transfer collision cell. Alternatively, the spectra can be generated by quadrupole precursor-ion selection, fragmentation in the trap collision cell, separation of the product ions by IMS, and subsequent fragmentation in the transfer region. Combining the latter approach with in-source fragmentation using a higher cone voltage results in MS⁴-like spectra.

Figure 3 illustrates all stages of CID fragmentation performed using the SYNAPT G2 HDMS in ESI positive-ion mode. Figure 3A shows the complete mass spectrum of midazolam N-glucuronide (m/z 502) at a cone voltage of 40 V. Increasing the cone voltage to 70 V induces in-source dissociation of the N-glucuronide, resulting in the formation of the aglycone at m/z 326 (Fig. 3B). Quadrupole selection of this aglycone ion and fragmentation in the trap, or in the transfer, T-Wave collision cells results in an MS^3 spectrum (Fig. 3C). MS^3 like spectra can alternatively be obtained via quadrupole precursor ion selection, fragmentation in the trap collision cell, separation of the product ions by ion mobility separation and subsequent fragmentation in the transfer region (spectra not shown). Combining the latter approach with in-source fragmentation using a higher cone voltage (as applied in this example to produce the aglycone ion) results in MS^4 -like spectra (Figs. 3D and E).

Because the MS^4 product ions are generated after the ion mobility separation, the precursor and product ions can be aligned on the basis of their drift time (Fig. 2). Such alignment allows the extraction of MS^4 -like spectra from the two-dimensional plot of m/z versus drift time or from the ion mobility “drift time” chromatograms (mobilograms).




Figure 3. Multi-stage mass spectrometry up to MS^4 on the SYNAPT G2 HDMS.

Mobilograms of the total ion current and main MS⁴ precursor ions appear in Figure 4. Spectra can be generated from this plot as they are for conventional LC-MS data. Where the precursor ions are not fully separated in IMS, a partial separation and drift-time alignment of the extracted ion chromatograms can suffice to indicate whether product ions are assignable to a specific precursor ion.




Figure 4. Overlay of the ion mobilograms of the total ion current (TIC) and MS⁴ precursor ions 326, 291, and 244.

The mobilograms in Figure 4 demonstrate that at the drift time of the m/z 244 precursor ion (3.86 ms), other ions are co-eluting. Therefore, a drift-time alignment of the extracted ion mobilograms at 3.86 ms and subtraction of the background was applied, producing the MS⁴ spectrum of m/z 244 on the Synapt G2 HDMS. Figure 5 shows corresponding midazolam N-glucuronide spectra obtained on an LTQ-Orbitrap mass spectrometer using standard precursor-ion selection. A comparison of the two data sets (Figs. 3 and 5) shows the spectra to be generally comparable.




Figure 5. Multi-stage mass spectrometry up to MS⁴, as performed on the LTQ-Orbitrap.

CONCLUSIONS

This example demonstrates the application of multi-stage mass spectrometry up to MS⁴ on a SYNAPT G2 HDMS System using quadrupole selection in combination with ion mobility separation.

The interpretation of the MSⁿ spectra obtained via ion mobility separation is assisted by using the drift-time relationship between an ion and its downstream fragments.

The product ion spectra obtained were comparable to those obtained on an LTQ-Orbitrap.

The demonstrated approach generates, in one step, product ion spectra and ion trees in a generic manner using a fast scan speed, which is compatible with modern chromatographic approaches.

References

1. Eckers C, Laures AM, Giles K, Major H, Pringle S. Evaluating the utility of ion mobility separation in combination with high-pressure liquid chromatography/mass spectrometry to facilitate detection of trace impurities in formulated drug products. *Rapid Commun Mass Spectrom.* 2007;21(7):1255–63.
2. Cuyckens F, Dillen L, Cools W, Bockx M, Vereyken L, de Vries R, Mortishire-Smith RJ. Identifying metabolite ions of peptide drugs in the presence of an in vivo matrix background. *Bioanalysis.* 2012 Mar;4(5):595–604.
3. Blech S, Laux R. Resolving the microcosmos of complex samples: UPLC/travelling wave ion mobility separation high resolution mass spectrometry for the analysis of in vivo drug metabolism studies. *Int. J. Ion Mobility.* 2013;16:5–17.
4. Williams JP, Bugarcic T, Habtemariam A, Giles K, Campuzano I, Rodger PM, Sadler PJ. Isomer separation and gas-phase configurations of organoruthenium anticancer complexes: ion mobility mass spectrometry and modeling. *J Am Soc Mass Spectrom.* 2009 Jun;20(6):1119–22.
5. Dear GJ, Munoz-Muriedas J, Beaumont C, Roberts A, Kirk J, Williams JP, Campuzano I. Sites of metabolic substitution: investigating metabolite structures utilising ion mobility and molecular modelling. *Rapid Commun Mass Spectrom.* 2010 Nov 15;24(21):3157–62.
6. Cuyckens F, Wassvik C, Mortishire-Smith RJ, Tresadern G, Campuzano I, Claereboudt J. Product ion mobility as a promising tool for assignment of positional isomers of drug metabolites. *Rapid Commun Mass Spectrom.* 2011 Dec 15;25(23):3497–503.

Waters

THE SCIENCE OF WHAT'S POSSIBLE.®

Waters, SYNAPT, HDMS, High Definition, and The Science of What's Possible are registered trademarks of Waters Corporation. All other trademarks are the property of their respective owners.

©2014 Waters Corporation. Produced in the U.S.A. September 2014 720005100EN AG-PDF

Waters Corporation
34 Maple Street
Milford, MA 01757 U.S.A.
T: 1 508 478 2000
F: 1 508 872 1990
www.waters.com