6.01053.4x0 Metrosep C Supp 2 - XXX/4.0

6.01053.410	Metrosep C Supp 2 - 100/4.0
6.01053.420	Metrosep C Supp 2 - 150/4.0
6.01053.430	Metrosep C Supp 2 - 250/4.0

DE**Säulenmaterial**

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EN

Column material

Polystyrene/divinylbenzene copolymer, with carboxyl groups, particle diameter 5 µm.

Dimensions

6.01053.410	100 x 4.0 mm
6.01053.420	150 x 4.0 mm
6.01053.430	250 x 4.0 mm

pH range

0 - 12 (eluent)

0 - 14 (sample)

Temperature range

10 °C - 60 °C

Recommended standard temperature: 40 °C

Maximum pressure

25.0 MPa (250 bar)

Flow rate

6.01053.410	max. 3.8 mL/min
6.01053.420	max. 3.1 mL/min
6.01053.430	max. 2.0 mL/min

Recommended flow rate: 1.0 mL/min (for all lengths)

Application

Determination of monovalent and divalent cations with suppression and conductivity detection. Optimized separation of sodium and ammonium with sizable differences in concentrations.

Eluent

- Standard eluent
5.0 mmol/L nitric acid and 50 µg/L rubidium (made of RbNO₃ salt).
- Permitted organic eluent additives
0 - 100% acetone and acetonitrile, no alcohol.
- Permitted organic additives for the sample matrix
0 - 100% acetone, acetonitrile and alcohol.

Preparation

Rinse the column with eluent for 3 h.

Guard column

Metrosep C Supp 2 Guard/4.0 (6.01053.500)

Storage

Store the column in standard eluent at ambient temperature.

Regeneration

NOTE

Ensure that the maximum pressure is never exceeded during regeneration.
If the pressure becomes too high, reduce the flow rate.

1. Disconnect the column outlet from the downstream functional units such as suppressor or detector. Instead, collect the flow of liquid in a beaker.
2. Depending on the type of contamination, regenerate the column as follows.
 - a. Organic contamination
Regenerate the column in the direction opposite to the flow at a flow rate of 1 mL/min with the following solutions in succession.
 - 1 h with ultrapure water
 - 1 h with acetonitrile-water mixture (40:60)
 - 1 h with ultrapure water
 - b. Inorganic contamination
Regenerate the column in the direction opposite to the flow at a flow rate of 1 mL/min for 1 h with 50 mmol/L nitric acid.

General notes

- Sample solutions must be microfiltered (0.45 µm).
- The column is not resistant to alcohol additives in the eluent.
- The eluent has only a limited shelf life due to bacterial and fungal growth. We therefore recommend preparing fresh eluent regularly.
- To protect the separation column, we recommend using the pulsation absorber (6.2620.150) to reduce the injector pressure surges.

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