Agilent BioHPLC and AdvanceBio Reversed-Phase Columns User Guide

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This booklet provides general information on Agilent BioHPLC and AdvanceBio reversed-phase columns. For more detailed information about your specific phase or family, see:

www.agilent.com/chem/biocolumnchoices

Getting started

A QC column performance report, including a test chromatogram, is enclosed with every Agilent column. The QC test system has been modified from a standard system to minimize system dead volume, so it may vary from the system used in your lab. This allows a better evaluation of the column efficiency and assures a more consistent product. An optimized LC system will generate similar results to the chromatogram on your QC performance report.

If you have specific questions, contact the Technical Support team at: **www.agilent.com/chem/columnsupport**.

Using Your Column

Installation

- The direction of flow is marked on the column.
- 1.8 μm columns (ZORBAX RRHT, ZORBAX RRHD) can only be operated in the direction flow marked on the column.

For removable, zero-dead volume column connections Agilent recommends the use of the InfinityLab Quick Connect and Quick Turn family. Choices are:

Maximum System Pressure	Recommended Fitting	Part Numbers
Up to 400 bar	InfinityLab Quick Turn fitting (finger-tight)	Fitting: 5067-5966
Up to 800 bar	InfinityLab Quick Turn fitting (with Mounting tool)	Fitting: 5067-5966 Mounting Tool: 5043-0915
Up to 1300 bar	InfinityLab Quick Connect fitting	Fitting: 5067-5965

For more information and part numbers, please see to the Agilent InfinityLab Fitting Brochure (5991-5164EN).

InfinityLab Quick Connect assembly, p/n 5067-5961

InfinityLab Quick Turn fitting, p/n 5067-5966

Learn more at: www.agilent.com/chem/infinitylabfittings

Column conditioning

Every column is tested before shipment. Before the first use, the shipping solvent must be replaced with eluent, taking care that all components are miscible and soluble. If mobile phase additives are used (such as buffers or ion-pair reagents), it is advisable to do an intermediate flush with a mobile phase of the correct composition, but without these additions. Flushing with 10 to 20 column volumes should help in transitioning to your mobile phase. Check that the column has been properly equilibrated before use. This will ensure reproducibility and help prevent retention time drifting. When using formic acid as a mobile phase additive, condition the column as recommended in the Column Conditioning Method Conditions for Formic Acid table.

Column Conditioning Method Conditions for Formic Acid					
Column id (mm)	Mobile Phase	Flow Rate (mL/min)	Column Temp. (°C)	Time (hrs)	After Conditioning
2.1	95/5 H ₂ O/CH ₃ CN + 0.1% formic acid	0.1	60	4	Flush and store in 100% CH ₃ CN
3.0	95/5 H ₂ O/CH ₃ CN + 0.1% formic acid	0.2	60	4	Flush and store in 100% CH ₃ CN
4.6	95/5 H ₂ O/CH ₃ CN + 0.1% formic acid	0.4	60	4	Flush and store in 100% CH ₃ CN

Note: Conditioning overnight, or up to 24 hours, may be beneficial, especially for longer columns or acidic analytes

Important safety considerations

- All connection points in liquid chromatographic systems are potential sources of leaks. Users should be aware of the toxicity or flammability of their mobile phases.
- Because of the small particle size, dry column packings are respirable. Agilent does not recommend removing the column end fittings and exposing the media. Columns should only be opened by trained personnel in a well-ventilated area.
- Please adhere to operating pressure limits noted for each column (see Table 1). Exceeding these limits will compromise chromatographic performance and column lifetime, and could be unsafe.

Table 1. Maximum operating parameters – columns up to 4.6 mm.

Column	Particle Size	Pressure Limit
AdvanceBio Peptide Mapping	2.7 μm	600 bar
AdvanceBio Peptide Plus		
AdvanceBio EC-C18		
AdvanceBio Oligonucleotide		
AdvanceBio Amino Acid Analysis		
AdvanceBio RP-mAb	3.5 μm	400 bar
ZORBAX 300SB	3.5 μm , 5 μm , 7 μm	
Poroshell 300	5 μm	1200 bar
ZORBAX RRHD	1.8 μm	
PLRP-S	3 μm	275 bar
	5 μm , 8 μm , 10 μm	207 bar
	10 to 15 μm , 15 to 20 μm , 30 μm	103 bar

Other operating tips

- Reverse flow will not usually harm the column, but should be avoided except to attempt removal of a clogged frit (see “column care”).
- It is recommended that the flow rate is started at a reduced rate and then gently increased to the desired operating flow rate.
- Always use high purity reagents and chromatography grade solvent to prepare your mobile phase. Degas and filter all mobile phase prior to use.
- Disassembling a column will degrade column performance.
- An inline filter or guard column may be used to protect your column and increase its lifetime.
- If the column is used outside of recommended pH ranges for column phase (see Table 2), a reduced lifetime will result.

- Columns should not be maintained at elevated pH or elevated temperature when not in use (see Table 2).
- New columns may contain a mixture of organic solvents and water, which may contain buffer salts. See Table 3 for details of the shipping solvents. Initially, care should be taken not to pass any mobile phase through the column that may cause a precipitate to form or may not be fully miscible.
- Agilent BioHPLC and AdvanceBio columns are compatible with all common eluents used for analysis of biomolecules.

Table 2. Column operating parameters – pH and temperature.

Column	Recommended pH Range	Maximum Operating Temperature
AdvanceBio Peptide Mapping	2.0 to 8.0	60 °C
AdvanceBio Peptide Plus		
AdvanceBio EC-C18		
ZORBAX RRHD	1.0 to 8.0	80 °C
ZORBAX 300SB-C8, C3, CN		
ZORBAX 300SB-C18	1.0 to 8.0	90 °C
AdvanceBio RP-mAb		
Poroshell 300		
Poroshell 300 Extend-C18	2.0 to 11.0	60 °C below pH 8, 40 °C above pH 8
AdvanceBio Oligonucleotide	3.0 to 11.0	65 °C
AdvanceBio Amino Acid Analysis		
PLRP-S	1.0 to 14.0	200 °C

Note: All silica-based packings have some solubility in pH >6 aqueous mobile phases. When using silicabased columns at pH >6, best column lifetime is obtained at lower temperatures (40 °C max) using low buffer concentrations in the range of 0.01 to 0.02 M.

Table 3. Shipping solvents.

Column	Shipping Solvent	Compatibility
AdvanceBio Peptide Plus	100% acetonitrile	Water and all common organic solvents. Avoid tetrahydrofuran (THF).
AdvanceBio Peptide Mapping	Mixture of acetonitrile and water	
AdvanceBio EC-C18		
ZORBAX RRHD		
ZORBAX 300SB		
AdvanceBio RP-mAb	Mixture of acetonitrile and water	Water and all organic solvents, including N,N-dimethylformamide and dimethyl sulfoxide
Poroshell 300 and 300 Extend	Mixture of methanol and water	
AdvanceBio Oligonucleotide	Mixture of acetonitrile and water	Water and all common organic solvents. Modifiers including HFIP and TEAA
AdvanceBio Amino Acid Analysis		
PLRP-S	7:1 acetonitrile/water	Aqueous organic solvents including N,N-dimethylformamide and dimethyl sulfoxide. 100% aqueous is not recommended as it will reduce column performance and lifetime

Mobile phase selection and operating temperatures

The bonded stationary phase is nonpolar in nature and is best used with polar mobile phases, such as methanol/water or acetonitrile/water mixtures. Increasing the amount of organic component reduces the retention time of the sample.

When the maximum temperature is used for prolonged periods this will reduce column lifetime.

Using 100% aqueous eluents with PLRP-S columns will significantly reduce the column lifetime and may result in a rapid deterioration in peak width and symmetry.

Recommended starting gradients

Separations of biological molecules typically use gradient conditions – increasing the amount of the organic component to achieve elution from the column. Separations of peptides, polypeptides, and proteins, are most commonly carried out using acidic eluents with trifluoroacetic acid (TFA) or formic acid (FA) used as the modifier for pH control and/or as an ion-pairing additive to achieve the desired retention and selectivity. Organics, such as acetonitrile, methanol, and ethanol, are used for elution, with acetonitrile being the most commonly used.

Additional information on peptide separations can be found in chapter 11, pages 497–508, *Introduction to Modern Liquid Chromatography, Second Edition*. L.R. Snyder and J. J. Kirkland, (John Wiley & Sons, 1979).

Column care

Columns with a particle size of 1.8 μm have a 0.5 μm inlet frit. Particulates will block the column inlet frits and should be removed before the sample is analyzed. Where this is not possible, an inline filter or guard column should be used to protect the analytical column, and increase the column lifetime. The AdvanceBio Peptide Mapping column has a 2 μm inlet frit. Samples that contain bio particulate matter larger than 2 μm will plug the inlet frit, so guard columns are recommended for use with these samples. It is best to filter your samples before injecting them onto any column. You should also check that all eluents are freshly prepared and filtered before use.

See www.agilent.com/chem/guards for more information about guard columns.

Cleaning your column/extending column life

For columns that can be backflushed (particles $>1.8\ \mu\text{m}$), clean the column in the reverse direction. Start with a stronger (less polar) solvent.

1. Disconnect column from detector and run wash solvents into a beaker.
2. Start with your mobile phase without buffer salts (water/organic). Run 10 to 20 column volumes through.
3. Next, use 100% organic (methanol or acetonitrile).
4. Check pressure to see if it has returned to normal. If not, then:
5. Discard column or consider stronger conditions, such as 75% acetonitrile: 25% isopropanol.
6. Increase to 100% isopropanol, 100% methylene chloride, or 100% hexane (if you use methylene chloride or hexane, you will need to flush the column with isopropanol prior to use and before returning to your mobile phase as it is not miscible with aqueous).

Aggressive cleanup cycles using 80% 1 M NaOH, or 1 M HCl with 20% organic, can be used with the PLRP-S columns.

For columns with $1.8\ \mu\text{m}$ particles, do not backflush the column – you can attempt the cleanup procedure as detailed above but maintain direction of flow.

Storage recommendations

In general, columns may be safely stored for short periods in most mobile phases. Long-term storage of silica-based, bonded phase columns should be in a pure organic solvent. If the column has previously been used with a buffered mobile phase, the buffer should first be removed by purging the column with 20 to 30 column volumes of a 50:50 mixture of methanol or acetonitrile and water, followed by 20 to 30 column volumes of the pure solvent. Before storing, end-fittings should be tightly capped with end-plugs to prevent packing from drying out.

To protect equipment, it is desirable to remove salts from the instrument and column by purging the column with the same mobile phase without the buffer (e.g. using 60:40 ACN/H₂O to remove a 60:40 ACN/0.02 M phosphate buffered mobile phase). Re-equilibration is rapid with the original mobile phase when using this approach and any possibility of corrosion from the salts is eliminated.

Tips for the best chromatographic results

- Optimize your instrument by minimizing tubing lengths between components. This will reduce extracolumn volume and band broadening. Use 0.12 mm id red tubing for Fast LC/high efficiency columns. Learn about capillary options at www.agilent.com/chem/lccapillaries.
- Ensure that the data collection rate is optimized for your column. Use a higher collection rate for Fast LC columns (ZORBAX RRHD).
- Use sample filtration or other sample preparation techniques appropriate for your sample. Learn more at www.agilent.com/chem/sampleprep.
- Use certified lamps in your LC instruments for best performance.

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
Agilent AdvanceBio columns

For faster, more consistent biopharmaceutical analysis

Agilent AdvanceBio is a family of state-of-the-art biocolumns designed to deliver consistent, exceptional performance for the separation and characterization of peptides and proteins. The science behind AdvanceBio columns helps advance accuracy, provide greater productivity, and eliminate interferences that can impede progress. AdvanceBio columns are rigorously tested by Agilent to ensure great results and are backed by Agilent's 60-day full satisfaction warranty.

Find the right column for your separation

Try the LC Column and Sample Prep Navigator

Get started now at www.agilent.com/chem/navigator.

Agilent ordering information

For more information on our products and services,
visit our web site at **www.agilent.com**

For technical support or to place an order,
visit **www.agilent.com/chem/columnsupport**

To place an order,
visit **www.agilent.com/en-us/contact-us/wheretobuy**

Agilent offers a complete line of sample preparation products
to support LC and LC/MS applications.

To find parts and supplies to support high-performance
bioseparation, such as stainless steel-coated PEEK capillaries
for high-pressure bio-inertness and robustness, visit
www.agilent.com/chem/lcsupplies

To receive a copy of the Infinity Supplies catalog, and other
key resources, visit **www.agilent.com/chem/getguides**

DE.3423263889

This information is subject to change without notice.

Agilent Technologies, Inc. 2018, 2020

Printed in USA, August 25, 2020

820000-997