



B-Adrenoceptor-blocking drugs

Analysis of esmolol and flumolol in whole blood

Application Note

Clinical Research

Authors

Agilent Technologies, Inc.

Introduction

According to K. Kylberg-Hanssen et al, sodium dodecyl sulfate (SDS) is an effective esterase inhibitor in the determination of labile β -blockers: esmolol and flumolol. The limit of detection (coefficient of variation: 10 - 15%) is 5 nmol/L of whole blood.



Agilent Technologies

Conditions


Technique : GC-capillary
Column : Agilent CP-Sil 8 CB, 0.32 mm x 25 m fused silica
WCOT CP-Sil 8 CB (0.12 μ m) (Part no. CP7741)
Temperature : esmolol determination:
100 °C (1 min) \rightarrow 185 °C 10 °C/min;
185 °C (15 min) \rightarrow 250 °C, 30 °C/min;
250 °C (7 min)
flumolol determination:
100 °C (1 min) \rightarrow 200 °C, 15 °C/min;
200 °C (10 min) \rightarrow 250 °C, 30 °C/min;
250 °C (7 min)
Carrier Gas : He, 60 kPa
Injector : Splitter, 20 mL/min
T = 270 °C
Detector : ECD
T = 350 °C

Courtesy : Reprinted with permission from Clinical Chemistry
Vol. 31, No. 6 1985. Copyright American Association
for Clinical Chemistry, Inc.


Peak identification

1. flumolol (100 nmol/L)
2. Int. Standard (H 163/37) (500 nmol/L)

Blood sample



Blood sample with added flumolol



Analytical procedure

Collect 5 mL of blood in a tube containing 0.5 mL SDS solution (200 g/L). After mixing transfer 1.0 g of the sample to a 15 mL centrifuge tube. Add 100 μ L of the internal standard solution (5 μ mol in phosphate buffer, pH = 6.0, 0.1 mol/L). Add 4 mL of toluene and 0.5 mL of phosphate buffer (pH = 11.8, 0.2 mol/L). Shake (15 min) and centrifuge. Immerse the tube in an acetone/solid CO₂ mixture to freeze the aqueous layer.


Pour the organic layer into a 5 mL centrifuge tube and evaporate under a stream of N₂ at 35 °C. Redissolve the residue in 100 μ L of ethylacetate and 25 μ L of pentafluoropropionic anhydride. After reacting for 20 min at room temperature, evaporate under nitrogen at 35 °C. Dissolve the residue in toluene (500 μ L) and inject 1.5 μ L.

Acknowledgement


Agilent Technologies thanks Kerstin Kylberg-Hanssen, Department of Bioanalytical Chemistry, AB Hässle, S-431 83 Mölndal, Sweden for all technical information and the chromatogram.

Peak identification


1. esmolol (300 nmol/L)
2. Int. Standard (H 163/37) (500 nmol/L)



Blood sample



Blood sample with added esmolol



www.agilent.com/chem

For Research Use Only. Not for use in diagnostic procedures.
This information is subject to change without notice.

© Agilent Technologies, Inc. 2015

Printed in the USA

Updated 18 May, 2015

First published prior to 11 May, 2010

A00320



Agilent Technologies