



Hydrocarbons, C₁ – C₄

Application Note

Energy & Fuels

Authors

Agilent Technologies, Inc.

Introduction

Retention times of C₁ - C₄ hydrocarbons on the Agilent CP-SilicaPLOT column are not influenced by polar impurities such as water, sulfur compounds or carbon dioxide. Only the Al₂O₃ PLOT column provides a higher selectivity for C₄-isomers.




Agilent Technologies

Conditions

Technique : GC-capillary
Column : Agilent CP-SilicaPLOT, 0.32 mm x 30 m, fused silica
PLOT (df = 4 µm) (Part no. CP8567)
Temperature : 60 °C
Carrier Gas : He, 150 kPa (1.5 bar, 21 psi)
Injector : Split,
T = 250 °C
Detector : FID
T = 250 °C
Concentration range : 1% in N₂

Peak identification

1. methane
2. ethane
3. ethylene
4. acetylene
5. propane
6. cyclopropane
7. propylene
8. isobutane
9. butane
10. 1-butene
11. propyne (methylacetylene)
12. 1,3-butadiene
13. trans-2-butene
14. isobutene
15. cis-2-butene



www.agilent.com/chem

This information is subject to change without notice.

© Agilent Technologies, Inc. 2011

Printed in the USA

31 October, 2011

First published prior to 11 May, 2010

A01294



Agilent Technologies