



Solvents

Application Note

Materials Testing & Research

Authors

Agilent Technologies, Inc.

Introduction

Gas chromatography with the Agilent PoraBOND Q column separates four chlorinated solvents in under 13 minutes.



Agilent Technologies


Conditions

Technique : GC
Column : Agilent PoraBOND Q, 0.53 mm x 25 m fused silica
(df = 10 μ L) (Part no. CP7354)
Temperature : 100 °C (2 min) \rightarrow 120°C, (0 min) 5 °C/min \rightarrow
240 °C, 20 °C/min, (2 min)
Carrier Gas : Helium, 65 kPa, 10 psi
Injector : Purge & Trap procedure
Detector : FID
Concentration Range : ca. 100 ppm level

Courtesy : David Black, Department of Civil Engineering,
North Carolina State University, USA

Peak identification

1. trans-1,2 dichloroethylene
2. cis-1,2-dichloroethylene
3. trichloroethylene
4. tetrachloroethylene



www.agilent.com/chem

This information is subject to change without notice.

© Agilent Technologies, Inc. 2011

Printed in the USA

31 October, 2011

First published prior to 11 May, 2010

A01905



Agilent Technologies