

Application Bulletin


Of interest to: General analytical laboratories, Metallurgical laboratories

A 1, 10

Determination of zinc by potentiometric titration with potassium hexacyanoferrate(II)

Summary

Zinc, such as that occurring as a constituent of light alloys, can be determined by precipitation titration with potentiometric end-point indication. The determination of zinc next to cadmium is also possible.



Apparatus and accessories

- Titrimo or Titrimo with Dosino or Dosimat
- Magnetic Swing-out Stirrer
- Exchange unit
- Pt-Titrode 6.0431.100 with electrode cable 6.2104.020
- Thermostated titration vessel, including titration vessel lid 6.1414.010 and titration vessel with thermostatic jacket 6.1418.250

Reagents

- Titrant, $\text{K}_4[\text{Fe}(\text{CN})_6]$ $c = 0.05 \text{ mol/L}$:
Dissolve 21.12 g $\text{K}_4[\text{Fe}(\text{CN})_6] \cdot 3 \text{H}_2\text{O}$ in dist. H_2O and fill up to 1 L.
- Sulfuric acid: H_2SO_4 conc. p.A.
- Zn standard:
Dissolve 22.0 g $\text{ZnSO}_4 \cdot 7 \text{H}_2\text{O}$ in dist. H_2O , add 1 mL conc. H_2SO_4 and fill up to 1 L with dist. H_2O .
1 mL = 5.0 mg Zn

Analysis

Pipet 5 mL sample, resp. standard solution, and approx. 50 mL dist. H_2O into the titration vessel. Add carefully 5 mL H_2SO_4 and heat the mixture up to 65 ... 70 °C, then titrate with $\text{K}_4[\text{Fe}(\text{CN})_6]$ in the MET mode of the titrator (parameters, see under appendix).

Calculation

1 mL $K_4[Fe(CN)_6]$ $c = 0.05$ mol/L = 4.904 mg Zn

g/L Zn = EP1 * C01 / C00

C00 = Sample weight in mL

C01 = 4.904

Remarks

- Do not allow temperature to fall below 60 °C during titration.
- The titrant must be standardized with Zn-standard because the stoichiometric factor does not quite agree.
- Cadmium also reacts with $K_4[Fe(CN)_6]$. Precipitation, however, is slow and takes place after that of Zn. It must be indicated potentiometrically. Thus, it is possible to determine Zn next to Cd.


Literature

- Vojlosnikova, M.T. / Kozlovkij, M.T. / Songina, A.O.
Eine amperometrische Methode zur Bestimmung sehr kleiner Zinkmengen
Zavodsk. Laborat. 23, (1957) 273-276 (Russ)
Ref: Fresenius, Z.Anal.Chem. 159, (1957) 59
- Agasjan, P.K.
Zur Bestimmung von Zink und Cadmium nebeneinander
Zavodsk. Laborat. 24, (1958) 532-534 (Russ)
Ref: Fresenius, Z.Anal.Chem. 167, (1959) 281
- Kao, S.S. / Chuang, W.T.
Simultantitration von Zink und Cadmium mit Hexacyanoferrat(II) - Lösung nach der Dead-Stop-Methode
Acta Chim. Sinica 24, (1958) 25-29 (Chin)
Ref: Fresenius, Z.Anal.Chem. 167, (1959) 282
- Malur, J. / Treptow, H.
Untersuchung über die Grenzen der potentiometrischen Simultanbestimmung von Zink und Cadmium
Fresenius, Z.Anal.Chem. 183, (1961) 426 - 432

Figures

```
'pa
736 GP Titrimo      04268 736.0011
date 99-09-06      time 10:00  2
MET U              AB11CdZn
parameters
>titration parameters
V step             0.05 ml
titr.rate          max. ml/min
signal drift       50 mV/min
equilibr.time      26 s
start V:           OFF
pause              10 s
dos.element:      internal D0
meas.input:        1
temperature        63.1 °C
>stop conditions
stop V:            abs.
stop V             20 ml
stop U             OFF mV
stop EP            9
filling rate       max. ml/min
>statistics
status:            OFF
>evaluation
EPC                50 mV
EP recognition:    all
fix EP1 at U      OFF mV
pK/HNP:           OFF
>preselections
req.ident:         OFF
req.smpl size:     all
activate pulse:    OFF
```


Fig. 1 Parameter report Titrimo



```
'fm
736 GP Titrimo      04268 736.0011
date 99-09-03      time 14:09  6
MET U              AB11
>calculations
Zn=EP1*C01/C00;2;g/l
C00=                2.5
C01=                4.904
=====
```

```
'fr
736 GP Titrimo      04268 736.0011
date 99-09-03      time 14:09  6
U(init)            8 mV MET U  AB11
smpl size          2.5 ml
EP1                2.614 ml    -57 mV
Zn                 5.13 g/l
stop V reached
=====
```

Fig. 2 Titration curve and results, Zn determination



```
'fm
736 GP Titrino      04268 736.0011
date 99-09-06     time 10:01  2
MET U             AB11CdZn
>calculations
Zn=EP1*C01/C00;2;g/l
Cd=(EP2-EP1)*C02/C00;2;g/l
C00=              2
C01=              5.62
C02=              4.904
```

```
-----
'fr
736 GP Titrino      04268 736.0011
date 99-09-06     time 10:00  2
U(init)          -29 mV MET U AB11CdZn
sml size         2 ml
EP1              2.041 ml      -40 mV
Zn               5.74 g/l
Cd               missing EP
```

Fig. 3 Titration curve and results, zinc next to cadmium