

Application Bulletin


Of interest to: Metal industry

A 10

Determination of nickel in gold and silver plating baths by potentiometric titration

Summary

A potentiometric method for the determination of nickel in gold and silver electroplating baths is described. The titration is done with KCN. Gold and silver are removed before titration by a reduction process. It is also possible to determine nickel in alloys, etc. (see the literature reference).



Apparatus and accessories

- Titrino or Titrande with Dosino or Dosimat
- Magnetic Swing-out Stirrer
- Exchange unit
- Ag Titrode with Ag₂S coating 6.0430.100 with electrode cable 6.2104.020

Reagents

- Titrant, $c(\text{KCN}) = 0.2 \text{ mol/L}$:
Dissolve 4 g KOH and 13.03 KCN in dist. H₂O and fill up to 1 L. Standardize by titrating against Ni standard.
- Ni standard, $c(\text{Ni}) = 0.01 \text{ mol/L}$:
Dissolve 2.809 g NiSO₄ * 7 H₂O in dist. H₂O and fill up to 1 L.
1 mL = 0.5869 mg Ni
- Hydrazine solution: $w(\text{hydrazine sulfate}) = 2\%$ in dist. H₂O
- Ammonia: $w(\text{NH}_3) = 25\%$
- Digestion acids: conc. HCl and conc. HNO₃.

Sample preparation

Heat 10.0 mL of the bath sample together with 10 mL HCl. Evaporate down to $\frac{1}{4}$ of the original volume. Work in a fume cupboard, since cyanide is given off!!! Add conc. HNO₃ drop by drop until the precipitate is redissolved and then evaporate down to $\frac{1}{4}$ of the volume again. After cooling, rinse with dist. H₂O into a 50 mL graduated flask, fill up to the mark and mix.

Analysis

Pipet out 10.0 mL of the sample solution as prepared above, (corresponding to 2 mL original bath) into a beaker, add 10 mL dist. H₂O, 10 mL ammonia and 5 mL hydrazine solution. Boil for a short time and (use a fume cupboard) and after cooling, titrate with c(KCN) = 0.2 mol/L.

Calculations

1 mL c(KCN) = 0.2 mol/L = 2.935 mg Ni

g/L Ni = EP1 * C01 * C02 / C00

C00 = Sample size in mL (2)

C01 = 2.935

C02 = Titer of the KCN solution

Literature


- Trepka-Bloch, E.
Potentiometrische Bestimmung von Nickel neben Ag, As, Bi, Co, Cu und Fe
Chemist Analyst **43**, (1954) 63-65
Ref: Fresenius, Z. Anal. Chem. **147**, (1955) 143
- Luke, C.L.
New rapid method for the determination of nickel in ferrous and ferromagnetic metals
Anal. Chem. **33**, (1961) 96-98

Figures

```
'pa
785 DMP Titrimo      02287  785.0010
□R□date 1999-07-15   time 08:54    2
MET U                *****
parameters
>titration parameters
  V step              0.10 ml
  dos.rate            max. ml/min
  signal drift        50 mV/min
  equilibr.time       26 s
  start V:            OFF
  pause               0 s
  meas.input:         1
  temperature         25.0 C
>stop conditions
  stop V:             abs.
  stop V              5 ml
  stop U              OFF mV
  stop EP             9
  filling rate        max. ml/min
>statistics
  status:             OFF
>evaluation
  EPC                 30 mV
  EP recognition:     all
  fix EP1 at U       OFF mV
  pK/HNP:            OFF
>preselections
  req.ident:          OFF
  req.smpl size:      OFF
  limit smpl size:   OFF
  activate pulse:     OFF
-----
```

Fig. 1 Parameters Titrimo

```
'fr
785 DMP Titrino      02287
785.0010
R□date 1999-07-15   time 08:22
1
card label:Appl.751
U(init)      94 mV MET U
*****
smpl size    10 ml
EP1          1.900 ml      262
mV
Ni           0.56 g/l
time        188 s
stop V reached
=====
```



```
'cu
785 DMP Titrino      02287
785.0010
R□date 1999-07-15   time 08:22
1
start V      0.000 ml MET U
*****
2.0 ml/div   dU=100.0 mV/div
```

```
'BMP-File:
C:\METROHM\DATASCAN\ASWNH4.bmp
□2      =====
```

Fig. 2 Results and titration curve