

# Application Bulletin

Of interest to: Electroplating industries, Metals

A 10

## Potentiometric analysis of tin plating baths

### Summary

Potentiometric titration methods for the analysis of acid and alkaline tin plating baths are presented. Following methods are described: tin(II) / tin(IV) / total tin, free fluoroboric acid or free sulfuric acid, chloride in acidic tin baths, free hydroxide and carbonate in alkaline tin baths.

### Apparatus and accessories

- Titrino or Titrande with Dosino or Dosimat
- Magnetic swing-out stirrer
- Exchange unit(s)
- Pt Titrode 6.0431.100 with electrode cable 6.2104.020
- Comb. pH glass electrode 6.0255.100
- Ag Titrode with Ag<sub>2</sub>S coating 6.0430.100

### Reagents

These are described under the individual analyses.

### 1. Iodometric tin determination

To increase the accuracy of this analysis, 10.0 mL bath sample is pipetted into a 100 mL graduated flask, filled up to the mark with distilled water and mixed well.

#### Reagents:

- $c(\text{Iodine solution}) = 0.05 \text{ mol/L}$
- $w(\text{HCl}) = 36\%$
- Iron powder p.a.

#### 1.1. Tin(II)

Add 15 mL HCl and 50 mL dist. H<sub>2</sub>O to 10.0 mL diluted sample (corresponding to 1 mL original bath) in a beaker and titrate with  $c(\text{iodine solution}) = 0.05 \text{ mol/L}$  against the Pt Titrode.

**1.2. Tin(IV) and total tin**

In compliance with the tin content, pipet between 10.0 - 50.0 mL of the diluted sample (1 - 5 mL original bath) into a wide-necked Erlenmeyer flask and add 50 mL HCl. Stirring well, add ca. 1 g iron powder in small portions and when reaction subsides, warm up until the iron powder is entirely dissolved. Cool immediately and titrate with c(iodine solution) = 0.05 mol/L against the Pt Titrode.

**Calculations:**

1 mL c(Iodine solution) = 0.05 mol/L = 5.9345 mg Sn

g/L Sn = EP1 x C01 / C00

C00 = sample size in mL original bath


C01 = 5.9345

**Figures:**

<pre>'pa 751 GPD Titrimo      05268  751.0011 date 2000-06-05    time 10:26    4 MET U              AB90 Sn. parameters &gt;titration parameters   V step           0.10 ml   dos.rate         max. ml/min   signal drift     30 mV/min   equilibr.time    32 s   start V:         OFF   pause            0 s   dos.element:     internal D0   meas.input:      1   temperature      25.0 °C</pre>	<pre>&gt;stop conditions   stop V:          abs.   stop V           6 ml   stop U           OFF mV   stop EP          9   filling rate     max. ml/min &gt;statistics   status:          OFF &gt;evaluation   EPC              30 mV   EP recognition:  greatest   fix EP1 at U    OFF mV   pK/HNP:         OFF &gt;preselections   req.ident:       OFF   req.smpl size:  OFF   activate pulse:  OFF   -----</pre>
--	---

**Fig. 1** Parameter report Titrimo, iodometric tin determination

```
'fr
751 GPD Titrimo      05268  751.0011
date 2000-06-05    time 10:26    4
U(init)            391 mV MET U  AB90 Sn.
smpl size          1.0 ml id#1    Sn(II)
EP1                3.030 ml      290 mV
Sn                 17.98 g/l
stop V reached
```



**Fig. 2** Titration curve iodometric tin determination

## 2. Free fluoroboric acid or free sulfuric acid

### Reagents:

- $c(\text{NaOH}) = 1 \text{ mol/L}$
- Sodium sulfate p.a.

### Analysis:

Dilute 10.0 mL bath sample to approx. 50 mL with dist.  $\text{H}_2\text{O}$  in a beaker. While stirring dissolve approx. 5 g  $\text{Na}_2\text{SO}_4$  and titrate afterwards with  $c(\text{NaOH}) = 1 \text{ mol/L}$  against the comb. pH glass electrode. The flat potential jump at  $\text{pH} = \text{approx. } 3.2$  is evaluated.

### Calculations:

1 mL  $c(\text{NaOH}) = 1 \text{ mol/L} = 87.81 \text{ mg HBF}_4$  or  $49.037 \text{ mg H}_2\text{SO}_4$

$\text{g/L HBF}_4 = \text{EP1} \times \text{C01} / \text{C00}$

$\text{g/L H}_2\text{SO}_4 = \text{EP1} \times \text{C02} / \text{C00}$


$\text{C00} = \text{Sample size in mL original sample (10)}$

$\text{C01} = 87.81$

$\text{C02} = 49.037$

### Figures:

```
'fr
751 GPD Titrimo           05268   751.0011
date 2000-06-06         time 08:23     3
pH(init)           1.82   DET pH   AB90 SO4
smpL size           5.0 ml
EP1+                9.932 ml           3.33
H2SO4               97.41 g/l
#EP's not corresponding
stop V reached
```



**Fig. 3** Titration curve free sulfuric acid

## 3. Chloride determination in acidic tin baths

### Reagents:

- $c(\text{AgNO}_3) = 0.1 \text{ mol/L}$
- $w(\text{HNO}_3) = 65\%$

**Analysis:**

Pipet 5.0 mL bath solution into a beaker and dilute to approx. 50 mL with dist. H<sub>2</sub>O. Add 2 mL HNO<sub>3</sub> and titrate with c(AgNO<sub>3</sub>) = 0.1 mol/L against the Ag Titrode (Ag<sub>2</sub>S-coating).

**Calculations:**

1 mL c(AgNO<sub>3</sub>) = 0.1 mol/L = 3.5453 mg Chloride


g/L Chloride = EP1 x C01 / C00

C00 = sample size in mL original bath (5)

C01 = 3.5453

**Figures:**

```
'fr
751 GPD Titrino           05268   751.0011
date 2000-06-06         time 09:02     5
U(init)                236 mV DET U   AB90 Cl
EP1                    3.473 ml          97 mV
Chlorid                12.31 g/l
stop V reached
=====
```



**Fig. 4** Titration curve chloride determination

**4. Free hydroxide and carbonate in alkaline baths**

**Reagents:**

- c(HCl) = 1 mol/L
- w(BaCl<sub>2</sub>) = 25%

**Analysis:**

Add 50 mL BaCl<sub>2</sub> to 10.0 mL bath sample in a wide-necked Erlenmeyer flask and boil for a short time. Allow to cool and slowly titrate the still warm solution with c(HCl) = 1 mol/L against the comb. pH glass electrode.

**Calculations:**

Two endpoints are obtained. The consumption up to EP1 corresponds to NaOH, between EP1 and EP2 to tin and between EP2 and EP3 to carbonate.

1 mL c(HCl) = 1 mol/L = 40.0 mg NaOH or 106.0 mg Na<sub>2</sub>CO<sub>3</sub>

g/L NaOH = EP1 x C01 / C00

g/L Na<sub>2</sub>CO<sub>3</sub> = (EP3 - EP2) x C02 / C00

C00 = Sample size in mL original sample (10)

C01 = 40


C02 = 106

**Figures:**

<pre>'pa 751 GPD Titrimo      05268  751.0011 date 2000-06-06    time 10:43    6 DET U              AB90 OH- parameters &gt;titration parameters   meas.pt.density      4   min.incr.           10.0 µl   dos.rate             max. ml/min   signal drift         25 mV/min   equilibr.time        34 s   start V:             OFF   pause                0 s   dos.element:         internal D0   meas.input:          1   temperature          25.0 °C</pre>	<pre>&gt;stop conditions   stop V:              abs.   stop V               20 ml   stop U               OFF mV   stop EP              9   filling rate         max. ml/min &gt;statistics   status:              OFF &gt;evaluation   EPC                  5   EP recognition:      all   fix EP1 at U        OFF mV   pK/HNP:              OFF &gt;preselections   req.ident:           OFF   req.smpl size:       OFF   activate pulse:      OFF</pre>
--	---

**Fig. 5** Parameter report Titrimo, free NaOH and carbonate

```
'fr
751 GPD Titrimo      05268  751.0011
date 2000-06-06    time 10:43    6
U(init)            -277 mV DET U    AB90 OH-
smpl size          10.0 ml
EP1                6.911 ml          -77 mV
EP2                11.657 ml         196 mV
EP3                14.329 ml         286 mV
NaOH                27.64 g/l
Na2CO3             28.32 g/l
stop V reached
```



**Fig. 6** Titration curve NaOH / (Sn) / Na<sub>2</sub>CO<sub>3</sub>

**Literature**

- Metrohm Ti Application Note No. T-5, T-21, T-23
- Wild, P.W.  
Moderne Analysen für die Galvanik  
Eugen G. Leuze Verlag, D-88348 Saulgau/Württ. 1972
- Jelinek, T.W.  
Prozessbegleitende Analytik in der Galvanotechnik  
Eugen G. Leuze Verlag, D-88348 Saulgau/Württ. 1999  
ISBN 3-87-480-135-7