

Potentiometric analysis of acidic lead and lead/tin plating baths

Of interest to:
Metals, electroplating
A 10

Summary

This Bulletin describes the potentiometric determination of lead, tin(II) and free fluoroboric acid.

Instruments and accessories

- Titrino or Titrando with Dosino or Dosimat
- Magnetic Stirrer
- Exchange Units (for NaOH possibly with 6.1608.040 PE reagent bottle)
- 6.0502.140 ion-selective copper electrode (Cu ISE)
- 6.0726.107 double-junction Ag/AgCl reference electrode [filled with $c(\text{KCl}) = 3 \text{ mol/L}$]
- 6.0255.100 Profitrode, combined double-junction pH glass electrode with separable ground-joint diaphragm
- 6.0431.100 Pt Titrode
- 6.2104.020 and 6.2106.020 electrode cables

1. Determination of lead

Reagents

- $c(\text{Na}_2\text{EDTA}) = 0.1 \text{ mol/L}$
- Buffer solution pH = 10; Dissolve 54 g NH_4Cl and 350 mL $w(\text{NH}_3) = 25\%$ in dist. water and fill up to 1 L.
- $\text{Cu}(\text{NH}_4)_2\text{EDTA}$, $c = 0.1 \text{ mol/L}$ (Merck no. 105217)

Analysis

Pipet a portion of the bath sample (should contain ca. 100 mg Pb) into a glass beaker and fill up to ca. 50 mL with dist. water. Add 5 mL each $\text{Cu}(\text{NH}_4)_2\text{EDTA}$ and buffer solution pH = 10 and, while stirring, allow to react for 30 s. Afterwards, titrate in the MET mode of the titrator with $c(\text{Na}_2\text{EDTA}) = 0.1 \text{ mol/L}$ against the Cu ISE. Work with a starting volume of 50% of the expected titrant consumption. The parameter «pause» is set at 30 s.

Calculation

$$1 \text{ mL } c(\text{Na}_2\text{EDTA}) = 0.1 \text{ mol/L} = 20.72 \text{ mg Pb}$$
$$\text{g/L Pb} = \text{EP1} * \text{C01} / \text{C00}$$
$$\text{C00} = \text{sample size in mL}$$
$$\text{C01} = 20.72$$

2. Determination of tin(II)

Reagents

- Iodine solution, $c(\text{I}_2) = 0.05 \text{ mol/L}$ (0.1 N); Dissolve 25 g potassium iodide in 40 mL dist. water. Add 12.7 g iodine and while shaking, bring it into solution. Finally fill up to 1 L with dist. water. Titer determination e.g. against As(III).
- $w(\text{H}_2\text{SO}_4) \approx 30\%$

Analysis

Mix approx. 50 mL dist. water and 20 mL H_2SO_4 in a beaker and purge with nitrogen. After addition of 2.0 mL bath sample, titrate immediately (under nitrogen) with $c(\text{I}_2) = 0.05 \text{ mol/L}$ (Pt Titrode).

Calculation

$$1 \text{ mL } c(\text{I}_2) = 0.05 \text{ mol/L} = 5.9345 \text{ mg Sn(II)}$$
$$\text{g/L Sn(II)} = \text{EP1} * \text{C01} / \text{C00}$$
$$\text{C00} = 2.0 \text{ (sample size in mL)}$$
$$\text{C01} = 5.9345$$

3. Determination of free fluoroboric acid

Reagents

- Sodium hydroxide solution, $c(\text{NaOH}) = 1 \text{ mol/L}$

Analysis

Dilute 10.0 mL bath sample with approx. 50 mL dist. water in a beaker and titrate (not too fast) with $c(\text{NaOH}) = 1 \text{ mol/L}$. The flat potential jump at $\text{pH} \approx 3.2$ is used for evaluation (combined pH glass electrode).

Calculation

$1 \text{ mL } c(\text{NaOH}) = 1 \text{ mol/L} = 87.81 \text{ mg HBF}_4$
 $\text{g/L free HBF}_4 = \text{EP1} * \text{C01} / \text{C00}$

$\text{C00} = 10.0$ (sample size in mL)

$\text{C01} = 87.81$

Literature

- Metrohm Application Bulletin No. 101
Complexometric titrations with the Cu-ISE
Metrohm Ltd., Herisau
- Metrohm Application Note T-21
Sn(II) and sulfuric acid in tin plating bath
Metrohm Ltd., Herisau
- Metrohm Application Note T-24
Metal contents of alkaline plating baths for cadmium, copper, lead or zinc
Metrohm Ltd., Herisau
- P. W. Wild
Moderne Analysen für die Galvanotechnik
Eugen G. Leuze Verlag, Saulgau, 1972
- T. W. Jelinek
Prozessbegleitende Analytik in der Galvanotechnik
Eugen G. Leuze Verlag, Saulgau, 1999
ISBN 3-87-480-135-7

Figures




Fig. 1: Parameter report Titrino for the Pb determination.

Fig. 2: Titration curve Pb.




Fig. 3: Titration curve Sn(II).

```
'pa
736 GP Titrino          04268  736.0011
date 15-01-18          time 09:53    10
DET pH                  HBF4
parameters
>titration parameters
  meas.pt.density      4
  min.incr.            10.0 µl
  titr.rate            max. ml/min
  signal drift         50 mV/min
  equilibr.time        20 s
  start V:             OFF
  pause                0 s
  dos.element:         internal D0
  meas.input:          1
  temperature          25.0 °C
>stop conditions
  stop V:              abs.
  stop V               8 ml
  stop pH              OFF
  stop EP              9
  filling rate         max. ml/min
>statistics
  status:              OFF
>evaluation
  EPC                  5
  EP recognition:      all
  fix EP1 at pH        3.2
  fix EP2 at pH        OFF
  pK/HNP:              OFF
>preselections
  req.ident:           OFF
  req.smpl size:       OFF
  activate pulse:      OFF
=====
```




Fig. 4: Parameter report Titrino for the determination of free HBF₄.

Fig. 5: Titration curve free HBF₄.