

# Determination of Mercury in Plastics/ Plastic Waste (Granulated or Flaked Particles)

LECO Corporation; Saint Joseph, Michigan USA

## Instrument: AMA254



### Sample Preparation

The sample particles can be weighed directly into the boat. To ensure a sufficient combustion do not exceed the sample weight range used in this application.

### Sample Weight

Maximum 100 to 150 mg; balance precision of 0.1 mg or better

### Accessories

614-822-114 Large Nickel Boats

### Calibration Samples

LECO 502-813 Fly Ash, LECO 502-499 (BCR 143r), LECO 502-649 Dry Sludge (NIST 2781), or other suitable reference material

### Analysis Time

~8 minutes

### Method Profile

Drying Time:	60 seconds
Decomposition Time:	200 seconds
Cuvette Clear Time:	45 seconds
Dosing Delay Time:	0 seconds
Cell Selection:	Auto Select
Metric for Calculations:	Peak Area

### Procedure

- Determine the blank as follows.
  - Enter "Blank" from the drop-down menu under the "Name" column.
  - Click "Analyze"; the door will open and the nickel loop will be presented.
  - Carefully place a 614-822-114 Large Nickel Boat into the nickel loop using clean tweezers.
  - Click "OK" in the "Load Sample" window; the door will close and the analysis sequence will start automatically.
  - Repeat steps 1a through 1d two more times. The system and boats will be purged of any interfering elements.
- Calibrate the instrument as defined in the instructional manual.
  - Analyze various sample weights of a relevant reference material in accordance to the absolute amount of mercury required to calibrate an appropriate dynamic range. The calibration samples are weighed into the 614-822-114 Large Nickel Boat.

- Enter each calibration sample with the appropriate ID code from the drop-down menu, and sample weight from an external balance measurement.
- Click "Analyze"; the door will open and the nickel loop will be presented.
- If there is a boat in the nickel loop, remove it and keep for later use.
- Carefully place the calibration sample boat into the nickel loop using clean tweezers.
- Click "OK" in the "Load Sample" window; the door will close and the analysis sequence will start automatically.
- Repeat steps 2a through 2f as per the calibration procedures.

*Note: The first analyzed sample after a long delay should be discarded. This sample should be considered a conditioner for the system, and not used for the actual calibration.*

- Complete a calibration by following the calibration procedure as outlined in the manual.
  - Verify the calibration by analyzing one of the calibration samples again. It should be within the expected tolerances. If not, repeat steps 2a through 2i.
- Analyze the samples as follows.
    - Weigh ~100 mg of the high concentration sample into a 614-822-114 Large Nickel Boat.  
*NOTE: Use ~200 mg for low concentration samples.*
    - Enter a sample identification in the Name column and the sample weight in the Mass column.
    - Click "Analyze"; the door will open and the nickel loop will be presented.
    - If there is a boat in the nickel loop, remove it and keep for later use.
    - Carefully place the sample boat into the nickel loop using clean tweezers.
    - Click "OK" in the "Load Sample" window; the door will close and the analysis sequence will start automatically.

## Typical Results

### Plastic Particles, Fine • Nominal Value 0.1 ppm

Sample Weight (mg)	ng	ppm
148.0	14.889	0.1006
123.2	12.159	0.0987
148.5	15.355	0.1034

Mean Value: 0.1009 ppm

SD: 0.0024 ppm

RSD: 2.38%

### Plastic Particles, Coarse • Nominal Value 0.1 ppm

Sample Weight (mg)	ng	ppm
107.7	10.328	0.0959
119.3	13.147	0.1102
98.7	10.245	0.1038

Mean Value: 0.1033 ppm

SD: 0.0072 ppm

RSD: 6.93%

### LECO Corporation

3000 Lakeview Avenue • St. Joseph, MI 49085 • Phone: 800-292-6141 • Fax: 269-982-8977  
info@leco.com • www.leco.com • ISO-9001:2000 • No. FM 24045 • LECO is a registered trademark of LECO Corporation.