

Ion Chromatography Assay for Ammonia in Adenosine

Lipika Basumallick and Jeffrey Rohrer
Thermo Fisher Scientific, Sunnyvale, CA USA

Keywords

Suppressed Conductivity, Fast Separation, Pharmaceuticals, Drug Products, Dionex IonPac CS12A-5 μm Column

Introduction

Adenosine (Figure 1) is a naturally occurring purine nucleoside that forms from the breakdown of adenosine triphosphate, the primary energy source for living cells. This nucleoside is composed of adenine attached to ribose in the furanose conformation via a β -N9-glycosidic bond. The major therapeutic uses of adenosine are for treating surgical and nerve pain, pulmonary hypertension, and irregular heartbeat; for controlling blood pressure during anesthesia/surgery; and for cardiac stress tests. For these indications, adenosine is administered either as a bolus intravenous injection or as an intravenous infusion.




Figure 1. The structure of adenosine.

Adenosine is produced by the condensation of ribose and adenine biochemically or by fermentation. Ammonia is an impurity in adenosine preparations and the U.S. Pharmacopeia (USP) adenosine monograph describes an



assay to determine the amount of ammonia in adenosine. That assay measures ammonia by comparing the color of an adenosine sample and a 0.4 mg/L ammonium chloride standard solution after the addition of an alkaline mercuric/potassium iodide solution to each. The acceptance criteria is defined as "The Sample solution does not exhibit a more intense yellow color than that of the Standard solution (NMT [not more than] 4 ppm of ammonia)."¹ However, that assay is subjective, potentially exposes the analyst to mercury, and generates hazardous waste.

Goal: To develop an ion chromatography (IC)-based method for the determination of ammonia in adenosine that can meet the assay requirements and replace the USP adenosine monograph's color-based method

Equipment and Software

- Thermo Scientific™ Dionex™ ICS-5000+ Reagent-Free™ HPIC™ system, capable of supporting high-pressure IC, including:
 - SP Single Pump or DP Dual Pump
 - EG Eluent Generator
 - DC Detector/Chromatography Compartment
 - CD Conductivity Detector
- Thermo Scientific Dionex AS-AP Autosampler with 10 µL PEEK Sample Loop (P/N 036104)
- Thermo Scientific Dionex EGC III Methanesulfonic Acid (MSA) Eluent Generator Cartridge (P/N 074535)
- Thermo Scientific Dionex CR-CTC Continuously Regenerated Cation Trap Column (P/N 066262)
- Thermo Scientific™ Dionex™ Chromeleon™ Chromatography Data System software

Reagents and Standards

- Deionized (DI) water, 18 MΩ-cm resistance or better
- Combined Six-Cation Standard-I (Fisher Scientific P/N DX040187)
- Adenosine 99+% (Fisher Scientific P/N AC16404-0050)
- Ammonia Standard, 1000 mg/mL in water (Fisher Scientific P/N US-ICC-101)

Conditions

Columns:	Thermo Scientific™ Dionex™ IonPac™ CG12A-5 µm Guard, 3 × 30 mm (P/N 057184) Dionex IonPac CS12A-5 µm Analytical, 3 × 150 mm (P/N 057185)
Eluent:	33 mM MSA
Flow Rate:	0.7 mL/min
Inj. Volume:	10 µL
Column Temp:	30 °C
Detector Temp:	30 °C
Detection:	Suppressed conductivity, Thermo Scientific™ Dionex™ CSRS™ 300 Cation Self-Regenerating Suppressor, 2 mm (P/N 064557) or Dionex SC-CSRS 300 Salt-Converter Cation Self-Regenerating Suppressor (P/N 067529), autosuppression recycle mode, power setting 68 mA
Back Pressure:	2600 psi
Background	
Conductance:	<0.250 µS
Noise:	~0.1–0.2 nS

Preparation of Solutions and Reagents

Ammonium Standards

Prepare standards gravimetrically by making appropriate dilutions of a commercial 1000 mg/L standard with DI water. Store standard solutions at 4 °C when not in use.

Sample Preparation

Adenosine sample solution (use sample within 24 h)

In the current USP monograph for ammonia assay in adenosine, sample preparation steps are: (i) suspend 0.5 g in 10 mL of water, (ii) stir for 30 s, (iii) pass through a coarse filter, (iv) dilute the filtrate to 15 mL and use the filtrate for the assay.¹

For this IC-based method, prepare adenosine samples according to the USP monograph steps described above, except do not dilute the filtrate to 15 mL; instead, directly inject the filtrate into the IC system.

General Design of Robustness Study

The robustness of an analytical procedure is a measure of its capacity to remain unaffected by small variations in the procedural parameters. The robustness of each method was studied by determining the response of a 1 mg/L ammonia standard solution under typical variations of analytical conditions. Variations in the following were tested:

- Mobile phase concentration (± 2 mM MSA*)
- Flow rate ($\pm 10\%$)
- Temperature of column (± 2 °C)
- Column-to-column (column sets from two different production batches, Column Sets 1 and 2)

* *Not typical for eluent generation but rather for manually prepared eluents*

Table 2. Precision for ammonia.

Day	Sample	Retention Time (min)	Retention Time RSD	Peak Area ($\mu\text{S} \cdot \text{min}$)	Peak Area RSD	Asymmetry	Resolution
1	1 mg/L Standard	2.144	0.09	0.165	0.19	1.3	—
	4 mg/L Standard	2.148	0.09	0.4385	0.08	1.4	—
	10 mg/L Standard	2.150	<0.01	0.9416	0.22	1.6	—
	Adenosine Solution	2.147	<0.01	0.0055	1.83	1.2	2.2
2	1 mg/L Standard	2.144	0.09	0.1651	0.88	1.4	—
	4 mg/L Standard	2.148	0.09	0.4403	0.01	1.4	—
	10 mg/L Standard	2.152	0.09	0.9433	0.24	1.6	—
	Adenosine Solution	2.144	<0.01	0.0051	2.27	1.2	2.2
3	1 mg/L Standard	2.149	0.09	0.1673	0.34	1.3	—
	4 mg/L Standard	2.152	0.09	0.4451	0.06	1.4	—
	10 mg/L Standard	2.158	0.09	0.9561	0.38	1.6	—
	Adenosine Solution	2.151	0.09	0.0052	1.79	1.2	2.2

The ICH guidelines recommend that repeatability be assessed using a minimum of nine determinations covering the specified range for the procedure (i.e., three concentrations and three replicates of each concentration).³ The retention time RSDs were <0.1 and the peak area RSDs ranged from 0.01 to 2.3 (Table 2).

Detection Limit

The limit for ammonia in adenosine in the current monograph is NMT 4 ppm (i.e., mg/L) of ammonia.

USP/ICH defines limit of detection (LOD) in terms of the signal-to-noise ratio (S/N) as 2:1 or 3:1 and limit of quantitation (LOQ) as 1:10. Currently, USP defines S/N as $2H/h$, where H is the peak height from the middle of the noise band to the top of the peak. The value of h is the “difference between the largest and smallest noise values observed equal to at least five times the width at the half height of the peak and, if possible, situated equally around the peak of interest” or the “range of the noise in a chromatogram obtained after injection or application of a blank, observed over a distance equal to at least five times the width at half height of the peak in the chromatogram obtained with the prescribed reference solution, and if possible, situated equally around the place where this peak is found.”^{2,3}

According to the new recommendation outlined in Hinshaw and Dolan’s article, S/N must be measured as S/N_{p-p} , with S defined as measurement of peak height from the middle of a noise band to the highest point of the peak (this is synonymous with H); N_{p-p} is defined as peak-to-peak baseline noise, or noise measured over $\geq 5 \times$ peak width measured at half of the peak height.⁴

Column: Dionex IonPac CS12A-5 μm
 Eluent: 33 mM MSA
 Flow Rate: 0.7 mL/min
 Inj. Volume: 10 μL
 Temperature: 30 °C
 Samples: A. 0.025 mg/L Ammonium Standard
 B. DI Water Blank




Figure 3. Chromatograms of (A) a 0.025 mg/L ammonium standard and (B) a DI water blank.

For this IC method, the LOD and LOQ for ammonia in adenosine—i.e., concentrations that resulted in peaks that were $3 \times$ and $10 \times$ of the noise (both as h and N_{p-p})—are 0.001 and 0.004 mg/L, respectively.

Figure 3 shows ammonium at 0.025 mg/L (Chromatogram A) and a DI water blank with a background level of ~ 0.008 mg/L ammonia (Chromatogram B).

Linearity

The USP/ICH recommendations for establishing linearity of an impurity in a drug substance or a finished product are a minimum of five concentrations ranging from 50 to 120% of the acceptance criteria.^{2,3}




Figure 4. Calibration curve for ammonium using a Dionex CSRS 300 cation self-regenerating suppressor.

For ammonium, linearity was investigated in the range of 0.025 to 10 mg/L (0.025, 0.05, 1.0, 2.0, 4.0, 6.0, and 10.0 mg/L). The highest concentration investigated was 2.5× the acceptance criteria (NMT 4 mg/L) for ammonia in adenosine. Weak bases—like ammonium—are partially dissociated and thus give a nonlinear response as the concentration increases. The coefficient of determination was 0.9992 using a quadratic curve-fitting function for ammonium (Figure 4) in the range of 0.025 to 10 mg/L.

For a narrow range (0.025–2 mg/L), response can be linear and the coefficient of determination for a linear fit is 0.9984. Linear response can be extended by converting




Figure 5. Calibration plot for ammonium using a Dionex SC-CSRS 300 suppressor.

the weak base to the fully ionized salt form using a Dionex SC-CSRS 300 suppressor. Figure 5 shows the linear calibration plot for ammonium with a Dionex SC-CSRS 300 suppressor over the 0.025–50 mg/L range (coefficient of determination 0.9989).

Robustness

The robustness of an analytical procedure is a measure of its capacity to remain unaffected by small variations in the procedural parameters.

An ammonium standard (1 mg/L) was used as the check standard. The data are summarized in Table 3. The peak asymmetry for ammonia ranged from 1.3 to 1.4 (USP) in the 1 mg/L standard and 1.2 to 1.25 in the adenosine solution. The resolution of ammonia from sodium in the adenosine solution ranged from 1.96 to 2.32 under the different conditions listed in Table 3. The retention time (RT) RSD was <0.2 and the peak area RSD ranged from 0.17 to 1.8.

Table 3. Robustness.

Condition	Sample	RT (min)	RT RSD	Peak Area (μS * min)	Peak Area RSD	Asymmetry (USP)	Resolution (USP)
Column Set 1							
30 °C, 0.7 mL/min, 33 mM MSA	Combined Six-Cation Standard-I	2.14	0.19	0.3226	0.17	1.36	1.96
30 °C, 0.7 mL/min, 33 mM MSA	1 mg/L Standard	2.144	0.09	0.165	0.19	1.33	—
30 °C, 0.7 mL/min, 33 mM MSA	Adenosine Solution	2.147	<0.01	0.0055	1.83	1.2	2.24
30 °C, 0.7 mL/min, 35 mM MSA	1 mg/L Standard	2.082	0.11	0.1675	0.53	1.4	—
30 °C, 0.7 mL/min, 35 mM MSA	Adenosine Solution	2.083	<0.01	0.0058	0.87	1.22	2.14
30 °C, 0.7 mL/min, 31 mM MSA	1 mg/L Standard	2.223	<0.01	0.1669	1.12	1.34	—
30 °C, 0.7 mL/min, 31 mM MSA	Adenosine Solution	2.225	0.11	0.006	0.75	1.2	2.32
32 °C, 0.7 mL/min, 33 mM MSA	1 mg/L Standard	2.115	0.11	0.1673	0.14	1.37	—
32 °C, 0.7 mL/min, 33 mM MSA	Adenosine Solution	2.117	<0.01	0.0061	0.39	1.22	2.19
28 °C, 0.7 mL/min, 33 mM MSA	1 mg/L Standard	2.17	0.22	0.1663	0.79	1.36	—
28 °C, 0.7 mL/min, 33 mM MSA	Adenosine Solution	2.17	<0.01	0.0061	0.73	1.22	2.23
30 °C, 0.65 mL/min, 33 mM MSA	1 mg/L Standard	2.287	<0.01	0.1763	0.26	1.35	—
30 °C, 0.65 mL/min, 33 mM MSA	Adenosine Solution	2.283	<0.01	0.0067	1.72	1.23	2.22
30 °C, 0.75 mL/min, 33 mM MSA	1 mg/L Standard	2.017	<0.01	0.1533	0.20	1.36	—
30 °C, 0.75 mL/min, 33 mM MSA	Adenosine Solution	2.02	<0.01	0.0057	1.23	1.25	2.18
Column Set 2							
30 °C, 0.7 mL/min, 33 mM MSA	Combined Six-Cation Standard-I	2.16	<0.01	0.3111	<0.01	1.23	1.98
30 °C, 0.7 mL/min, 33 mM MSA	1 mg/L Standard	2.162	0.09	0.1667	0.87	1.29	—
30 °C, 0.7 mL/min, 33 mM MSA	Adenosine Solution	2.168	0.09	0.0052	1.49	1.2	2.17

Conclusion

The Dionex IonPac CS12A-5 μm column with its small-diameter resin enables a fast and accurate isocratic method for the determination of ammonia in adenosine. A run time of 5 min yields high sample throughput. The electrolytically generated eluent and the self-regenerating suppressor contribute to the robustness and sensitivity of the method. In addition, only water is used for eluent generation; no hazardous chemicals are required.

The precision (RT RSD <0.1, peak area RSD <2), accuracy (average recovery 88–116%), linearity, detection and quantitation limits, and robustness for the assays presented in this study show that the proposed IC method for ammonia determination in adenosine meets the analytical performance characteristics outlined in USP General Chapter <1225>.

References

1. Adenosine 2079. *U.S. Pharmacopeial Convention (USP) Erata for USP-NF*; Rockville, MD, USP35-NF30.
2. International Conference on Harmonisation of Technical Requirements for Registration of Pharmaceuticals for Human Use (ICH). [Online] www.ich.org/fileadmin/Public_Web_Site/ICH_Products/Guidelines/Quality/Q2_R1/Step4/Q2_R1__Guideline.pdf (accessed Nov 20, 2013).
3. Validation of Compendial Methods, General Chapter <1225>. *U.S. Pharmacopeia/National Formulary*; Rockville, MD, 2012; pp 877–881.
4. Hinshaw, J.V.; Dolan, J.W. Stimuli to the Revision Process: Signal-to-Noise Measurements from Chromatographic Data. *Pharmaceutical Forum* 2012, 38 (3) [Online] www.usppf.com/pf/pub/index.html (accessed Oct 18, 2013).

www.thermofisher.com/dionex

©2016 Thermo Fisher Scientific Inc. All rights reserved. All trademarks are the property of Thermo Fisher Scientific Inc. and its subsidiaries. This information is presented as an example of the capabilities of Thermo Fisher Scientific Inc. products. It is not intended to encourage use of these products in any manners that might infringe the intellectual property rights of others. Specifications, terms and pricing are subject to change. Not all products are available in all countries. Please consult your local sales representative for details.

Africa +43 1 333 50 34 0
Australia +61 3 9757 4300
Austria +43 810 282 206
Belgium +32 53 73 42 41
Brazil +55 11 3731 5140
Canada +1 800 530 8447
China 800 810 5118 (free call domestic)
 400 650 5118

Denmark +45 70 23 62 60
Europe-Other +43 1 333 50 34 0
Finland +358 9 3291 0200
France +33 1 60 92 48 00
Germany +49 6103 408 1014
India +91 22 6742 9494
Italy +39 02 950 591

Japan +81 6 6885 1213
Korea +82 2 3420 8600
Latin America +1 561 688 8700
Middle East +43 1 333 50 34 0
Netherlands +31 76 579 55 55
New Zealand +64 9 980 6700
Norway +46 8 556 468 00

Russia/CIS +43 1 333 50 34 0
Singapore +65 6289 1190
Sweden +46 8 556 468 00
Switzerland +41 61 716 77 00
Taiwan +886 2 8751 6655
UK/Ireland +44 1442 233555
USA +1 800 532 4752

AN70794_E 08/16S

Thermo
 S C I E N T I F I C

Part of Thermo Fisher Scientific