

Screening Method for 30 Pesticides in Green Tea Extract Using Automated Online Sample Preparation with LC-MS/MS

Yang Shi, Catherine Lafontaine, Thermo Fisher Scientific, Franklin, MA, USA

Fangting Ye, Zheng Jiang, Le Ma, Ting Liu, Haijian Wang, Thermo Fisher Scientific, Shanghai, China

Key Words

- Transcend TLX-1 System
- TurboFlow Technology
- TSQ Access MAX
- Food Safety

Introduction

Analysis of pesticide residues has been one of the most important tasks of food safety laboratories. Mass spectrometers (MS), with liquid chromatography coupled to triple stage quadrupole mass spectrometers (LC-MS/MS), have been the main tools used in pesticide residue analysis. There is a consensus that sample preparation is becoming the bottleneck to the entire workflow. Traditional sample preparation methods, usually involving liquid-liquid extraction (LLE) or solid phase extraction (SPE), can be time-consuming and labor-intensive. In addition, low recovery, matrix interference and poor reproducibility are among other major concerns. In recent years, a rapid processing method, QuEChERS, has gained popularity. The QuEChERS method makes it easier and less expensive for analytical chemists to examine pesticide residues in various food matrices¹. However, some reports show matrix interference tends to be severe after QuEChERS, and the mass spectrometer is more vulnerable to contamination by highly complex food matrices².

In this study, we describe an easy, comprehensive, on-line screening LC method using a Thermo Scientific Transcend TLX-1 system powered by Thermo Scientific TurboFlow technology to analyze multiple pesticide residues in green tea extract. Figure 1 illustrates a typical Transcend™ TLX-1 system with the Thermo Scientific TSQ Access MAX triple stage quadrupole mass spectrometer.

Goal

Develop a rapid and sensitive automated online sample preparation LC-MS/MS method to screen for multiple pesticides in green tea extract.

Experimental

The matrix standard curve

One gram of Chinese green tea was extracted using 10 mL HPLC grade acetonitrile followed by 15 minutes of ultra-sonication. The extract was then filtered through a 0.45 µm membrane filter. The resultant solution was used to prepare the matrix calibrators and QC samples. The matrix calibrant concentrations are 6.25 µg/L, 12.5 µg/L, 25 µg/L, 50 µg/L and 100 µg/L, respectively. The matrix QC sample concentration is 10 µg/L.

TurboFlow™ Method Parameters

System:	Transcend TLX-1 system controlled by Thermo Scientific Aria OS 1.6.3 software
Column:	TurboFlow Cyclone 0.5 x 50 mm
Injection Volume:	10 µL
Loading Solvent:	0.1% formic acid in water
Loading Flow Rate:	1.5 mL/min
Eluting Solvent:	0.1% formic acid in methanol



Figure 1. Typical layout of a Transcend TLX-1 system with a TSQ Access MAX™ triple stage quadrupole mass spectrometer.

HPLC Method Parameters

Analytical Column: Thermo Scientific Hypersil GOLD
2.1 x 100 mm, 3 μ m
Solvent A: 0.1% formic acid in water
Solvent B: 0.1% formic acid in methanol

Mass Spectrometer Parameters

MS: TSQ Quantum Access MAX
MS Ionization Source: Heated Electrospray Ionization (H-ESI)
Ion Polarity: Positive ion mode
Spray Voltage: 2 KV
Sheath Gas Pressure (N_2): 30 arbitrary units
Auxiliary Gas Pressure (N_2): 15 arbitrary units
Vaporizer Temperature: 300 $^{\circ}$ C
Capillary Temperature: 300 $^{\circ}$ C
Collision Gas Pressure: 1.5 mTorr




Figure 2. TurboFlow method schematic as viewed in the Aria OS software.

Results and Discussion

The Extraction and Separation of 30 Pesticide Residues

In 2006, Japan released the most stringent pesticide-related regulation in history entitled “Positive List System for Agricultural Chemical Residues in Foods”³. Since Japan is China’s major tea importer, the limits discussed in the current study follow this regulation. As described in the Experimental section, the tea matrix standard samples are 6.25 μ g/L, 12.5 μ g/L, 25 μ g/L, 50 μ g/L and 100 μ g/L, respectively. The matrix QC samples are 10 μ g/L. Figure 3 shows the representative chromatograms at 6.25 μ g/L, which has been determined as the lower limit of quantitation (LLOQ). The data demonstrate that 30 pesticides were well separated with good peak shape. The peaks’ signal to noise ratios are far greater than the required 10:1 at the LLOQ. Table 1 shows the linear curve for these 30 analytes. All R^2 values are between 0.993-0.999. The relative standard deviation (RSD) for 6 consecutive injections of 6.25 μ g/L calibrator was in the range of 2.85% -7.48%.

Background Reduction Effects using TurboFlow Technology

By using the Transcend TLX system with TurboFlow technology, the background noise and interference peaks are reduced significantly. Figure 4 compares chromatograms of Clomazone at 6.25 μ g/L in tea extract using standard HPLC (top) and the TurboFlow method (bottom). The left panel (A-1 and B-1) shows the primary transition of m/z 240 > 125. The right panel (A-2 and B-2) shows the secondary transition of m/z 240 > 89. It clearly shows the effectiveness of background reduction using TurboFlow technology while the signal to noise ratio increased by 3 and 4 times for m/z 125 and 89 transitions, respectively. The area responses of both peaks also increase by more than 50% due to the minimization of ion suppression incurred by matrix. We also noticed the mass spectrometry response become more stable across the entire tested concentration range, thus improving the method reliability.

A Simple Method Optimization Process

During TurboFlow method development, the sample loading condition, elution solvents and many other parameters may need to be optimized. Aria™ OS 1.6.3 operation software for Transcend systems offers a method variable function. By utilizing this unique tool, different parameters can be easily tried using the same method in a single batch. For example, in this study, one of the critical steps was to find the optimal solvent content in the transfer loop to elute the target analytes completely from the TurboFlow column without introducing unnecessarily high organic solvent into the analytical column. We compared 5 different concentration ratios of 0.1% formic acid in acetonitrile to 0.1% formic acid in water (10:90, 30:70, 50:50, 70:30 and 90:10). The results indicated that with the increase of organic content, the target compounds were more completely washed off from TurboFlow column. However, once the organic concentration reached 50%, the elution strength was approaching a balance. Therefore, we chose 50:50 as the optimal elution ratio of organic to aqueous solvent in the transfer loop. Another example of method optimization appears in Figure 5, showing the effects of the loading flow rate on Dimethametryn’s elution peak shape. All these tests were done in just one sample batch without writing multiple methods, which simplified the method development process and improved method reliability.

The Comparison of TurboFlow Technology with Two of the Most Popular Pesticides Sample Preparation Methods

As shown in Figure 6, we compared a TurboFlow method and two currently popular methods for pesticide residue sample preparation, SPE and QuEChERs. A typical SPE method involves equilibrating the cartridge, loading, washing and eluting analytes. It usually takes about 1 week to process 100 samples. Although QuEChERs was designed to simplify sample preparation, it still requires two-step centrifugation and concentration. A few days are typically required to prepare 100 samples with QuEChERs. TurboFlow technology minimizes preparation of 100 samples to less than 3 hours, dramatically improving the efficiency and throughput of this routine lab test.

Table 1: Standard curve linearity and QC results for the 30 pesticides in tea extract.

Compound	RT (min)	Parent ion (m/z)	Product ion (m/z)	Collision Energy (V)	Linear Curve	R ²	CV% (n = 6) QC = 10 µg/L
Prometon	4.82	226.0	184.0 142.1	20 27	Y=167343+396533X	0.999	2.99%
Ametryn	5.07	228.0	186.0 96.0	26 34	Y=83264.1+194461X	0.999	2.85%
Dimethametryn	5.68	256.1	186.1 158.1	21 27	Y=166875+605055X	0.999	3.29%
Mefenoxam	5.79	280.0	220.0 192.0	17 20	Y=460109+272420X	0.998	4.23%
Monolinuron	5.85	215.0	126.0 99.0	17 36	Y=-10985.6+18335.3X	0.998	6.51%
Isoprocarb	5.94	194.0	95.0 137.0	16 11	Y=-18662+12428.2X	0.999	6.43%
Dimethachlor	6.01	256.0	224.0 148.0	15 28	Y=-23531.9+96341.3X	0.997	5.53%
Clomazone	6.05	240.0	125.0 89.0	20 37	Y=-43447.5+42181.6X	0.998	6.37%
Furalaxyl	6.21	302.0	242.0 270.0	15 10	Y=358101+267257X	0.998	4.85%
Azoxystrobin	6.33	404.0	372.0 329.0	15 33	Y=538988+377945X	0.997	4.00%
Triadimefon	6.39	294.0	197.0 225.0	19 19	Y=-20167.4+16685.8X	0.997	7.31%
Ethoprophos	6.41	243.0	131.0 97.0	21 33	Y=-13814+14313.8X	0.997	7.48%
Iprobenfos	6.52	289.0	205.0 91.0	12 23	Y=53008.6+137376X	0.999	6.28%
Isoprothiolane	6.57	291.0	189.0 231.0	22 12	Y=123106+87284X	0.998	6.00%
Flutolanil	6.60	324.0	242.0 262.0	26 18	Y=6077+47866X	0.998	6.56%
Propiconazole	6.65	342.0	159.0 69.0	30 31	Y=-17113.2+36428.7X	0.997	6.74%
Benalaxyl	6.78	326.0	148.0 208.0	25 20	Y=172291+126493X	0.997	5.92%
Pirimiphos-methyl	6.81	306.0	164.0 108.0	22 33	Y=227752+204491X	0.994	4.73%
Picoxystrobin	6.82	368.0	145.1 205.0	22 7	Y=320093+78661.3X	0.993	4.03%
Diazinon	6.90	305.0	169.0 153.0	24 26	Y=182248+386247X	0.998	4.96%
Thiazopyr	6.95	397.0	335.0 275.0	30 40	Y=-5052.12+18434.8X	0.997	6.47%
Piperophos	7.09	354.0	171.0 143.0	25 33	Y=142671+143459X	0.996	4.68%
Trifloxystrobin	7.13	409.0	186.0 206.0	21 16	Y=-18755.2+43150.6X	0.998	6.22%
Tebufenpyrad	7.16	334.0	145.0 117.0	28 36	Y=-3267.09+9390.51X	0.998	7.01%
Piperonyl butoxide	7.25	356.0	177.0 119.0	13 33	Y=-300922+175066X	0.996	4.16%
Pyriproxyfen	7.34	322.0	96.0 185.2	16 27	Y=-19160.9+56881.6X	0.999	4.73%
Tralkoxydim	7.39	330.0	284.0 138.0	15 20	Y=-8119.47+46536.4X	0.997	5.01%
Fenazaquin	7.55	307.0	161.0 57.0	18 23	Y=-56587.3+97365.3X	0.998	2.76%
Butralin	7.58	296.0	240.0 222.0	15 20	Y=-2485.92+25777.9X	0.998	4.72%
DEF	7.85	315.0	169.0 113.0	15 25	Y=-8658.91+7992.12X	0.998	3.89%




Figure 3. Selected ion chromatograms at LLOQ of 6.25 µg/L for all 30 analytes (same as the order in Table 1).




Figure 4: Comparison of chromatograms of Clomazone at 6.25 µg/L in tea extract using standard HPLC (top) and the TurboFlow method (bottom).





Figure 5. Effect of the loading flow rate on Dimethametryn's elution peak shape.



In addition to these offices, Thermo Fisher Scientific maintains a network of representative organizations throughout the world.

Africa-Other
+27 11 570 1840

Australia
+61 3 9757 4300

Austria
+43 1 333 50 34 0

Belgium
+32 53 73 42 41

Canada
+1 800 530 8447

China
+86 10 8419 3588

Denmark
+45 70 23 62 60

Europe-Other
+43 1 333 50 34 0

Finland/Norway/Sweden
+46 8 556 468 00

France
+33 1 60 92 48 00

Germany
+49 6103 408 1014

India
+91 22 6742 9434

Italy
+39 02 950 591

Japan
+81 45 453 9100

Latin America
+1 561 688 8700

Middle East
+43 1 333 50 34 0

Netherlands
+31 76 579 55 55

New Zealand
+64 9 980 6700

Russia/CIS
+43 1 333 50 34 0

South Africa
+27 11 570 1840

Spain
+34 914 845 965

Switzerland
+41 61 716 77 00

UK
+44 1442 233555

USA
+1 800 532 4752

Figure 6. Comparison of the TurboFlow method to SPE and QuEChERS.

Conclusion

A quick, automated online sample preparation LC-MS/MS method has been developed that is sensitive enough to screen the tested pesticides in tea extracts. The method detection and quantitation limits are significantly lower than the strictest limits set by the Japanese government. TurboFlow technology eliminates the need for time-consuming sample preparation procedures such as SPE and QuEChERS. By using Aria OS software, the method development and optimization process is greatly simplified.

References

1. Steven J. Lehotay, Quick, Easy, Cheap, Effective, Rugged and Safe (QuEChERS Approach for Determining Pesticide Residues. *Methods in Biotechnology*, 2006, 19, 239-261
2. Steven J. Lehotay, Michelangelo Anastassiades, Ronald E. Majors, QuEChERS, a Sample Preparation Technique that Is "Catching On": An Up-to-Date Interview with the Inventors, LCGC North America, July 2010
3. The Japanese Positive List System for Agricultural Chemical Residues in Foods (enforcement on May 29, 2006)
4. GB/T 19648-2005 (466 Pesticides analysis method in Fruit and Vegetable GC-MS and LC-MS-MS)

www.thermoscientific.com

Legal Notices: ©2011 Thermo Fisher Scientific Inc. All rights reserved. All trademarks are the property of Thermo Fisher Scientific Inc. and its subsidiaries. This information is presented as an example of the capabilities of Thermo Fisher Scientific Inc. products. It is not intended to encourage use of these products in any manners that might infringe the intellectual property rights of others. Specifications, terms and pricing are subject to change. Not all products are available in all countries. Please consult your local sales representative for details.



Thermo Fisher Scientific, San Jose, CA USA is ISO Certified.

AN63386_E 01/11S