

IC Application Note No. C-117


Title: MiPT – Metrohm intelligent Partial Loop Technique

Summary: Calibration of lithium, sodium, ammonium, zinc, potassium, magnesium and calcium applying partial loop technique using cation chromatography with direct conductivity detection.

This technique allows a calibration range of 1:100 (e.g. 1 µg/L to 100 µg/L corresponding to 2 µL to 200 µL injected volume) out of 1 calibration solution. Applying the full range of partial loop injection to the samples one calibration covers a sample concentration range of 1 to 10'000 e.g. 2 µL of a 10 mg/L solution corresponds to the highest calibration level (100 µg/L) while 200 µL of a 1 µg/L solution corresponds to the lowest calibration level.


Sample: 100 µg/L multi-cation standard solution
Sample Preparation: –

Column: 6.1050.420 Metrosep C 4 – 150
Eluent: 2.5 mmol/L oxalic acid
Flow: 1.0 mL/min
Injection Volume: 250 µL loop; variable volumes 2...200 µL



Calibration: 1...100 µg/L; 2...200 µL	Li ⁺	Na ⁺	NH ₄ ⁺	Zn ²⁺	K ⁺	Mg ²⁺	Ca ²⁺
Correlation coefficient	0.99999	0.99999	0.99999	0.99994	0.99999	0.99999	0.99990
Percentage standard deviation (%RSD)	0.41	0.42	0.45	1.58	0.30	0.58	2.04

Example calibration curve:



Function: $A = 4.68855\text{E-}5 + 5.37919\text{E-}6 \times Q$
 Relative standard deviation 0.296085 %
 Correlation coefficient 0.999998

Carryover test:

200 μL	Sample $\mu\text{g}/\text{L}$	Blank (ultrapure water) $\mu\text{g}/\text{L}$	Carryover %
Lithium	10'000	0.1	0.001

Carryover was evaluated by injection of a blank (ultrapure water) immediately after injection of a 10'000 $\mu\text{g}/\text{L}$ -Lithium standard.

Precision:

Lithium	Volume injected μL	RSD (n=6) %
500 $\mu\text{g}/\text{L}$	2	1.38
500 $\mu\text{g}/\text{L}$	10	0.51
500 $\mu\text{g}/\text{L}$	20	0.16
500 $\mu\text{g}/\text{L}$	400	0.30
500 $\mu\text{g}/\text{L}$	200	0.02

Principle of MiPT: Dosino and 2-mL buffer tubing is required.

