


Thermo. Titr. Application Note No. H-086

Title: Determination of Sodium in Potato Chips

Scope: Determination of sodium in potato chips (crisps) and similar fried snack foods

Principle: Titration of an unfiltered suspension of the sample with a standardized solution of aluminium containing a stoichiometric excess of potassium ions in the presence of ammonium hydrogen difluoride at ~ pH 3 to give an exothermic reaction, forming insoluble NaK_2AlF_6 .



The titrant is standardized against a solution prepared from anhydrous sodium sulfate or sodium carbonate

Reagents: *Titrant:* Mixed 0.5 mol/L $\text{Al}(\text{NO}_3)_3$, 1.1 mol/L KNO_3 solution.
Conditioning reagent: 300 g/L $\text{NH}_4\text{F}\cdot\text{HF}$

Method: *Basic Experimental Parameters:*

Titrant delivery rate (mL/min.)	4
Titrant pre-dose (mL)	0.2
Delay before titration commences (sec)	60
No. of exothermic endpoints	1
Data smoothing factor (DSF)	45
Stirring speed (802 stirrer)	15

Sample preparation: Crush the potato chips to coarse crumbs in the bag. Weigh accurately approximately 10 g of crushed chips into a clean, dry 200 mL blender container. Add 150 mL D.I. water by bulb pipette. Place on a blender for 60 secs., until a stable suspension with minimum settling is produced.

Titration program: The titration program is set up to pre-dose 0.2 mL of titrant before the actual titration starts. A delay of 60 seconds is programmed to allow an initial precipitate seed surface upon which subsequent precipitate may grow rapidly. This helps to improve the precision of the method. The pre-dose volume is added to the titrant volume at the endpoint to obtain the actual titrant consumption.

Titration: Weigh accurately 25-30 g of homogenized suspension into a titration vessel. Add 20 mL D.I. water and 5 mL 300 g/L NH₄F·HF solution.

Examples: *Analysis of flavoured potato chips (crisps) – Woolworths Home Brand Variety pack*

	Chip Flavour	Sodium, mg/100g	
		As per label	Analyzed here
	Regular	570	875 ± 15 (n=5)
	Chicken	874	857 ± 13 (n=5)
	Barbecue	898	1132 ± 13 (n=5)
	Salt & Vinegar	1100	1545 ± 20 (n=5)
	Sour Cream & Onion	810	826 ± 14 (n=5)

Calculations:

1. Calculation of mass of chips in aliquot.

Mass of chips in blender = c

Mass of water in blender = w

Mass of aliquot taken for analysis = s


Mass of chips in aliquot (m) = $s \times c / (w+c)$

2. mg Na/100mg product:

$$= \frac{((\text{titre, mL} - \text{blank, mL}) \times M_{Al} \times AW_{Na} \times 100)}{m}$$

Thermometric Titration Plot:

Legend:
Red = solution
temperature curve
Black = second
derivative curve (for
endpoints)



Titration of sodium in barbecue flavoured chips