

# Thermo. Titr. Application Note No. H-091

**Title:** Standardization of tetrasodium EDTA solutions with standard magnesium solution

**Scope:** Standardization of 1mol/L tetrasodium EDTA ( $\text{Na}_4\text{EDTA}$ ) solutions by titration with standard magnesium solution.

**Principle:** Aliquots of standard  $\text{Mg}^{2+}$  solution are titrated with  $\text{Na}_4\text{EDTA}$  solution to an endothermic endpoint. From a plot of mmol  $\text{Mg}^{2+}$  (x-axis) against mL  $\text{Na}_4\text{EDTA}$  (y-axis), the gradient of the linear regression is computed. The molarity of the  $\text{Na}_4\text{EDTA}$  solution is equal to the reciprocal of the gradient.

**Reagents:** *Titrant:* 1mol/L  $\text{Na}_4\text{EDTA}$   
*Standard:* 0.2mol/L  $\text{Mg}^{2+}$   
 $\text{NH}_3/\text{NH}_4\text{Cl}$  buffer. Dissolve 87.5g  $\text{NH}_4\text{Cl}$  in 568mL 28% w/v  $\text{NH}_3$  solution and dilute to 1000mL with DI water

**Method:** *Basic Experimental Parameters:*

Titration delivery rate (mL/min.)	4
No. of exothermic endpoints	1
Data smoothing factor	55
Stirring speed (802 stirrer)	8
Delay before start (secs.)	10

*Preparation of standard  $\text{Mg}^{2+}$  solution*

Prepare the Mg metal by scraping A.R. Mg ribbon free of all oxide with a sharp blade until uniformly shiny. Cut into ~20mm lengths. Weigh 2.4305g, and transfer to a 500mL volumetric flask. Add 200mL DI water to the flask. Holding the flask at an angle, slowly add 25mL concentrated HCl down the neck. The flask should then be placed at an angle of approximately  $45^\circ$  so that the effervescence is directed against the walls of the flask. When all effervescence has ceased and all metal has dissolved, cool the flask contents to room temperature, and make to volume with DI water.

### Titration Program

A titration program is set up to pre-dose 10mL buffer solution, after which the titration commences automatically.

### Determination of molarity

Into individual titration vessels, pipette 10, 15, 20, 25 and 30mL standard 0.2mol/L  $Mg^{2+}$  solution, equivalent to 2, 3, 4, 5 and 6mmol  $Mg^{2+}$  respectively. Add DI water to bring the total volume in each vessel to approximately 30mL. Either titrate individually or place in the rack of a 814 Sample Processor for automated titration.


Plot mmol  $Mg^{2+}$  (x-axis) against mL  $Na_4EDTA$  titrated (y-axis) and determine the gradient and y-intercept by linear regression. Determine the molarity of the  $Na_4EDTA$  titrant as the reciprocal of the gradient.

### Example:

Standard  $Mg^{2+}$  solution: weighed 2.4285g Mg, made to 500mL = 0.19984mol/L  $Mg^{2+}$

mL $Mg^{2+}$ soln.	mmol $Mg^{2+}$	mL $Na_4EDTA$ titrated
30.00	5.9951	5.456
25.00	4.9959	4.565
20.00	3.9967	3.652
15.00	2.9975	2.751
10.00	1.9984	1.864

### Linear Regression Plot:



**Calculation:**

$1 \text{ mol Mg}^{2+} \equiv 1 \text{ mol Na}_4\text{EDTA}$ ,  
thus  $M \text{ Na}_4\text{EDTA} = 1/\text{gradient}$   
 $= 1/0.90061 = 1.1104 \text{ mol/L Na}_4\text{EDTA}$

**Thermometric Titration Plot:****Legend:**

Red = solution  
temperature curve  
Black = second  
derivative curve (for  
endpoints)

