

KF Application Note No. K- 18

Title: Water in cyclopropyl methyl ketone

Summary: The water content of cyclopropyl methyl ketone is determined according to Karl Fischer by coulometric titration using special reagents for aldehydes and ketones.

Sample: Cyclopropyl methyl ketone

Sample Preparation: none

Instruments and Accessories: 737 KF Coulometer, cell without diaphragm, 728 Magnetic Stirrer, printer

Analysis: Fill Hydranal Coulomat AK into the measuring cell and condition it until the drift is steady and below 10 ug/min. Inject ca. 300 mg sample with a syringe (the exact sample mass is determined by difference weighing) and start the automatic determination.

Reagents:
Hydranal Coulomat AK (Riedel-de Haën)

Results: Sample 1: AVG(5) = 0.166 +/- 0.002 % water
Sample 2: AVG(5) = 532 +/- 6 ppm water

Settings: 737 KF Coulometer

smpl.req:	on
d.start	20 ug/min
extr.	0 s
stop drift:	auto
delay time	3 s
report:	full