


Trace analysis of iodide in sodium chloride applying amperometric detection



Pure sodium chloride contains much less iodide than e.g. table salt which usually is fortified with it. Trace determination of iodide is easily performed applying ion chromatography with amperometric detection. This detection mode is specially selective and sensitive. The actual separation is achieved using a Metrosep A Supp 5 - 250/4.0 column. The detection happens at a silver working electrode. LOQ is at approximately 1.0 µg/L (in solution) and 50 µg/kg in the sample. The use of a shorter column might further improve the LOQ.

Results

| Anion | Conc. [µg/L] | RSD [% , n = 3] | Spike 1 [µg/L] | Recovery [%] | Spike 2 [µg/L] | Recovery [%] |
|--------|---------------|-----------------|----------------|--------------|----------------|--------------|
| Iodide | 1.55 | 4.2 | 6.75 | 104.1 | 12.03 | 104.8 |
| | Conc. [µg/kg] | | | | | |
| Iodide | 77.5 | | | | | |

Sample

Sodium chloride

Sample preparation

20 g dissolved and diluted to 1 L with ultrapure water

Columns

| | |
|-----------------------------|------------|
| Metrosep A Supp 5 - 250/4.0 | 6.1006.530 |
| Metrosep A Supp 5 Guard/4.0 | 6.1006.500 |

Solutions

| | |
|--------|---|
| Eluent | 6.4 mmol/L sodium carbonate 2.0 mmol/L sodium hydrogen carbonate |
|--------|---|

Parameters

| | |
|--------------------|------------|
| Flow rate | 0.8 mL/min |
| Injection volume | 50 µL |
| Pmax | 15 MPa |
| Recording time | 20 min |
| Column temperature | 40 °C |

PAD Parameters


| | |
|---------------------|---------------|
| Cell | Wall-Jet cell |
| Working electrode | Silver |
| Reference electrode | Ag/AgCl |
| Spacer | 50 µm |
| Measuring potential | 0.15 V |
| Range | 200 nA |
| Temperature | 35 °C |
| Mode | DC |

Analysis

Amperometric detection

Instrumentation

| | |
|--|------------|
| 930 Compact IC Flex Oven/Deg | 2.930.2160 |
| IC Amperometric Detector | 2.850.9110 |
| 858 Professional Sample Processor | 2.858.0020 |
| IC equipment Wall-Jet cell: without electrodes | 6.5337.000 |
| Ag working electrode | 6.1257.240 |
| Ag/AgCl reference electrode | 6.1257.720 |



www.metrohm.com

 **Metrohm**