

VA Application Note No. V - 137

Title: Suppressor «Cupraspeed» in acid copper baths (Atotech)

Summary: Determination of suppressor «Cupraspeed» in acid copper baths by dilution titration (DT) using cyclic voltammetric stripping (CVS).

Sample: Acid copper electroplating bath

Sample preparation: None

Analysis of suppressor «Cupraspeed»

Electrolyte Virgin make-up solution (VMS)
CuSO₄, H₂SO₄ and NaCl concentrations according to the supplier specifications.

Measuring solution 100 mL VMS

Working electrode (WE) **Pt-RDE:**
Drive shaft6.1246.000
+ Pt tip for CVS6.1204.160

Auxiliary electrode (AE) **Pt**.....6.0343.000

Reference electrode (RE) Reference system: Ag/AgCl/KCl (3 mol/L)6.0728.020
Intermediate electrolyte: KNO₃ sat.:H₂O (3:1) ..6.1245.010

Parameters		
	Working electrode	RDE (hydrodynamic measurement)
Stirrer speed	2000 rpm	
Mode	CVS	
Calibration technique	DT	
Start potential	1.575 V	
First vertex potential	-0.25 V	
Second vertex potential	1.575 V	
Voltage step	0.006 V	
Sweep rate	0.1 V/s	
Peak potential (Cu)	0.2 V ± 0.2 V	
Evaluation ratio (Q/Q(0))	0.7	

Determination of suppressor «Cupraspeed»

