

High Throughput Screening and Characterization of Bispecifics Using Native Ion Mobility Mass Spectrometry

Caroline S. Chu¹; Te-Wei Chu¹; Gregory O Staples¹; Patrick D Perkins¹; Andy Gieschen²; Christian Klein¹; Carol H. Ball³; Ning Tang¹
¹Agilent Technologies, Santa Clara, CA; ²Agilent Technologies, La Jolla, CA; ³Agilent Technologies, Wilmington, DE

ASMS 2016
ThP #560



Introduction

Bispecific (bsAbs) are an emerging group of biotherapeutics. While there are over 30 bsAbs in development, only two bsAbs are approved for therapy.¹ Bispecific antibodies are a unique class of antibodies combining the specificities of two monoclonal antibodies (mAbs) to bind two different targets at the same time.^{1,2} Generating bsAbs creates different combinations of the heavy and light chains from the original two mAbs.^{1,2} The heterogeneity of the resulting bsAbs poses an analytical challenge for high throughput screening for determining the optimal antibody pairs. Here we present an LC/MS method using native and denaturing ion mobility mass spectrometry the screening and characterization of bsAbs.

Experimental

A method for generating a bsAb was adapted from Debaene et al.³ using an in-house IgG1 (mAb A) and a commercial IgG1 standard (mAb STD). The two IgG1s were mixed with 5mM glutathione-reduced (GSH) for 24 hours at 37°C. The resulting bsAb and both starting IgG1s were buffer exchanged into 100 mM ammonium acetate buffer using Micro Bio-spin 30 columns. The samples were introduced into the IM-QTOF using a nanoLC interface and sample introduction was performed isocratically at a flowrate of 0.4 µL/min with 200 mM ammonium acetate, pH 7. Drift times and collision cross sectional areas were determined using the IM-MS Browser. Protein deconvolution was performed using the maximum entropy algorithm found in BioConfirm for MassHunter Qualitative Analysis.

Native Ion Mobility Mass Spectrometry for Screening of IgGs




Figure 2: Native IM-MS profiling of the different subclasses of IgG1 standard (above left), IgG2-κ (above right), IgG3-κ (bottom left), and IgG4-λ (bottom right) using nanoLC with the G1992A interface with 200 mM Ammonium Acetate, pH 7.0 as the mobile phase.




Figure 3: Native IM-MS profiling of the two IgG1s selected for generating a bispecific. The charge envelope, drift spectra and drift scope of the mAb STD (above left) and mAb A (above right) were profiled.

Preliminary Results for Bispecific




Figure 4: Native MS of the bispecific including the charge envelope and corresponding drift spectra (above). Zoom in of the drift scope of the charge envelope with the putative bispecific annotated and the charge envelope from mAb A (top right). Zoom in of the m/z values corresponding to the bispecific with the drift spectra, indicating the presence of two conformations.




Figure 5: Overlay of the drift spectra at the different charge states of mAb STD (left), mAb A (middle), and the resulting bispecific (right).

Agilent 6560 Ion Mobility Q-TOF

- Allows use of drift gases other than N₂
- Automatically controls pressure in drift tube of pure gas
- Available as an upgrade kit

Alternate Drift Gas Upgrade Kit

- Small molecule tune with fragile molecule options.
- Tunes entire system, including IM
- Capability to schedule and run in a worklist

SWARM Autotune

- IM Molecular feature extraction (iMFE) for compound extraction based on m/z , retention time, and drift time
- Annotations on drift scope

Ion Mobility Molecular Feature Extraction, iMFE

- Calculate collision cross section (CCS) for compounds
- Determine the CCS for compounds in a mixture identified with iMFE

Single Field Collision Cross Section, CCS

- Alternating frame acquisition with high/low collision energy
- Ramped collision energy per drift cycle
- Exportable MS/MS to Spectrum Mill and MassHunter Qualitative Analysis

Alternating High/Low Collision Energy



Conclusion

- High throughput intact mass analysis of bispecific antibodies using native and denaturing ion mobility mass spectrometry
- ¹Spiess, Christoph, Qianting Zhai, and Paul J. Carter. "Alternative molecular formats and therapeutic applications for bispecific antibodies." *Molecular immunology* 67.2 (2015): 95-106.
- ²Kontermann, Roland E., and Ulrich Brinkmann. "Bispecific antibodies." *Drug discovery today* 20.7 (2015): 838-847.
- ³Debaene, François, et al. "Time resolved native ion-mobility mass spectrometry to monitor dynamics of IgG4 Fab arm exchange and "bispecific" monoclonal antibody formation." *Analytical chemistry* 85.20 (2013): 9785-9792.