BHIMADZU

Multi-residue veterinary drug analysis of >200 compounds using MRM Spectrum mode by LC-MS/MS

David Baker¹, Laetitia Fages², Eric Capodanno², Neil Loftus¹ ¹Shimadzu, Manchester, UK ² Phytocontrol, Nimes, France

I. Introduction

Veterinary drugs are used for therapeutic, metaphylactic, prophylactic and growth promotion purposes. To provide an assurance that food from animals is safe with regards to residues of veterinary medicines, regulatory authorities have established Maximum Residue Limits (MRL's) for certain drugs in target tissues and animal species and has also identified pharmacologically active compounds that are prohibited and considered a hazard at any level (EU regulation EC 37/2010; Commission Decision 2003/181/EC; 21CFR Part 556 Tolerances for Residues of New Animal Drugs in Food). In this work, we describe a method that delivers highly sensitive and selective triple quadrupole detection together with MRM Spectrum mode to reduce false positive and false negative reporting. MRM Spectrum mode acquires a high number of fragment ion transitions for each target compound generating a fragmentation spectra that could be used in routine library searching and compound verification using reference library match scores.

2. Materials and Methods

Samples of beef, egg, honey, milk and salmon were extracted and spiked in the calibration range 0.001 – 0.1 mg/kg. Repeatability was assessed at low and high concentrations. Samples were measured using a Nexera UHPLC and the LCMS-8060 triple quadrupole detector (Table 1). Over 200 veterinary drugs were targeted, with more than 2,000 MRM transitions in both ESI +/- during the 12minute gradient.

| Liquid chromatography | | | | |
|---|--|---------------------|---|----------------------------|
| UHPLC Analytical column Column temperature / Flow rate Solvent A Solvent B Binary Gradient | 40°C ; 0.4mL/i 0.1% formic ad | nyl (100 x 2.1, 2.7 | • | olution %B 2 Stop |
| Mass spectrometry | | | | |
| Mass spectrometer Pause time/dwell time Polarity switching time Scope | Shimadzu LCMS-8060 1 msec/3 msec Pos/neg switching time set to 5 msec 218 drugs in positive ion mode (including internal standards) 11 drugs in negative ion mode Structure Analytics (in house development tool) | | | |
| Source temperatures (interface; heat block; DL) Gas flows (nebulising; heating; drying) | 350°C; 300°C; 150°C 3L/min; 10 L/min; 10L/min | | | |

Table 1. LC and MS/MS acquisition parameters used to create the LC-MS/MS method.

2.1 Advantages of MRM Spectrum mode

The method is straightforward to set-up using conventional MRM optimization procedures and acquisition windows (scheduled MRM) resulting in high data densities and a high data sampling rate across a peak. This approach generates a consistent loop time and sampling rate producing reliable quantitation and peak integration. It also enables greater flexibility in routine veterinary drug monitoring programs by supporting a change in the qualifier and quantifier ion selection. The number of precursor-fragment ion transitions used to generate a product ion spectrum is only limited by the chemical structure of the veterinary drug.

Figure 2. MRM reference spectrum for chlortetracycline with assigned fragment structures. MRM Spectrum mode combines MRM with the generation of a product ion spectrum. The product ion spectrum can be used for compound identification by searching a library. As the collision energy was optimized for each fragment ion to generate a product ion spectrum, the library spectrum is highly specific and selective. (Each precursor-fragment ion transitions structure may be verified by using Structure Analytics to show commonly described losses and charge migration; the hydrogen deficit is shown in brackets).






Figure 1. MRM Spectrum mode was used to acquire a high number of fragment ions for each veterinary drug target. For chlortetracycline, 11 precursor-fragment ion transitions were acquired using optimized collision energies. By acquiring a high number of fragment ion transitions, each target veterinary drug had a corresponding fragmentation spectra which could be used in routine library searching and compound verification using reference library match scores. (Chlortetratcycline is a tetracycline class of antimicrobial. According to the Sixth ESVAC report published in 2016, of the overall sales of antimicrobials in the 29 EU countries in 2014, the largest amounts, expressed as a proportion of mg/PCU, were accounted for by tetracyclines (33.4 %), penicillins (25.5%) and sulfonamides (11.0 %). Chlortetracycline was selected as a representative target).





177.00 162.10

Figure 4. Virginiamycin S1 was spiked into an extract of salmon at a concentration of 10pg/uL and MRM Spectrum mode was used to quantify and identify 212 drugs (the method included 2,009 repeatedly injected (n=10). The library score was above 99 in all injections (injection 1, 5 and 10 shown MRM transitions) without compromising limits of detection, linearity or repeatability compared to above) using 11 MRM's to generate a product ion spectrum (Structure Analytics was used to propose a conventional 2MRM method. fragment structures; the hydrogen deficit is shown in brackets). %RSD for oxytetracycline, sulfadimethoxine, ormetoprim, and virginiamycin spiked into salmon (n=10; 10pg/uL) acquired using a conventional 2 MRM method compared to a MRM spectrum method with a higher number of fragment ions. [Compound selection based on antibiotics used in aquaculture;



https://www.ncbi.nlm.nih.gov/pmc/articles/PMC4254636/].

pesticide was dependent upon the chemical structure with typically more than 10 fragment ions searching and compound verification and identification.