

# Analysis of the California list of Pesticides and Mycotoxins in Edibles

Nathaly Reyes-Garcés; Colton Myers; Ashlee Gerardi

Restek Corporation


## Introduction

- The use of cannabis for medicinal and/or recreational purposes has become legal in several states. Regulations that permit the use of different forms of cannabis demand effective and reliable analytical strategies to ensure the safety of cannabis users.
- Pesticide content is one of the main parameters tested in cannabis and cannabis-derived products due to the risks that these compounds pose for human health.
- The main challenges associated with pesticide testing rely on the broad range of physicochemical properties of these compounds, the low action levels requested by the regulations, and the complexity and diversity of matrices to be analyzed.
- The purpose of this work is to present sample preparation and instrumental strategies for the accurate quantitation of the California list of pesticides and mycotoxins in cannabis products.

## Goal

To provide an effective workflow for the simultaneous analysis of the California list of pesticides and mycotoxins in brownies using both LC-MS/MS and GC-MS/MS for instrumental analysis.

## Method Development: Sample Preparation



**Figure 1.** Sample preparation workflow

**For LC-MS/MS analysis:** take 750  $\mu$ L of extract and mix it with 250  $\mu$ L of water. Inject.

**For GC-MS/MS analysis:**

- Transfer the remaining extract to a tube with dSPE sorbents (magnesium sulfate + PSA) (cat.# 26215).
- Vortex and centrifuge. Dilute 500  $\mu$ L of this extract with 500  $\mu$ L of 1% acetic acid in acetonitrile. Inject.

## Method Development: LC/GC-MS/MS


**Table 1.** LC-MS/MS conditions (ionization: ESI)

<b>Column</b>	Raptor ARC-18 2.7 $\mu$ m, 100 mm x 2.1 mm (Restek Cat.# 9314A12)			
<b>Guard Column</b>	Raptor ARC-18 EXP Guard Column Cartridge 2.7 $\mu$ m, 5 x 2.1 mm (cat.# 9314A0252)			
<b>Mobile Phase A</b>	Water, 2 mM ammonium formate, 0.1% formic acid			
<b>Mobile Phase B</b>	Methanol, 2 mM ammonium formate, 0.1% formic acid			
<b>Time Program</b>	Time (min.)	%B	Time (min.)	%B
	0	5	10.5	100
	1.5	65	10.6	5
	8.5	95	12.0	5
	9.5	100		
<b>Other parameters</b>	Column T: 40 $^{\circ}$ C; autosampler T: 10 $^{\circ}$ C; injection volume: 2 $\mu$ L			
<b>Instrument</b>	Shimadzu LCMS-8060			


**Table 2.** GC-MS/MS conditions (ionization: EI)

<b>Column</b>	Rxi-5ms (cat.# 13423)
<b>Injection</b>	Splitless, 1 $\mu$ L (0.5 min splitless time, 14 mL/min split flow)
<b>Liner</b>	Topaz 4.0 mm ID Single Taper Inlet Liner w/ Wool (cat.# 23447)
<b>Inj. T</b>	250 $^{\circ}$ C
<b>Purge Flow</b>	5 mL/min
<b>Oven</b>	90 $^{\circ}$ C (hold 1 min) to 310 $^{\circ}$ C (hold 10 min) by 25 $^{\circ}$ C/min
<b>Carrier Gas</b>	He, at a constant flow of 1.4 mL/min
<b>Transfer line T</b>	290 $^{\circ}$ C
<b>Source T</b>	330 $^{\circ}$ C
<b>Instrument</b>	Thermo Trace 1310-TSQ 8000

## Results and Discussion



**Figure 2.** LC-MS/MS amenable pesticides



**Figure 3.** GC-MS/MS amenable pesticides

As seen in Table 3, the proposed methodology showed satisfactory results in the quantification of all the pesticides and mycotoxins regulated by the state of California. The LOQ values obtained were significantly below the action levels established by the state of CA in cannabis goods (not inhalable).

**Table 3.** Figures of merit corresponding to pesticides and mycotoxins analyzed in brownies.

Compound	LOQ ng/g	R <sup>2</sup>	100 ng/g Acc. Precision %	Compound	LOQ ng/g	R <sup>2</sup>	100 ng/g Acc. Precision %
Daminozide	25	0.9954	102	Spinosyn A (71%)	3.5	0.9994	102
Acephate	10	0.9944	104	Diazinon	5	0.9995	101
Thiamethoxam	5	0.9979	106	Coumaphos	5	0.9997	102
Methomyl	5	0.9996	104	Clofentezine	5	0.9997	103
Oxamyl	5	0.9986	104	Spinosyn D (29%)	1.5	0.999	101
Imidacloprid	10	0.9979	103	Spinosyn J (80%)	4	0.9991	105
Dimethoate	5	0.9994	101	Spinosyn L (20%)	1	0.9991	100
Acetamiprid	5	0.9991	103	Trifloxystrobin	5	0.9997	102
Thiacloprid	5	0.9993	106	Prallethrin	25	0.9996	98
Aldicarb	5	0.9988	99	Hexythiazox	5	0.9996	104
Naled	25	0.9962	105	Cyfluthrin	50	0.9988	107
Mevinphos I (79%)	4	0.9991	104	Pyrethrin I (54%)	5.4	0.9998	104
Mevinphos II (21%)	2	0.9981	106	Pyrethrin II (34%)	26	0.999	100
Carbofuran	5	0.9994	105	Etoazole	5	0.9998	102
Carbaryl	5	0.9997	103	Piperonyl Butoxide	5	0.9998	103
Dichlorvos	5	0.9949	101	Chlorpyrifos	5	0.9994	101
Propoxur	5	0.9993	106	Permethrin-cis (41%)	4.1	0.9994	105
Chloanthraniliprole	10	0.9992	105	Permethrin-trans (59%)	5.9	0.9999	106
Imazalil	5	0.9993	100	Fenpyroximate	5	0.9998	103
Metalaxyl	5	0.9996	103	Bifenthrin	5	0.9994	103
Azoxystrobin	5	0.9998	103	AbamectinB1a	10	0.9999	105
Myclobutanil	5	0.9997	104	Cypermethrin	25	0.9991	93
Phosmet	5	0.9997	103	Etofenprox	5	0.9995	107
Spiroxamine	5	0.9987	102	Pyridaben	10	0.9989	101
Fenoxycarb	5	0.9995	103	Acequinocyl	5	0.9987	103
Methiocarb	5	0.9997	104	Fonicamid	10	0.9993	101
Spiromesifen	25	0.9994	103	Fipronil	10	0.9993	102
Boscalid	5	0.9998	106	Fludioxonil	5	0.9995	104
Paclobutrazol	5	0.9996	103	Captan (GC)	10	0.9941	103
Malathion	5	0.9995	102	Chlordane (GC)	25	0.9939	106
Dimethomorph I (39%)	4	0.9994	101	Chlorfenapyr (GC)	25	0.9953	102
Dimethomorph II (61%)	3	0.9994	103	Methyl parathion (GC)	5	0.9976	103
Tebuconazole	5	0.9996	104	PCNB (GC)	5	0.9975	103
Bifenazate	5	0.9999	104	Cyfluthrin (GC)	5	0.9983	103
Fenhexamid	10	0.9992	103	Cypermethrin (GC)	10	0.9986	103
Propiconazole	5	0.9997	105	Aflatoxin G2	5	0.9987	104
Spirotetramat	5	0.9990	104	Aflatoxin G1	5	0.9984	96
Ethoprophos	5	0.9997	103	Aflatoxin B2	5	0.9996	99
Kresoxym-methyl	5	0.9993	104	Ochratoxin A	10	0.9943	112
Category I pesticides, LOQ < or = to 100 ng/g				Aflatoxin B1	5	0.999	96

Since this methodology only uses 3.5 mL of extraction solvent per sample, a significant reduction in solvent usage/waste is also possible.

## Conclusions

A simple analytical workflow involving SPE and dSPE in combination with LC-MS/MS and GC-MS/MS analysis was proven to be effective in the quantitation of CA pesticides and mycotoxins regulated in cannabis. Satisfactory results in terms of linearity, accuracy, and precision were attained for all the target compounds at the three evaluated concentration levels (10, 100 and 500 ng/g).

## References

- N. Reyes and C. Myers. Analysis of pesticides and mycotoxins in cannabis brownies. Technical article: <https://www.restek.com/pdfs/FFAN3149-UNV.pdf>