

Introduction

Per/Polyfluoroalkyl substances (PFAS), are compounds that have uniquely desirable properties for use in various industries. However, their wide-ranging use leads to emission into the environment, and as PFAS are persistent and bioaccumulate in the environment and wildlife, they are contaminants of concern. Monitoring PFAS precursors present in an environmental sample may impact decisions in treatment processes at remediation sites and help deduce possible degradation products that could exist in the environment. Consequently, scientists are contributing newly identified PFAS structures and spectra to various publicly available databases: growing the list of precursors and degradation by-products, some listing thousands of PFAS.

There are estimated to be >4,000 PFAS in the environment while analytical standards for only around ~100 are readily available. Thus, we need strategies for screening and identifying a larger array of PFAS using high-resolution accurate mass (HRAM) mass spectrometry like quadrupole time of flight (Q/TOF) and sophisticated yet user friendly software tools. This study summarizes some approaches to screen and identify new and emerging PFAS on an LC-Q/TOF instrument without the need for every analytical standard.

Quantification with HRMS

Quantification of a target list of PFAS using a LC-Q/TOF instrument is possible due to its improved sensitivity and robustness. In this study, a 5 mL sample of water was diluted with 5 mL of methanol, acidified with acetic acid, filtered and centrifuged. The a 30 µL injection of the 50/50 methanol/water was injected to the Agilent 1290 Infinity II LC coupled to a 6546 Q-TOF (Santa Clara, USA) to quantify 22 PFAS across multiple classes. Fig 1 denoted the chromatogram of the 22 PFAS.




Fig 1. Chromatogram of 22 PFAS analyzed on LC-Q/TOF




Fig 2. Low level EICs and calibration curves for PFBA, PFOS and PFOA

Screening Strategies for PFAS

Using Libraries & Databases

One of the most common ways to use a LC-Q/TOF and screen for compounds including PFAS is to use personal compound databases (PCD) or libraries (PCDL). There are several different forms of this but typically a PCD is defined as a group of compounds (PFAS in this case) where information including molecular mass, formula and optionally retention time is provided. A PCDL includes the information included in a PCD but in addition also has collected MS/MS spectra to provide fragment information thus adding additional specificity and confidence in identification. These can be commercially curated and available or user-collected.




Fig 3. Information available in the Agilent PCDL for PFHxA

Increasing scope with Metabolite and Transformation Products

PFAS are stable in the environment generally but it is possible for them to degrade and transform in the environment especially under harsh conditions or during certain weathering and water treatment processes. To account for this, it can be useful to increase the scope of screening strategies by adding metabolites and transformation products of PFAS. Here, we use a software available in the Agilent MassHunter package that can automatically calculate biotransformation products of a compound given a chemical formula.




Fig 4. Agilent Biotransformation software that calculated multiple possible PFOS metabolites for increasing the scope of screening

Identification of 6:2 FTS in a wastewater effluent sample using the database. The mass accuracy was <2 ppm but also, the isotopic fidelity match (shown in inset) was excellent compared to theoretical spectra. Subsequently, an analytical standard was procured and the RT was an exact match too thus adding multiple layers for increased confidence in ID.




Fig 5. Identification of 6:2 FTS in a wastewater effluent sample with good match to the library for mass and isotope fidelity along with RT confirmation with the standard

Retention Time Prediction for Increased Confidence

Adding retention time information is another parameter that increases confidence in ID. However, PFAS standards are very limited and so make RT collection for many very difficult. A model to predict RTs for a wide range of PFAS was developed by first projecting RTs from a validated method analyzing a chemically diverse set of PFAS (Figure 6A), to increase the model training set size. Measured and projected RTs were then regressed on predicted physicochemical properties, including LogP (Figure 6C) and LogS (Figure 6D) as well as the number of -CF₂- subunits (Figure 6B) in the chemical structure. A weighted average of predicted RTs from individual regressions was calculated.




Fig 7. Suspect screening results for Perfluoro(2-ethoxyethane) sulfonic acid using a PCD and RT prediction to increase confidence of ID without need for a standard.

A Personal Compound Database (PCD) containing predicted RTs was curated from US EPA's PFAS Inventory List. These suspect PFAS were appended to the MassHunter Quantitative 10.0 data analysis method, and screened for. The Screening Summary Report is shown in Figure 5 for a selected sample.

Perfluoro (2-ethoxyethane) sulfonic acid had a predicted retention time of 4.56 minutes and was putatively detected at 3.78 minutes in this sample, which is within an expected window given error in the prediction model. This suspect must be subject to further inspection to determine if the putative identification is real. However, the predicted retention time, accurate mass results, stable isotope pattern provides much data for this inspection.




Fig 6. Correlation of other PFAS RTs and chemical properties to RTs measured in this study

Identifying Homologous series using Repeating Units & Kendrick Mass Defect

Many PFAS are separated by a -CF₂ unit and this is a critical tool in identifying whole series of PFAS in a sample by knowing just one homologue. Since these compounds exhibit a negative mass defect due to the F atom, and making use of the repeating unit of 44 Da; software like the Agilent Mass Profiler Professional can identify repeating series with specified mass units assisting in ID of series of PFAS in a sample.




Fig 8. Agilent Mass Profiler Professional identifying series of PFAS using repeating units of -CF₂ in the sample.




Fig 9. Agilent MassHunter workflow utilizing Kendrick Mass Defect to identify specific features in the sample followed by formula generation

In Silico Fragment Predictors for Structure Identification of Unknown PFAS

Presently, there are >4000 PFAS thought to be created and many thousands of them are still unknown. Further, The availability of PFAS standards of a vast majority of them is non-existent. So, without standards, putative identifications rely on knowledge of fragment annotation and physicochemical properties. Software, such as Agilent's Molecular Structure Correlator can be used to correlate MS/MS spectra with chemical structures.