



# Agilent Sample Preparation The Pesticide Analysis Workflow

---

*Insights into Preparing your Cannabis Sample and Triple  
Quad Mass Spec Analysis*

Joan Stevens, Ph.D.

Sample Preparation Application Specialist

# Outline

- Introduction to Plant Material and Edibles
- Existing Sample Prep Techniques
  - Extraction, QuEChERS
  - SPE, dispersive SPE
- Evaluation of Extraction
  - Data from Initial Studies
  - Insights and observations
- Optimization Sample Preparation Techniques
  - Sample extraction, cleanup and dilution
  - EMR-Lipid dispersive SPE
- Future Investigations, Aspects and Considerations

***Agilent products and solutions are intended to be used for cannabis quality control and safety testing in laboratories where such use is permitted under state/country law.***

# Cannabis and Cannabis-Based Products: Pesticide Analysis

No tolerances have been established for marijuana, because of its illegal federal status and because the pesticide companies have yet to embark on the lengthy and expensive process of testing their products on cannabis..... *Cannabis Now Issue 19*



# Pesticide Analytes and their action levels in Oregon

Range 0.1-1 ppm or 100-1000 ppb




Analyte	Chemical Abstract Services (CAS) Registry number	Action level ppm	Analyte	Chemical Abstract Services (CAS) Registry number	Action level ppm
Abamectin	71751-41-2	0.5	Imazalil	35554-44-0	0.2
Acephate	30560-19-1	0.4	Imidacloprid	138261-41-3	0.4
Acequinocyl	57960-19-7	2	Kresoxim-methyl	143390-89-0	0.4
Acetamiprid	135410-20-7	0.2	Malathion	121-75-5	0.2
Aldicarb	116-06-3	0.4	Metalaxyl	57837-19-1	0.2
Azoxystrobin	131860-33-8	0.2	Methiocarb	2032-65-7	0.2
Bifenazate	149877-41-8	0.2	Methomyl	16752-77-5	0.4
Bifenthrin	82657-04-3	0.2	Methyl parathion	298-00-0	0.2
Boscalid	188425-85-6	0.4	MGK-264	113-48-4	0.2
Carbaryl	63-25-2	0.2	Myclobutanil	88671-89-0	0.2
Carbofuran	1563-66-2	0.2	Naled	300-76-5	0.5
Chlorantraniliprole	500008-45-7	0.2	Oxamyl	23135-22-0	1
Chlorfenapyr	122453-73-0	1	Paclobutrazol	76738-62-0	0.4
Chlorpyrifos	2921-88-2	0.2	Permethrins*	52645-53-1	0.2
Clofentezine	74115-24-5	0.2	Phosmet	732-11-6	0.2
Cyfluthrin	68359-37-5	1	Piperonyl_butoxide	51-03-6	2
Cypermethrin	52315-07-8	1	Prallethrin	23031-36-9	0.2
Daminozide	1596-84-5	1	Propiconazole	60207-90-1	0.4
DDVP (Dichlorvos)	62-73-7	0.1	Propoxur	114-26-1	0.2
Diazinon	333-41-5	0.2	Pyrethrins†	8003-34-7	1
Dimethoate	60-51-5	0.2	Pyridaben	96489-71-3	0.2
Ethoprophos	13194-48-4	0.2	Spinosad	168316-95-8	0.2
Etofenprox	80844-07-1	0.4	Spiromesifen	283594-90-1	0.2
Etoxazole	153233-91-1	0.2	Spirotetramat	203313-25-1	0.2
Fenoxycarb	72490-01-8	0.2	Spiroxamine	118134-30-8	0.4
Fenpyroximate	134098-61-6	0.4	Tebuconazole	80443-41-0	0.4
Fipronil	120068-37-3	0.4	Thiacloprid	111988-49-9	0.2
Flonicamid	158062-67-0	1	Thiamethoxam	153719-23-4	0.2
Fludioxonil	131341-86-1	0.4	Trifloxystrobin	141517-21-7	0.2
Hexythiazox	78587-05-0	1			

\* Permethrins should be measured as cumulative residue of cis- and trans-permethrin isomers (CAS numbers 54774-45-7 and 51877-74-8).




† Pyrethrins should be measured as the cumulative residues of pyrethrin 1, cinerin 1 and jasmolin 1 (CAS numbers 121-21-1, 25402-06-6, and 4466-14-2 respectively).

# Constituents of Cannabis Plants: Complex

- **Nitrogenous compounds (27 known)**
  - Amino acids (18),
  - Proteins (3)
  - Glycoproteins (6)
  - Enzymes (2)
- **Sugars and related compounds (34)**
  - Hydrocarbons (50)
  - Simple alcohols (7)
  - Aldehydes (13)
  - Ketones (13)
  - Simple acids (21)
- **Fatty acids (22)**
  - Simple esters (12)
  - Lactones (1)
  - Steroids (11)
- **Terpenes (120)**
- **Non-cannabinoid phenols (25)**
- **Cannabinoids (66)**
- **Flavonoids (21)**
  - Vitamins (1) [Vitamin A]
  - Pigments (2)
  - Elements (9).




## Structure of Terpenes



# Sample Preparation: Pesticide Analysis

- Complex matrix associated with cannabis plant material and edibles needs to be addressed
  - CBDs are in large amounts (10-20%), THCA can interfere with analysis with broad interference (100,000-200,000 ppm)
  - Terpenes and other non-cannabinoid compounds are also in large ppm quantities (10-5000 ppm)
  - Pesticides in 500 ppb amounts (0.00005% )



LC/MSMS analysis after QuEChERS and dSPE cleanup (universal)

# Sample Preparation Techniques: Pesticide Analysis

- Methanol extraction with C18 SPE
- Acetonitrile extraction with C18 SPE
- Acetonitrile extraction with both NH<sub>2</sub> and C18 SPE
- QuEChERS extraction with SPE cleanup
- QuEChERS extraction with dispersive SPE cleanup
  - Plant material and edibles considered dry matrix; < 60% water
  - Addition of water required for QuEChERS extraction/partitioning: step 1
  - Use of Salts: Na Acetate, Citrates, NaCl with MgSO<sub>4</sub>
  - Super-saturates the water with the salts allowing separation of water from the ACN
  - Analytes of interest extract/transfer into the ACN layer
  - Clean-up of co-extractive matrix: step 2




# Basic Protocol: Pesticide Analysis



Grind flower or plant material (weigh 1.5 g)  
SPEX GenoGrinder, 2 Ceramic Homogenizers, 5 min at 1500 rpm

Individual grinding in 50 mL disposable centrifuge tubes eliminates washing, carryover and cross contamination

2010  
Geno/Grinder



Dispersive SPE: Universal or All Food Types: (2 mL) 50 mg PSA, 50 mg C18, 7.5 mg GCB, 150 mg MgSO<sub>4</sub> (5982-0028)  
(15 mL) 400mg PSA, 400 mg C18, 45 mg GCB, 1200 MgSO<sub>4</sub> (5982-0029)




# Basic Protocol: Pesticide Analysis



Grind flower or plant material (weigh 1.5 g)  
SPEX GenoGrinder, 2 Ceramic Homogenizers, 5 min at 1500 rpm

Individual grinding in 50 mL disposable centrifuge tubes eliminates washing, carryover and cross contamination

2010  
**Geno/Grinder**



Dispersive SPE: Universal or All Food Types: (2 mL) 50 mg PSA, 50 mg C18, 7.5 mg GCB, 150 mg MgSO<sub>4</sub> (5982-0028)  
(15 mL) 400mg PSA, 400 mg C18, 45 mg GCB, 1200 MgSO<sub>4</sub> (5982-0029)


# Basic Protocol: Pesticide Analysis



Grind flower or plant material (weigh 1.5 g)  
SPEX GenoGrinder, 2 Ceramic Homogenizers, 5 min at 1500 rpm

Individual grinding in 50 mL disposable centrifuge tubes eliminates washing, carryover and cross contamination



2010  
Geno/Grinder



Dispersive SPE: Universal or All Food Types: (2 mL) 50 mg PSA, 50 mg C18, 7.5 mg GCB, 150 mg MgSO<sub>4</sub> (5982-0028)  
(15 mL) 400mg PSA, 400 mg C18, 45 mg GCB, 1200 MgSO<sub>4</sub> (5982-0029)

# Optimization and Custom Dispersive SPE: Why is it Important

Step 1:



QuEChERS Extraction/Partitioning  
or  
ACN (1% acetic acid) extraction

- Advantages or Disadvantages
- Addition of water
- Addition of 1% acetic acid



Step 2:

Clean-up: Dispersive SPE

- Advantage over SPE, less steps, no manifold, vacuum
- Enhanced (custom) clean-up or additional clean-up

# Optimization and Custom Dispersive SPE: Why is it Important

Step 1:



QuEChERS Extraction/Partitioning  
or  
ACN (1% acetic acid) extraction

- Advantages or Disadvantages
- Addition of water
- Addition of 1% acetic acid



Step 2:

Clean-up: Dispersive SPE

- Advantage over SPE, less steps, no manifold, vacuum
- Enhanced (custom) clean-up or additional clean-up

# Optimization and Custom Dispersive SPE: Why is it Important


Step 1:



QuEChERS Extraction/Partitioning  
or  
ACN (1% acetic acid) extraction

- Advantages or Disadvantages
- Addition of water
- Addition of 1% acetic acid

Step 2:





Clean-up: Dispersive SPE


- Advantage over SPE, less steps, no manifold, vacuum
- Enhanced (custom) clean-up or additional clean-up

# Optimization and Custom Dispersive SPE: Why is it Important

Step 1:



Step 2:



QuEChERS Extraction/Partitioning  
or  
ACN (1% acetic acid) extraction


- Advantages or Disadvantages
- Addition of water
- Addition of 1% acetic acid

Clean-up: Dispersive SPE


- Advantage over SPE, less steps, no manifold, vacuum
- Enhanced (custom) clean-up or additional clean-up

# Optimization and Custom Dispersive SPE: Why is it Important

Step 1:



Step 2:



QuEChERS Extraction/Partitioning  
or  
ACN (1% acetic acid) extraction


- Advantages or Disadvantages
- Addition of water
- Addition of 1% acetic acid

Clean-up: Dispersive SPE


- Advantage over SPE, less steps, no manifold, vacuum
- Enhanced (custom) clean-up or additional clean-up

# Optimization and Custom Dispersive SPE: Why is it Important

Step 1:



Step 2:



QuEChERS Extraction/Partitioning  
or  
ACN (1% acetic acid) extraction

- Advantages or Disadvantages
- Addition of water
- Addition of 1% acetic acid



Clean-up: Dispersive SPE

- Advantage over SPE, less steps, no manifold, vacuum
- Enhanced (custom) clean-up or additional clean-up




# Optimization and Custom Dispersive SPE: Why is it Important

Step 1:



QuEChERS Extraction/Partitioning  
or  
ACN (1% acetic acid) extraction

- Advantages or Disadvantages
- Addition of water
- Addition of 1% acetic acid



Step 2:





Clean-up: Dispersive SPE


- Advantage over SPE, less steps, no manifold, vacuum
- Enhanced (custom) clean-up or additional clean-up

# Optimization and Custom Dispersive SPE: Why is it Important

Step 1:



Step 2:



QuEChERS Extraction/Partitioning  
or  
ACN (1% acetic acid) extraction


- Advantages or Disadvantages
- Addition of water
- Addition of 1% acetic acid

Clean-up: Dispersive SPE

- Advantage over SPE, less steps, no manifold, vacuum
- Enhanced (custom) clean-up or additional clean-up


# Optimization and Custom Dispersive SPE: Why is it Important

Step 1:




QuEChERS Extraction/Partitioning  
or  
ACN (1% acetic acid) extraction

- Advantages or Disadvantages
- Addition of water
- Addition of 1% acetic acid



Step 2:




Clean-up: Dispersive SPE

- Advantage over SPE, less steps, no manifold, vacuum
- Enhanced (custom) clean-up or additional clean-up


# Optimization and Custom Dispersive SPE: Why is it Important

Step 1:




QuEChERS Extraction/Partitioning  
or  
ACN (1% acetic acid) extraction

- Advantages or Disadvantages
- Addition of water
- Addition of 1% acetic acid



Step 2:





Clean-up: Dispersive SPE


- Advantage over SPE, less steps, no manifold, vacuum
- Enhanced (custom) clean-up or additional clean-up

# Optimization and Custom Dispersive SPE: Why is it Important

Step 1:



Step 2:




QuEChERS Extraction/Partitioning  
or  
ACN (1% acetic acid) extraction

- Advantages or Disadvantages
- Addition of water
- Addition of 1% acetic acid

Clean-up: Dispersive SPE

- Advantage over SPE, less steps, no manifold, vacuum
- Enhanced (custom) clean-up or additional clean-up

# Cannabis mixture cleanup – QuEChERS with 1x dSPE vs. 2 x dSPE custom dispersive steps



GC MS 7010 Triple Quad TIC Scan, 0.5  $\mu$ L injection, Gain 10  
Custom dispersive: Mixture of PSA, C18-EC, GCB, **additional sorbent**, and MgSO<sub>4</sub>

# Sample Preparation: Insights and Observations

## Step 1: QuEChERS or ACN Extraction


- Addition of Water, vortex 15-30 minutes



5 min



10 min



20 min with indicator pH range




5 min with indicator pH range

Water 50 mM Na  
Acetate pH 3.5

Water/ACN (1:1 v/v)


Addition of water to cannabis flower/plant turns basic almost immediately

Why is this a concern.....degradation of base sensitive pesticides (fungicides)  
Captan, Folpet, Dichlofluanid, Chlorothalonil, Dicofol, Tolyfluanid



# Recovery Comparison using Custom dSPE with and without Water

Custom dSPE Recovery Experimental Recovery Data



EXP 1a: AOAC QuEChERS: ACN (1%AA) 13.75 mL, dSPE custom, 1xdSPE


EXP 1b: AOAC QuEChERS: ACN (1%AA) 13.75 mL, dSPE custom, 2xdSPE

EXP 2a: ACN (1%AA) 13.75 mL only no QuEChERS salts and no H2O, dSPE custom, 1xdSPE


EXP 2b: ACN (1%AA) 13.75 mL only no QuEChERS salts and no H2O, dSPE custom, 2xdSPE




# QuEChERS AOAC Extraction/Partitioning with Custom dSPE versus Universal dSPE: Hops



# QuEChERS AOAC Extraction/Partitioning with Custom F dSPE versus Universal dSPE: Hops



# Pesticide Recovery after QuEChERS AOAC with Custom F dSPE or Universal dSPE: Hops




Concentration of spiked pesticides: 100 ppb

Daminozide is not extracted into the ACN layer after QuEChERS AOAC salt partitioning

# Pesticide Recovery after QuEChERS AOAC with Custom F dSPE or Universal dSPE: Hops

QuEChERS AOAC with Custom F dSPE versus Universal dSPE


Custom Dispersive SPE Did Not Improve Matrix Removal in Cannabis



Concentration of spiked pesticides: 100 ppb


Daminozide is not extracted into the ACN layer after QuEChERS AOAC salt partitioning

# Pesticide Recovery after QuEChERS AOAC with Custom F dSPE or Universal dSPE: Hops




Concentration of spiked pesticides: 100 ppb

# Pesticide Recovery after QuEChERS AOAC with Custom F dSPE or Universal dSPE: Hops




Concentration of spiked pesticides: 100 ppb

# Pesticide Recovery after QuEChERS AOAC with Custom F dSPE or Universal dSPE: Hops




# Comparison of Cannabis Extract versus Extract with 20 x Dilution

## Full Scan GC/MSMS






# Comparison of Cannabis Extract versus Extract 20 x Dilution with Cleanup for GCMSMS



# Overlay of Cannabis Extracts Full Scan GC/MSMS



# Dilution better for Quantitation & Recovery?


Analyte	No dilution	1:5 Dilution	1:10 Dilution	1:20 Dilution	1:50 Dilution	1:100 Dilution
Acetamidiprid	47.1 ± 3.7	88.5 ± 2.3	94.6 ± 3	99.2 ± 6.8	104.5 ± 6.4	116.7 ± 5.9
Aflatoxin B1	31.1 ± 0.9	80 ± 1.3	90 ± 3.8	92.1 ± 7.7	110.2 ± 8.8	112 ± 3.4
Buprofezin	4.1 ± 0.3	16.6 ± 1.1	32.1 ± 1.2	47.3 ± 5.8	84.6 ± 9	109.8 ± 9.9
Chlorantraniliprole	8.2 ± 0.7	31.1 ± 1.1	50.9 ± 5.2	65.5 ± 10.9	94.2 ± 10.2	101.9 ± 20.4
Lenacil	16.9 ± 1.8	57.7 ± 5.9	75.7 ± 5.4	82.4 ± 3.6	97.5 ± 9.2	106.9 ± 11.5
Methomyl	24.9 ± 1.9	51.3 ± 3.5	65.3 ± 6.5	77.1 ± 7.5	107.3 ± 10.7	113.3 ± 14.6
Profenofos	2.6 ± 0.3	12.5 ± 0.6	20.4 ± 2.1	34.2 ± 1.7	59 ± 8.3	88.5 ± 15.5
Propamocarb	105.9 ± 8	102.4 ± 3.4	101 ± 3	100.2 ± 2.8	114.5 ± 3.4	111.5 ± 6.8
Proquinazid	11.4 ± 0.6	42.8 ± 0.9	57.9 ± 1	68.3 ± 3.9	96 ± 11.6	111.3 ± 6.8
Sudan red 7B	19.8 ± 2.5	57.4 ± 4.2	61 ± 6.9	70.5 ± 2.4	84 ± 5	81.5 ± 7.2

Compound recoveries in black pepper extremely complex matrix

# Dilution better for Quantitation & Recovery?

Analyte	No dilution	1:5 Dilution	1:10 Dilution	1:20 Dilution	1:50 Dilution	1:100 Dilution
Acetamiprid	47.1 ± 3.7	88.5 ± 2.3	94.6 ± 3	99.2 ± 6.8	104.5 ± 6.4	116.7 ± 5.9
Aflatoxin B1	31.1 ± 0.9	80 ± 1.3	90 ± 3.8	92.1 ± 7.7	110.2 ± 8.8	112 ± 3.4
Buprofezin	4.1 ± 0.3	16.6 ± 1.1	32.1 ± 1.2	47.3 ± 5.8	84.6 ± 9	109.8 ± 9.9
Chlorantraniliprole	8.2 ± 0.7	31.1 ± 1.1	50.9 ± 5.2	65.5 ± 10.9	94.2 ± 10.2	101.9 ± 20.4
Lenacil	16.9 ± 1.8	57.7 ± 5.9	75.7 ± 5.4	82.4 ± 3.6	97.5 ± 5.5	100.5 ± 11.7
Methomyl	24.9 ± 1.9	51.3 ± 3.5	65.3 ± 6.5	77.1 ± 7.5	107.3 ± 11.5	114.5 ± 11.5
Profenofos	2.6 ± 0.3	12.5 ± 0.6	20.4 ± 2.1	34.2 ± 1.7	59 ± 8.5	77.1 ± 11.5
Propamocarb	105.9 ± 8	102.4 ± 3.4	101 ± 3	100.2 ± 2.8	114.5 ± 11.5	114.5 ± 11.5
Proquinazid	11.4 ± 0.6	42.8 ± 0.9	57.9 ± 1	68.3 ± 3.9	96 ± 11	101.9 ± 20.4
Sudan red 7B	19.8 ± 2.5	57.4 ± 4.2	61 ± 6.9	70.5 ± 2.4	84 ± 5	101.9 ± 20.4

Compound recoveries in black pepper extremely complex matrix



Agilent 6495B: LC/MSMS

1:100 dilution required for acceptable recoveries >90% of the compounds

# Recovery of Problematic Pesticides from Cannabis Extraction with 20 x Dilution LC/MSMS

Compound	Recovery%	Average ng/mL
Abamectin	79	0.63
Acephate	86	0.69
Acequinocyl-Hydroxy	87	0.69
Azoxystrobin	88	0.70
Bifenazate	87	0.69
Daminozide	95	0.76
Hexythiazox	78	0.62
Piperonyl butoxide	89	0.72
Prallethrin	92	0.74
Pyrethrin I	88	3.50
Pyrethrin II	91	3.64
Spinosyn A	82	3.64
Spinosyn D	74	0.66
Spirotetramat	82	0.66
Spiroxamine	94	0.75



N=5

# THC and CBD High Fat Products



# Cannabis-Based Products: Unique Dispersive SPE

## EMR-Lipid Innovative Sorbent!



When “activated” by water EMR-Lipid Sorbent Selectively traps lipids.




**Size Exclusion:** Unbranched hydrocarbon chains (lipids) enter the sorbent; bulky analytes do not.

**Sorbent Chemistry:** Lipid chains that enter the sorbent are trapped by hydrophobic interactions.



# Cannabis-Based Products: Unique Dispersive SPE

## EMR-Lipid Innovative Sorbent!



When “activated” by water EMR-Lipid Sorbent Selectively traps lipids.

**Size Exclusion:** Unbranched hydrocarbon chains (lipids) enter the sorbent; bulky analytes do not.

**Sorbent Chemistry:** Lipid chains that enter the sorbent are trapped by hydrophobic interactions.



**EMR-Lipid Mechanism – Size exclusion and hydrophobic interaction**



... and what does it do?





# ... and what does it do?

**EMR sorbent removes LIPIDS it interacts with the long aliphatic chain NOT the functional groups**

What are Lipids?

A class of naturally occurring hydrocarbon containing compounds commonly known as fats and oils







# ... and what does it do?

**EMR sorbent removes LIPIDS it interacts with the long aliphatic chain NOT the functional groups**


What are Lipids?

A class of naturally occurring hydrocarbon containing compounds commonly known as fats and oils




# What Does EMR *NOT* Interact With?

EMR does **NOT** remove analytes of interest




Insecticides and anthelmintics




Imidazole pesticides

Organophosphate Pesticides



Organochlorine Pesticides




Fumonisin B2

Aflatoxins and Mycotoxins



# Cannabis-Based Product: THC-Butter


- 1.5 g of THC-butter was added to a 50 mL centrifuge tube
- Add 15 mL of ACN (1%acetic acid)
- Vortex 30 minutes, centrifuge 5000 rpm, 5 min
- Add 5 mL water to EMR-Lipid, vortex immediately 1 min
- Add 5 mL of ACN (1% acetic acid) extract, vortex immediately 2 min
- Centrifuge 5000 rpm, 5 min
- Analyze by LCMSMS



Pesticide Screened:

Abamectin B1-a
Abamectin B1-b
Azoxystrobin
Bifenazate
Etoxazole
Imazalil
Imidacloprid
Malathion
Myclobutanil
Permethrin-1
Permethrin-2
Spinosyn A
Spinosyn D
Spiromesifen
Spirotetramat
Tebuconazole

# Extraction of Butter: MS1 Scan



Blue Chromatogram: 1.5 g butter, AOAC fatty dSPE (PSA and C18EC)  
Green Chromatogram: 1.5 g butter, AOAC salts, EMR-Lipid dispersive SPE

# EMR-Lipid Cartridge Formate

- SPE format with gravity flow through design
- Accommodate mLs of extract
- Very applicable to LC/MSMS work flow
- Continuing the work with high fat products for pesticide analysis
- Oils, cookies, chocolate based products
- Application work available at product release early fall

# Insights and Observations :



# Insights and Observations :

- The active ingredients can cause interference and instrument issues: cannabinoids and terpenes

# Insights and Observations :

- The active ingredients can cause interference and instrument issues: cannabinoids and terpenes
- If QuEChERS AOAC method is being used you should add the ACN (1% acetic acid immediately after the addition of water, before vortexing for 20-30 minutes

# Insights and Observations :

- The active ingredients can cause interference and instrument issues: cannabinoids and terpenes
- If QuEChERS AOAC method is being used you should add the ACN (1% acetic acid immediately after the addition of water, before vortexing for 20-30 minutes
- Do not use QuEChERS EN method, water and ACN mix will be too basic for labile compounds based on the weaker buffering capacity of citrate salts used in the QuEChERS EN

# Insights and Observations :

- The active ingredients can cause interference and instrument issues: cannabinoids and terpenes
- If QuEChERS AOAC method is being used you should add the ACN (1% acetic acid immediately after the addition of water, before vortexing for 20-30 minutes
- Do not use QuEChERS EN method, water and ACN mix will be too basic for labile compounds based on the weaker buffering capacity of citrate salts used in the QuEChERS EN
- QuEChERS methodology: Salting/Partitioning and dSPE may negatively affect recovery; e.g. Daminozide

# Insights and Observations :

- The active ingredients can cause interference and instrument issues: cannabinoids and terpenes
- If QuEChERS AOAC method is being used you should add the ACN (1% acetic acid immediately after the addition of water, before vortexing for 20-30 minutes
- Do not use QuEChERS EN method, water and ACN mix will be too basic for labile compounds based on the weaker buffering capacity of citrate salts used in the QuEChERS EN
- QuEChERS methodology: Salting/Partitioning and dSPE may negatively affect recovery; e.g. Daminozide
- Promising results with ACN (1% acetic acid) extraction, no QuEChERS salts

# Insights and Observations :

- The active ingredients can cause interference and instrument issues: cannabinoids and terpenes
- If QuEChERS AOAC method is being used you should add the ACN (1% acetic acid immediately after the addition of water, before vortexing for 20-30 minutes
- Do not use QuEChERS EN method, water and ACN mix will be too basic for labile compounds based on the weaker buffering capacity of citrate salts used in the QuEChERS EN
- QuEChERS methodology: Salting/Partitioning and dSPE may negatively affect recovery; e.g. Daminozide
- Promising results with ACN (1% acetic acid) extraction, no QuEChERS salts
- Dispersive SPE or SPE will remove some of the matrix co-extractives

# Insights and Observations :

- The active ingredients can cause interference and instrument issues: cannabinoids and terpenes
- If QuEChERS AOAC method is being used you should add the ACN (1% acetic acid immediately after the addition of water, before vortexing for 20-30 minutes
- Do not use QuEChERS EN method, water and ACN mix will be too basic for labile compounds based on the weaker buffering capacity of citrate salts used in the QuEChERS EN
- QuEChERS methodology: Salting/Partitioning and dSPE may negatively affect recovery; e.g. Daminozide
- Promising results with ACN (1% acetic acid) extraction, no QuEChERS salts
- Dispersive SPE or SPE will remove some of the matrix co-extractives
- Dilution had the greatest effect on removing matrix interferences, consider incorporating in your protocol, require very sensitive triple quadrupole

# Insights and Observations :

- The active ingredients can cause interference and instrument issues: cannabinoids and terpenes
- If QuEChERS AOAC method is being used you should add the ACN (1% acetic acid immediately after the addition of water, before vortexing for 20-30 minutes
- Do not use QuEChERS EN method, water and ACN mix will be too basic for labile compounds based on the weaker buffering capacity of citrate salts used in the QuEChERS EN
- QuEChERS methodology: Salting/Partitioning and dSPE may negatively affect recovery; e.g. Daminozide
- Promising results with ACN (1% acetic acid) extraction, no QuEChERS salts
- Dispersive SPE or SPE will remove some of the matrix co-extractives
- Dilution had the greatest effect on removing matrix interferences, consider incorporating in your protocol, require very sensitive triple quadrupole
- Is the existing clean-up techniques enough for long term analysis, issues with analysis and instrument maintenance, check you sample prep methodology by full scan both LC and GC/MSMS