


Agilent Technologies

Investigating Miniaturization in GPC/SEC



Graham Cleaver
GPC Business Development

What are the advantages of Minaturiztion?



How to Achieve Higher Resolution



- High efficiency GPC/SEC columns packed with small diameter particles
- High pore volume GPC/SEC columns packed with “multipore” particles, with a near-linear molecular weight range appropriate for your samples, that will increase resolution
- Shorter lengths and wider diameters that allow high linear velocity to be used in order to greatly reduce run time and still maintain acceptable resolution over a wide molecular weight range

PLgel 3um / 5um


PL aquagel-OH 5um

PlusPore



RL Rapide

PL Rapide Aqua

High Resolution with Smaller Particles



PlusPore – High Performance GPC Columns




PlusPore calibration curves


PlusPore selection guide

Column	MW range (g/mol) (PS)	Nominal particle size (µm)	Typical efficiency (p/m)	Recommended calibrants	Frit porosity (µm)
PolyPore	200 to 2,000,000	5	>60,000	EasiCal PS-1or EasiVial PS-H	2
ResiPore	up to 500,000	3	>80,000	EasiCal PS-2 or EasiVial PS-M	2
MesoPore	up to 25,000	3	>80,000	Polystyrene S-L-10 Kit,	2
OligoPore	up to 3,300	6	>55,000	Polystyrene S-L2-10 Kit	2

Avoid Dislocations



Comparison of PolyPore with conventional individual pore size GPC columns



- Individual pore size columns can exhibit dislocations where the Mw resolution ranges overlap
- Mixed bed columns have a wide linear range that prevent dislocations

High Resolution OligoPore Separations


Polystyrene Standards

Columns: 2xOligoPore, 300x7.5mm (PL1113-6520)

Eluent: THF


Flow Rate: 1.0ml/min

Detector: RI



- The OligoPore column is filled with 6 μ m packing material, but due to a very large pore volume it gives increased resolution compared to a 3 μ m Plgel column
- As a result of the larger particle size, OligoPore is very resilient to extra-column dispersion

Polystyrene Mw 580 – Oligopore 250x4.6mm




MW Range: up to 3,300 (g/mol)

Nominal Particle Size: 6 μ m

Typical Efficiency: >55,000 p/m


Resolution maintained as speed is increased



Polystyrene Mw 580 resolved using Oligopore 4.6x250mm


Different flow rates overlaid to show that faster doesn't sacrifice resolution. The chromatograms have been normalised to better illustrate the differences

High Speed MesoPore Columns



Conditions

Column: 2 x MesoPore, 4.6 x 250 mm (PL1513-5325)
Sample: Epoxy resin
Eluent: THF
Flow rate: 0.35 and 1.2 mL/min
Inj vol: 4 μ L
System: 1260 Infinity GPC/SEC System, UV, 254 nm



Easy Method Transfer from Standard to rapid GPC on MesoPore 250x4.6mm GPC columns

MW Range: up to 25,000 (g/mol)

Nominal Particle Size: 3 μ m


Typical Efficiency: >80,000 p/m

Very High Speed PL-Rapide Separations

PL Rapide columns reduce analysis times while maintaining the excellent solvent compatibility and mechanical stability of all GPC columns from Agilent

Conditions

Column: PL Rapide L, 100 x 10 mm
Sample: Epoxy resin
Eluent: THF
Flow Rate: 1.0, 2.0 and 3.0 mL/min
Detector: UV, 254 nm




Description	MW range (g/mol)
PL Rapide H, 150 x 7.5 mm	500 to 10,000,000
PL Rapide H, 100 x 10mm	500 to 10,000,000
PL Rapide M, 150 x 7.5 mm	200 to 2,000,000
PL Rapide M, 100 x 10 mm	200 to 2,000,000
PL Rapide L, 150 x 7.5 mm	200 to 500,000
PL Rapide L, 100 x 10 mm	200 to 500,000
PL Rapide F, 150 x 7.5 mm	up to 3,300
PL Rapide F, 100 x 10 mm	up to 3,300



Reproducible Results

Samples	Columns
Agilent EasiVial PS-H and PS-M	Agilent PLgel 10 μ m 10 ⁶ 7.5 x 300 mm
	Agilent PLgel 5 μ m 10 ⁵ 7.5 x 300 mm
	Agilent PLgel 5 μ m 10 ⁴ 7.5 x 300 mm
	Agilent PLgel 5 μ m 10 ³ 7.5 x 300 mm
Kraton	3 x Agilent PL Rapide L, 10 x 100 mm Linear MW range up to 500,000 g/mol
	2 x Agilent ResiPore, 4.6 x 250 mm Linear MW range up to 500,000 g/mol

Columns:	4 x PL gel	3 x PL Rapide L	2 x ResiPore
Mobile phase:	THF		
Flow rate:	1.0 mL/min	1.5 mL/min	0.6 mL/min
Sample conc.:	1.5 mg/mL		
Inj. volume:	100 μ L	5 μ L	2 μ L
Temperature:	30 °C	50 °C	50 °C
Run time:	50 minutes	15 minutes	15 minutes



Consistent Results

methods can be transferred with confidence without the risk of absorbance effects or other interactions between the analytes and stationary phase.

Chromatography and MW results run on PL gel, PL Rapide L and Resipore columns


Columns	Tr (Peak 2)	Rs	N/m (Peak 2)	a	Area (%)	Height (%)
4 x PLgel	28.46	1.22	5653	1.05	8	7
3 x PL Rapide L	7.41	1.13	23727	1.06	7	7
2 x ResiPore	6.66	1.10	14510	1.05	8	8

Average and precision of molecular weight results run on PL gel, PL Rapide L and Resipore columns

	Mp	Mn	Mw	PD	Mp	Mn	Mw	PD
	Peak 1				Peak 2			
Mean	110262	107964	111181	1.03	56041	51015	53208	1.05
%RSD	3	1	1	0	5	8	5	4

Potential for 2DLC/GPC


PS standards UV Chromatograms – 10x50 Mix-E at 6ml/min



Comparison of molecular weight distributions derived from 10 x 50 mm and 7.5 x 600 mm versions of the Agilent PLgel MIXED-E column.


Agilent PLgel dimensions	Mn (g/mole)	Mw (g/mole)	Mz (g/mole)	Mz+1 (g/mole)
10 x 50 mm	1.369	2.938	4.957	6.517
7.5 x 600 mm	1.509	2.912	4.562	5.732

Solvent Usage – Cost Saving



**70% saving in analysis time
and 55% saving in solvent
usage**

**Use a narrower id
column for an 82%
saving on solvent usage**



Summary

Shorter run times

High speed, high resolution separations can be achieved using Multiporous Polystyrene/Divinyl Benzene particles

Save cost with reduced solvent consumption

Smaller i.d. columns and faster flow rates reduces solvent consumption by as much as 80% compared to conventional columns.

No column dislocations or absorption

Linear column calibrations improve peak shape and accuracy of results

Increased performance on an analytical GPC system

Miniaturized high pore volume columns packed with PS/DVB particles operate at low pressure and elevated flow rates

Any Questions?

