

Instrument: CN628

Carbon/Nitrogen in Soil and Plant Tissue

LECO Corporation; Saint Joseph, Michigan USA

Sample Preparation

Samples must be of uniform consistency to produce suitable results. Samples should be ground to pass a 0.5 mm screen.

Carbon and nitrogen results for soil and plant tissue materials are normally reported on a dry basis. The materials can either be dried prior to analysis or the moisture content determined and the values corrected. Please see the note addendum at the end of this document for details on drying these materials.

Accessories

502-186 Tin Foil Cup

Calibration Samples

502-092 EDTA, 502-642 Phenylalanine, 501-050 Nicotinic Acid, or other suitable pure compounds.

Analysis Parameters*

Furnace Temperature	950°C
Afterburner Temperature	850°C

Element Parameters

	Carbon	Nitrogen
Analyze	Yes	Yes
Baseline Delay Time	0 seconds	10 seconds
Minimum Analysis Time	20 seconds	40 seconds
Comparator Level	100.00	100.00
Endline Time	1 seconds	2 seconds
Conversion Factor	1.00	1.00
Significant Digits	5	5
IR Baseline Time	1 second	—
TC Baseline Time	—	10 seconds

Burn Profile

Burn Steps	Time (seconds)	Furnace Flow
1	20 seconds	High
2	150 seconds	Medium
3	30 seconds	High

Ballast Parameters

Equilibrate Time	30 seconds
Not Filled Timeout	300 seconds

Aliquot Loop

Fill Pressure Drop	200 mm Hg
Equilibrate Pressure Time	8 seconds

*Refer to CN628 Operator's Instruction Manual for Method Parameter definitions.

Procedure

1. Prepare instrument for operation as outlined in the operator's instruction manual.
2. Determine blank.
 - a. Enter 1.0000 g mass into Sample Login (F3) using Blank as the sample name.
 - b. Select 10 replicates.
 - c. Initiate the analysis sequence (F5).
 - d. Set the blank using at least 5 results following the procedure outlined in the operator's instruction manual.
 - e. The standard deviation of the last 5 blanks should be less than or equal to 0.001% (10 ppm) for nitrogen. Additional blanks beyond the recommended 10 may need to be analyzed in order to achieve the recommended precision.
3. Calibrate.
 - a. Weigh ~0.15 g of EDTA into a 502-186 Tin Foil Cup and seal.
 - b. Enter sample mass and identification into Sample Login (F3).
 - c. Transfer sample to the appropriate position in the sample carousel.
 - d. Repeat steps 3a through 3c a minimum of five times.
 - e. Initiate the analysis sequence (F5).
 - f. Calibrate the instrument using single standard calibration (fixed at origin) following the procedure outlined in the operator's instruction manual.

Note: A CN628 can be calibrated using several replicates of a single mass range (nominal 0.15 g) of EDTA utilizing a single standard calibration. The calibration can be verified by analyzing a pure compound that is different than the material used for calibration, such as phenylalanine (~0.1 g) or nicotinic acid (~0.1 g). Multi-point (fractional weight or multiple calibration samples) may also be used to calibrate if desired.

4. Analyze Samples.
 - a. Weigh ~0.25 g of the soil or plant tissue sample into a 502-186 Tin Foil Cup and seal.
 - b. Enter mass and sample identification into Sample Login (F3).
 - c. Transfer sample to the appropriate position on the sample carousel.
 - d. Repeat steps 4a through 4c for each sample to be analyzed.
 - e. Initiate the analysis sequence (F5).

5. Atmospheric Blank.

Some atmosphere will be trapped with the sample when it is encapsulated in the tin foil cup. This will cause biased nitrogen results at low nitrogen concentrations (particularly with soil samples), therefore an atmospheric blank should be determined and entered using the following procedure: Analyze an inert material such as LECO 501-427 Com-Aid several times using similar weights of the Com-Aid to the weight of samples being analyzed. Enter the actual weight of the Com-Aid (Com-Aid should be baked-off in a muffle furnace at ~1000°C for 15 minutes, allowed to cool,

and stored for up to 24 hours in a desiccator until used). The nitrogen value obtained is considered the atmospheric blank and can be automatically compensated using the CN628 software. Refer to the operator's instruction manual for details regarding the setting of the atmospheric blank.

Notes

Sample Drying Instructions

- Soils—Samples should be dried at 105°C for one hour prior to analysis.
- Plant Tissues—Samples should be dried between 80°C and 85°C for two hours prior to analysis.

Typical Results

(Based on a single standard calibration with 0.15 g of 502-092 EDTA)

Sample	Mass g	% Carbon	% Nitrogen
Soil 502-308*	0.2504	2.42	0.186
Lot: 1016	0.2529	2.40	0.185
2.35 ±0.07% C	0.2497	2.41	0.185
0.183 ±0.010% N	0.2501	2.40	0.182
	0.2510	2.40	0.185
	0.2502	2.40	0.187
	0.2497	2.40	0.186
	0.2498	2.40	0.185
	0.2502	2.40	0.187
	0.2520	2.39	0.186
	X =	2.40	0.185
	s =	0.01	0.001

Soil 502-309*	0.2511	13.06	1.12
Lot: 1009	0.2478	13.04	1.11
13.01 ±0.27% C	0.2493	13.14	1.12
1.12 ±0.04% N	0.2512	13.01	1.12
	0.2518	13.12	1.12
	0.2498	13.07	1.12
	0.2525	13.17	1.12
	0.2523	12.97	1.11
	0.2518	12.95	1.11
	0.2496	12.99	1.11
	X =	13.05	1.12
	s =	0.07	0.01

Sample	Mass g	% Carbon	% Nitrogen
Orchard Leaves	0.2522	50.2	2.27
502-055**	0.2511	50.4	2.28
Lot: 1032	0.2471	50.3	2.27
50.4 ±0.4% C	0.2497	50.3	2.27
2.28 ±0.04% N	0.2480	50.4	2.27
	0.2540	50.3	2.26
	0.2480	50.2	2.26
	0.2506	50.5	2.27
	0.2509	50.2	2.27
	0.2494	50.2	2.26
	X =	50.3	2.27
	s =	0.1	0.01

Tobacco 502-082**	0.2499	46.7	2.56
Lot: 1015	0.2502	46.6	2.55
46.3 ±0.4% C	0.2513	46.6	2.56
2.56 ±0.04% N	0.2519	46.6	2.56
	0.2495	46.5	2.55
	0.2501	46.5	2.55
	0.2498	46.5	2.56
	0.2493	46.5	2.55
	0.2511	46.6	2.55
	0.2509	46.6	2.55
	X =	46.6	2.56
	s =	0.1	0.01



LECO Corporation | 3000 Lakeview Avenue | St. Joseph, MI 49085 | Phone: 800-292-6141 | 269-985-5496

info@leco.com • www.leco.com | ISO-9001:2015 Q-994 | LECO is a registered trademark of LECO Corporation.

LCRM = LECO Certified Reference Material; LRM = LECO Reference Material and are registered trademarks of LECO Corporation.