

Carbon in Refractory and Reactive Metals and Their Alloys[†]

LECO Corporation; Saint Joseph, Michigan USA

Instrument: CS600

Sampling and Sample Preparation

Surface contamination on the sample can cause significant errors in the analytical data; therefore, care must be taken to ensure a clean, representative sample is analyzed. Solid samples should be abraded with a clean file, rinsed in acetone, and dried with warm air prior to analysis. Samples that cannot be abraded due to irregular shapes should be rinsed in a suitable solvent such as acetone and dried with warm air. Care must be taken to remove all traces of the solvent. If a sample is porous, it is advisable to avoid using solvents, as it will be difficult to remove all traces of the solvent by drying.

Method Reference

ASTM E1941

Accessories

528-018 or 528-018HP Ceramic Crucibles and 619-880 Ceramic Crucible Covers (preheated)*; 502-492 High Purity Copper Accelerator and 502-231 High Purity Iron Chip Accelerator

**Ceramic crucibles and covers are baked in a muffle or tube furnace (LECO TF10) at 1250°C for a minimum of 15 minutes, or at 1000°C for 40 minutes. The crucibles/covers are removed from the furnace, allowed to cool for 1 to 2 minutes, and transferred to a desiccator for storage. If the crucibles/covers are not used within four hours, they should be re-baked.*


Calibration Samples

NIST or other suitable reference materials.

Note: Due to limited availability of refractory/reactive reference materials, steel calibration samples can be used to calibrate.

Method Parameters

Purge Time (seconds)	10
Delay Time (seconds)	20
Furnace Low Power (%)	100
Furnace High Power (%)	100
Furnace Ramp Rate	0
	Carbon
Minimum Timeout (seconds)	40
Comparator Level	1.00
Significant Digits	4 or 5
Integration Delay	0



Procedure

1. Prepare instrument for operation as outlined in the operator's instruction manual.
2. Determine blank.
 - a. Enter 1.0000 g mass into Sample Login (F3) using Blank as the sample name.
 - b. Add ~1.0 g of iron chip and ~1.5 g of copper accelerator to a prepared crucible.
 - c. Place a prepared crucible lid on the crucible.
 - d. Place the crucible on the furnace pedestal (or appropriate autoloader position if so equipped) and initiate Analyze (F5).
 - e. Repeat steps 2a through 2d a minimum of three times.
 - f. Set blank following procedure outlined in operator's instruction manual.
3. Calibrate/Drift Correct.
 - a. Weigh ~0.5 g calibration/drift sample into a prepared crucible, enter mass and sample identification into Sample Login (F3).
 - b. Add ~1.0 g of iron chip accelerator and ~1.5 g of copper accelerator on top of sample.
 - c. Place a prepared crucible lid on the crucible.
 - d. Place the crucible on the furnace pedestal (or appropriate autoloader position if so equipped) and initiate Analyze (F5).
 - e. Repeat steps 3a through 3d a minimum of three times for each calibration/drift sample intended for calibration/drift.
 - f. Calibrate/drift correct using the procedure outlined in the operator's instruction manual.
4. Analyze Samples.
 - a. Weigh ~0.5 g sample into prepared crucible, enter mass and sample identification into Sample Login (F3).
 - b. Add ~1.0 g of iron chip accelerator ~1.5 g of copper accelerator on top of sample.
 - c. Place a prepared crucible lid on the crucible.
 - d. Place crucible on furnace pedestal (or appropriate autoloader position if so equipped) and initiate Analyze (F5).

[†]Ti, Zr, W, Ta, Mo, Nb, Hf, and Re

Typical Results*

Sample	Mass g	C %
Titanium	0.5353	0.0247
Wire	0.5215	0.0246
	0.5106	0.0252
	0.5377	0.0245
	0.4677	0.0244
	0.4969	0.0246
	0.4590	0.0246
	0.4981	0.0245
	0.5281	0.0245
	0.4285	0.0246
	X =	0.0246
	s =	0.0002

NIST	0.5005	0.0250
SRM 173a	0.4989	0.0251
Titanium alloy	0.5031	0.0251
Chip	0.4973	0.0249
@ 0.025% C	0.4988	0.0251
	0.5041	0.0249
	0.5010	0.0249
	0.4986	0.0249
	0.4996	0.0250
	0.5023	0.0251
	X =	0.0250
	s =	0.0001

Sample	Mass g	C %
NIST	0.4983	0.0156
SRM 360a	0.5032	0.0158
Zircaloy-2	0.5056	0.0158
Chip	0.5061	0.0157
(C not certified)	0.4983	0.0154
	0.5094	0.0156
	0.5015	0.0156
	0.5032	0.0156
	0.5099	0.0154
	0.5091	0.0160
	X =	0.0157
	s =	0.0002

Tantalum	0.5029	0.0022
Powder	0.4962	0.0020
	0.5024	0.0022
	0.5069	0.0021
	0.4974	0.0021
	0.4977	0.0022
	0.5060	0.0024
	0.5092	0.0022
	0.5066	0.0022
	0.5004	0.0022
	X =	0.0022
	s =	0.0001

*Based on single-standard calibration using NIST SRM 173a



LECO Corporation

3000 Lakeview Avenue • St. Joseph, MI 49085 • Phone: 800-292-6141 • Fax: 269-982-8977
 info@leco.com • www.leco.com • ISO-9001:2008 • No. FM 24045 • LECO is a registered trademark of LECO Corporation.