



Can I Really Enjoy the Benefits of the Latest LC Column Technology?

Stephen Luke

LC Columns Product Manager

Agilent Technologies

Organization and Laboratory Needs

Organization needs

- Increased capacity
- Shorter time to market
- Increased profitability



Laboratory needs


- Better use of resources
- Increase productivity
- Reduce costs




Laboratory Needs

Translated to liquid chromatography

Better use of resources
 Increase productivity
 Reduce costs




Use all the instruments in the lab
 Run fast with high resolution
 Get long column lifetimes



Laboratory Needs

What's stopping you achieving them?


Better use of resources
Increase productivity
Reduce costs



Use all the instruments in the lab
Run fast with high resolution
Get long column lifetimes


Yes, but...

- I don't all have UHPLC instruments
- I can't / don't want to change my methods much
- My "dirty" samples clog my columns
- My columns fail in the mobile phases we use




Superficially Porous Particles for LC

Make every LC and LC/MS in your lab work harder



Poroshell 120 2.7 µm




Poroshell 120 4 µm


www.agilent.com/chem/discoverporoshell

Laboratory Needs Translated to liquid chromatography

Better use of resources
Increase productivity
Reduce costs



Use all the instruments in the lab
Run fast with high resolution
Get long column lifetimes



Equation for Pressure Drop Across an HPLC Column

$$\Delta P = \frac{\eta \cdot L \cdot v}{\theta \cdot d_p^2}$$

ΔP = Pressure Drop

η = Fluid Viscosity

L = Column Length

v = Flow Velocity

d_p = Particle Diameter

θ = Dimensionless Structural Constant
~ 600 For Packed Beds in LC


- Pressure and delay/dispersion volume are key components of instrument compatibility
- Shorter column lengths and larger particle diameters reduce column pressure



Use All the Instruments in the Lab

4 μm pressures suitable for any instrument

Flow Rate and Pressure



Column Dimensions: 4.6 x 100 mm

Mobile phase A: 0.1 % formic acid in water

Mobile phase B: 0.1 % formic acid in acetonitrile

Temperature: 35°C


Pressure reading at 5% B

× Eclipse Plus C18 1.8 μm ▲ Poroshell 120 EC-C18 2.7 μm
◆ Poroshell 120 EC-C18 4 μm ● Eclipse Plus C18 5 μm




Application note [5991-5510EN](#)

Volume Considerations for LC Systems

- The time taken for solvent to fill the delay volume represents an isocratic hold at the start of a run. For fast methods low delay volumes are required
- Sharp (efficient) peaks have a low peak volume. Low dispersion volumes are needed to avoid broadening of low volume peaks
- Both effects are most pronounced with narrow (2.1 mm) ID columns




Use All the Instruments in the Lab


| Instrument / Preference | Column | Details |
|--|---|---|
|  <p>UHPLC Low delay/dispersion Up to 1300 bar</p> | <p>Poroshell 120 2.7 μm</p> | <p>Pressure 50% of sub-2 μm TPP</p> |
|  <p>HPLC/UHPLC Mid delay/dispersion Up to 600 bar</p> | | |
|  <p>HPLC High delay/dispersion Up to 400 bar</p> | <p>Poroshell 120 4 μm</p> | <p>Typical pressure < 200 bar</p> |

Laboratory Needs Translated to liquid chromatography

Better use of resources
Increase productivity
Reduce costs




Use all the instruments in the lab
Run fast with high resolution
Get long column lifetimes



Fundamental Resolution Equation

$$R_s = \frac{\sqrt{N}}{4} \frac{(\alpha - 1)}{\alpha} \frac{k}{(k + 1)}$$




To increase resolution:

- Increase retention (k)
- Change selectivity (α)
- Increase efficiency (N)

Fundamental Resolution Equation

$$R_s = \frac{\sqrt{N} (\alpha - 1)}{4 \alpha} \frac{k}{(k + 1)}$$




- Retention and selectivity are impacted by changes to chemistry (mobile phase and stationary phase)
- Using a shorter column with high efficiency will maintain resolution in a shorter run time

Run Fast with High Resolution

All 3 van Deemter terms are reduced with SPP



van Deemter equation:

$$h = A + B/v + C \cdot v$$



Lower h = higher efficiency!


- **A term** – eddy diffusion
 - Particle size & packing quality
 - Narrow particle size distribution
- **B term** – longitudinal diffusion
 - Less mobile phase in the column
 - Reduced diffusion
- **C term** – mass transfer
 - Shorter diffusion paths
 - More effect on large molecules



Increase Resolution

4 μm gives good resolution at low pressure

Resolution vs Flow Rate
4-methyl phenol / 2-methyl phenol



▲ Eclipse Plus 1.8um 6,7 ■ Poroshell 120 2.7um 6,7
 ✕ Poroshell 120 4um 6,7 ◆ Eclipse Plus 5um 6,7

Column Dimensions: 4.6 x 100 mm

Mobile phase A: 0.1 % formic acid in water

Mobile phase B: 0.1 % formic acid in acetonitrile

Temperature: 35°C

| % B | Time (min) | | | | | | |
|--------------------|------------|----|-------|------|-----|------|------|
| 5 | 4 | 2 | 1.33 | 1 | 0.8 | 0.67 | 0.34 |
| 40 | 34 | 17 | 11.33 | 8.5 | 6.8 | 5.67 | 2.84 |
| 40 | 40 | 20 | 13.33 | 10 | 8 | 6.67 | 3.34 |
| 5 | 42 | 21 | 14 | 10.5 | 8.4 | 7 | 3.5 |
| 5 | 50 | 25 | 16.67 | 12.5 | 10 | 8.34 | 4.17 |
| Flow rate (mL/min) | 0.5 | 1 | 1.5 | 2 | 2.5 | 3 | 3.5 |

| Results at 2mL/min | R_s | P (bar) |
|---------------------------------|-------|---------|
| Eclipse Plus 1.8 μm | 3.1 | 422 |
| Poroshell 120 2.7 μm | 2.7 | 293 |
| Poroshell 120 4 μm | 2.5 | 169 |
| Eclipse Plus 5 μm | 2.0 | 103 |




Application note [5991-5510EN](#)

Other SPP Column Performance

Poor retention and peak shape from non-Agilent columns



Mobile phase A: 20 mM K_2HPO_4/KH_2PO_4 pH 7
Mobile phase B: Methanol
Flow rate: 1.2 ml/min

1. Uracil
2. Propranolol
3. Butyl Paraben
4. Dipropylphthalate
5. Amitriptyline



Increase Resolution

4 μm column increases resolution






Notoginseng, China Pharmacopoeia, edition 2015

Application note [5991-5554EN](#)

Run Fast with High Resolution

Shorter 4 μm column gives fast, high resolution runs

System Suitability Method Requirement: $N > 4000$, $R_s > 11.5$



Mobile phase: 50:49:1 MeCN:H₂O:Acetic acid

Flow rate: 1.2 mL/min

Peak ID

1. Naproxen

2. Butyrophenone

Application note [5991-5408EN](#)

Return on Investment

\$43K savings annually with 160% ROI / year

| | | |
|--------------|-----------------------|---|
| Scope | Number of instruments | 1 |
| | Number of years | 5 |


| | | Conventional | Fast |
|---------------|-----------------------------------|--------------|-------|
| Method | Analysese per year (1 instrument) | 5,000 | 5,000 |
| | Run time (minutes) | 10 | 6 |
| | Flow rate (mL/min) | 1.2 | 1.2 |
| | Column lifetime (injections) | 1,000 | 1,000 |
| | Development time (hours) | 0 | 80 |

| | | Conventional | Fast | Conventional | Fast | Saving |
|--------------|--------------------------------------|------------------------|------------------------|------------------|------------------|------------------|
| Value | | <i>Amount/analysis</i> | <i>Amount/analysis</i> | <i>Cost/year</i> | <i>Cost/year</i> | <i>Cost/year</i> |
| | Overhead per instrument (\$/hour) | \$20.83 | \$12.50 | \$104,167 | \$62,500 | <u>\$41,667</u> |
| | Cost of mobile phase purchase (\$/L) | \$0.75 | \$0.45 | \$3,750 | \$2,250 | <u>\$1,500</u> |
| | Cost of mobile phase disposal (\$/L) | \$0.06 | \$0.04 | \$300 | \$180 | <u>\$120</u> |
| | Columns (\$) | \$0.55 | \$0.55 | \$2,750 | \$2,750 | <u>\$0</u> |
| | <u>Total</u> | <u>\$22.19</u> | <u>\$13.54</u> | <u>\$110,967</u> | <u>\$67,680</u> | |
| | Saving | | \$8.66 | | \$43,287 | |

| | | Conventional | Fast | Cost |
|--------------|-----------------------------------|--------------------|--------------------|--------------------|
| Costs | | <i>Cost/method</i> | <i>Cost/method</i> | <i>Cost/method</i> |
| | Overhead per instrument (\$/hour) | \$0 | \$10,000 | <u>\$10,000</u> |
| | Method development (\$/hour) | \$0 | \$14,000 | <u>\$14,000</u> |
| | Instrument | \$0 | \$0 | <u>\$0</u> |
| | <u>Total</u> | <u>\$0</u> | <u>\$24,000</u> | |

| | | Conventional | Fast |
|-------------|-------------|--------------|-------------|
| ROI% | Value | \$0.00 | \$216,433 |
| | Cost | \$0.00 | \$24,000 |
| | ROI% / year | | 160% |

Increase Resolution



Selectivity Impacts Resolution Most

- Change bonded phase
- Change mobile phase

$$R_s = \frac{\sqrt{N} (\alpha - 1)}{4 \alpha} \frac{k}{(k + 1)}$$

Increase Resolution

Choice of bonded phase selectivity - 12 chemistries

| Best all around | Best for low pH mobile phases | Best for high pH mobile phases | Best for alternative selectivity | Best for more polar compounds |
|--|--|---|---|---|
| Poroshell 120 EC-C18 2.7 µm, 4 µm | Poroshell 120 SB-C18 2.7 µm | Poroshell HPH-C18 2.7 µm, 4 µm | Poroshell 120 Bonus-RP 2.7 µm | Poroshell 120 SB-Aq 2.7 µm |
| Poroshell 120 EC-C8 2.7 µm, 4 µm | Poroshell 120 SB-C8 2.7 µm | Poroshell HPH-C8 2.7 µm, 4 µm | Poroshell 120 PFP 2.7 µm, 4 µm | Poroshell 120 EC-CN 2.7 µm |
| Poroshell 120 Phenyl-Hexyl 2.7 µm, 4 µm | | | | Poroshell 120 HILIC 2.7 µm, 4 µm |



Wall Chart [5991-6240EN](#)

Increase Resolution

The right selectivity

| Column | H | S* | A | B | C _(2.8) | C _(7.0) | F _s (pH 2.8) |
|----------------------------|------|-------|-------|------|--------------------|--------------------|-------------------------|
| Poroshell 120 EC-C18 | 1.02 | 0.01 | -0.13 | 0.00 | 0.16 | 0.12 | 0.0 |
| Poroshell 120 EC-C8 | 0.88 | 0.01 | -0.23 | 0.02 | 0.13 | 0.09 | 6.0 |
| Poroshell 120 Phenyl-Hexyl | 0.75 | -0.08 | -0.39 | 0.02 | 0.14 | 0.14 | 13.1 |
| Poroshell 120 PFP | 0.63 | -0.52 | -0.52 | 0.43 | -0.11 | N/A | 85.5 |

Hydrophobic Subtraction Model (HSM) Data provided by Dwight Stoll

F_s factor describes the similarity of two column selectivities. A small F_s indicates that two columns are very similar, while a large factor indicates that two columns are very different. Calculated according to the following equation:

| Hydrophobic Subtraction Model | |
|-------------------------------|----|
| Hydrophobicity | H |
| Steric interaction | S* |
| Hydrogen-bond acidity | A |
| Hydrogen-bond basicity | B |
| Ion-exchange capacity | C |

$$F_s = \left\{ [12.5(H_2 - H_1)]^2 + [100(S_2^* - S_1^*)]^2 + [30(A_2 - A_1)]^2 + [143(B_2 - B_1)]^2 + [83(C_2 - C_1)]^2 \right\}^{\frac{1}{2}}$$


Further details at:

<http://www.usp.org/app/USPNF/columnsDB.html>

Reversed-phase only
so Poroshell 120 HILIC
not included

Increase Resolution


The right selectivity



Application note [5991-6283EN](#)

Increase Resolution

Very similar selectivity - Poroshell 120 HPH-C18 vs EC-C18




MP-A 10mM Ammonium formate/water
adj. to pH3 using Formic Acid
MP-B ACN
Flow – 0.42 mL/min
Column Temp. Ambient
1 uL injection
Detection 254 nm
Gradient:

| time, min | %B |
|-----------|----------|
| 0 | 5 |
| 4 | 95 |
| 5 | 95 |
| 6 | 5 |
| 7 | stop run |

Increase Resolution


Very different selectivity – pH10 vs pH3



Increase Resolution

Change in retention with mobile phase pH


More retention for non-charged analytes
(i.e. acids at low pH and bases at high pH)



Increase Resolution

Very different selectivity at high and low pH

Mixture of acidic and basic analytes




Application note [5991-6525EN](#)

Advantages of High pH with LC-MS

Improved retention, resolution and response at high pH



Drugs of abuse analysis by LC-MS



Application note [5991-6523EN](#)

Other High-pH SPP Column Performance

Poor retention and peak shape from non-Agilent column



| Compound | HPH Column T _f | Column Y T _f |
|---------------|---------------------------|-------------------------|
| Doxepin | 1.68 | 2.97 |
| Nortriptyline | 1.36 | 1.96 |
| Amitriptyline | 1.25 | 2.41 |
| Trimipramine | 1.07 | 1.60 |

Application note [5991-6512EN](#)

Instrument Considerations at High pH



- In general, Agilent Infinity modules operate over a pH range of 1 to 12.5
- pH < 2.3
 - Solvents must not contain acids that attack stainless steel
- pH > 9.5
 - Replace standard (Vespel) rotor seals in all valves with either Tefzel or PEEK seals*
 - Replace standard glass solvent inlet filters with stainless steel filters
 - Be aware that quartz flow cells in detectors will etch slowly. Do not leave high pH solvents stationary in the cells for extended periods

* Examples:

5068-0171 Rotor Seal PEEK FL for 1290 Infinity Binary Pump


5068-0170 Rotor Seal PEEK FL for 1290 Infinity High Performance Autosampler

5068-0172 Rotor Seal PEEK FL for 1290 Infinity Quaternary Pump


5068-0223 Rotor Seal PEEK for 6-Column Selector Valve

Laboratory Needs Translated to liquid chromatography

Better use of resources
Increase productivity
Reduce costs




Use all the instruments in the lab
Run fast with high resolution
Get long column lifetimes



Costs of Short Column Lifetime

- More columns to purchase
- Disruption to workflow
- Re-work required



Column Lifetime



Get Long Column Lifetimes


Causes of Column Failure

- Blocked frits
- Contamination/blockage
- Voiding of packed bed
- Dissolution of particles
- Loss of bonded phase





Can be Avoided Using

- Wide frit porosity
- Guard columns
- Robust packed bed
- High stability particles
- Robust bonding



Get Long Column Lifetimes

>1800 Injections - no performance change



A: 0.2% HCOOH in H₂O, B: 0.2% HCOOH in CH₃CN

0.833 mL/min

| | | | |
|------|------|------|------|
| Time | 0.00 | 1.25 | 2.50 |
|------|------|------|------|

| | | | |
|----|----|----|----|
| %B | 10 | 15 | 27 |
|----|----|----|----|

40 °C P_{max} = 550 bar

Agilent Poroshell 120 SB-C18, 2.1 x 100 mm, 2.7 μm

Sig=210,4nm, Ref=Off

2-μL, 3-mm micro flow cell (PN G1315-60024)

Sample: 2 μL of freshly brewed green tea

(brewed from a commercial tea bag in 6 oz of initially boiling water for six minutes)

2 μm inlet frit prevents clogging and extends column lifetime


Get Long Column Lifetimes

Advantages of guard columns



Accelerated Lifetime Test

Similac sample (milk substitute diluted 300:1) containing 2 sulfa drugs
Peak width change indicating column failure



No Guard

- Column failure @ inj. 70; new column required

With Guard



- Guard failure @ inj. 80; guard replaced
- Same column used throughout analysis

By installing a guard column when using dirtier samples, one can extend the life of their column, and utilize more inexpensive guard columns rather than analytical column replacements

Get Long Column Lifetimes

Approaches for longer lifetime at high pH

| Approach | Comments |
|---|--|
| Totally porous silica-hybrid particles | Do not have the efficiency and moderate pressure of superficially porous particles |
| Bonding chemistry on superficially porous silica particles | Do not have the lifetime of silica-hybrid particles |
| Chemical modification of the outer layer of superficially porous silica particles | Combine the advantages of silica-hybrid and superficially porous particles |



Get Long Column Lifetimes


Stability at high pH

Column: Poroshell HPH-C18, 2.1 x 50 mm, 4 um

A: 10 mM, pH 10 ammonium bicarbonate
 B: Acetonitrile
 25 °C, 0.4 mL/min, 2.1 x 50 mm column
 220 nm, 4 nm ref off
 Agilent 1290 Infinity, PEEK rotor seals

| Time | % B |
|------|-----|
| 0.0 | 5 |
| 5.0 | 95 |
| 5.1 | 5 |
| 7.0 | 5 |

1. 2-Hydroxy-5-methyl-benzaldehyde
2. 4-Chloro-cinnamic acid
3. Acetophenone
4. Quinine sulfate
5. Nortriptyline
6. Amitriptyline
7. Hexanophenone



Application note [5991-6525EN](#)

Other High-pH SPP Column Performance

Drifting retention times

Column: Column Y C18, 2.1 x 50 mm, 2.6 μ m


Mobile phase A: 10 mM, pH 10 Ammonium Bicarbonate

Mobile phase B: Acetonitrile

Flow rate: 0.4 ml/min


Temperature: 25 C

1. 2-Hydroxy-5 methyl-benzaldehyde
2. 4-Chloro-cinnamic acid
3. Acetophenone
4. Quinine Sulfate
5. Nortriptyline
6. Hexanophenone
7. Amitriptyline




Poroshell 120

Make every LC and LC/MS in your lab work harder



Poroshell 120 2.7 µm

Efficiency 90% of < 2 µm TPP
Pressure 50% of < 2 µm TPP
12 chemistries, including HPH
2 µm inlet frit



Poroshell 120 4 µm


Efficiency 2x 5 µm TPP
Pressure often below 200 bar
7 chemistries, including HPH
2 µm inlet frit

www.agilent.com/chem/discoverporoshell

TPP = totally porous particle

Summary

- The pressure and peak widths of Poroshell 4 μm particles allow you to use all the instruments in the lab
- High efficiency Poroshell particles and the right bonded phase selectivity enable fast runs with high resolution
- 2 μm inlet frits and innovative chemistries like HPH provide long column lifetimes – even at high pH
- You really can enjoy the benefits of the latest LC column technology



stephen.luke@agilent.com
www.agilent.com/chem/discoverporoshell