

Agilent RapidFire 365 High-throughput Mass Spectrometry System

Data Analysis Guide

Data Analysis with RapidFire Integrator or MassHunter Quantitative Analysis	2
Data analysis with MassHunter Quantitative Analysis	2
To convert data for use with MassHunter Quantitative Analysis	3
Installation	5
To install Agilent RapidFire Integrator software	5
Location, Format, and Preparation of Data for RapidFire Integrator	6
Location of data	6
Data files	7
Exact mass platemap	9
Loading RapidFire-QQQ Data in RapidFire Integrator	10
To prepare Agilent MassHunter data	10
To prepare AB Sciex Analyst data	11
To prepare Thermo Fisher Xcalibur data	13
To open RapidFire-QQQ data	16
Loading RapidFire-TOF/QTOF Data in RapidFire Integrator	17
To open RapidFire-TOF/QTOF data	17
Analyzing Data in RapidFire Integrator	20



Data Analysis with RapidFire Integrator or MassHunter Quantitative Analysis

RapidFire 365 data can be acquired from triple-quadrupole (QQQ) or time-of-flight (TOF/QTOF) LC/MS instruments controlled by the Agilent MassHunter Data Acquisition software or other software, such as Analyst or Xcalibur.

After data acquisition is completed (see *G9530-90000 RapidFire 365 User Guide*), you can process the data using RapidFire Integrator as described in this guide.

- The Agilent RapidFire Integrator data analysis software can be used to integrate all of the data points collected in a single Agilent G9530AA RapidFire 365 High-throughput Mass Spectrometry System run.
- Each sequence of a RapidFire-MS batch is processed individually. There is no limit on the number of plates.

For data that was acquired with Agilent MassHunter Data Acquisition software, RapidFire 365 data can also be processed using Agilent MassHunter Workstation Software Quantitative Analysis. For more information, see:

- “Data analysis with MassHunter Quantitative Analysis” below, and
- *G3335-90152 Agilent MassHunter Workstation Software Quantitative Analysis Familiarization Guide*.

Data analysis with MassHunter Quantitative Analysis


During the RapidFire-MassHunter data acquisition, data files (.D) are generated for each sequence.

The total ion chromatogram (TIC) can be viewed in real time in MassHunter Data Acquisition (MHAcq), and can be inspected further in MassHunter Qualitative Analysis.


After a run has completed, the original data files (long sequence **X.D** files that are generated by MassHunter Data Acquisition and **RFDDataBase.XML** files that are generated by RapidFire) must be split into individual injection files for the standard data analysis by MassHunter Quantitative Analysis software.

To convert data for use with MassHunter Quantitative Analysis

- 1 In the RapidFire user interface, click **File > Convert MS Data**.




- 2 On the Convert MS Data dialog box:
 - a Select the whole-experiment data directory (**Data Path**) you want to convert.
 - b Select where to save the resulting single-injection data files (**Save To**). The data path of the latest run is displayed as a default selection.
 - c Click **Convert**.



Data Analysis with RapidFire Integrator or MassHunter Quantitative Analysis

To convert data for use with MassHunter Quantitative Analysis

- 3 View the list of converted files as shown in the following example. The single-peak data files are automatically named **InjXXXXX-barcode-well.d**, where:
- "XXXXX" is the injection order in the RapidFire batch.
 - "barcode" is the plate name.
 - "well" is the plate position of the sample on the RapidFire platemap.



- 4 Proceed with viewing and processing the single-injection data files, as described in MassHunter Quantitative Analysis online Help and manual.

Installation

To install Agilent RapidFire Integrator software

The RapidFire Integrator data analysis software can be installed on any personal computer.

Before you begin

Contact your Agilent field support engineer to request a copy of the latest RapidFire Integrator program, as well as the following prerequisite Windows components:

- Microsoft .NET 4.0 framework
- Visual C++ runtime libraries

To install RapidFire Integrator

- 1 Copy the entire **RapidFire Integrator** folder and paste it into **C:\Agilent**.
- 2 Create a shortcut to the **RapidFire Integrator.exe** file and place it on the computer desktop.

Location, Format, and Preparation of Data for RapidFire Integrator

The following reference material applies to RapidFire Integrator data analysis software.

Location of data

RapidFire data is stored in the folder **D:\ MassHunter\ Data\ RapidFire** on the MS computer.

- A new data folder is created each day with the date as the folder name.
- A new folder is created when a new experiment (in Plates mode) or a new batch (in Sequences mode) is started in RapidFire acquisition software.

Example The first run performed on April 5, 2013 is saved in the following folder on the MS computer:

D:\ MassHunter\ Data\ RapidFire \ 2013\ April\ 5\ 1

NOTE

The data file folders that are created by MS computer settings may be different than the scheme described above. Refer to the **BASE_DATA_DIR** variable in **RFMassHunterS.cfg** on your MS computer, which is located in **C:\ Agilent\ RapidFire Communicator\ cfgs**.

NOTE


If the RapidFire front end is connected to an AB Sciex MS, then the data are stored in **D:\ Analyst Data\ RapidFire** and **BASE_DATA_DIR = "D:\ Analyst Data\ RapidFire"** in **RFAAOS.cfg** on your MS computer.

Data files

A data set for analysis with the RapidFire Integrator software includes the following types of files. *These files must be present in the same folder.*

- **batch.rftime**
- **sequenceX.d** mass spectrometer (MS) data files in MassHunter, *or* **sequenceX-xic-YZ.rfxic** in AB Sciex Analyst or Thermo Fisher Xcalibur
- A **platemap.tofmap.txt** tab-delimited text file, for RapidFire-TOF/QTOF systems. The original **platemap.tofmap.txt** file is first empty. Modify it as described on [page 9](#).

The following example data folder of a RapidFire qTOF shows the three required file types:



If the data is acquired within a single experiment in Plates mode, then multiple 96-well or 384-well plates are analyzed within a single folder.

The following columns of data are saved in the **batch.rftime** text file:

- **Plate identity**
- **Injection number**
- **Sequence number**
- **Row**
- **Column**
- **Sip Time**, the time stamp for the actuation of valve 2 from inject to load position, as well as the start of elution to MS.
- **Sip sensor value**. The displayed value is **1** if the optical sip sensor detected the presence of liquid, or **0** if it did not.

Example

	A	B	C	D	E	F	G	H
1	plate	sip	seq	row	col	siptime	sip sensor	
2	-----							
3	DEMO PLATE	1	1	2001	2	3.304	0	WASH STATION
4	DEMO PLATE	2	1	1	1	12.529	1	
5	DEMO PLATE	3	1	1	3	21.341	1	
6	DEMO PLATE	4	1	1	5	30.238	1	
7	DEMO PLATE	5	1	2001	2	38.486	1	
8	DEMO PLATE	6	1	2	1	47.582	1	B1
9	DEMO PLATE	7	1	2	3	56.322	1	B3
10	DEMO PLATE	8	1	2	5	65.164	1	B5
11	DEMO PLATE	9	1	2001	2	73.49	1	
12	DEMO PLATE	10	1	3	1	82.726	1	
13	DEMO PLATE	11	1	3	3	91.705	1	
14	DEMO PLATE	12	1	3	5	100.552	1	
15	DEMO PLATE	13	1	2001	2	108.996	1	
16	DEMO PLATE	14	1	4	1	118.225	1	
17	DEMO PLATE	15	1	4	3	127.228	1	
18	DEMO PLATE	16	1	4	5	136.047	1	

NOTE

- The aqueous and organic wash stations, WASH1 and WASH2, have (row, column) addresses of (2001, 2) and (2001, 3), respectively.
- The matrix vials MAT1, MAT2, MAT3, and MAT4 have (row, column) addresses of (1001, 2), (1001, 3), (1001, 4), and (1001, 5), respectively.
- If data were acquired in Plates mode, make sure to name or rename the MS data files to **sequence1.d**.
- If the RapidFire front end is connected to an AB Sciex QQQ, then the **.WIFF** MS data files are automatically processed by the RapidFire software and transformed into text files called **sequenceX-xic-Q1---Q3.rfxic**.

Exact mass platemap

For data that were acquired via RapidFire-QTOF in *MS/MS mode*, you don't need to modify the **platemap.tofmap.txt** file.

For RapidFire-TOF/QTOF experiments in *MS mode*, modify **platemap.tofmap.txt** to specify the exact masses of the compounds of interest for each well that was injected. The RapidFire Integrator uses these values to extract the appropriate ion chromatograms from the TOF data of each well. For data that are acquired in MS mode, the columns of data in the **platemap.tofmap.txt** text file are:

- Sequence number
- Row
- Column
- Exact mass 1 (or chemical formula 1)
- Exact mass 2 (or chemical formula 2), and so on

More columns can be appended as needed.

Example



	A	B	C	D	E	F
1	Sequence	Row	Column	EXACT MASS	Internal standard	
2						
3	1	2001	2	507.2344911	493.14382	WASH STATION
4	1	1	1	507.2344911	493.14382	
5	1	1	3	507.2344911	493.14382	
6	1	1	5	507.2344911	493.14382	
7	1	2001	2	327.1946771	493.14382	
8	1	2	1	327.1946771	493.14382	B1
9	1	2	3	327.1946771	493.14382	B3
10	1	2	5	327.1946771	493.14382	B5
11	1	2001	2	427.0921388	493.14382	
12	1	3	1	427.0921388	493.14382	
13	1	3	3	427.0921388	493.14382	
14	1	3	5	427.0921388	493.14382	
15	1	2001	2	318.1038131	493.14382	
16	1	4	1	318.1038131	493.14382	
17	1	4	3	318.1038131	493.14382	
18	1	4	5	318.1038131	493.14382	

Loading RapidFire-QQQ Data in RapidFire Integrator

A run in Plates mode is defined as the data that are collected without setting up a batch. Unlike Sequences mode, data acquired in Plates mode are named and saved in folders possibly with *different* names and locations on the RapidFire and MS computers.

If the data are collected in Plates mode with MassHunter, Analyst, or Xcalibur, then prepare the data as described in this section before you analyze it with RapidFire Integrator.

If the data are collected in Sequences mode, then proceed with “[To open RapidFire-QQQ data](#)” on page 16.


To prepare Agilent MassHunter data

Repeat this process for each RapidFire-MS run.

- 1 Rename the MS data files to **sequence1.d**.
- 2 Copy the following files to a unique folder, such as **\data sequenceX**:
 - **sequence1.d**.
 - **batch.rftime** (from the RapidFire computer, see “[Data files](#)” on page 7).

To prepare AB Sciex Analyst data

- 1 Open the Analyst program.
 - a Click **Explore > Open Data File**.




- b Select the appropriate **.WIFF** file from the list.
- c Click to select the sample of interest, then click **OK**.

The MS chromatogram traces for the selected sample appear in the main window.

- 2 Create extracted ion chromatograms (XIC) as follows:
 - a Click on the **XIC** button on the Analyst toolbar.
 - b Click to select the masses of interest, which are the individual MS chromatograms for each of the MRM channels that were monitored during the experiment.

Loading RapidFire-QQQ Data in RapidFire Integrator


To prepare AB Sciex Analyst data



- 3 Save each XIC as a text file as follows:
 - a Right-click on the panel for the XIC of interest and select **Save to Text File**.
 - b Select the folder where you want to save the file. Use the same directory where the RapidFire **batch.rftime** data file is already saved.
 - c Name the MS data files **sequence1-xic-XXX.txt**, where XXX is any unique string of alphanumeric characters, such as **sequence1-xic-MSprod.txt**.


To prepare Thermo Fisher Xcalibur data

- 1 Open the Xcalibur program.
 - a Click **File > Open**.




- b Select the appropriate **.RAW** file from the list.
 - c Click **Open**.
- The MS chromatogram traces for the selected sample appear in the main window.
- 2 Create extracted ion chromatograms (XIC) as follows:
 - a Right-click right on the total ion chromatogram (TIC) and select **AutoFilter**.
 - b Delete the TIC panel at the top of the window.

Loading RapidFire-QQQ Data in RapidFire Integrator To prepare Thermo Fisher Xcalibur data



- 3 Save each XIC as a text file as follows:
 - a Right-click on the panel for the XIC of interest and select **Export to Clipboard (Chromatograms)**.




- b Paste the selection into Excel.
- c Save it as a text file (.TXT) in the same directory where the RapidFire **batch.rftime** data file is already saved.
- d Name the MS data files **sequence1-xic-XXX.txt**, where XXX is any unique string of alphanumeric characters, such as **sequence1-xic-MSprod.txt**.

To open RapidFire-QQQ data

- 1 To start the Agilent RapidFire Integrator software, double-click the desktop shortcut



- 2 Click **File > Load RF300-QQQ Data**.



- 3 When the Browse For Folder window appears:
 - a Select the folder that contains the **.RFTIME** and **.RFXIC** or **.D** files.
 - b Click **OK**.
- 4 Continue with “[Analyzing Data in RapidFire Integrator](#)” on page 20.


Loading RapidFire-TOF/QTOF Data in RapidFire Integrator

To open RapidFire-TOF/QTOF data

- 1 To start the Agilent RapidFire Integrator software, double-click the desktop shortcut.




- 2 Click **File > Load RF360-TOF/QTOF Data**.



- 3 When the Browse For Folder window appears:
 - a Select the folder that contains the **.D**, **.RFTIME** and **platemap.TOFMAP.TXT** files.
 - b Click **OK**.
- 4 Click **Data > Extract TOF/QTOF XICs**.


- 5 For data that were acquired with a TOF or with a QTOF in *MS mode*:
 - a Set the following parameters on the TOF Options Form dialog box:
 - **Display Time Range** (set to all by default)
 - **Polarity**
 - **Error (ppm)**, which determines the accuracy
 - Whether or not to **Use Internal Standard**, and if marked, then enter the mass of the internal standard. To extract multiple internal standards in every well of the RapidFire-TOF run, enter their masses separated by commas (,).




- b Click **OK** to close the TOF Options Form dialog box.
- 6 For data that were acquired with a QTOF in *MS/MS mode*:
 - a Click **Add** on the QTOF Options dialog box.
 - b In the Add Product Ion dialog box:
 - Select **Sequence** and **Precursor Ion** from the drop-down fields.
 - Type in the mass or chemical formula of one **Product Ion** of interest and the **Error** (ppm).
 - Click **OK**.
 - c Repeat for all Precursor Ions and all Product Ions of interest.

Loading RapidFire-TOF/QTOF Data in RapidFire Integrator

To open RapidFire-TOF/QTOF data




- d When you are finished, click **OK** to close the QTOF Options dialog box.
- 7 Select **Minutes** for the MS data files time units.




Analyzing Data in RapidFire Integrator

- 1 Click to select the sequence of interest in the left pane, as shown for Sequence 1 in the following example.



- 2 Review the data using the following features:

- To display the list of m/z , click the + sign next to the sequence.
- To display the following menu options, right-click in the left pane:



- Click to select the number of graphs to display in the upper right area of the window.
- To select an area of the graph to expand, drag the mouse. Repeat to expand the area further.
- To zoom out each successive level, right-click in the graph and select **Un-zoom**.


3 Integrate peaks as follows:

a Zoom in on a known landmark peak within the sequence.

- Any peak within the experiment can be selected as a landmark.
- Data analyses can be based on any mass.
- The title of the plot identifies monitored analytes by their m/z .

To navigate from one m/z to another, use the scroll bar on the right side of the window.


b Select the identity of the peak from the **Sample** list.



c Define the peak to integrate for area-under-the-curve (AUC) calculations in *either* of the following ways:

- Drag the mouse in a thin slab rectangle from the front edge of the peak to its tail end.
- Use the keyboard to type in the T_1 and T_2 times.

The defined time range is displayed in the T_1 and T_2 boxes above the graph. The peak area is also emphasized by teal coloring as shown in the following example.




d To process the data for the current sequence, click the summation button.



- All masses within the same sequence are processed together.
- AUCs of injections turn blue, but background signals remain red.

e Repeat Steps 3a. through 3d. for each sequence of the run.



- 4 To save partial progress in data analysis before all of the sequences of the run are analyzed, click **File > Save As**.
- The results are saved in the specified **.RFD** file.
 - Click **File > Open** to reopen experiment files that have been partially processed.

5 To export data that were acquired with a QQQ:

- Click **File > Export Peak Data List**.

A **.RFPKS** text file is created, which contains a list of injections.

You can analyze this file with a spreadsheet program. The following is an example file opened in Microsoft Excel.

	A	B	C	D	E	F
1	Sequence 1					
2	PlateName	Row	Column	xic-342---203.1	xic-345---206.1	Injection Time (sec)
3	Plate1	1	1	18800	3346420	8.42
4	Plate1	1	2	67890	6215320	18.42
5	Plate1	1	3	352140	6442985	28.23
6	Plate1	1	4	1707265	5712430	37.58
7	Plate1	1	5	4602310	6090260	47.2
8	Plate1	1	6	5283070	7169040	56.89
9	Plate1	1	7	4763010	6029000	66.08
10	Plate1	1	8	5072060	6568350	75.56
11	Plate1	1	9	4681235	6430260	84.95
12	Plate1	1	10	4934730	6827290	94.59
13	Plate1	1	11	23465	49370	103.96
14	Plate1	1	12	10939460	14423820	113.34
15	Sequence 2					
16	PlateName	Row	Column	xic-342---203.1	xic-345---206.1	Injection Time (sec)
17	Plate1	2	1	18230	3743510	7.72
18	Plate1	3	2	81380	6190690	17.83
19	Plate1	2	3	410120	7120610	27.55
20	Plate1	3	4	2192220	7106175	37.47

6 To export data that were acquired with a TOF or QTOF in *MS mode*:

- Click **File >Export Peak Data TOF**.

A **.RFPKS** text file is created, which contains a list of injections.

Only the masses of interest are reported for each well.

You can further analyze the results with a spreadsheet program. The following is an example file opened in Microsoft Excel.

	A	B	C	D	E	F	G	H	I	J
1	Sip	Sequence	PlateName	Row	Column	Injection time (sec)	Mass 1	Abundance 1	Mass 2	Abundance 2
2	1	1	DEMO PLATE	2001	2	9.23	507.2345	10648	493.1438	2747
3	2	1	DEMO PLATE	1	1	18.46	507.2345	2992271	493.1438	297438
4	3	1	DEMO PLATE	1	3	27.27	507.2345	3061670	493.1438	490474
5	4	1	DEMO PLATE	1	5	36.17	507.2345	2900104	493.1438	475369.5
6	5	1	DEMO PLATE	2001	2	44.41	327.1947	0	493.1438	4295
7	6	1	DEMO PLATE	2	1	53.51	327.1947	1819135	493.1438	574327
8	7	1	DEMO PLATE	2	3	62.25	327.1947	1888354	493.1438	552930
9	8	1	DEMO PLATE	2	5	71.09	327.1947	1403209	493.1438	517027
10	9	1	DEMO PLATE	2001	2	79.42	427.0921	41	493.1438	4003.5
11	10	1	DEMO PLATE	3	1	88.65	427.0921	1135663	493.1438	527903
12	11	1	DEMO PLATE	3	3	97.63	427.0921	195475	493.1438	503490
13	12	1	DEMO PLATE	3	5	106.48	427.0921	292988	493.1438	492337.5
14	13	1	DEMO PLATE	2001	2	114.92	318.1038	311	493.1438	2913
15	14	1	DEMO PLATE	4	1	124.15	318.1038	31407	493.1438	545323
16	15	1	DEMO PLATE	4	3	133.16	318.1038	29637	493.1438	505078.5
17	16	1	DEMO PLATE	4	5	141.98	318.1038	34226	493.1438	490114

7 To export data that were acquired with a QTOF in *MS/MS mode*:

- Click **File > Export Peak Data List**.

A **.RFPKS** text file is created, which contains a list of injections.

Only the **precursor ion : product ion** pairs of interest are reported for each sequence.

You can further analyze the results with a spreadsheet program. The following is an example file opened in Microsoft Excel.

	A	B	C	D	E	F	G
1	Sequence 1						
2	PlateName	Row	Column	208.2379 : 160.0000	211.2564 : 163.0000	276.3509 : 205.0000	Injection Time (sec)
3	Std Curve	1	3	27	7826	0	10.47
4	Std Curve	2	3	0	18471	202	27.1
5	Std Curve	3	3	395	16137	768	43.8
6	Std Curve	4	3	145	13475	406	60.38
7	Std Curve	5	3	0	19239	530	76.86
8	Std Curve	6	3	254	14069	1586	93.48
9	Std Curve	7	3	1019	20907	3064	109.94
10	Std Curve	8	3	1263	17357	5734	126.44
11	Std Curve	1	4	47	30059	106	142.87
12	Std Curve	2	4	85	34318	0	159.34
13	Std Curve	3	4	0	33395	0	176.12
14	Std Curve	4	4	120	37298	0	192.44
15	Sequence 2						
16	PlateName	Row	Column	208.2379 : 160.0000	211.2564 : 163.0000	276.3509 : 205.0000	Injection Time (sec)
17	Std Curve	1	3	68	8270	0	11.24
18	Std Curve	2	3	134	19820	257	27.8
19	Std Curve	3	3	77	15757	196	44.26
20	Std Curve	4	3	221	17702	632	60.83

8 To export RapidFire-QQQ data *by plates*:

- Click **File > Export Peak Data Plate-wise**.

Keep in mind the following caveats:

- Data for MAT# and WASH# injections are *not* displayed.
- This option is only offered if each sample well within a plate is visited only once.

You can analyze this file with a spreadsheet program, as shown in the following example in Microsoft Excel.

	A	B	C	D	E	F	G	H	I	J	K	L	M
1	XIC = xic-342---203.1												
2		1	2	3	4	5	6	7	8	9	10	11	12
3	A	17060	69240	351195	1675510	4485350	5145550	4675030	5006445	4630220	4912200	5900	10379525
4	B	22620	101930	409910	2519030	5382330	6074400	5543170	6310075	5374720	5889960	5441500	9831590
5	C	15530	82915	487280	2191820	5966590	5331980	6205310	5569820	6738300	5162370	6531670	7544630
6	D	35660	116230	635105	2999890	7014370	6974840	7239930	7294740	6773310	6273410	6765560	12194250
7	E	28270	123070	619530	3241380	7240720	7195915	7231920	6897140	6673390	6537440	6863870	14057730
8	F	16730	113440	648320	3263800	3135335	7091210	6938345	5974780	6590740	6263830	6233490	13339040
9	G	9220	138810	613570	3266685	7388100	6034660	7144950	6912480	5846570	5996250	6499790	7195990
10	H	5000	126890	702220	3093510	7312870	7475130	5911630	5634800	6832050	6414380	6990580	9278240
11													
12	XIC = xic-345---206.1												
13		1	2	3	4	5	6	7	8	9	10	11	12
14	A	3225680	6069120	6334215	5616810	5918680	6971700	5893470	6458745	6351660	6765690	13070	13552035
15	B	3762920	7471300	7134560	8075850	7261610	7952300	7361700	8446565	7231980	8255770	7480010	12304715
16	C	4614840	6195005	7830980	7109350	7752745	7066970	8158290	7669750	9088970	7261570	8688110	10037330
17	D	5137720	8598850	9067560	9328840	9354920	9053750	9416945	9146200	9264970	8556380	9120250	15759350
18	E	5040400	8489620	8636100	9187590	9358630	9060130	9115380	9003460	9253780	9332390	9280810	18907170
19	F	4970310	8287770	9052850	8977820	4413770	8956190	8779035	8012300	8917535	8807140	8532365	16720920
20	G	5260530	8750170	8687950	9301155	9386290	7878040	9091835	8854535	8031880	8598980	8799570	9218770
21	H	5259840	8857620	9088550	8940810	9437290	9680240	7971770	7499840	9280300	9121100	9257000	11785570

9 To export RapidFire TOF or QTOF data acquired with *MS mode*:

- Click **File > Export Peak Data TOF**.

A **.RFPKS** text file is created, which contains a list of injections.

Only the masses of interest are reported for each well.

You can further analyze the results with a spreadsheet program. The following is an example file opened in Microsoft Excel.

	A	B	C	D	E	F	G	H	I	J	K	L	M	N	O	P
1	Sequence 1												INTERNAL STANDARD			Injection time (sec)
2	PlateName	Row	Column	260.1572	319.1038	328.1947	340.1383	365.0324	368.1564	378.2103	423.1413	428.0921	494.1438	497.2452	508.2345	
3	DEMO PLATE	2001	2	4277	0	6	323	681	1738	9332	750	1085	2747	12245	10648	9.23
4	DEMO PLATE	1	1	0	139	0	126	8994	23	30920	675	224	297438	4205	2992271	18.46
5	DEMO PLATE	1	3	0	265	0	0	24523	73.5	28588	1319.5	819	490474	8052	3061670	27.27
6	DEMO PLATE	1	5	0	77.5	0	0	13564	100	0	836.5	327	475369.5	6326	2900104	36.17
7	DEMO PLATE	2001	2	0	54	0	0	299	0	0	233	0	4295	2197	4893	44.41
8	DEMO PLATE	2	1	0	479	1819135	0	16647	32469	629	698	424	574327	8337	82630	53.51
9	DEMO PLATE	2	3	0	1079	1888354	0	1880	59220	228	851	482	552930	8459	23723	62.25
10	DEMO PLATE	2	5	0	708	1403209	0	2532	58457	21	486	224	517027	9908	9391	71.09
11	DEMO PLATE	2001	2	0	638	1295.5	0	109.5	0	0	470	41	4003.5	2713.5	3307	79.42
12	DEMO PLATE	3	1	0	960	2248	86	10454	52	277	412	1135663	527903	8916	25129	88.65
13	DEMO PLATE	3	3	0	487	0	747	991	8.5	0	350	195475	503490	7489.5	11622	97.63
14	DEMO PLATE	3	5	0	429	0	880.5	632.5	167	0	309	292988	492337.5	7719	7239.5	106.48
15	DEMO PLATE	2001	2	0	311	0	0	28	5	0	339	0	2913	2301	2620	114.92
16	DEMO PLATE	4	1	0	31407	0	0	838	0	0	430	14744	545323	39815	14282	124.15
17	DEMO PLATE	4	3	0	29637	0	0	300.5	0	4.5	212.5	743	505078.5	340302	6610	133.16
18	DEMO PLATE	4	5	0	34226	0	0	554	0	0	206	315	490114	40181	4813	141.98

Analyzing Data in RapidFire Integrator
To open RapidFire-TOF/QTOF data

www.agilent.com

In This Book

This guide has instructions for installing and using the Agilent RapidFire Integrator.

Windows is a registered trademarks of Microsoft Corporation in the United States and/or other countries.

© Agilent Technologies, Inc. 2013

Printed in USA
Revision A, December 2013



G9530-90001



Agilent Technologies