

# Welcome to the webinar

## Sample purification and GC – MS/MS for dioxin analysis

Waldemar Weber  
Product manager GCMS  
Shimadzu Europa GmbH



# Background

## ➤ **Market trend**


- In the past, Dioxins in food had been analysed by GC HRMS
- Recently, GC MS/MS was also confirmed as official method

## ➤ **Requirement of customer**

- Start analysis immediately without adjusting the analysis conditions
- Create reports showing items required by EU regulations
- Compare the respective quantitative capabilities of GC MS/MS and GC HRMS
- Fast and simple sample preparation

# Manual vs. Automated Sample Prep.

## Manual sample preparation



New fully automatic sample preparation system for Dioxin analysis

# Features of the Miura GO-EHT System

- High performance clean-up

- 4 columns
- Heated purification and elution

Better purification

Faster elution

- User-friendly design

- Easy to apply sample
- Easy connections
- Disposable parts

No washing

No cross-contamination

- Unique flow switching




Perfect separation of PCBs and dioxins  
in one fraction each with excellent quality

# New Miura GO-EHT Features

## GO-EHT series

Touch panel




7.5 inch

Power supply

100-240 V, 50/60 Hz

Method programs



Flexible

Basic method time


80 mins (possible to further optimize)

Functions

Dry-air leak check before start

Data log

P, T, Operation/Alarm










# Dioxin Analyzer

## ➤ 1. Method Files Registered with the Optimal Conditions

- Method files for DXNs, PCBs, and BFRs



- Automatic adjustment of retention times (AART)





# Dioxin Analyzer

- **2. Report creation tool, capable of outputting items required by regulation**
  
- **3. Method file performance confirmed by the analysis of 44 types/201 samples of foods and feeds**

- Animal feed product
- BétaïneHCL 95%
- Bovine milk
- Compound fish food
- Diverse nature
- Fish
- Game liver
- Goat liver
- Milk
- Mussels
- Ovine fat
- Oysters
- Poultry eggs
- Powder
- Salmon
- Scallops
- Shrimp
- Veal fat
- Yellow Pigment
- Bétaïne anhydre
- Bovine fat
- Bovine muscle
- Dairy product
- Eels
- Fresh product
- Goat fat
- Grasses
- Molluscs
- Oilcake
- Ovine liver
- Pork fat
- Poultry muscle
- QC
- Sardine
- Shellfish
- Thréronine
- Vitamine K4
- other


# Steps from the Preparation for Analysis

**Step 1**



**Preparation for analysis**

**Step 2**



**Creation of method files  
AART, Calibration curve**

**Step 3**



**Analysis of samples**

**Step 4**

ID	Compound Name	LOQ	TRQ	TRQ	TRQ
			Lower limit	Method limit	Upper limit
1	1,2,3,7,8-Tetrachlorodibenzo-p-dioxin	0.0020	0.00050	0.00142	0.00284
2	1,2,3,7,8-Pentachlorodibenzo-p-dioxin	0.0001	0.00050	0.00250	0.00500
3	1,2,3,4,7,8-Hexachlorodibenzo-p-dioxin	0.0005	0.00070	0.00070	0.00210
4	1,2,3,6,7,8-Hexachlorodibenzo-p-dioxin	0.0014	0.00080	0.00080	0.00080
5	1,2,3,4,7,8-Hexachlorodibenzo-p-dioxin	0.0005	0.00110	0.00170	0.00210
6	1,2,3,4,6,7,8-Heptachlorodibenzo-p-dioxin	0.0005	0.00205	0.00205	0.00205
7	Octachlorodibenzo-p-dioxin	0.0008	0.00019	0.00019	0.00019
<b>Sum of PCDDs</b>					
8	1,2,3,7,8-Tetrachlorodibenzofuran	0.0014	0.00050	0.00050	0.00050
9	1,2,3,7,8-Pentachlorodibenzofuran	0.0010	0.00054	0.00054	0.00054
10	1,2,3,4,7,8-Hexachlorodibenzofuran	0.0009	0.00040	0.00040	0.00040
11	1,2,3,4,6,7,8-Hexachlorodibenzofuran	0.0013	0.00170	0.00170	0.00170
12	1,2,3,6,7,8-Hexachlorodibenzofuran	0.0012	0.00150	0.00150	0.00150
13	1,2,3,4,6,7,8-Hexachlorodibenzofuran	0.0007	0.00099	0.00099	0.00099
14	1,2,3,7,8,9-Hexachlorodibenzofuran	0.0021	0.00210	0.00210	0.00210
15	1,2,3,4,6,7,8-Hexachlorodibenzofuran	0.0006	0.00041	0.00041	0.00041
16	1,2,3,4,7,8-Hexachlorodibenzofuran	0.0004	0.00009	0.00009	0.00009
17	Octachlorodibenzofuran	0.0001	0.00001	0.00001	0.00001
<b>Sum of PCDFs</b>					
			<b>0.02125</b>	<b>0.02125</b>	<b>0.02125</b>

**Creation of reports**




# Analysis result

I.D.	Compound Name	Calibration Point Concentration						Average RRF	RRF (level 1)	Dev(%) *Criteria <30%
		Level 1 (pg/uL)	Level 2	Level 3	Level 4	Level 5	Level 6			
1	2,3,7,8-Tetrachlorodibenzo-p-dioxin	0.025	0.050	0.100	0.250	0.500	1.000	1.212	1.144	5.60
2	1,2,3,7,8-Pentachlorodibenzo-p-dioxin	0.025	0.050	0.100	0.250	0.500	1.000	1.089	0.990	9.11
3	1,2,3,4,7,8-Hexachlorodibenzo-p-dioxin	0.025	0.050	0.100	0.250	0.500	1.000	1.106	1.157	-4.62
4	1,2,3,6,7,8-Hexachlorodibenzo-p-dioxin	0.025	0.050	0.100	0.250	0.500	1.000	1.043	1.043	-0.06
5	1,2,3,7,8,9-Hexachlorodibenzo-p-dioxin	0.025	0.050	0.100	0.250	0.500	1.000	1.039	0.936	9.95
6	1,2,3,4,6,7,8-Heptachlorodibenzo-p-dioxin	0.025	0.050	0.100	0.250	0.500	1.000	1.033	1.115	-7.86
7	Octachlorodibenzo-p-dioxin	0.050	0.100	0.200	0.500	1.000	2.000	1.180	1.261	-6.84
8	2,3,7,8-Tetrachlorodibenzofuran	0.025	0.050	0.100	0.250	0.500	1.000	1.159	1.213	-4.65
9	1,2,3,7,8-Pentachlorodibenzofuran	0.025	0.050	0.100	0.250	0.500	1.000	1.047	0.974	6.94
10	2,3,4,7,8-Pentachlorodibenzofuran	0.025	0.050	0.100	0.250	0.500	1.000	1.038	0.962	7.35
11	1,2,3,4,7,8-Hexachlorodibenzofuran	0.025	0.050	0.100	0.250	0.500	1.000	1.106	1.358	-22.81
12	1,2,3,6,7,8-Hexachlorodibenzofuran	0.025	0.050	0.100	0.250	0.500	1.000	1.052	1.134	-7.82
13	2,3,4,6,7,8-Hexachlorodibenzofuran	0.025	0.050	0.100	0.250	0.500	1.000	1.000	0.923	7.67
14	1,2,3,7,8,9-Hexachlorodibenzofuran	0.025	0.050	0.100	0.250	0.500	1.000	1.021	1.205	-18.09
15	1,2,3,4,6,7,8-Heptachlorodibenzofuran	0.025	0.050	0.100	0.250	0.500	1.000	1.097	1.157	-5.46
16	1,2,3,4,7,8,9-Heptachlorodibenzofuran	0.025	0.050	0.100	0.250	0.500	1.000	1.056	1.080	-2.27
17	Octachlorodibenzofuran	0.050	0.100	0.200	0.500	1.000	2.000	0.981	0.975	0.66


# Result of the robustness test

- Transition of the sensitivity of STD (0.05 pg/uL)



Lifetime was over than 500 injection

# Method validation



# Dioxins S<sup>3</sup> Smart Solution Systems



Extraction



Clean-up



Analysis

**Thank You for your Attention**