

High Speed vs. High Resolution Analyses in HPLC: A Critical Performance Comparison of Column Options Using Poppe and Kinetic Plots

Jason C. Link, William E. Barber, William Long, Anne Mack
Agilent Technologies


PittCon 2010

Improvements to HPLC


Approach	Specialized Instrument Required	Key Advantages	Major Limitations
Sub-2 μm Particle	Yes	Very high plate counts in short analysis times	Extra-column broadening, frictional heating
Superficially Porous Particle	No	High plate counts at relatively lower pressure	Limited commercial phases
High Temperature	Yes (above 100°C)	High efficiency maintained at high mobile phase velocity	Solutes degradation, Limited number of stable stationary phases
Monolith	No	High column permeability	Batch to batch reproducibility, Limited column dimensions

Poroshell 120 for Small Molecules

1.8 μm totally porous



2.7 μm Poroshell 120



- High efficiency of Poroshell 120 particles similar to sub-2 μm totally porous particles
- Gains due to shorter mass transfer and narrower particle size distribution
- Low back pressure (40-60% of sub-2 μm totally porous particles)


Efficiency improvement with superficially porous particles

(van Deemter equation) $h = A + B / v + C \cdot v$

- **A term – eddy diffusion and flow distribution**
 - Column particle size and column packing quality impact this
 - **Tight particle size distribution improves the A term**
- **B term – longitudinal diffusion**
 - Impact in superficially porous particle not yet determined
- **C term – mass transfer component**
 - **Mass transfer is improved by using shorter diffusion paths, improving the C term**
 - This is improved with a superficially porous particle
 - The C term has more effect on large molecules than on small molecules


Poroshell 120 pore and particles size distribution

BJH Adsorption $dV/d\log(w)$ Pore Volume




Poroshell 120 particles have an average pore size of 120 Å.

Particle Size Distribution Comparison with Totally Porous Particles




Normalized Particle Size Distribution

Poroshell 120 particles size distribution




Particle Size Distribution of Poroshell 120



SEM Picture of Poroshell 120


Poroshell 120 particles also show spherical and smooth surface.

Backpressure comparison of 2.7 μm Poroshell 120 with 1.8 μm totally porous particles



- Superficially porous particles have 40%-50% back pressure of 1.8 μm totally porous particles.


Which Yields Better Performance



- Which one performs better under *optimized* conditions
- Relative performance as a function of separation goal
 - Fast separation (low efficiency with short column)
 - Slow separation (high efficiency with long column)

Theoretical Poppe plot

- Efficiency (N) – Higher N means higher resolving power
- Separation speed (t_0/N) – Lower t_0/N equates to faster separation



3.5µm

$P_{max} = 400 \text{ bar}$, $T = 40 \text{ }^\circ\text{C}$
van Deemter equation

$A = 1.0$, $B = 5.0$, $C = 0.05$

$D_m = 1 \times 10^{-5} \text{ cm}^2/\text{sec}$

$\eta = 6.6 \times 10^{-4} \text{ Pa/sec}$

$\phi = 500$, $\varepsilon_e = 0.38$, $\varepsilon_i = 0.30$

$t_0 = 1 \text{ sec}$

$L = 3.41 \text{ cm}$

$F = 4.36 \text{ mL/min}$

$N = 1,000$

$t_0 = 30 \text{ sec}$

$L = 18.7 \text{ cm}$

$F = 0.796 \text{ mL/min}$

$N = 19248$

$t_0 = 1000 \text{ sec}$

$L = 108 \text{ cm}$

$F = 0.138 \text{ mL/min}$

$N = 141960$


<http://homepages.gac.edu/~dstoll/calculators/optimize.html>



Agilent Technologies

Why use Poppe plots?

Poppe optimum vs H optimum



Poppe plot analysis allows for a fair comparison under optimized conditions

	van Deemter optimum		Poppe optimum	
	1.8 μm	3.5 μm	1.8 μm	3.5 μm
Point	1	2	3	4
L (cm)	5	5	5.96	16.2
n	10	10	11.9	32.2
F (ml/min)	0.44	0.23	0.52	0.73
ΔP (bar)	282	38	400	400
N	13900	7100	16,400	16,700
t_0 (sec)	13.4	26.1	13.4	26.1
N/t_0	1037	272	1224	640
% N		51%		102%
% N/t_0		26%		52%

$P_{\text{max}} = \mathbf{400 \text{ bar}}$, $T = 40 \text{ }^\circ\text{C}$ van Deemter equation
 $A = 1.0$, $B = 5.0$, $C = 0.05$ $D_m = 1 \times 10^{-5} \text{ cm}^2/\text{sec}$
 $\eta = 6.6 \times 10^{-4} \text{ Pa/sec}$ $\phi = 500$, $\varepsilon_e = 0.38$, $\varepsilon_i = 0.30$

Experimental Design

Sub-2 μm fully porous particle

ZORBAX Rapid Resolution
High Definition (RRHD) Eclipse
Plus C18

VS

2.7 μm superficially porous particle

Poroshell 120 EC-C18

1. Perform flow study on each column type with alkylphenones

$$h = A + B/v + C \cdot v$$

2. Transform van Deemter analysis into Poppe Plot
3. Test accuracy of method experimentally with columns connected in series to mimic a variety of column lengths

van Deemter metrics:


- Plate height (h)
- Mobile phase velocity (v)



Practically meaningful metrics:

- Plate count (N)
- Plate count per unit time (N/t₀)

Reduced van Deemter analysis




	Poroshell	RRHD
A	.70	0.81
B	5.6	4.7
C	0.06	0.13
h_{\min}	1.90	2.38

Agilent 1290 Infinity LC System
2.1mm x 50mm columns
Analyte: Hexanophenone
Temp: 40 °C
%ACN adjusted to achieve $k'=6$
for each column

In the reduced van Deemter plot:

- Lower C term on the Poroshell 120 column (better mass transfer)
- Lower h_{\min} on the Poroshell 120 Column

Theoretical Poppe plot




	A	B	C
1.8µm ZORBAX RRHD C18 @ 1000 bar	0.81	4.66	0.13
1.8µm ZORBAX RRHD C18 @ 550 bar	0.81	4.66	0.13

Higher pressures result in more efficient separations in longer analyses

<http://homepages.gac.edu/~dstoll/calculators/optimize.html>

Theoretical Poppe plot

Faster Separation ↓



	A	B	C
1.8 μ m ZORBAX RRHD C18 @ 1000 bar	0.81	4.66	0.13
2.7 μ m ZORBAX C18 @ 550 bar (hypothetical)	0.81	4.66	0.13
2.7 μ m Poroshell 120 @ 550 bar	0.71	5.6	0.06



- Poroshell 120 gives similar performance to ZORBAX RRHD **at 50% backpressure**

<http://homepages.gac.edu/~dstoll/calculators/optimize.html>

Fast Analysis (5cm)


2.7um Poroshell 120
EC-C18
2.1x5 cm
F = 1.94 mL/min, 40°C
56% ACN, 44% Water
549 bar

1.8um ZORBAX RRHD
Eclipse Plus C18
2.1x5cm
F = 1.82 mL/min, 40°C
59% ACN, 41% Water
1003 bar




Long Analysis (55cm)


2.7um Poroshell 120
EC-C18
2.1x55cm (3x15cm, 10cm)
F = 0.202 mL/min, 40°C
56% ACN, 44% Water
547 bar



1.8um ZORBAX RRHD
Eclipse Plus C18
2.1x55cm (3x15cm, 10cm)
F = 0.23 mL/min, 40°C
59% ACN, 41% Water
1001 bar



Theoretical Poppe plot



RRHD (40°C)				Poroshell 120 (40°C)		
L (mm)	N (Calc)	N (Exp)	% Err	N (Calc)	N (Exp)	% Err
50	4805	6119	27.35%	4333	5390	24.39%
300	69695	71643	2.80%	61220	60264	-1.56%
550	114361	102189	10.64%	101312	97170	-4.09%

Conclusions

**Sub-2 μm fully porous particle
ZORBAX RRHD Eclipse Plus
C18**

VS

**2.7 μm superficially porous particle
Poroshell 120 EC-C18**

Comparison under optimized conditions

1. Similar performance between Poroshell 120 at 550 bar and 1.8 μm ZORBAX RRHD at 1000 bar pressures
2. Increasing temperature tends to shift analysis toward slightly higher efficiencies, with a greater impact observed with Poroshell 120
3. More experimentation with other solvents, such as methanol, need to be evaluated to verify performance in a more viscous mobile phase environment

Conclusions

Which column option should be chosen?

Sub 2 μ m totally porous or 2.7 μ m Poroshell 120

For an HPLC (<600 bar), consider Poroshell 120

- Unless the chromatography requires scale-up in particle sizes

For a UHPLC, both choices are possible

- Consider the analysis goals: scalability, selectivity, loadability
- Sub 2 μ m will be an excellent choice for the highest pressure UHPLC systems

Acknowledgements

- Jia Liu¹
- Hongli Wang¹
- Sami Chaana¹
- Charles Lofton¹
- Xiaoli Wang²

¹ Agilent Technologies

² AstraZeneca Pharmaceuticals