

Using the JetClean Self-Cleaning Ion Source to Extend Maintenance Free Operation

Bruce Quimby¹, Michael Szelewski², Elizabeth Almasi¹, Ken Brady¹

¹Agilent Technologies Inc., Wilmington, DE, USA ²M Szelewski LLC, Hockessin, DE, USA

Introduction

Phthalate Analysis

Matrix and column bleed deposits degrade instrument response over time, necessitating routine source cleaning. The process requires removal and abrasive cleaning of the source, reducing productivity. The patented JetClean self-cleaning ion source (JetClean) was developed to maintain consistent MS response in difficult matrices for extended periods of time via controlled introduction of hydrogen into the MS source. This poster describes the application of JetClean to the analysis of phthalate esters (phthalates). Due to the phthalates' adverse health effects (suspected endocrine disruptors) their use is limited by international regulations.

Phthalates can exhibit several undesirable characteristics in GC/MS analysis:

- **Non-linearity:** some analytes have reduced response at the lower cal levels.
- **Peak Tailing:** compounds stick to the source with some ions tailing more than others.
- **Sensitivity:** Higher source temperatures are often used to improve linearity, but some compounds lose sensitivity.
- **Dropping Response:** Raw area response for replicate injections can exhibit significant loss of response with time.

These problems were investigated and addressed with a number of hardware and method changes, resulting in significantly improved results.


Experimental

Instrument Conditions and Test Mixture

GC				MSD				Agilent 5977 Xtr					
Ramp Initial	°C/min	°C	Hold min	Inlet Liner	Ultra Inert Single Taper w/ Wool (PN. 5190-2293)	Solvent Delay	5.0 min	SIM, TID ON					
Ramp 1	20	220	1.0	Column	HP-5MSUI (19091S-433UI)	Acquisition Mode	Atune						
Ramp 2	5	290	4.0	Column Flow	30m x 0.25 mm id x 0.25 um 1.0 mL/min, const flow	Tune	EMV						
Runtime	40 min			Splitter Pressure	2-way splitter 4.2 psig (initial)	Quad Temp	150 °C						
Inlet Mode	Split/Splitless (or MMI)			MSD Restictor	0.68 m x 0.12 mm id 1.2 mL/min const flow	Source Temp	300 °C						
Temp	290 °C			FID Restictor	0.37 m x 0.12 mm id	Transfer Line	280 °C						
Mode	Hot Pulsed Splitless, 15.0 psi			Split ratio	~1:1 MSD:FID	Drawout Lens	9 mm						
Init. Pressure	35 psi until 1.0 min					JetClean Flow	JetClean Mode						
Pulse Pressure	50 mL/min at 1.0 min												
Purge Flow	3 mL/min, switched												
Septum Purge	0.5 uL												
Injection volume													


Hardware Configuration

The system has a CFT splitter between the MSD and FID. The FID signal makes a good reference for peak shape, linearity, and precision. Use of the FID helps distinguish between detector effects and inlet, liner, and column effects.



Phthalate Test Mixture

The phthalate test mixture was diluted in isooctane at various levels to make the calibration standards for linearity and parameter optimization.




Results

Peak Shape Improvement

The first parameters investigated were source temperature and drawout lens diameter. A diameter of 9 mm instead of the typical 3 mm and a temperature of 300°C were found to be optimum. However, without hydrogen cleaning gas added, the peaks still showed tailing. The degree of tailing was not the same on all ions, EICs for multiple ions from the same peak had differing amounts of tailing.

The figure to the right shows that this differential tailing is greatly reduced with the continuous addition of 0.3 mL/min of H₂ to the source, resulting in TIC and FID peak shapes to be very comparable.



Linearity Improvement


The calibration parameters were studied, examining source temperature, H₂ addition, and tune type. In all cases, the 9 mm drawout lens was used. Calibrations were ISTD with 3 runs at each level.

Optimal results were found with the source at 300°C, 0.3 mL/min H₂ added, and using ATUNE.

Cal Range, pg	2.5 - 1250		2.5 - 1250		2.5 - 1250		2.5 - 1250	
	Stat	%RSD	r ²	%RSD	r ²	%RSD	r ²	%RSD
Source Temp °C	300	300	230	230	230	230	230	230
H ₂ , mL/min	0.3	0.3	0.3	0.3	0.3	0.3	no H ₂	no H ₂
	Tune	Atune	Atune	Atune	Etune	Etune	Etune	Etune
Tetradecane	0.9999	3.0	0.9999	5.9	0.9999	4.2	0.9979	22.9
Dimethyl	0.9998	2.4	0.9999	6.2	0.9993	6.9	0.9993	15.0
Diethyl	0.9998	1.9	0.9999	6.9	0.9990	8.8	0.9987	17.5
Disobutyl	0.9998	4.4	0.9999	8.8	0.9987	9.4	0.9990	15.9
Dibutyl	0.9999	10.9	0.9998	14.7	0.9983	14.0	0.9982	21.2
Bis(2-methoxyethyl)	0.9998	21.1	0.9990	21.9	0.9980	26.3	0.9980	35.8
Bis(4-methylpentyl)	0.9996	4.2	0.9999	9.5	0.9990	12.5	0.9991	20.0
Bis(2-ethoxyethyl)	0.9997	16.7	0.9994	20.8	0.9982	24.6	0.9983	32.0
Dipentyl	0.9998	3.8	0.9998	7.8	0.9978	15.7	0.9969	27.9
Dihexyl	0.9998	4.1	0.9999	6.8	0.9984	13.3	0.9988	20.0
Benzyl butyl	0.9997	5.9	0.9997	13.0	0.9982	19.7	0.9975	32.6
Bis(2-butoxyethyl)	0.9998	13.8	0.9997	20.3	0.9982	23.0	0.9990	30.3
Dicyclohexyl	0.9979	5.4	0.9989	23.6	0.9976	30.7	0.9971	42.3
Bis(2-ethylhexyl)	0.9997	3.8	0.9999	7.9	0.9981	13.7	0.9982	22.7
Di n-octyl	0.9997	2.9	0.9999	7.7	0.9984	14.1	0.9992	19.3
Dinonyl	0.9996	3.1	0.9998	9.8	0.9985	14.4	0.9993	20.0

Area Precision

One significant problem with phthalates is dropping response with replicate injections. This problem may be present even if no matrix is injected. Employing JetClean significantly increased precision, which was further improved by employing a Merlin Microseal instead of the standard spectrum.



Normalized area response of 60 consecutive injections of the phthalate mixture at 125pg level. All compounds show a remarkable replicate area precision. The Merlin Microseal contribution to the precision improvement is currently under investigation.

Summary and Conclusions

The GC/MS analysis of phthalates can be improved in terms of peakshape, linearity, and repeatability by incorporating these changes to the analysis method:

- Run the source temperature at 300°C. Values lower than that may result in tailing and dropping response. Higher values result in problems with some of the more thermally labile phthalates.
- Reduce drawout lens interactions with the phthalates by changing to a larger diameter, like 9 mm.
- Using ATUNE.U instead of ETUNE.U improves linearity and peakshape for phthalates.
- Addition of continuous hydrogen to the source during analysis with JetClean improves peakshape, linearity, and replicate precision.
- Use of a Merlin Microseal further improved response stability.

Cal Range, pg	2.5 - 1250		2.5 - 1250		2.5 - 1250		2.5 - 1250		2.5 - 1250
Stat	r ²	%RSD	r ²	%RSD	r ²	%RSD	r ²	%RSD	r ²
Source Temp C	300	300	230	230	230	230	230	230	300
H2, mL/min	0.3	0.3	0.3	0.3	0.3	0.3	0.3	0.3	100
	Tune	Atune	Atune	Atune	Atune	Etune	Etune	Etune	Atune
Tetradecane	0.9999	3.0	0.9999	5.9	0.9999	4.2	0.9999	4.2	0.9999
Dimethyl	0.9998	2.4	0.9999	6.2	0.9993	6.9	0.9993	6.9	0.9999
Diethyl	0.9998	1.9	0.9999	6.9	0.9990	8.8	0.9990	8.8	0.9999
Disobutyl	0.9998	4.4	0.9999	8.8	0.9987	9.4	0.9987	9.4	0.9998
Dibutyl	0.9999	10.9	0.9998	14.7	0.9983	14.0	0.9983	14.0	0.9997
Bis(2-methoxyethy	0.9998	21.1	0.9990	21.9	0.9980	26.3	0.9980	26.3	0.9994
Bis(4-methylpentyl	0.9996	4.2	0.9999	9.5	0.9990	12.5	0.9990	12.5	0.9998
Bis(2-ethoxyethyl)	0.9997	16.7	0.9994	20.8	0.9982	24.6	0.9982	24.6	0.9996
Dipentyl	0.9998	3.8	0.9998	7.8	0.9978	15.7	0.9978	15.7	0.9996
Diethyl	0.9998	4.1	0.9999	6.8	0.9984	13.3	0.9984	13.3	0.9997
Benzyl butyl	0.9997	5.9	0.9997	13.0	0.9982	19.7	0.9982	19.7	0.9995
Bis(2-butoxyethyl)	0.9998	13.8	0.9997	20.3	0.9982	23.0	0.9982	23.0	0.9995
Dicyclohexyl	0.9979	5.4	0.9989	23.6	0.9976	30.7	0.9976	30.7	0.9986
Bis(2-ethylhexyl)	0.9997	3.8	0.9999	7.9	0.9981	13.7	0.9981	13.7	0.9993
Di n-octyl	0.9997	2.9	0.9999	7.7	0.9984	14.1	0.9984	14.1	0.9995
Dinonyl	0.9996	3.1	0.9998	9.8	0.9985	14.4	0.9985	14.4	0.9993