

Extraction of Oil & Grease from Water Samples Using ISOLUTE® O&G Solid Phase Extraction Columns

This method is applicable to the determination of Oil and Grease as defined by EPA Method 1664 as Hexane Extractable Material (HEM, also referred to as Total Oil and Grease) as well as Silica Gel Treated-Hexane Extractable Material (SGT-HEM, also referred to as Total Petroleum Hydrocarbons or TPH).

Oil and Grease samples contain a range of polar and non-polar compounds. Using ISOLUTE® O&G SPE columns, both polar and non-polar compounds are extracted from water samples.

The compounds can be fractionated into the non-polar total petroleum hydrocarbon (TPH or SGT-HEM) fraction and remaining polar fraction using a dual elution solvent strategy.

If fractionation is not required, a single elution, eluting total Oil & Grease (or HEM) is used.

O&G or TPH levels are determined gravimetrically after gentle evaporation of the elution solvent.

ISOLUTE® O&G solid phase extraction columns are available in two sizes, and contain an integral depth filter which prevents blocking by highly particulate laden or contaminated wastewater samples.

Extraction Procedure

Format

ISOLUTE® O & G ISOLUTE® O&G 1 g/6 mL Column (Depth Filter)
(Part Number 753-0100-CD or

ISOLUTE® O & G ISOLUTE® O&G 3 g/70 mL Column (Depth Filter)
Part Number 753-0300-FD



Acidify the sample to pH of 1.9 to 2.1 with 6 M HCl to neutralize fatty acids. Add 10 mL of methanol to 1 L of sample. Chill to 4 °C for samples with high concentrations of fatty acids.

Solvation

Rinse the extraction column with 10 mL of methanol at 10 mL/min.

Equilibration

Rinse the extraction column with 10 mL of reagent water, acidified to pH between 1.9 and 2.1 with 6 M HCl, at 10 mL/min.

Sample Application

The sample may be loaded at rates not exceeding 10 mL/min. After the method has been optimized, increased loading rates should be tested. Loading rates between 50 and 100 mL/min are realistic.

Interference Elution

Rinse the sample bottle with 20 mL of reagent water acidified to pH ~2 with 6 M HCl. If high concentrations of fatty acids are present, repeat this rinse step. Load onto the extraction column. Add 10 mL acetone to the sample bottle. Shake well, making sure that the sides come into good contact with the acetone. **Dilute the acetone with 40 mL of reagent water acidified to pH ~2 with 6 M HCl.** Swirl. Load onto the extraction column. Repeat this acetone rinse until the bottle washings appear clear. Dry the extraction column for approximately 30 minutes on a vacuum manifold or with gas (nitrogen or carbon dioxide at 4 L/min). When dry, the column will no longer feel cool to the touch.

Analyte Elution

1. For determination of the non-polar fraction (TPH, or SGT-HEM), elute with 2 volumes of 4 mL hexane into a tared collection tube. Include a 2 minute soak step with each elution volume. Follow the instructions for evaporation in step 3.
2. Elute the polar fraction into a separate tared collection tube using 2 volumes of 4 mL hexane/THF (1:1, v/v). Include a 2 minute soak step with each elution volume. If separate determination of polar and non-polar fractions is not required, skip step 1.
3. Concentrate both fractions to near dryness under a gentle stream of nitrogen at a maximum temperature of 35 °C. When the solvent is almost gone, weigh the tubes at 1 minute intervals until weight loss is less than 1 mg. Total Oil and Grease is the combined weight of residue from both elutions (or the THF/hexane elutions alone if step 1 is skipped).

Ordering Information

Part Number	Description	Quantity
753-0100-CD	ISOLUTE® O&G 1 g/6 mL Column (Depth Filter)	30
753-0300-FD	ISOLUTE® O&G 3 g/70 mL Column (Depth Filter)	16

General Comments

1. Due to the heterogeneous nature of oil and grease samples, the choice of column size used should be determined by sample type. The 6 mL configuration (p/n 753-0100-CD) should be used for samples containing SGT-HEM concentration < 50 ppm, and HEM concentration < 50 ppm. For higher concentrations, the ISOLUTE® O&G column is available in a higher capacity, wider column format (p/n 753-0300-FD). This format facilitates the rapid extraction of samples containing higher concentrations of Oil & Grease, or high levels of particulate material. Samples can be processed without chilling using this format. Alternatively, the sample volume may be reduced.
2. Due to the nature of the analytes, the bottle washing steps after sample loading are very important, as analytes do stick to the walls of the sample bottle. For this reason, sample splitting is not recommended for Oil and Grease samples.

EUROPE

Main Office: +46 18 565900
 Toll Free: +800 18 565710
 Fax: +46 18 591922
 Order Tel: +46 18 565710
 Order Fax: +46 18 565705
 order@biotage.com
 Support Tel: +46 18 56 59 11
 Support Fax: + 46 18 56 57 11
 eu-1-pointsupport@biotage.com

NORTH & LATIN AMERICA

Main Office: +1 704 654 4900
 Toll Free: +1 800 446 4752
 Fax: +1 704 654 4917
 Order Tel: +1 704 654 4900
 Order Fax: +1 434 296 8217
 ordermailbox@biotage.com
 Support Tel: +1 800 446 4752
 Outside US: +1 704 654 4900
 us-1-pointsupport@biotage.com

JAPAN

Tel: +81 3 5627 3123
 Fax: +81 3 5627 3121
 jp_order@biotage.com
 jp-1-pointsupport@biotage.com

CHINA

Tel: +86 21 68162810
 Fax: +86 21 68162829
 cn_order@biotage.com
 cn-1-pointsupport@biotage.com

KOREA

Tel: +82 31 706 8500
 Fax: +82 31 706 8510
 korea_info@biotage.com
 kr-1-pointsupport@biotage.com

INDIA

Tel: +91 11 45653772
 india@biotage.com

Distributors in other regions
 are listed on www.biotage.com

Literature Number:IST1005A.V.1

© 2020 Biotage. All rights reserved. No material may be reproduced or published without the written permission of Biotage. Information in this document is subject to change without notice and does not represent any commitment from Biotage. E&OE. A list of all trademarks owned by Biotage AB is available at www.biotage.com/legal. Other product and company names mentioned herein may be trademarks or registered trademarks and/or service marks of their respective owners, and are used only for explanation and to the owners' benefit, without intent to infringe.