

Where are my PFAS coming from?

Ruth Marfil-Vega, PhD

Senior Market Manager - Environmental

In today's presentation



1. Why PFAS?
2. Experimental plan
3. Results
4. Conclusions
5. Q&A

Why?

Plan

Results

Conclusions

Q&A

Why PFAS?



Why PFAS?

A systematic evaluation of the sources of PFAS in the laboratory



Why?

Plan

Results


Conclusions

Q&A

Experimental Plan

Tubing Type	With Degasser		Without Degasser	
	Without Delay Column	With Delay Column	Without Delay Column	With Delay Column
PEEK			✓	✓
FEP	✓	✓		
LLDPE	✓	✓	✓	


FEP: Fluorinated Ethylene Propylene



Experimental Plan

- ❑ **Targets – EPA 533.**
- ❑ **LC and MS conditions – listed in Table.**
 - ✓ Initial LC conditions: 5%, 10%, 20%, 30% B.
- ❑ **Equilibration times/sequence:**
 - ✓ Standard operation.
 - ✓ 30 min delay.
 - ✓ 120 min delay.



Shimadzu LCMS 8060	
Mobile Phase A	20 mM Ammonium Acetate in Water
Mobile Phase B	Methanol
Ionization Mode	ESI Negative
Analytical Column	Shim-pack Velox C18 (50 x 2.1, 1.8 μM)
Delay Column	Shimadzu Nexcol C18 (50 x 3.0, 5 μM)
Column Temperature	40 °C
Injection Volume	2 μL
Sample Temperature	8 °C
Interface Temperature	100 °C
Desolvation Temperature	160 °C



Specific results

- ✓ Blanks (80:20 MeOH:H2O):
 - ✓ 120 min equilibration time.
- ✓ LC conditions:
 - 5 %B, 10 %B, 20 %B, 30 %B.



Tubing Type	With Degasser		Without Degasser	
	Without Delay Column	With Delay Column	Without Delay Column	With Delay Column
PEEK			✓	✓
FEP	✓	✓		
LLDPE	✓	✓	✓	



Specific results

- ✓ Blanks (80:20 MeOH:H2O) with 120 min equilibration time.
- ✓ LC conditions:
5 %B, 10 %B, 30 %B.



Tubing Type	With Degasser		Without Degasser	
	Without Delay Column	With Delay Column	Without Delay Column	With Delay Column
PEEK			✓	✓
FEP	✓	✓		
LLDPE	✓	✓	✓	




Specific results

- ✓ Blanks (80:20 MeOH:H2O)
 - 120 min equilibration time
 - 30 min equilibration time
 - No extended equilibration time
- ✓ LC condition: 5%B

Tubing Type	With Degasser		Without Degasser	
	Without Delay Column	With Delay Column	Without Delay Column	With Delay Column
PEEK			✓	✓
FEP	✓	✓		
LLDPE	✓	✓	✓	



Overall results



Contamination from PFBA comes from other sources different from LCMS or mobile phase.

Why?

Plan




Results

Conclusions

Q&A

Overall results

- Different tubing, no degasser, no delay column; blanks injected after 120 minutes of equilibration.



Contamination in PFPeA MRM in blank injections with three types of tubing tested.

Why?

Plan


Results

Conclusions

Q&A

Overall results


❑ Effect of LC equilibration time.



- ✓ Blanks (80:20 MeOH:H₂O).
- ✓ FEP tubing, degasser, no delay column.

Area of background peaks increases with longer equilibration time before injection.

Overall results



- ✓ Standard (10 ppb in vial) with no equilibration time.
- ✓ FEP tubing, degasser, no delay column.
- ✓ LC conditions: **5 %B, 10 %B, 30 %B.**

Percent of organic reduces PFAS background without sensitivity loss.

- ✓ Blanks (80:20 MeOH:H₂O) with 120 min equilibration time.
- ✓ FEP tubing, degasser, no delay column.
- ✓ LC conditions: **5 %B, 10 %B, 30 %B.**

Summary results

**Compounds meeting identification criteria^(*)
in blanks after 120 min equilibration time, without delay column.**

Tubing	PEEK (no degasser)				FEP (degasser)			LLDPE (no degasser)		
	5 %B	10 %B	20 %B	30 %B	5 %B	10 %B	30 %B	5 %B	20 %B	30 %B
LC initial										
PFHpA					X	X		X		
6-2 FTS	X	X	X		X	X		X	X	
PFOA	X	X			X	X		X	X	
PFNA	X	X	X		X	X	X	X	X	X
PFOS	X	X	X	X	X	X	X	X	X	X
PFDA	X	X	X	X	X	X	X	X	X	X
PFUnA	X				X	X	X	X	X	X
PFDoA	X				X	X		X	X	X


(*) Presence of primary MRM and RT within established window.

Targets included in EPA 533 NOT listed in this table: do not meet both identification criteria.



Future work

- ❑ Systematic evaluation of vials and caps, in combination with different LC configurations.



Comparisons of a polypropylene cap and vial (blue) to PTFE cap with silanized glass vial (orange).

Chromatogram from 30th consecutive injection of blank.



Conclusions



- ✓ All types of tested tubing (PEEK, FEP, and LLDPE) displayed PFAS compounds (targets from EPA 533).
- ✓ Presence or absence of in-line degasser did not significantly contribute to PFAS background in samples.
- ✓ Contamination was more severe with longer equilibration times before injection.
- ✓ Initial LC conditions can reduce PFAS background without significant loss in signal.
- ✓ Installation of delay column eliminates measurable PFAS background from all types of tubing tested (Non fluorinated tubing is not essential!).

Acknowledgments



Elizabeth Barnabe
& Christopher Gilles

Why?

Plan

Results

Conclusions

Q&A



Visit OneLabOneEarth.com



Ruth Marfil-Vega, PhD
rmmarfilvega@shimadzu.com