

# Analysis of Trace Persistent Organic Pollutants (POPs) in Dust Using a Novel Multi-Mode Ionization Source and High Resolution Time-of-Flight Mass Spectrometry

David E. Alonso, Todd Richards and Joseph E. Binkley | LECO Corporation, Saint Joseph, Michigan, USA


## Introduction

- Household dust is complex: Numerous, chemically diverse constituents in a wide concentration range
- Dust can be used to estimate human exposure to contaminants (e.g., Pesticides)
- Analysis of dust has primarily been conducted using targeted methods and different analytical instrumentation

## Objectives

- To use a novel multi-mode source with high resolution time-of-flight mass spectrometry to analyze dust
- To use enhanced chromatography and HRTOFMS to identify trace POPs in the complex samples

## Analytical Platform



Mass Accuracy: 1 PPM  
Acquisition Speed: 200 sps  
Enhanced Chromatographic Resolution: GCxGC  
Resolution: Up to 50,000  
Multimode Source (MMS): EI, PCI and ECNI

Figure 1. Pegasus® HRT® 4D and Multi-Mode Source™ (MMS)

## Dust Sample Preparation

- Dust: 1) NIST SRM 2585, 2) Office Sample, and 3) Household Sample
- Extraction:





Figure 2. Dust General Extraction Procedure

## Instrument Acquisition Parameters

Gas Chromatograph	Agilent 7890B with LECO Dual Stage Quad Jet Modulator
Injection	2µL liquid injection, Splitless, 70°C to 300°C at 500°C/min
Carrier Gas	He @ 1.0 mL/min, Corrected Constant Flow
Primary Column	HP-5MS UI, 30 m x 0.25 mm i.d. x 0.25 µm
Secondary Column	BPX-50, 0.60 m x 0.10 mm x 0.10 µm
Temperature Program	70°C (1 min) to 140°C then ramp 8°C/min to 300 °C, hold for 4 min
Modulation Period	Secondary oven maintained +40 °C relative to primary oven
Transfer Line	800 seconds
Mass Spectrometer	LECO Pegasus® HRT® 4D
Ion Source Temperature	EI, 250°C; PCI, ECNI, 165°C
Acquisition Mode	High Resolution, R ≥ 25,000 for m/z 219, Mass Accuracy ≤ 1 ppm
Ionization	EI, PCI/ECNI (Reagent Gas = CH <sub>4</sub> )
Mass Range (m/z)	EI 50-1000; PCI 60-1000; ECNI 30-1000
Acquisition Rate	125 sps

Table 1. Pegasus HRT® 4D Instrument Acquisition Parameters

## NIST SRM: EI, ECNI and PCI-HRT+ 4D Data



- Hard ionization (Fragmentation = More Information)
- Non-selective (Universal)
- Good sensitivity

- Soft ionization (Less Fragmentation)
- Selective (Electron Rich Compounds)

Figure 3. EI, ECNI, and PCI NIST SRM 2585 Data Collected using the Multi-Mode Source

## Advantages of GCxGC-HRT+ 4D Data




Figure 4. Contour Plot Displaying Some of the Major Compounds in NIST SRM 2585

## NIST SRM 2585: Representative Compounds

Name	Formula	R.T. (s)	Similarity	Mass Accuracy (ppm)
Benzyl alcohol	C <sub>8</sub> H <sub>10</sub> O	496.2364	880	
Benzene	C <sub>6</sub> H <sub>6</sub>	126.1212	620	
Benzoic acid	C <sub>7</sub> H <sub>6</sub> O <sub>2</sub>	122.0200	940	
Benzofuran	C <sub>8</sub> H <sub>6</sub> O	126.0200	940	
Benzonitrile	C <sub>7</sub> H <sub>5</sub> N	104.0400	560	
Benzophenone	C <sub>14</sub> H <sub>10</sub> O	206.0700	510	
Benzothiazole	C <sub>7</sub> H <sub>5</sub> N <sub>2</sub> S	155.0400	510	
Benzotriazole	C <sub>7</sub> H <sub>5</sub> N <sub>3</sub>	155.0400	510	
Benzyl chloride	C <sub>7</sub> H <sub>7</sub> Cl	126.5212	620	
Benzyl bromide	C <sub>7</sub> H <sub>7</sub> Br	206.0700	510	
Benzyl iodide	C <sub>7</sub> H <sub>7</sub> I	270.0300	510	
Benzyl acetate	C <sub>9</sub> H <sub>10</sub> O <sub>2</sub>	134.0400	560	
Benzyl benzoate	C <sub>14</sub> H <sub>14</sub> O <sub>2</sub>	226.0900	510	
Benzyl propyl ether	C <sub>10</sub> H <sub>12</sub> O	150.0600	560	
Benzyl phenyl ether	C <sub>12</sub> H <sub>14</sub> O	178.0800	560	
Benzyl phenyl sulfide	C <sub>12</sub> H <sub>14</sub> S	178.0800	560	
Benzyl phenyl sulfone	C <sub>12</sub> H <sub>12</sub> S <sub>2</sub> O	226.0900	510	
Benzyl phenyl selenide	C <sub>12</sub> H <sub>14</sub> Se	226.0900	510	
Benzyl phenyl telluride	C <sub>12</sub> H <sub>14</sub> Te	270.0300	510	
Benzyl phenyl tin	C <sub>12</sub> H <sub>14</sub> Sn	270.0300	510	
Benzyl phenyl lead	C <sub>12</sub> H <sub>14</sub> Pb	270.0300	510	
Benzyl phenyl borane	C <sub>12</sub> H <sub>14</sub> B	206.0700	510	
Benzyl phenyl fluoride	C <sub>12</sub> H <sub>13</sub> F	206.0700	510	
Benzyl phenyl chloride	C <sub>12</sub> H <sub>13</sub> Cl	206.0700	510	
Benzyl phenyl bromide	C <sub>12</sub> H <sub>13</sub> Br	270.0300	510	
Benzyl phenyl iodide	C <sub>12</sub> H <sub>13</sub> I	334.0300	510	
Benzyl phenyl selenide	C <sub>12</sub> H <sub>13</sub> Se	334.0300	510	
Benzyl phenyl telluride	C <sub>12</sub> H <sub>13</sub> Te	378.0300	510	
Benzyl phenyl tin	C <sub>12</sub> H <sub>13</sub> Sn	378.0300	510	
Benzyl phenyl lead	C <sub>12</sub> H <sub>13</sub> Pb	378.0300	510	
Benzyl phenyl borane	C <sub>12</sub> H <sub>13</sub> B	206.0700	510	
Benzyl phenyl fluoride	C <sub>12</sub> H <sub>12</sub> F	206.0700	510	
Benzyl phenyl chloride	C <sub>12</sub> H <sub>12</sub> Cl	206.0700	510	
Benzyl phenyl bromide	C <sub>12</sub> H <sub>12</sub> Br	270.0300	510	
Benzyl phenyl iodide	C <sub>12</sub> H <sub>12</sub> I	334.0300	510	
Benzyl phenyl selenide	C <sub>12</sub> H <sub>12</sub> Se	334.0300	510	
Benzyl phenyl telluride	C <sub>12</sub> H <sub>12</sub> Te	378.0300	510	
Benzyl phenyl tin	C <sub>12</sub> H <sub>12</sub> Sn	378.0300	510	
Benzyl phenyl lead	C <sub>12</sub> H <sub>12</sub> Pb	378.0300	510	
Benzyl phenyl borane	C <sub>12</sub> H <sub>11</sub> B	206.0700	510	
Benzyl phenyl fluoride	C <sub>12</sub> H <sub>11</sub> F	206.0700	510	
Benzyl phenyl chloride	C <sub>12</sub> H <sub>11</sub> Cl	206.0700	510	
Benzyl phenyl bromide	C <sub>12</sub> H <sub>11</sub> Br	270.0300	510	
Benzyl phenyl iodide	C <sub>12</sub> H <sub>11</sub> I	334.0300	510	
Benzyl phenyl selenide	C <sub>12</sub> H <sub>11</sub> Se	334.0300	510	
Benzyl phenyl telluride	C <sub>12</sub> H <sub>11</sub> Te	378.0300	510	
Benzyl phenyl tin	C <sub>12</sub> H <sub>11</sub> Sn	378.0300	510	
Benzyl phenyl lead	C <sub>12</sub> H <sub>11</sub> Pb	378.0300	510	
Benzyl phenyl borane	C <sub>12</sub> H <sub>10</sub> B	206.0700	510	
Benzyl phenyl fluoride	C <sub>12</sub> H <sub>10</sub> F	206.0700	510	
Benzyl phenyl chloride	C <sub>12</sub> H <sub>10</sub> Cl	206.0700	510	
Benzyl phenyl bromide	C <sub>12</sub> H <sub>10</sub> Br	270.0300	510	
Benzyl phenyl iodide	C <sub>12</sub> H <sub>10</sub> I	334.0300	510	
Benzyl phenyl selenide	C <sub>12</sub> H <sub>10</sub> Se	334.0300	510	
Benzyl phenyl telluride	C <sub>12</sub> H <sub>10</sub> Te	378.0300	510	
Benzyl phenyl tin	C <sub>12</sub> H <sub>10</sub> Sn	378.0300	510	
Benzyl phenyl lead	C <sub>12</sub> H <sub>10</sub> Pb	378.0300	510	
Benzyl phenyl borane	C <sub>12</sub> H <sub>9</sub> B	206.0700	510	
Benzyl phenyl fluoride	C <sub>12</sub> H <sub>9</sub> F	206.0700	510	
Benzyl phenyl chloride	C <sub>12</sub> H <sub>9</sub> Cl	206.0700	510	
Benzyl phenyl bromide	C <sub>12</sub> H <sub>9</sub> Br	270.0300	510	
Benzyl phenyl iodide	C <sub>12</sub> H <sub>9</sub> I	334.0300	510	
Benzyl phenyl selenide	C <sub>12</sub> H <sub>9</sub> Se	334.0300	510	
Benzyl phenyl telluride	C <sub>12</sub> H <sub>9</sub> Te	378.0300	510	
Benzyl phenyl tin	C <sub>12</sub> H <sub>9</sub> Sn	378.0300	510	
Benzyl phenyl lead	C <sub>12</sub> H <sub>9</sub> Pb	378.0300	510	
Benzyl phenyl borane	C <sub>12</sub> H <sub>8</sub> B	206.0700	510	
Benzyl phenyl fluoride	C <sub>12</sub> H <sub>8</sub> F	206.0700	510	
Benzyl phenyl chloride	C <sub>12</sub> H <sub>8</sub> Cl	206.0700	510	
Benzyl phenyl bromide	C <sub>12</sub> H <sub>8</sub> Br	270.0300	510	
Benzyl phenyl iodide	C <sub>12</sub> H <sub>8</sub> I	334.0300	510	
Benzyl phenyl selenide	C <sub>12</sub> H <sub>8</sub> Se	334.0300	510	
Benzyl phenyl telluride	C <sub>12</sub> H <sub>8</sub> Te	378.0300	510	
Benzyl phenyl tin	C <sub>12</sub> H <sub>8</sub> Sn	378.0300	510	
Benzyl phenyl lead	C <sub>12</sub> H <sub>8</sub> Pb	378.0300	510	
Benzyl phenyl borane	C <sub>12</sub> H <sub>7</sub> B	206.0700	510	
Benzyl phenyl fluoride	C <sub>12</sub> H <sub>7</sub> F	206.0700	510	
Benzyl phenyl chloride	C <sub>12</sub> H <sub>7</sub> Cl	206.0700	510	
Benzyl phenyl bromide	C <sub>12</sub> H <sub>7</sub> Br	270.0300	510	
Benzyl phenyl iodide	C <sub>12</sub> H <sub>7</sub> I	334.0300	510	
Benzyl phenyl selenide	C <sub>12</sub> H <sub>7</sub> Se	334.0300	510	
Benzyl phenyl telluride	C <sub>12</sub> H <sub>7</sub> Te	378.0300	510	
Benzyl phenyl tin	C <sub>12</sub> H <sub>7</sub> Sn	378.0300	510	
Benzyl phenyl lead	C <sub>12</sub> H <sub>7</sub> Pb	378.0300	510	
Benzyl phenyl borane	C <sub>12</sub> H <sub>6</sub> B	206.0700	510	
Benzyl phenyl fluoride	C <sub>12</sub> H <sub>6</sub> F	206.0700	510	
Benzyl phenyl chloride	C <sub>12</sub> H <sub>6</sub> Cl	206.0700	510	
Benzyl phenyl bromide	C <sub>12</sub> H <sub>6</sub> Br	270.0300	510	
Benzyl phenyl iodide	C <sub>12</sub> H <sub>6</sub> I	334.0300	510	
Benzyl phenyl selenide	C <sub>12</sub> H <sub>6</sub> Se	334.0300	510	
Benzyl phenyl telluride	C <sub>12</sub> H <sub>6</sub> Te	378.0300	510	
Benzyl phenyl tin	C <sub>12</sub> H <sub>6</sub> Sn	378.0300	510	
Benzyl phenyl lead	C <sub>12</sub> H <sub>6</sub> Pb	378.0300	510	
Benzyl phenyl borane	C <sub>12</sub> H <sub>5</sub> B	206.0700	510	
Benzyl phenyl fluoride	C <sub>12</sub> H <sub>5</sub> F	206.0700	510	
Benzyl phenyl chloride	C <sub>12</sub> H <sub>5</sub> Cl	206.0700	510	
Benzyl phenyl bromide	C <sub>12</sub> H <sub>5</sub> Br	270.0300	510	
Benzyl phenyl iodide	C <sub>12</sub> H <sub>5</sub> I	334.0300	510	
Benzyl phenyl selenide	C <sub>12</sub> H <sub>5</sub> Se	334.0300	510	
Benzyl phenyl telluride	C <sub>12</sub> H <sub>5</sub> Te	378.0300	510	
Benzyl phenyl tin	C <sub>12</sub> H <sub>5</sub> Sn	378.0300	510	
Benzyl phenyl lead	C <sub>12</sub> H <sub>5</sub> Pb	378.0300	510	
Benzyl phenyl borane	C <sub>12</sub> H <sub>4</sub> B	206.0700	510	
Benzyl phenyl fluoride	C <sub>12</sub> H <sub>4</sub> F	206.0700	510	
Benzyl phenyl chloride	C <sub>12</sub> H <sub>4</sub> Cl	206.0700	510	
Benzyl phenyl bromide	C <sub>12</sub> H <sub>4</sub> Br	270.0300	510	
Benzyl phenyl iodide	C <sub>12</sub> H <sub>4</sub> I	334.0300	510	
Benzyl phenyl selenide	C <sub>12</sub> H <sub>4</sub> Se	334.0300	510	
Benzyl phenyl telluride	C <sub>12</sub> H <sub>4</sub> Te	378.0300	510	
Benzyl phenyl tin	C <sub>12</sub> H <sub>4</sub> Sn	378.0300	510	
Benzyl phenyl lead	C <sub>12</sub> H <sub>4</sub> Pb	378.0300	510	

Ave. 880/1000

Table 2. Representative Compounds in NIST SRM 2585

## NIST SRM 2585: Halogenated Organic Compounds




Figure 5. Selected Halogenated Compounds in SRM 2585

Name	Formula	R.T. (s)	Similarity	Mass Accuracy (ppm)
Benzyl chloride	C <sub>7</sub> H <sub>7</sub> Cl	488.1952	873	0.87
Benzyl bromide	C <sub>7</sub> H <sub>7</sub> Br	560.2263	866	0.54
2-Chlorophenyl isocyanate	C <sub>7</sub> H <sub>5</sub> ClNO	600.2264	937	0.31
3-Chlorophenol	C <sub>6</sub> H <sub>4</sub> ClO	664.2784	896	1.32
4-Chloroaniline	C <sub>6</sub> H <sub>4</sub> ClN	696.3160	933	1.52
p-Ethylbenzyl chloride	C <sub>9</sub> H <sub>9</sub> Cl	704.2624	884	0.30
4-Chloroaniline	C <sub>6</sub> H <sub>4</sub> ClN	704.3184	816	1.13
p-Ethylbenzyl chloride	C <sub>9</sub> H <sub>9</sub> Cl	720.2704	931	0.21
α-Bromostyrene	C <sub>8</sub> H <sub>7</sub> Br	912.3228	861	-0.23
3,4-Dichlorophenyl isocyanate	C <sub>7</sub> H <sub>3</sub> Cl <sub>2</sub> NO	944.3358	820	1.93
2,5-Dichloroaniline	C <sub>6</sub> H <sub>3</sub> Cl <sub>2</sub> N	1352.7124	933	-0.25
1-Chlorododecane	C <sub>12</sub> H <sub>25</sub> Cl	1504.3656	949	-0.80
Dodecane, 2-bromo-	C <sub>12</sub> H <sub>25</sub> Br	2206.3788	873	N/A
2,2',4,4'-Tetrachlorodiphenyl ether	C <sub>12</sub> H <sub>4</sub> Cl <sub>4</sub> O	2488.2976	781	N/A
TCEP	C <sub>10</sub> H <sub>10</sub> Cl <sub>3</sub> P	2632.3392	893	N/A
TCPP	C <sub>10</sub> H <sub>10</sub> Cl <sub>3</sub> P	2704.3008	918	N/A
TCDF	C <sub>10</sub> H <sub>8</sub> Cl <sub>4</sub> P	2776.3024	816	N/A
Chlorophene	C <sub>6</sub> H <sub>4</sub> Cl <sub>2</sub> O	2904.4120	824	1.22
TDCPP	C <sub>10</sub> H <sub>8</sub> Cl <sub>3</sub> P	3912.4904	902	N/A
2,2',4,4'-Tetrabromodiphenyl ether	C <sub>12</sub> H <sub>4</sub> Br <sub>4</sub> O	4384.6528	831	0.34
2,2',4,4'-Tetrabromodiphenyl ether	C <sub>12</sub> H <sub>4</sub> Br <sub>4</sub> O	4556.7296	892	-0.45

Table 3. SRM 2585 Halogenated Compounds

## NIST SRM 2585: EI, Pesticides




Figure 7. Pesticides in SRM 2585

