

STREAMLINED METHOD DEVELOPMENT FOR ACTIVE PHARMACEUTICAL INGREDIENTS IN COLD AND FLU MEDICATION USING A SYSTEMATIC PROTOCOL

Authors; Margaret Maziarz(1), Paul Rainville(1) Yukari Haramaki(2)
Affiliations (1)Waters Corporation, (2)Waters Asia Pacific

INTRODUCTION

In this work, we present the development of a MS-compatible UPLC method for analysis of eight active pharmaceutical ingredients (APIs) found in common over-the-counter (OTC) cold and flu medication¹.

Active Pharmaceutical Ingredient (API)	Formula	Monoisotopic mass (m/z)
Acetaminophen	C ₉ H ₉ NO	151.06
Dextromethorphan HBr	C ₁₈ H ₂₅ BrNO	351.11
Phenylephrine HCl	C ₉ H ₁₀ ClNO ₂	203.07
Chlorpheniramine maleate	C ₂₀ H ₂₃ ClN ₂ O ₄	390.13
Ibuprofen	C ₁₃ H ₁₈ O ₂	206.13
Pseudoephedrine HCl	C ₁₀ H ₁₆ ClNO	201.09
Guaifenesin	C ₁₀ H ₁₄ O ₄	198.09
Doxylamine succinate	C ₂₁ H ₂₈ N ₂ O ₅	388.20

Table 1. Compounds for method development.

METHODS

Sample Preparation

Separate stock solutions were prepared in methanol at 1.0 mg/mL and subsequently diluted with 90:10 water/ methanol to make a mixture with 100 µg/mL of each analyte.

UPLC Method

LC System	ACQUITY UPLC® H-Class PLUS with PDA & ACQUITY QDa™ Detectors			
Column	ACQUITY UPLC BEH C18, 2.1 x 50, 1.7-µm			
Mobile Phase	A: 10 mM Ammonium acetate in water with 0.2% of ammonium hydroxide B: Methanol with 0.2% ammonium hydroxide			
Flow Rate	0.6 mL/min			
Column Temp.	40 °C			
Sample Temp.	15 °C			
Gradient	Step	Time (minutes)	Solvent A (%)	Solvent B (%)
	1	Initial	95.0	5.0
	2	5.0	10.0	90.0
	3	5.5	10.0	90.0
	4	5.6	95.0	5.0
5	7.5	95.0	5.0	


Table 2. UPLC conditions of the final method.

RESULTS AND DISCUSSION

Method development

A systematic protocol² that includes scouting, screening, and optimization steps was employed to develop a UPLC method for analysis of eight APIs (Table 1).

In the rapid scouting, low and high pH experiments were performed to quickly identify condition for best retention of the active ingredients. Mass spectral data from an ACQUITY QDa Detector was used to track the elution order of each analyte over the pH experiments (Figure 1).



SampleName	152.0 (min)	166.1 (min)	168.1 (min)	199.1 (min)	205.2 (min)	271.2 (min)	272.2 (min)	275.1 (min)
1 APIs, Low pH	0.658	0.243	0.204	1.581	3.683	0.576	1.806	1.381
2 APIs, High pH	0.474	1.988	0.765	1.613	1.753	3.148	3.924	3.367

Figure 1. Rapid scouting with low and high pH. Separation with mass-to-charge (m/z) ratio for each peak (A). Empower MS Peak Tracking report displays retention time of each peak (B).

pH	Column	Strong Solvent	Total Peaks	Total Peaks Rs >=2.0	Total Peaks Tailing <=1.5	Min k'	Lowest Rs	RT of First Peak	RT of Last Peak
1 High pH	CSH C18	ACN	8	7	4	0.98	4.817	0.474	3.92
2 Low pH	CSH C18	ACN	8	5	8	-0.15	1.355	0.204	3.68

Figure 2. Empower 3 custom scoring report for rapid scouting. Conditions with best separation were ranked based on the numbers of peaks that met the performance goals. High pH ranked highest.

The high pH condition from the scouting step was screened with an ACQUITY UPLC CSH C18 and ACQUI- TY UPLC BEH C18 columns using methanol and acetonitrile solvents, respectively. The scoring report was used to analyze the chromatographic data and showed that the ACQUITY UPLC BEH C18 with methanol provided best separation with the highest number of peaks for the USP resolution and peak tailing (Figure 3).

Sample	Column	Strong Solvent	pH	Total Peaks	Total Peaks Rs >=2.0	Total Peaks Tailing <=1.5	Lowest Rs	Max Peak Tailing	Min k'	RT of Last Peak
1 APIs mx	BEH C18	MeOH	High pH	8	7	7	2.353	2.3	1.38	4.89
2 APIs mx	CSH C18	ACN	High pH	8	7	5	4.562	3.3	1.36	3.93
3 APIs mx	BEH C18	ACN	High pH	8	7	5	4.675	2.0	1.38	3.71
4 APIs mx	CSH C18	MeOH	High pH	8	6	7	1.212	1.5	2.27	4.84

Figure 3. Screening with columns and solvents. Empower 3 scoring report shows that the ACQUITY UPLC BEH C18 column and methanol provided best separation.

The ACQUITY UPLC BEH C18 with methanol method was optimized by studying gradient slope, column temperature, pH and wavelength. In addition, the use of MS-compatible buffers was investigated to further improve separation and peak symmetry for the analytes. Addition of ammonium acetate to the mobile phase with 0.2% of ammonium hydroxide improved chromatographic separation and reduced peak tailing (Figure 4).




Figure 4. Mobile phase optimization. 0.1% ammonium hydroxide in water and methanol (A). 0.2% ammonium hydroxide in 10 mM ammonium acetate and in methanol (B). UV at 215 nm.

System Suitability of the final method

Performance was measured by evaluating repeatability of five replicate injections of the sample according to the specifications defined in the USP General Chapter <621> Chromatography³. System suitability results were well below the USP acceptance criteria (Figure 5).

Name	Peak Label	# of Inj.	Ave RT	Ave k'	%RSD RT	%RSD Peak Areas	Ave USP Resolution	Ave USP Tailing
1 Acetaminophen	ACE	5	0.706	2.2	0.21	0.35		1.1
2 Phenylephrine	PHE	5	1.114	4.1	0.11	0.10	11.9	1.2
3 Guaifenesin	GUA	5	2.241	9.2	0.04	0.33	31.7	1.1
4 Pseudoephedrine	PSE	5	2.543	10.6	0.02	0.35	7.4	1.2
5 Ibuprofen	IBU	5	3.731	15.9	0.06	0.33	21.3	1.4
6 Doxylamine	DOX	5	4.100	17.6	0.01	0.29	6.8	1.2
7 Chlorpheniramine	CHL	5	4.332	18.7	0.01	0.23	5.7	1.1
8 Dextromethorphan	DEX	5	4.804	20.8	0.01	0.33	11.0	1.2

Analysis of cold and flu medication

Over-the-counter (OTC) cold and flu samples were prepared in 90:10 water/methanol diluent to the working concentration (Table 3) and filtered prior analysis

Medication	API concentration in sample solution (µg/mL)
Mucinex® Cold, Flu and Sore Throat Maximum Strength Syrup	Acetaminophen: 325 Dextromethorphan HBr: 10 Guaifenesin: 200 Phenylephrine HCl: 5
Vicks™ NyQuil Severe Caplets	Acetaminophen: 325 Dextromethorphan HBr: 10 Phenylephrine HCl: 5 Doxylamine succinate: 6.25
CVS Health Sinus PE + Allergy tablets	Chlorpheniramine maleate: 80 Phenylephrine HCl: 200 Acetaminophen: 325
Tylenol® Cold & Flu Severe Caplets	Dextromethorphan HBr: 10 Guaifenesin: 200 Phenylephrine HCl: 5

Table 3. OTC cold and flu medication for analysis by UPLC

Spectral purity or homogeneity of each active ingredient was checked using peak purity tools in the Empower Software (Figure 6). The UV peak purity plot showed that the phenylephrine purity angle was below the threshold angle, confirming spectral homogeneity (Figure 6B). The Empower 3 Mass Analysis Window showed one mass (m/z) across the entire peak, specific for phenylephrine (Figure 6C). Using both the UV and MS spectral data enabled spectral homogeneity confirmation to ensure that analytes are not subject to interference with any excipients of the sample formulations.




Figure 6. Peak purity determination of phenylephrine API. Mucinex syrup sample, UV at 215 nm (A). UV peak purity plot (B). Empower 3 Mass Analysis window with peak purity spectrum (C).




Figure 7. Analysis of cold and flu medication samples confirms that each API peak is spectrally homogenous

CONCLUSION

Quick and efficient development of reproducible and robust methods

- Reduce development time with a systematic approach
- Automate columns and solvents screening using the UPLC H-Class PLUS System
- Easily identify and track peaks by using the ACQUITY QDa Mass Detector
- Perform quick data analysis and method selection with Empower scoring reports

References

1. Maziarz M, Rainville P. Robust and Rapid Method for Analysis of Active Pharmaceutical Ingredients in Multi-Component Cold and Flu Medication, Waters Application Note. 720006523EN. 2019.
2. Maziarz M, McCarthy SM, Wrona M. Improving Effectiveness in Method Development by Using a Systematic Screening Protocol for a USP Method Waters Application Note 720005026EN. 2014
3. USP General Chapter, <621>, USP45-NF36, Chromatography, The United States Pharmacopeia Convention, Official August 2017