

**RIVA 2016**  
**Poster number**

Title: **Using the JetClean  
Self-Cleaning Ion Source  
to Extend Maintenance  
Free Operation**

Authors: **Bruce Quimby<sup>1</sup>, Michael  
Szelewski<sup>2</sup>, Elizabeth Almasi<sup>1</sup>, Ken  
Brady<sup>1</sup>**

**<sup>1</sup>Agilent Technologies Inc.,  
Wilmington, DE, USA <sup>2</sup>M Szelewski  
LLC , Hockessin, DE, USA**



# Introduction

## Phthalate Analysis

Matrix and column bleed deposits over time degrade instrument response, necessitating routine source cleaning. The process requires removal and abrasive cleaning of the source, leading to lost productivity. The patented JetClean self-cleaning ion source (JetClean) was developed to maintain consistent MS response in difficult matrices for extended periods of time through the introduction of carefully controlled hydrogen into the MS source. This poster describes the application of JetClean to the analysis of phthalate esters (phthalates). Due to the phthalates' adverse health effects (plausible endocrine disruptors) their use is limited by international regulations.

Phthalates can exhibit several undesirable characteristics in GC/MS analysis:

- **Non-linearity:** some analytes have reduced response at the lower cal levels.
- **Peak Tailing:** compounds stick to the source with some ions tailing more than others.
- **Sensitivity:** Higher source temperatures are often used to improve linearity, but some compounds lose sensitivity.
- **Dropping Response:** Raw area response for replicate injections can exhibit significant loss of response with time.

These problems were investigated and addressed with a number of hardware and method changes, resulting in significantly improved results.


# Experimental

## Instrument Conditions and Test Mixture

| GC               |   |     |                | MSD  |                           |
|------------------|---|-----|----------------|--|---------------------------|
| Agilent 7890B    |   |     |                | Agilent 5977 Xtr                             |                           |
| Ramp Initial     | °C/min  | °C  | Hold min       | Solvent Delay                                | 5.0 min                   |
| Ramp 1           | 20  | 60  | 1.0            | Acquisition Mode                             | SIM, TID ON               |
| Ramp 2           | 5   | 220 | 1.0            | Tune   | Atune                     |
| Runtime          | 40 min  | 290 | 4.0            | EMV  | Gain 1                    |
| Inlet Temp       | Split/Splitless (or MMI)                              |     |                | Quad Temp                                    | 150 °C                    |
| Mode             | 290 °C  |     |                | Source Temp                                  | 300 °C                    |
| Init. Pressure   | Hot Pulsed Splitless, 15.0 psi                        |     |                | Transfer Line                                | 280 °C                    |
| Pulse Pressure   | 35 psi until 1.0 min                                  |     |                | Drawout Lens                                 | 9 mm                      |
| Purge Flow       | 50 mL/min at 1.0 min                                  |     |                | JetClean Flow                                | 0.3 mL/min H <sub>2</sub> |
| Septum Purge     | 3 mL/min, switched                                    |     |                | JetClean Mode                                | continuous                |
| Injection volume | 0.5 µL  |     |                |  |                           |
| Inlet Liner      | Ultra Inert Single Taper w/ Wool (PN. 5190-2293)      |     | Column Flow    | 1.0 mL/min, const flow                       |                           |
| Column           | HP-5MSUI (19091S-433UI)<br>30m x 0.25 mm id x 0.25 µm |     | Splitter       | 2-way splitter                               |                           |
| Column Flow      | 1.0 mL/min, const flow                                |     | Pressure       | 4.2 psig (initial)                           |                           |
|                  |   |     | MSD Restrictor | 0.68 m x 0.12 mm id<br>1.2 mL/min const flow |                           |
|                  |   |     | FID Restrictor | 0.37 m x 0.12 mm id                          |                           |
|                  |   |     | Split ratio    | ~1:1 MSD:FID                                 |                           |


## Hardware Configuration

The system has a CFT splitter between the MSD and FID. The FID signal makes a good reference for peak shape, linearity, and precision. Use of the FID helps distinguish between detector effects and inlet, liner, and column effects.



## Phthalate Test Mixture

The phthalate test mixture was diluted in isooctane at various levels to make the calibration standards for linearity and parameter optimization.




| Peak Number | Compound Name       | Tgt Ion |
|-------------|---------------------|---------|
| 1           | Diisobutyl          | 149     |
| 2           | Dibutyl             | 149     |
| 3           | Bis(2-methoxyethyl) | 59      |
| 4           | Bis(4-methylpentyl) | 149     |
| 5           | Bis(2-ethoxyethyl)  | 149     |
| 6           | Dipentyl            | 149     |
| 7           | Dihexyl             | 149     |
| 8           | Benzyl butyl        | 149     |
| 9           | Bis(2-butoxyethyl)  | 149     |
| 10          | Dicyclohexyl        | 149     |
| 11          | Bis(2-ethylhexyl)   | 149     |
| 12          | Di n-octyl          | 149     |
| 13          | Dinonyl             | 149     |

# Results

## Peak Shape Improvement

The first parameters investigated were source temperature and drawout lens diameter. A diameter of 9 mm instead of the typical 3 mm and a temperature of 300°C were found to be optimum. However, without hydrogen cleaning gas added, the peaks still showed tailing. The degree of tailing was not the same on all ions, EICs for multiple ions from the same peak had differing amounts of tailing.

The figure to the right shows that this differential tailing is greatly reduced with the continuous addition of 0.3 mL/min of H<sub>2</sub> to the source, resulting in TIC and FID peak shapes to be very comparable.



## Linearity Improvement

The calibration parameters were studied, examining source temperature, H<sub>2</sub> addition, and tune type. In all cases, the 9 mm drawout lens was used. Calibrations were ISTD with 3 runs at each level.


Optimal results were found with the source at 300°C, 0.3 mL/min H<sub>2</sub> added, and using ATUNE.

< 0.999  
< 30%  
< 20%  
< 10%

| Stat                    | 2.5 - 1250     |       | 2.5 - 1250     |       | 2.5 - 1250     |       | 2.5 - 1250     |       | 2.5 - 1250        |                   |
|-------------------------|----------------|-------|----------------|-------|----------------|-------|----------------|-------|-------------------|-------------------|
|                         | r <sup>2</sup> | %RSD  | r <sup>2</sup> | %RSD  | r <sup>2</sup> | %RSD  | r <sup>2</sup> | %RSD  | r <sup>2</sup>    | %RSD              |
| Source Temp °C          | 300            | 300   | 230            | 230   | 230            | 230   | 230            | 230   | no H <sub>2</sub> | no H <sub>2</sub> |
| H <sub>2</sub> , mL/min | 0.3            | 0.3   | 0.3            | 0.3   | 0.3            | 0.3   | 0.3            | 0.3   | no H <sub>2</sub> | no H <sub>2</sub> |
| Tune                    | Atune          | Atune | Atune          | Atune | Etune          | Etune | Etune          | Etune | Etune             | Etune             |
| Tetradecane             | 0.9999         | 3.0   | 0.9999         | 5.9   | 0.9999         | 4.2   | 0.9999         | 4.2   | 0.9979            | 22.9              |
| Dimethyl                | 0.9998         | 2.4   | 0.9999         | 6.2   | 0.9993         | 6.9   | 0.9993         | 6.9   | 0.9993            | 15.0              |
| Diethyl                 | 0.9998         | 1.9   | 0.9999         | 6.9   | 0.9990         | 8.8   | 0.9990         | 8.8   | 0.9987            | 17.5              |
| Disobutyl               | 0.9998         | 4.4   | 0.9999         | 8.8   | 0.9987         | 9.4   | 0.9987         | 9.4   | 0.9990            | 15.9              |
| Dibutyl                 | 0.9999         | 10.9  | 0.9998         | 14.7  | 0.9983         | 14.0  | 0.9983         | 14.0  | 0.9982            | 21.2              |
| Bis(2-methoxyethyl)     | 0.9998         | 21.1  | 0.9990         | 21.9  | 0.9980         | 26.3  | 0.9980         | 26.3  | 0.9980            | 35.8              |
| Bis(4-methylpentyl)     | 0.9996         | 4.2   | 0.9999         | 9.5   | 0.9990         | 12.5  | 0.9990         | 12.5  | 0.9991            | 20.0              |
| Bis(2-ethoxyethyl)      | 0.9997         | 16.7  | 0.9994         | 20.8  | 0.9982         | 24.6  | 0.9982         | 24.6  | 0.9983            | 32.0              |
| Dipentyl                | 0.9998         | 3.8   | 0.9998         | 7.8   | 0.9978         | 15.7  | 0.9978         | 15.7  | 0.9969            | 27.9              |
| Dihexyl                 | 0.9998         | 4.1   | 0.9999         | 6.8   | 0.9984         | 13.3  | 0.9984         | 13.3  | 0.9988            | 20.0              |
| Benzyl butyl            | 0.9997         | 5.9   | 0.9997         | 13.0  | 0.9982         | 19.7  | 0.9982         | 19.7  | 0.9975            | 32.6              |
| Bis(2-butoxyethyl)      | 0.9998         | 13.8  | 0.9997         | 20.3  | 0.9982         | 23.0  | 0.9982         | 23.0  | 0.9990            | 30.3              |
| Dicyclohexyl            | 0.9979         | 5.4   | 0.9989         | 23.6  | 0.9976         | 30.7  | 0.9976         | 30.7  | 0.9971            | 42.3              |
| Bis(2-ethylhexyl)       | 0.9997         | 3.8   | 0.9999         | 7.9   | 0.9981         | 13.7  | 0.9981         | 13.7  | 0.9982            | 22.7              |
| Di n-octyl              | 0.9997         | 2.9   | 0.9999         | 7.7   | 0.9984         | 14.1  | 0.9984         | 14.1  | 0.9992            | 19.3              |
| Dinonyl                 | 0.9996         | 3.1   | 0.9998         | 9.8   | 0.9985         | 14.4  | 0.9985         | 14.4  | 0.9993            | 20.0              |

## Area Precision

One significant problem with phthalates is dropping response with replicate injections. This problem may be present even if no matrix is injected. Employing JetClean significantly increased precision, which was further improved by employing a Merlin Microseal instead of the standard septum.



Normalized area response of 60 consecutive injections of the phthalate mixture at 125pg level. All compounds show a remarkable replicate area precision. The Merlin Microseal contribution to the precision improvement is currently under investigation.

## Summary and Conclusions

- The GC/MS analysis of phthalates can be improved in terms of peak shape, linearity, and repeatability by incorporating these changes to the analysis method:
- Run the source temperature at 300°C. Values lower than that may result in tailing and dropping response. Higher values result in problems with some of the more thermally labile phthalates.
- Reduce drawout lens interactions with the phthalates by changing to a larger diameter, like 9 mm.
- Using ATUNE.U instead of ETUNE.U improves linearity and peak shape for phthalates.
- Addition of continuous hydrogen to the source during analysis with JetClean improves peak shape, linearity, and replicate precision.
- Use of a Merlin Microseal further improved response stability.