

Assay Analysis of Tetracycline Hydrochloride (USP-38 method):

SAMPLE PREPARATION:

Diluent: Mix 680ml of 0.1M ammonium oxalate and 270 ml of Dimethylformamide.

Standard solution: 0.5mg/ml of Tetracycline HCl in Diluent.

Resolution solution: 100µg/ml of Tetracycline HCl and 25µg/ml of USP 4-Epianhydrotetracycline HCl in Diluent.

CHROMATOGRAPHIC CONDITIONS:

Instrument: UltiMate 3000 LC

Column: Acclaim 120-C8 (4.6*250mm, 5µm, p/n: 059141, lot no.: 018-01-158)

Mobile phase: Dimethylformamide, 0.1M ammonium oxalate and 0.2M dibasic ammonium phosphate (27:68:5), adjust if necessary with 3N ammonium hydroxide or 3N phosphoric acid to a pH of 7.6 to 7.7

Separation Mode: Isocratic

Column temperature: 25°C

Flow rate: 2.0 mL/min

Injection Volume: 20 µl

Detector wavelength: UV280 nm


Run Time: 25min

System Suitability Results:

Sr. No.	Parameters	USP Criteria	Obtained Results
1	Resolution between 4-Epianhydrotetracycline HCl and Tetracycline	NLT 1.2	3.46
2	RRT of 4-Epianhydrotetracycline HCl with respect to Tetracycline HCl peak	About 0.9	0.82
3	The %RSD for replicate injection of standard solution	NMT 2.0%	0.4%

CHROMATOGRAMS:

System Suitability:



Standard Solution:

