

Impurities test for Caffeine (USP 38 method):

SAMPLE PREPARATION:

System Suitability Solution: 0.02 mg/ml of Theophylline in Mobile phase.

Standard solution: Transfer 5.0mg of Caffeine standard to 25 ml volumetric flask, add 5.0 ml of System suitability solution and 10 ml of Mobile phase, make up volume with mobile phase.

Sample solution: 0.2 mg/ml of Caffeine in Mobile phase.

CHROMATOGRAPHIC CONDITIONS:

Instrument: UltiMate 3000 LC

Column: Acclaim 120 C18 (4.6*150mm, 5um, p/n 059148, lot no.:018-01-171)

Mobile phase: Acetonitrile: Tetrahydrofuran:Buffer (25:20:955) adjust pH to 4.5 with Glacial acetic acid.

Buffer: 0.82gm/lit of anhydrous sodium acetate.

Separation Mode: Isocratic

Column temperature: 25°C

Flow rate: 1.0 mL/min

Injection Volume: 10 µl

Detector wavelength: UV 275nm


Run Time: 15 minutes

System Suitability Results:

Sr. No.	Parameters	USP Criteria	Obtained Results
1	Resolution b/w Theophylline and Caffeine Peak	NLT 6.0	12.3
2	Tailing Factor for Theophylline and Caffeine peaks	NMT 2.0	For Both peaks 1.1

CHROMATOGRAMS:

System Suitability:



Impurity Mix:

