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



TORION T-9


TORION T-9

Rapid Identification of Illicit Drug Substances Using Thermal Desorption Coupled with a Portable Toroidal Trap GC/MS System


- Methamphetamine
- Tapentadol
- Methadone
- Cannabidiol
- Hydrocodone
- Morphine
- Hydromorphone
- Codeine
- Diazepam
- Oxycodone
- Cannabinol
- Oxymorphone
- Fentanyl
- Clonazepam
- Alprazolam
- Buprenorphine

Rapid Identification of Illicit Drug Substances Using Thermal Desorption Coupled with a Portable Toroidal Trap GC/MS System

Introduction

The misuse of drug substances creates the need for advancements in the ability to rapidly determine the type of chemical compound present. Gas chromatography-mass spectrometry (GC/MS) has regularly been used by scientists to detect the presence of pain medication in the laboratory. Fortunately, novel innovation and miniaturization of GC/MS instrumentation are making this technique available for in-field investigators and researchers, which offers a simple, prompt, and onsite identification of a drug substance. This study describes the injection, separation, and identification of 16 drugs compounds in less than 10 minutes using portable gas chromatograph-toroidal ion trap mass spectrometry (PerkinElmer, Torion® T-9 Portable GC/MS) combined with a coiled microextraction (CME) sampling injector to provide an effective tool for onsite analysis of illicit drugs substances.

For research use only. Not for use in diagnostic procedures.

TORION T-9

MTBE and BTEX On-Site Rapid Screening of Contaminated Ground Water by Portable SPME-GC/MS Methodology

APPLICATION NOTE

Gas Chromatography/ Mass Spectrometry

Authors:
Portable GC/MS Team
PerkinElmer, Inc.
Shelton, CT


MTBE and BTEX On-Site Rapid Screening of Contaminated Ground Water by Portable SPME-GC/MS Methodology

Introduction

Methyl-tert-butyl ether (MTBE) and the combination of benzene, toluene, ethylbenzene, and xylenes (BTEX), are routinely monitored in soil and ground water to determine whether petroleum contamination has occurred due to

industrial runoff, or waste water effluent, in order to protect local tributaries and surrounding property^{1,2}. For industrial environmental monitoring programs, the ability to quickly analyze samples taken from nearby ground and water sources is key to environmental protection. The ability to analyze samples at the source, without the need for fixed laboratory analysis, can result in significant time and cost savings, and in the case of a spill, can lead to improved response times.

Sampling using the Custodion® SPME syringe, with separation and analysis using the Torion® T-9 portable GC/MS, allows rapid and on-site sample extraction and analysis, and provides actionable results within minutes. Using the portable Torion T-9 GC/MS, informed decisions can be made in critical situations, and can optimize follow-up activities by providing the user with the appropriate information to determine which, and if, samples should be sent to a laboratory for confirmatory analysis. In addition, on-site sample collection, extraction and analysis, using SPME and the Torion T-9 GC/MS, minimizes target analyte losses, compared to traditional techniques that risk analyte loss during complex storage and transport conditions, therefore ensuring the integrity of critical and valuable samples.

Sampling:	Solid phase microextraction (SPME)
SPME Phase	Polydimethylsiloxane/divinylbenzene/carboxen (PDMS/DVB/CAR)
GC Inj. Temp	270 °C
GC Column	MXT-5, 5 m x 0.1 mm, 0.4 µm df
GC Carrier Gas	Helium
GC Column Temp	50-200 °C at 2 °C/s
Transfer Line	270 °C
Injection Split	30:1 after 4 sec splitless
Mass Analyzer	Toroidal ion trap (TMS)
TMS Mass Range	41-500 Da
Ionization Mode	In-trap electron impact (EI)
Detector	Electron multiplier
Vacuum	Roughing and turbo molecular pumps
MS Resolution	Less than unit mass to 230 amu, nominal unit mass to 500 amu

ANALYTIKA

HYPHENATION

TORION T-9

Děkuji.

WWW.PESYSTEMS.CZ