

# Summary results using G7000A GC-QQQ for OCPs, Pesticides and PCBs in Marine Sediment and Biota



**Paul Zavitsanos**  
**Agilent Technologies**

**PITTCON 2009**  
**Chicago**

# **G7000A GC-QQQ**

## **OCPs and PCBs in Marine Sediment and Biota**



**Chris Sandy**

**GC-MS Applications Chemist**

**Agilent Technologies UK**

**Melissa Churley**

**GC-MS Applications Chemist**

**Agilent Technologies USA**

**Chicago 2009**



**Agilent Technologies**


# EA Evaluation Protocol

## 15 Target Analytes, 2 ISTDs


TS	Start Time	Parent Ion		Ret Time
1	2.0	219	a-HCH	6.19
		219	b-HCH	6.76
		219	Lindane	6.89
		219	d-HCH	7.45
2	7.80	256	PCB 28	8.25
3	8.80	292	PCB 52	9.15
4	9.30	263	Aldrin	9.47
5	9.90	193	Isodrin	10.24
6	10.90	360	PCB 155 (ISTD)	11.28
7	11.40	326	PCB 101	11.51
8	11.90	263	Dieldrin	12.18
		263	Endrin	12.56
9	12.70	326	PCB 118	12.84
10	13.00	360	PCB 153	13.25
		360	PCB 138	13.72
11	14.20	394	PCB 180	14.75
12	14.90	272	Mirex (ISTD)	15.09




# Example EI Spectrum – PCB 101



# Example EI Spectrum – g-HCH



# Example EI Spectrum - Endrin




# Samples

- **NIST SRM 2977 Mussel tissue**
- **NIST SRM 1944 New York / New Jersey Waterway Sediment**
- **The samples (approx 2-5g) are extracted using accelerated solvent extraction (Dionex ASE system) and then purified using size exclusion chromatography followed by alumina clean up. Biota samples are cleaned up using sulphuric acid for PCB analysis .**




## 2ul Splitless Injection Scan Mode – 5975 MSD




Agilent Technologies




## 2ul Splitless Injection Scan Mode – 5975 MSD



## 2ul Splitless Injection Scan Mode – 5975 MSD




## 2ul Splitless Injection SIM Mode – 5975 MSD



# OCPs / PCBs in Marine Sediments

## Agilent RTL Pesticide Method



# GC Method Agilent RTL Pest 2x Translated for AUX 4.0 psig

GC Method Translation - RTPSTX2.MXD

Criterion:  Translate Only  Best Efficiency  Fast Analysis  None **Speed gain: 2.00000**

	Original Method	Translated Method																																				
<b>Column</b>																																						
Length, m	30.00	<input type="checkbox"/> 15.00																																				
Internal Diameter, $\mu\text{m}$	250.0	<input type="checkbox"/> 250.0																																				
<b>Film</b>																																						
Thickness, $\mu\text{m}$	0.250	<input type="radio"/> Unlock																																				
Phase Ratio	250.0	<input type="radio"/> 0.250																																				
		<input checked="" type="radio"/> 250.0																																				
<b>Carrier Gas</b>	Helium	<input type="checkbox"/> Helium																																				
Enter one Setpoint		<input checked="" type="radio"/> Unlock																																				
Head Pressure, psi	19.664	<input type="radio"/> 17.177																																				
Flow Rate, mL/min	2.1000	<input type="radio"/> 2.3704																																				
Outlet Velocity, cm/sec	Very large	73.10																																				
Average Velocity, cm/sec	52.86	<input type="radio"/> 52.86																																				
Hold-up Time, min	0.945907	<input type="radio"/> 0.472954																																				
Outlet Pressure (absolute), psi	0.000	<input type="checkbox"/> 18.696																																				
Ambient Pressure (absolute), psi	14.696	<input type="checkbox"/> 14.696																																				
<b>Oven Temperature</b> 3-ramp Program																																						
	<table border="1"> <thead> <tr> <th>Ramp Rate</th> <th>Final Temp.</th> <th>Final Time</th> </tr> <tr> <th><math>^{\circ}\text{C}/\text{min}</math></th> <th><math>^{\circ}\text{C}</math></th> <th>min</th> </tr> </thead> <tbody> <tr> <td></td> <td>70.00</td> <td>2.000</td> </tr> <tr> <td>Ramp 1</td> <td>25.000</td> <td>150.00</td> </tr> <tr> <td>Ramp 2</td> <td>3.000</td> <td>200.00</td> </tr> <tr> <td>Ramp 3</td> <td>8.000</td> <td>280.00</td> </tr> </tbody> </table>	Ramp Rate	Final Temp.	Final Time	$^{\circ}\text{C}/\text{min}$	$^{\circ}\text{C}$	min		70.00	2.000	Ramp 1	25.000	150.00	Ramp 2	3.000	200.00	Ramp 3	8.000	280.00	<table border="1"> <thead> <tr> <th>Ramp Rate</th> <th>Final Temp.</th> <th>Final Time</th> </tr> <tr> <th><math>^{\circ}\text{C}/\text{min}</math></th> <th><math>^{\circ}\text{C}</math></th> <th>min</th> </tr> </thead> <tbody> <tr> <td></td> <td>70.00</td> <td>1.000</td> </tr> <tr> <td></td> <td>50.000</td> <td>150.00</td> </tr> <tr> <td></td> <td>6.000</td> <td>200.00</td> </tr> <tr> <td></td> <td>16.000</td> <td>280.00</td> </tr> </tbody> </table>	Ramp Rate	Final Temp.	Final Time	$^{\circ}\text{C}/\text{min}$	$^{\circ}\text{C}$	min		70.00	1.000		50.000	150.00		6.000	200.00		16.000	280.00
Ramp Rate	Final Temp.	Final Time																																				
$^{\circ}\text{C}/\text{min}$	$^{\circ}\text{C}$	min																																				
	70.00	2.000																																				
Ramp 1	25.000	150.00																																				
Ramp 2	3.000	200.00																																				
Ramp 3	8.000	280.00																																				
Ramp Rate	Final Temp.	Final Time																																				
$^{\circ}\text{C}/\text{min}$	$^{\circ}\text{C}$	min																																				
	70.00	1.000																																				
	50.000	150.00																																				
	6.000	200.00																																				
	16.000	280.00																																				
<b>Sample Information</b> None																																						


# Target RTs and Transitions

TS	Target	Ret Time	Parent Ion m/z	Quant m/z	Qual m/z	CV
1	a-HCH	6.19	219	147	183	20
	b-HCH	6.76	219	147	183	20
	Lindane	6.89	219	147	183	20
	d-HCH	7.45	219	147	183	20
2	PCB 28	8.25	256	186	151	20
3	PCB 52	9.15	292	220	257	20
4	Aldrin	9.47	263	193	228	30
5	Isodrin	10.24	193	123	157	30
6	PCB 155 (ISTD)	11.28	360	290		20
7	PCB 101	11.51	326	256	291	20
8	Dieldrin	12.18	263	193	228	30
	Endrin	12.56	263	193	228	30
9	PCB 118	12.84	326	256		20
10	PCB 153	13.25	360	290	325	20
	PCB 138	13.72	360	290	325	20
11	PCB 180	14.75	394	324	359	20
12	Mirex (ISTD)	15.09	272	237		20

# Calibration Range

Std Ref	Conc ppb	uL Injected	pg injected
Cal 5	200	2	400
Cal 4	100	2	100
Cal 3	40	2	80
Cal 2	20	2	40
Cal 1	8	2	16
Low Cal A	1.6	2	3.2
Low Cal B	0.8	2	1.6
Low Cal C	0.2	2	0.4
Low Cal D	0.1	2	0.2


# OCs / PCBs TIC MRM





Agilent Technologies




# Lindanes : MRM Quant Transition 219 → 147




# Cal Curves Isodrin, PCB 101





# MRM Quant Transition 263 → 193





# Lindanes : MRM 219 → 147, RMS S/N 200 fg injected





# SIM (5973) vs MRM : a-HCH 3.2 pg injected




# SIM (5973) vs MRM : Isodrin 3.2 pg injected



# SIM (5973) vs MRM : Dieldrin 3.2 pg injected




# PCB 52 : 0.63 ppb in Biota Extract






# PCB 180 : 19.3 ppb in Sediment Extract



# 'Results at a Glance' – PCB 153

Sample			PCB 153...	PCB 153 Results						Qu...	PCB 155 (ISTD) (IST...	
Name	Type	Acq. Date-Time	Exp. Conc.	RT	Resp.	Calc. Conc.	Final Conc.	Accuracy	Ratio	RT	Resp.	
EA Low Cal C	Cal	8/5/2008 3:06 PM	0.4000	13.256	392	0.2905	0.2905	72.6	41.5	11.280	1407548	
EA Low Cal B	Cal	8/5/2008 3:23 PM	0.8000	13.256	1153	0.8234	0.8234	102.9	35.3	11.277	1240816	
EA Low Cal A	Cal	8/5/2008 3:42 PM	1.6000	13.257	3195	1.7957	1.7957	112.2	34.0	11.277	1394530	
EA Cal 01	Cal	8/5/2008 4:01 PM	8.0000	13.257	13772	6.5645	6.5645	82.1	35.1	11.276	1341342	
EA Cal 02	Cal	8/5/2008 4:20 PM	20.0000	13.258	44311	16.9417	16.9417	84.7	35.0	11.275	1440849	
EA Cal 03	Cal	8/5/2008 4:39 PM	40.0000	13.259	93432	36.7860	36.7860	92.0	34.9	11.277	1238753	
EA Cal 04	Cal	8/5/2008 4:58 PM	100.0000	13.259	330146	103.6555	103.6555	103.7	35.9	11.276	1320111	
EA Cal 05	Cal	8/5/2008 5:17 PM	200.0000	13.258	781324	199.0773	199.0773	99.5	36.7	11.276	1468188	
EA Extract Biota Blank	Sample	8/5/2008 6:53 PM		13.256	409	0.0805	0.0805		69.7	11.277	6486277	
EA Biota CRM 2977	Sample	8/5/2008 7:31 PM		13.258	10272	1.5128	1.5128		42.2	11.280	5467023	
EA Biota CRM 2977	Sample	8/5/2008 8:09 PM		13.258	9834	1.4803	1.4803		35.1	11.280	5366849	
EA Biota CRM 2977	Sample	8/5/2008 8:47 PM		13.261	10925	1.5501	1.5501		41.1	11.284	5652983	
EA Extract Sediment Blank	Sample	8/5/2008 9:25 PM		13.261	9597	2.7890	2.7890		35.3	11.283	2516798	
EA Sediment CRM 1944	Sample	8/5/2008 10:03 PM		13.275	132267	25.1141	25.1141		43.7	11.295	2727390	
EA Sediment CRM 1944	Sample	8/5/2008 10:41 PM		13.278	135023	24.4716	24.4716		43.8	11.300	2868966	
EA Sediment CRM 1944	Sample	8/5/2008 11:19 PM		13.279	119553	23.8858	23.8858		28.8	11.301	2612484	



## Quant Reproducibility– ppb

	PCB 28	PCB 52	PCB 101	PCB 118	PCB 153	PCB 138	PCB 180
EA Extract Biota Blank	0.03	< 0.01	<0.01	0.04	ND	ND	ND
EA CRM 2977 Biota Extract	0.72	0.63	0.83	1.47	1.51	1.47	1.02
EA CRM 2977 Biota Extract	0.74	0.66	0.84	1.47	1.48	1.46	1.03
EA CRM 2977 Biota Extract	0.76	0.64	0.85	1.55	1.55	1.48	1.14
EA Extract Sediment Blank	0.94	2.93	1.59	2.49	2.79	1.57	4.36
EA 1944 Sediment Extract	46.69	23.29	18.81	23.03	25.11	25.10	19.36
EA 1944 Sediment Extract	44.87	21.83	18.60	23.65	24.47	25.05	19.29
EA 1944 Sediment Extract	45.28	22.28	18.29	22.25	23.88	24.41	18.87


# Pesticide/PCB results from the G7000A QQQ-GCMS



Melissa Churley, Santa Clara  
Charles Thomson, Santa Clara  
Bernhard Rothweiler, WAD  
Pasi Vuorinen, Finland  
Max Ruemler, WAD

# Pesticide STD

## 5 levels, 10 – 200 ppb



# Pesticide STD




## Results table (ppb)

	QQQ	Q	QQQ	Q	QQQ	Q
	strawberry	strawberry	carrot	carrot	parsley	parsley
propoxur					4.0	
chlorpropham			0.4		3.7	
triflualin					22.9	23.6
Carbofuran	1.1		0.5		1.5	
BHC Beta gama					0.2	
diazinon			122.8	112.4	0.1	
vinclozolin			212.1	187.5		
chlorpyriphos methyl			16.6	21.7		
malathion			146.5	139.2	0.3	
fenpropimorph			0.3			
pendimethalin			0.3		26.7	
endosulfan alfa					0.3	
mepanipyrim	64.4	65.2				
kresoxim methyl			10.4	10.8		
Endosulfan beta					0.2	
oxadixyl					4.7	
endosulfan sulfate			20	21.9		
quinoxifen			65	63.7		
triphenyl phosphate	ISTD	ISTD	ISTD	ISTD	ISTD	ISTD
phosmet			48.4	46.1	5	
cyhalothrin	54.4	49.1			125.8	125.5






# Pesticides in the sample

## Range detected, 0.1 to 200 ppb






**0.1 ppb Diazinon  
in parsley matrix**






**0.4 ppb  
Chlorpropham  
in carrot matrix**

# Pesticides in the sample

## Range detected, 0.1 to 200 ppb





**3.7 ppb Chlorpropham  
in parsley matrix**



**125.8 ppb Cyhalothrin  
in parsley matrix**




# PCB STDs 1 to 1000 ppb



## List:

- PCB18
- PCB28+31
- PCB52
- PCB44
- PCB101
- PCB149
- PCB118
- PCB153
- PCB138
- PCB180
- PCB170
- PCB194
- PCB209

# MRM data from 1 ppb PCB STD



# PCBs within fish matrix

	QQQ	Q	QQQ	Q
	eel	eel	trout	trout
PCB18	0.2		0.3	15.2
PCB28+31	1.1		31.4	
PCB52	6.7	3.8	2.2	0.7
PCB44	4.7	1.3	2.7	0.6
PCB101	18.4	10.7	16.1	7.1
PCB149	54.5	14.3	36.9	10.4
PCB118	34.4	22.6	14.9	9.5
PCB153	148	118.5	52.3	40.8
PCB138	154.8	108.7	54.6	34.1
PCB180	100.9	58.5	46.4	22.5
PCB170	97.1	30.7	43	10.6
PCB194	13.5	3.8	9.1	0.9
PCB209	82.5	56.9	136	88.7



# Results summary

## Pesticide Standard

- ❑ Concentrated on 29 from ~120 Pesticides found
- ❑ Compared STD results between single Q and QQQ
  - Scan spectra and SIM data comparable (Collision gas ON!)
  - Over 5 concentration levels (10 – 200 ppb) comparable good results

## QQQ benefits

- MRM data are more selective and sensitive (below 1 ppb)
- No matrix effects visible at the QQQ



# Results summary

## Pesticides in the sample

- Found pesticides in each of the three samples (strawberry, carrot, and parsley) detected compounds in the range of 0.1 to 200 ppb
- Required detection limit is 10 ppb

## PCBs in the samples and STD

- STD range 1 – 1000 ppb
- Required detection limit is 1 ppb
- All 14 PCBs detected in matrix (fish)



# Results summary

## Backflush

- Not running backflush at the Single Quad already effected the chromatographic results after about 20 injections
- Not visible at the QQQ GC/MS system due to backflush



# Pesticide STD

## Analytical QQQ benefits

- MRM data are more selective and sensitive
  - Found 0.1 to 200 ppb in sample
  - Detected all Pesticide in spiked customer sample
- No matrix effects visible at the QQQ
  - Lower possibility for false positives and negatives

